
Quantu 
PureFlow

Cardioplegia heat exchanger


Spectrum  *Medical*


EN	ENGLISH	Instructions for use	1
IT	ITALIANO	Istruzioni per l'uso	8
FR	FRANÇAIS	Mode d'emploi	15
DE	DEUTSCH	Gebrauchsanleitung	22
ES	ESPAÑOL	Instrucciones de uso	29
NL	NEDERLANDS	Gebruiksaanwijzing	36
SV	SVENSKA	Bruksanvisning	43
DA	DANSK	Brugsanvisning	50
FI	SUOMI	Käyttöohjeet	57
NO	NORSK	Bruksanvisning	64
CS	ČEŠTINA	Návod k použití	71
PL	POLSKI	Instrukcje dotyczące użytkowania	78

EN - Cardioplegia Heat Exchanger Performance	IT - Performance dello scambiatore di calore per cardioplegia
FR - Performances de l'échangeur thermique pour cardioplegie	DE - Leistung des Kardioplegiawärmetauschers
ES - Rendimiento del intercambiador de calor para cardioplegia	NL - Prestaties van de warmtewisselaar voor cardioplegie
SV - Prestanda för värmväxlare för kardioplegilösning	DA - Kardioplegi-varmevekslerens ydeevne
FI - Kardioplegiassa käytettävän lämmönvaihtimen suorituskyky	NO - Ytelse til kardioplegisk varmeveksler
CS - Výkon tepelného výměníku pro kardioplegii	PL - Wydajność wymiennika ciepła do kardioplegii


HX32V-C0 - HX32V-C4


HX22V-C0 - HX22V-C1 - HX22V-C4 - HX22V-C5

HX32V-C0 - HX32V-C4

HX22V-C0 - HX22V-C1 - HX22V-C4 - HX22V-C5

Quantum PureFlow Cardioplegia Heat Exchangers

		HX-C	
REF.	DESCRIPTION	CODE	DESCRIPTION
HX32V-C0	Cardioplegia Heat Exchanger High Flow 1/4	HX32V-C4	Cardioplegia Heat Exchanger High Flow 1/4 - PL
HX22V-C0	Cardioplegia Heat Exchanger Low Flow 3/16	HX22V-C4	Cardioplegia Heat Exchanger Low Flow 3/16 - PL
HX22V-C1	Cardioplegia Heat Exchanger Low Flow 1/4	HX22V-C5	Cardioplegia Heat Exchanger Low Flow 1/4 - PL

Cardioplegia Heat Exchanger

CONTENTS

- A. WARNINGS
- B. CONTRAINDICATIONS AND ADVERSE REACTIONS
- C. DESCRIPTION
- D. TECHNICAL CHARACTERISTICS
- E. INTENDED USE
- F. DEVICE PREPARATION AND SET-UP
- G. DEVICE PRIMING PROCEDURE
- H. DEVICE USE AND REMOVAL
- I. DEVICE REPLACEMENT
- J. MEDICAL DEVICES TO BE USED WITH THE DEVICE
- K. RETURN OF USED DEVICES AND NOTIFICATIONS
- L. LIMITED WARRANTY
- M. SYMBOLS USED

A. WARNINGS

CAUTION

- Carefully read the instructions before use.
- The device must be used in accordance with these instructions for use. QURA cannot be held responsible for damages resulting from incorrect use; the user is fully responsible for improper use or use different from that indicated.
- The device may only be used by adequately trained and qualified personnel.
- The device is sterilized by ethylene oxide and is sold STERILE. Sterility is guaranteed only when the package is intact. Do not use the device if the package is damaged, tampered with, open or wet.
- Visually inspect and carefully check the device before use. Inappropriate transport and/or storage conditions not compliant with the instructions may have damaged the product.
- Fragile, handle with care.
- Keep dry. Store in a dry place between 0°C/32°F and 60°C/140°F.
- Single-use product to be used only once and on one single patient. Reuse (on the same person or on other patients) may result in cross-contamination, infection and sepsis after the device comes into contact with blood/body fluids or any substances infused into the patient during the procedure. Reusing the product could also compromise its characteristics in terms of performance and safety.
- Do not resterilize, subject to further processing or clean for reuse.
- The device is not made with natural rubber latex and does not contain phthalates.
- The device is non-toxic and non-pyrogenic.
- Check the expiry date printed on the label. Do not use the device beyond this date.
- In order to avoid product contamination, the device must be used immediately after opening the sterile packaging and handled aseptically.
- The device must be used only in combination with the medical devices listed in section J "Medical devices to be used with the device" in these instructions for use. Carefully read the user manuals of the medical devices listed in section J before using the device.
- Do not allow alcohol, alcohol-based liquids, ethers, halogenated liquids (such as Halothane), or corrosive solvents (such as acetone) to come into contact with the device. This would cause damage such as to compromise the integrity and functionality of the device.
- The connection lines must be properly connected to prevent tube kinking or occlusion which could lead to a reduction in flow and/or an increase in circuit pressure.
- It is advisable to externally secure the tubes on the herringbone connectors of the device using cable-ties.
- Continuously and carefully monitor perfusion.
- A spare device must always be available during the procedure.
- Follow an anticoagulation protocol and monitor anticoagulation throughout the procedure. The benefit from using an extracorporeal support must however be evaluated against the risk of systemic anticoagulation and must be supported by a medical prescription.
- Frequently monitor both the patient and the device; do not leave the device unattended while it is running.
- After use, dispose of the device in accordance with the applicable regulations in force in the country of use.
- In the event of a malfunction, always keep the defective device and follow the instructions given in section J "Return of used devices and notifications" if it needs to be returned.
- For further information and/or in case of complaint, contact Qura S.r.l. or the distributor of the Qura product.
- For US market - Caution: Federal law restricts this device to sale by or on the order of a physician.

B. CONTRAINDICATIONS AND ADVERSE REACTIONS

The device may not be used for any purpose other than its intended use. The user is responsible for any use different from that intended. The following adverse reactions are possible: allergic reaction, haemorrhage, blood toxicity, embolism, excessive activation of blood components, bleeding, fever, hemolysis, hypercapnia and acidosis, hypotension, hypovolemia, hypoxia, compromised patient temperature control, infection, neurological dysfunction and stroke, organ dysfunction, renal failure, reproductive harm in male neonatal and paediatric patients, respiratory and pulmonary problems, sepsis, sterility, thrombosis, toxemia, carcinogenicity, mutagenicity, sensitization, cytotoxicity, vascular damage, death.

C. DESCRIPTION

The device (Fig. 1) has been designed to heat/cool cardioplegic solution during extracorporeal circulation. The device is equipped with an integrated bubble trap filter (ref. 7 in Fig.1). The device is available in the versions listed in Table 1 – section D. The device is sold sterile in a blister pack in a multiple box of 4 units. All the device surfaces in contact with blood are treated with a phosphorylcholine-based biocompatible coating. Contact with cardioplegic solution for a period of more than 6 hours is not permitted. The device has been designed to be powered by heaters coolers that use a glycol-based heat transfer fluid (HTF). The device must be used in combination with the medical devices listed in section J in these instructions for use.


Fig. 1

References:

- | | |
|---|------------------------------|
| 1. HTF inlet connector | 5. Purge Luer-lock connector |
| 2. HTF outlet connector | 6. Overpressure valve |
| 3. Cardioplegic solution inlet connector | 7. Bubble trap filter |
| 4. Cardioplegic solution outlet connector | 8. Temperature probe |

D. TECHNICAL CHARACTERISTICS

Characteristics	CODE						
	HX32V-C0	HX22V-C0	HX22V-C1	HX32V-C4	HX22V-C4	HX22V-C5	
	High Flow	Low Flow	Low Flow	High Flow	Low Flow	Low Flow	
Blood contact materials	Polycarbonate Polyethylene Polyester AISI 304 stainless steel sheet Epoxy resin Silicone Coating: phosphorylcholine						
Maximum operating pressure - cardioplegic solution side	750 mmHg / 100 kPa / 1 bar / 14.5 psi						
Overpressure valve opening	500-750 mmHg / 66-100 kPa / 0.6-1 bar / 9.7-14.5 psi						
Maximum operating pressure - HTF side	1500 mmHg / 200 kPa / 2 bar / 29 psi						
Maximum HTF compartment supply temperature	42°C						
Static priming volume (ml)	32	22	22	32	22	22	
Exchange surface area (m ²)	0.055	0.028	0.028	0.055	0.028	0.028	
Maximum cardioplegic solution flow rate (l/min)	1	0.5	0.5	1	0.5	0.5	
Filter cutoff (µm)	120						
Temperature probe	YSI 400						
Connections - cardioplegic solution side							
	IN	1/4" (6.35 mm)	3/16" (4.76 mm)	1/4" (6.35 mm)	1/4" (6.35 mm)	3/16" (4.76 mm)	1/4" (6.35 mm)
	OUT	1/4" (6.35 mm)	3/16" (4.76 mm)	1/4" (6.35 mm)	Pos Lock	Pos Lock	Pos Lock
Connections - HTF side (inlet/outlet)		CPC HFC profile	CPC APC profile	CPC APC profile	CPC HFC profile	CPC APC profile	CPC APC profile

Table 1

E. INTENDED USE

The Cardioplegia Heat Exchanger has been designed to heat/cool cardioplegic solution and physiological fluids, remove air and monitor the temperature during surgical procedures that require extracorporeal circulation and cardiopulmonary bypass for periods of maximum 6 hours. Any blood to be treated must contain anticoagulant. Contact with cardioplegic solution for longer than this period is not permitted.

F. DEVICE PREPARATION AND SET-UP

CAUTION

- A spare device must always be available during extracorporeal circulation.
- Check that the O-Ring is present on the HTF connectors (ref. 3c and 3d in Fig.2).
- The device has been designed to be used only with the dedicated holder provided by Spectrum Medical.
- During circuit set-up and before starting extracorporeal circulation, check that the spare medical devices and accessories indicated in the manual of the medical devices listed in section J are available.
- Be careful not to damage the connectors when connecting the tubes.
- Follow a suitable blood anticoagulation protocol before and during the procedure.
- The ACT (activated coagulation time) value must always be such as to guarantee anticoagulation suited to the procedure.

F.1. PREPARING THE HEATER COOLER

CAUTION

- The pressure of the HTF solution inside the device must not exceed the value indicated in Table 1 - section D.
- The heat exchanger temperature must not exceed the value indicated in Table 1 - section D.

1. Prepare the quick-coupling tubes provided with the Spectrum Medical Quantum Heater-Cooler System.

F.2. INSTALLING THE HEAT EXCHANGER

1. Remove the device from the sterile packaging.
2. Fit the device on the dedicated holder provided by Spectrum Medical. For more details, refer to the user manual of the Spectrum Medical Quantum HLM heart-lung machine (Quantum Pump Console User Guide).
3. Make sure the device is firmly attached to the holder.
4. Remove the protective caps from the connectors on the HTF side of the device (ref. 3c and 3d in Fig.2).
5. Connect the tubes of the Spectrum Medical Quantum Heater-Cooler System (ref. 4 in Fig. 2) to the device (ref. 3c and 3d in Fig. 2) using the female quick-coupling connectors.

F.3. CHECKING THE INTEGRITY OF THE DEVICE

1. Start priming of the HTF circuit following the instructions given in the user manual of the Spectrum Medical Quantum Heater-Cooler System.
2. The integrity of the device is guaranteed if:

- there are no HTF solution leaks from its compartment to the outside environment;
- there are no HTF solution leaks in the cardioplegia compartment.

F.4. CONNECTING THE HEAT EXCHANGER TO THE EXTRACORPOREAL CIRCUIT

CAUTION

- The maximum pressure of the cardioplegia compartment must not exceed the value indicated in Table 1 - section D.
- Only use tubes with dimensions compatible with the cardioplegia inlet/outlet connectors of the device (see Table 1 - section D).
- The overpressure valve activation pressure is indicated in Table 1 – section D.

1. Connect the device to the cardioplegia extracorporeal circuit as shown in Figure 2.
2. Remove the protective cap from the blood inlet connector (ref. 3a in Fig. 2).
3. Connect the appropriate tube to the inlet port of the device using an aseptic technique in order to maintain sterility.
4. Remove the protective cap from the blood outlet connector (ref. 3b in Fig. 2).
5. Connect the appropriate tube to the outlet port of the device using an aseptic technique in order to maintain sterility.
6. The purge Luer-lock connector (ref. 5 in Fig. 2) must be connected to a purge line equipped with a stopcock. The purge line must be equipped with a one-way valve to prevent backflow into the device.
7. Connect the free end of the purge line to a filtered Luer-lock connector of the cardiotomy/venous reservoir.
8. Connect an additional purge line to the overpressure valve (ref. 6 in Fig. 2).
9. Connect the other end of the overpressure valve purge line to the cardiotomy/venous reservoir.
10. Insert the YSI compatible temperature probe into the probe coupling of the device (ref. 7 in Fig. 2).


Fig. 2

References:

- | | | |
|--|--|--|
| 1. Cardioplegic solution and/or oxygenator cardioplegia outlet | 3b. Cardioplegic solution outlet connector | 5. Purge line |
| 2. Cardioplegia pump | 3c. HTF inlet connector | 6. Overpressure valve |
| 3. Cardioplegia heat exchanger | 3d. HTF outlet connector | 7. Temperature probe |
| 3a. Cardioplegic solution inlet connector | 4. Quantum Heater-Cooler System | 8. Connection line with cardioplegia cannula |

G. **DEVICE PRIMING PROCEDURE**

G.1. PRIMING THE CIRCUIT

CAUTION

1. Open the purge line (ref. 5 in Fig. 2).
2. Clamp the device outlet line and start the cardioplegia pump (ref. 2 in Fig. 2) at a low flow (30-40 ml/min).
3. Wait until the air in the device has been removed. If necessary, lightly tap the device to facilitate air removal.
4. Close the purge line (ref. 5 in Fig. 2).
5. Remove the clamp from the outlet line.
6. Completely fill the circuit.
7. Stop the pump and clamp the device outlet line.

G.2. REMOVING AIR BUBBLES

CAUTION

- Do not use inappropriate tools to remove any air bubbles, as they may damage the device and alter its operation.
1. Open the purge line.
 2. Remove the air by circulating solution inside the circuit at a flow of 100-150 ml/min.
 3. Verify the air removal.
 4. Close the purge line.
 5. Verify that all the air bubbles have been removed from the circuit.
 6. Stop the pump and clamp the device outlet line.

H. DEVICE USE AND REMOVAL

CAUTION

- To reduce the risk of air reaching the patient, make sure that the device and the circuit are free of air and check that they have been filled properly before starting the procedure.
- Visually check that all the air has been removed from the device.
- Carefully check for the absence of potential occlusions in the circuit before starting perfusion with cardioplegic solution.
- In case of leaks or other faults, remove the device in accordance with what is indicated in these instructions for use.
- Prolonged use of the device with stationary flow may cause the formation of clots or agglutination.
- If the device is used for blood cardioplegia and it is fitted downstream of the membrane oxygenator, make sure that the cardioplegic perfusion flow is always lower than that of systemic perfusion.
- Check that all the Luer-lock connections are secure.
- All accessory lines connected to the device must be perfectly hermetically sealed.

H.1. USING THE DEVICE DURING THE PROCEDURE

1. Let the HTF circulate at the desired temperature.
2. Remove the clamp from the device outlet line.
3. Start the cardioplegia pump at the desired flow rate.
4. Check the temperature and pressure of the cardioplegic solution by means of the sensor integrated in the device.
5. Repeat the operations (see steps 2, 3 and 4) each time you need to administer cardioplegic solution.

H.2. REMOVING THE DEVICE

CAUTION

- Remove and dispose of the device after each procedure in accordance with what is indicated in section A. Follow the instructions given in section K if used devices need to be returned.
1. Once perfusion is complete, stop the cardioplegia pump.
 2. Empty the HTF circuit following the instructions reported in the user manual of the Spectrum Medical Quantum Heater-Cooler System.
 3. Disconnect the tubes of the Spectrum Medical Quantum Heater-Cooler System heat transfer circuit.
 4. Disconnect the temperature probe.
 5. Remove the device from the cardioplegia circuit.

I. DEVICE REPLACEMENT

CAUTION

After six hours of use with perfusion fluid or if situations occur that, based on the judgement of the person responsible for perfusion, compromise the safety conditions for the patient (insufficient performance, loss of integrity of the HTF compartment or the cardioplegic solution compartment, etc.), proceed as follows to replace the device.

- If the device needs to be replaced during circulatory support, evaluate the potential risks for the patient before replacing the device.
 - Never disconnect the lines if the circuit is pressurized.
 - Visually inspect and carefully check the device before replacing it.
 - Handle the device aseptically throughout the replacement procedure.
 - Follow the instructions given in section K if the device needs to be returned.
1. Empty the HTF circuit following the instructions reported in the user manual of the Spectrum Medical Quantum Heater-Cooler System.
 2. Disconnect the tubes of the Spectrum Medical Quantum Heater-Cooler System.
 3. Remove the new device from the outer packaging and from the sterile packaging checking its integrity.
 4. Stop the cardioplegia pump following the instructions in the user manual of the Spectrum Medical Quantum HLM heart-lung machine (Quantum Pump Console User Guide).
 5. Place two clamps on both the cardioplegic solution inlet and outlet line of the device.
 6. Disconnect the temperature probe.
 7. Clamp the overpressure valve purge line and remove it.

8. Clamp the purge line and remove it.
9. Remove the device from the dedicated Spectrum Medical holder.
10. Cut the inlet and outlet lines at the point between the two clamps taking care to keep adequate lengths for reconnection.
11. Discard the device to be replaced.
12. Position the new device on the dedicated Spectrum Medical holder.
13. Fill the HTF compartment of the heat exchanger as described in section F.2 and check its integrity as described in section F.3.
14. Start the Spectrum Medical Quantum Heater-Cooler System.
15. Connect the cardioplegic solution inlet line to the cardioplegic solution inlet connector and the cardioplegic solution outlet line to the cardioplegic solution outlet connector of the new device.
16. Use cable-ties to externally lock the inlet and outlet tubes on the respective connectors of the device.
17. Connect the purge line to the new device.
18. Connect the overpressure valve purge line to the new device.
19. Connect the temperature probe.
20. Fill the device carrying out the operations described in section G.
21. Check that there is no air in the device.
22. If there is air in the device, follow the instructions given in section G.2.
23. Open the outlet line only after having checked that there is no air in the device.
24. Open the outlet line.
25. Resume perfusion.

J. MEDICAL DEVICES TO BE USED WITH THE DEVICE

CAUTION

- The user must respect the warnings and precautions and follow the instructions for use provided with the device.
- The connections of the circuit in which the cardioplegic solution circulates must be made using tubes of a diameter compatible with the dimensions of the connectors on the device in accordance with Table 1 - section D.
- The temperature must be monitored with a YSI 400 compatible probe.

The device is qualified to be used exclusively with:

- Spectrum Medical Quantum Heater-Cooler System.
- Spectrum Medical Quantum HLM heart-lung machine and accessories.

K. RETURN OF USED DEVICES AND NOTIFICATIONS

CAUTION

- Always keep any defective medical devices. Follow the instructions below if the devices need to be sent for investigation.
- Healthcare institutions are responsible for preparing and identifying the products for return shipments. If the product to be returned has been in contact with the patient's blood, it must be treated, packed and handled in compliance with the provisions of the legislation in force in the country where the product has been used. Do not send products that have been contaminated by infectious pathogenic agents.
- Qura stresses the importance of reporting any incident that occurred regarding the device to Qura and to the competent authority of the Member State where the product was used. Should the user be dissatisfied with anything related to the quality of the product, the user may notify Qura or the distributor of the Qura product. All notifications considered critical by the user must be reported with particular care and urgency. Below is the minimum information to be provided:
 - thorough description of the event and, if pertinent, the conditions of the patient;
 - identification of the product involved;
 - batch number of the product involved;
 - availability of the product involved;
 - all the indications the user considers useful in order to understand the origin of the dissatisfaction.
- Qura reserves the right to authorize the return of the product involved in the notification for assessments.

L. LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the purchaser pursuant to applicable law. Qura guarantees that all reasonable precautions have been taken during the manufacture of this medical device, as required by the nature of the device and its intended use. Qura guarantees that its medical device is capable of functioning as indicated in these Instructions For Use when used in accordance with that specified in these Instructions by qualified personnel and within the expiry date indicated on the packaging. However, Qura cannot guarantee that the user will use the device correctly, nor that incorrect diagnosis or therapy and/or the particular physical and biological characteristics of an individual patient do not affect the performance and effectiveness of the device or cause damaging consequences for the patient, even though the specified Instructions for Use have been respected. Whilst emphasising the need to adhere strictly to the Instructions for Use and to adopt all the precautions necessary for correct use of the device, Qura therefore cannot assume any responsibility for any incidents, harm, loss, expense or consequences arising directly or indirectly from improper use of, tampering with or modifications made to this device. Qura will replace the medical device if it is defective at the time it is put on the market or while being shipped by Qura up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties, explicit or implicit, written or verbal, including warranties of marketability and/or functionality. No person, including any representative, agent, dealer, distributor or intermediary of Qura or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. Qura disclaims any variations to this Limited Warranty and to the information/Instructions for Use as explicitly contained herein. The purchaser pledges to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with Qura, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The relationship between the parties regarding the contract (even if not stipulated in writing) to whom this warranty is issued, as well as any dispute related to it or in any way connected with it, as well as any relationship or dispute concerning this warranty, its interpretation and execution, nothing excluded and/or reserved, are governed exclusively by Italian law and jurisdiction. The court appointed in the event of the above-said disputes shall exclusively be the Court of Modena (Italy).

Qura S.r.l.
Via di Mezzo, 23
41037 Mirandola (MO), Italy
Tel. +39 0535-1803050
www.quramed.com

M. SYMBOLS USED

	Manufacturer		Temperature limit
	Use-by Date		Non-pyrogenic
	Batch code		For US Audiences Only
	Catalogue number		Not made with natural rubber latex
	Sterilized using ethylene oxide		Serial number
	Do not resterilize		This way up
	Do not use if the packaging is damaged		Fragile, handle with care
	Do not reuse		Keep dry
	Consult the instructions for use		Keep away from sunlight
	Caution		Open here
	Compliance with Directive 93/42/EC		Packaging unit

Quantum PureFlow Cardioplegia Heat Exchangers

HX-C

REF.	DESCRIPTION	CODE	DESCRIPTION
HX32V-C0	Cardioplegia Heat Exchanger High Flow 1/4	HX32V-C4	Cardioplegia Heat Exchanger High Flow 1/4 - PL
HX22V-C0	Cardioplegia Heat Exchanger Low Flow 3/16	HX22V-C4	Cardioplegia Heat Exchanger Low Flow 3/16 - PL
HX22V-C1	Cardioplegia Heat Exchanger Low Flow 1/4	HX22V-C5	Cardioplegia Heat Exchanger Low Flow 1/4 - PL

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Medtronic Australasia Pty Ltd
Qura S.r.l. - Via di Mezzo, 23 - 41037 Mirandola (MO) Italy