ACTIVATED CLOTTING TIME CARTRIDGES  
 CARTOUCHES POUR TEMPS DE COAGULATION ACTIVÉE  
 KARTUSCHEN ZUR BESTIMMUNG DER AKTIVIERTEN  
 GERINNUNGSZEIT  
 CARTUCHOS DE TIEMPO DE COAGULACIÓN ACTIVADA  
 PATRONEN VOOR GEACTIVEERDE STOLLINGSTIJD  
 CARTUCCE PER IL TEMPO DI COAGULAZIONE ATTIVATA  
 AKTIVOIDUN HYYTYMISAJAN TESTIPUTKET  
 TESTRÖR FÖR BESTÄMNING AV AKTIVERAD  
 KOAGULATIONSTID  
 ΦΙΑΛΙΔΙΑ ΓΙΑ ΤΟΝ ΕΝΕΡΓΟΠΟΙΗΜΕΝΟ ΧΡΟΝΟ ΠΗΞΗΣ  
 AKTIVE KOAGULATIONSTIDS KASSETTER  
 CARTUCHOS DE TEMPO DE COAGULAÇÃO ACTIVADA


Low Range (LR-ACT) / Valeurs faibles (LR-ACT) / Low Range (LR-ACT) / Bajo rango (LR-ACT) / Laag bereik (LR-ACT) / Basso range (LR-ACT) / Low Range (LR-ACT) / Lågt värde (LR-ACT) / Χαμηλού Εύρους (LR-ACT) / Low Range (LR-ACT) / Valores baixos (LR-ACT)

Recalcified (RACT) / Recalcifié (RACT) / Recalcified (RACT) / Recalcificado (RACT) / Recalcified (RACT) / Ricalcificata (RACT) / Recalcified (RACT) / Rekalcifierad (RACT) / Επανασβεστοποιημένα (RACT) / Recalcificeret (RACT) / Recalcificado (RACT)




High Range (HR-ACT) / Valeurs élevées (HR-ACT) / High Range (HR-ACT) / Alto rango (HR-ACT) / Hoog bereik (HR-ACT) / Ampio range (HR-ACT) / High Range (HR-ACT) / Högt värde (HR-ACT) / Υψηλού Εύρους (HR-ACT) / High Range (HR-ACT) / Valores altos (HR-ACT)

Instructions for Use ■ Mode d'emploi ■ Gebrauchsanweisung ■ Instrucciones de uso ■ Gebruiksaanwijzing  
 Istruzioni per l'uso ■ Käyttöohjeet ■ Bruksanvisning ■ Οδηγίες χρήσης ■ Brugsanvisning ■ Instruções de utilização

Medtronic Confidential

**Explanation of symbols on package labeling**  
**Explication des symboles des étiquettes sur l'emballage**  
**Erklärung der Symbole auf dem Verpackungsetikett**  
**Explicación de los símbolos en el etiquetado del envase**  
**Verklaring van de symbolen op de verpakkingslabels**  
**Spiegazione dei simboli sulle etichette della confezione**  
**Pakkausen merkintöjen selitykset**  
**Förklaring av symboler på förpackningsetiketten**  
**Επεξήγηση των συμβόλων στην ετικέτα της συσκευασίας**  
**Forklaring af symboler på emballagens mærkater**  
**Explicação dos símbolos nas etiquetas da embalagem**

Conformité Européenne (European Conformity). This symbol means that the device fully complies with the European Council Directive 98/79/EC (27 October 1998, on in vitro diagnostic medical devices)

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## INTENDED USE

### 1 INTENDED USE

For use in determining the activated clotting time of whole blood.

For In Vitro Diagnostic Use.

### 2 SUMMARY

The activated clotting time test was originally described by Hattersley in 1966 (Hattersley, 1966). The test and its variations have been widely used to monitor the anticoagulant effect of heparin.<sup>1,2,3,5,7,12</sup> Medtronic Activated Clotting Time Cartridges utilize either fresh whole blood or citrated whole blood. Three types of cartridges are available (Table I). The Low Range Activated Clotting Time Cartridge, LR-ACT (catalog # 402-01) and the High Range Activated Clotting Time Cartridge, HR-ACT (catalog # 402-03), are for use with fresh drawn whole blood, and the Recalcified Activated Clotting Time Cartridge, RACT (catalog # 402-02) is for use with citrated whole blood. The two fresh draw ACT cartridges are distinguished by their respective sensitivities to the anticoagulant heparin. The HR-ACT cartridge has a decreased response to heparin which favors its use in cardiovascular surgery. The LR-ACT cartridge has increased response to heparin, which favors its use in therapeutic heparin administration.<sup>14</sup>

### 3 PRINCIPLES OF THE PROCEDURE

The activated clotting time test measures the clotting time of fresh whole blood activated by surface contact. Hattersley's original procedure used diatomaceous earth as the activating agent. Diatomaceous earth exhibits differing degrees of contact activation depending on its purity, the source of the deposit, and the particle size and shape. Activated Clotting Time Cartridges use a specially prepared and standardized kaolin as the activating agent. This kaolin exhibits superior characteristics in contact activation of blood coagulation and minimizes lot to lot variation in activation characteristics.

The end point of the test is the detection of clot (fibrin) formation. Medtronic Automated Coagulation Timers detect clot formation by measuring the rate of fall of the plunger-flag mechanism contained in each cartridge channel. The plunger assembly falls rapidly through an unclotted sample, but the fibrin web formed during clotting impedes the rate of descent. This is detected by a photo optical system located in the cover assembly of the instrument. Automated Coagulation Timers perform simultaneous duplicate tests using the two channel cartridge. The clotting times for each channel are displayed in seconds following completion of the test. By depressing the front panel display switch on the ACT II® instrument the average and difference of the duplicate tests can be displayed.

### 4 REAGENTS

Table 1 describes the cartridge reagents (0.1 ml per channel) and maximum prewarm times for the cartridges and cartridges + sample:

Table 1.

Cartridge Type	Reagent	Cartridge Max. Prewarm Limits	Cartridge + Sample Prewarm Limits
HR-ACT	12% Kaolin 0.05M CaCl <sub>2</sub> HEPES Buffer <sup>a</sup> Sodium Azide <sup>b</sup>	12 hours	Non-applicable. Perform test immediately.
LR-ACT	0.75% Kaolin <sup>c</sup> 0.0025M CaCl <sub>2</sub> HEPES Buffer Sodium Azide	12 hours	Non-applicable. Perform test immediately.
RACT	2.2% Kaolin <sup>c</sup> 0.05M CaCl <sub>2</sub> HEPES Buffer Sodium Azide	12 hours	5 min

<sup>a</sup> HEPES (hydroxyethyl-piperazine-ethanesulfonic acid) buffer.

<sup>b</sup> Sodium Azide is a bacteriostatic agent.

<sup>c</sup> Kaolin concentration may vary slightly.

### 5 STORAGE

Store the Activated Clotting Time Cartridges from 2 °C to 25 °C (36 °F to 77 °F) for the expiration dating period indicated on the cartridge container. Do not use if cartridge appears discolored, cracked or disfigured, or if the reagent appears to have evaporated or be contaminated.

Do not freeze the cartridges.

## INSTRUMENT

Store cartridges in box until time of use.

## 6 INSTRUMENT

Medtronic Activated Clotting Time Cartridges are specifically designed for use in Medtronic Automated Coagulation Timers. For instrument instructions consult the operator's manual.

If blood should get into the light paths, it is critical the instrument be cleaned as soon as possible. In order to ensure proper instrument performance, it is important to clean the instrument a minimum of every thirty (30) days or more. An Actuator Cleaning Kit (catalog # 303-58) is provided with each instrument. See details in Cleaning Kit package insert.

## 7 CARTRIDGES

The following Medtronic Activated Clotting Time Cartridges are available for use with the Medtronic Automated Coagulation Timer (Table 2).

Table 2.

Catalog #	Product	Sample Type	Channel Sample Volume	Intended Use
402-01	LR-ACT	Fresh drawn	0.2 ml	Therapeutic (0-1.5 Units Heparin/ml)
402-02	RACT	Citrated	0.2 ml	Therapeutic (0-1.5 Units Heparin/ml)
402-03	HR-ACT	Fresh drawn	0.4 ml	Cardiovascular Cath Lab (1 Unit Heparin/ml or more)

## 8 SPECIMEN COLLECTION AND PREPARATION

Improperly handled samples can yield erroneous coagulation test results; therefore, samples must be collected carefully and handling techniques must be consistent.

If blood is obtained by venipuncture, the venipuncture should be clean and the first several milliliters of blood discarded to avoid contamination of the sample with tissue thromboplastin. Blood should not be taken from heparinized indwelling catheters or other anticoagulated lines without thoroughly flushing the line.

If testing fresh whole blood, samples should be run as quickly as possible after being drawn. Baseline samples, which contain no anticoagulant, should be run within one (1) minute. Heparinized samples should be run within two (2) minutes.

If testing citrated whole blood, test the sample immediately or store the sample at room temperature and test samples from any one patient in the same time frame. For example, if a patient sample is tested 15 minutes after drawing then subsequent samples for that patient should be tested at 15 minutes.

Citrated samples should be tested within 6 hours when stored at room temperature. Ensure that the sample has the proper amount of citrate (a 1:9 volume ratio of 3.2% or 3.8% sodium citrate to blood).

## 9 PROCEDURE

This product is intended for use only by trained individuals.

A. Materials Supplied:

High Range Activated Clotting Time Cartridge (HR-ACT), catalog # 402-03

or

Recalcified Activated Clotting Time Cartridge (RACT), catalog # 402-02

or

Low Range Activated Clotting Time Cartridge (LR-ACT), catalog # 402-01

B. Materials required, not supplied:

Automated Coagulation Timer with heat block at  $37^{\circ} \pm 0.5^{\circ}\text{C}$  ( $98.6^{\circ} \pm 0.9^{\circ}\text{F}$ ).

Syringe with needle (preferably blunt 18 or 19 gauge, Syringe catalog # 300-01 and Blunt Needle catalog # 300-02) or device that delivers 0.200 ml or 0.400 ml, depending on the test performed.

Fresh whole blood or citrated whole blood sample, depending on the test performed.

C. Sequence of Operation:

ACT II®

(1) Set Incubate off for the HR-ACT and LR-ACT Cartridges with fresh whole blood.

## RESULTS

- (2) Set Incubate on for the RACT Cartridges with citrated whole blood. For instrument instructions consult the operator's manual.
- (3) Tap the cartridge to resuspend the activator. Make sure any activator clinging to the cartridge wall is dislodged. Tapping the bottom of the cartridge against a solid surface or flicking with a finger will facilitate this dislodging.
- (4) HR-ACT and LR-ACT Cartridges
  - a. Warm the cartridge in the heat block to  $37^{\circ} \pm 0.5^{\circ}\text{C}$  ( $98.6^{\circ} \pm 0.9^{\circ}\text{F}$ ) (allow 3-5 minutes) then add the required amount of blood sample to each cartridge channel (HR-ACT: 0.4 ml/ channel; LR-ACT: 0.2 ml/channel). Each channel is marked with fill lines that indicate upper and lower limits of the sample volume for each channel.
  - b. Start the test by rotating the actuator heat block to the closed position.
- (5) RACT Cartridges
  - a. Add the required amount of blood sample (0.2 ml/channel) to each channel of the cartridge. Each channel is marked with fill lines that indicate the upper and lower limits of the sample volume.
  - b. Simultaneously warm the cartridge and blood sample to  $37^{\circ} \pm 0.5^{\circ}\text{C}$  ( $98.6^{\circ} \pm 0.9^{\circ}\text{F}$ ) by setting Incubate on and starting the test by rotating the actuator heat block to the closed position. The test will automatically begin after the incubation cycle is complete.


Figure 1. Filling the Cartridge

- (6) Results are displayed in seconds in two formats: a) the test time for each channel, or b) the average and difference of the channels. Pressing the front panel switch controls the display format.
- (7) Following termination of the test, dispose of the cartridge, syringe, and needle in approved biohazard disposable containers.

### ACT Plus™

- (1) From the Main Menu select [Cartridge Type], press the key until the correct cartridge is displayed. Press **Enter** to confirm selection.  
**Note:** Lot Numbers and Expiration Dates for cartridges and controls must be entered prior to running a test.
- (2) From the Main Menu select **Enter ID**. Enter the **Patient ID** and **User ID** numbers. Press **Enter** to confirm.
- (3) Prewarm the cartridge for 3 - 5 minutes before collecting the test sample.
- (4) Tap the cartridge to resuspend the reagent.
- (5) Fill each cartridge chamber with the appropriate sample volume:

Cartridge	Sample Type	Sample Volume per Channel
HR-ACT	Fresh Whole Blood	0.4 ml
LR-ACT	Fresh Whole Blood	0.2 ml
RACT	Citrated Whole Blood	0.2 ml

- (6) Insert the cartridge into the ACT Plus™ and close the actuator heat block to start test.
- (7) Clot formation is signaled by an audible tone, the actuator heat block opens and the results are displayed.

## 10 RESULTS

The response of the activated clotting time to heparin varies considerably from individual to individual.<sup>1,2,4</sup> Various drugs affect the activated clotting time; in particular, drugs which inhibit platelet activation.<sup>15</sup> A number of other factors can affect the response of an individual's activated clotting time to heparin (e.g. antithrombin III levels, heparin potency, coagulation

#### LIMITATIONS OF THE PROCEDURE

factor deficiencies, sample activation, consumptive coagulopathies, excessive sample dilution, and sample temperature) and should be taken into account when interpreting the results of activated clotting time tests.<sup>8,9,11</sup>

Duplicate channel results should fall within 10% of each other for baseline (unheparinized) samples and within 12% of each other for extended or heparinized samples.

Sample Calculation:

Ch 1 clotting time = 210 seconds  
Ch 2 clotting time = 200 seconds  
Mean clotting time = 205 seconds  
Difference = 10 seconds  
12% of the mean = 25 seconds

The difference of 10 seconds is less than 12% of the mean, 25 seconds. These results are acceptable.

#### 11 LIMITATIONS OF THE PROCEDURE

To obtain accurate and consistent activated clotting times, the following conditions must be observed:

1. Both the instrument heat block and cartridge temperature should be at  $37^{\circ} \pm 0.5^{\circ}\text{C}$  ( $98.6^{\circ} \pm 0.9^{\circ}\text{F}$ ). (The cartridge will come to temperature in approximately 3 minutes in the instrument heat block.)
2. Activator must be thoroughly resuspended.
3. Blood must be free of tissue thromboplastin and run as soon as possible after drawing (fresh drawn), or stored in a non-activating, capped container (citrated whole blood).
4. Patient diagnosis and medications should be noted. Medications can alter clotting times.<sup>15</sup> Samples which give unexplained abnormal values should be redrawn and retested. If the results are still abnormal and no cause is apparent, a sample should be citrated and sent to your clinical lab for special coagulation analysis.
5. To optimize precision, all technique variables should be held constant from test to test.
6. Methodologies used in cardiopulmonary bypass surgery continually change as does the population demographics of bypass patients. Patient sensitivity to heparin, dilution of clotting factors by the extracorporeal circuit, the use of citrated blood products, the use of antiplatelet drugs, hypothermia, fluctuating calcium levels, a change in platelet function or count, unknown coagulopathies, etc., are influenced by these changes.<sup>13</sup> Excessive blood dilution during bypass, frequently due to a desire to limit the use of heterologous blood products, can lead to postsurgically hemodiluted patients. This dilution can be great enough that clotting times, as reflected by the activated clotting time, are extended due to the dilution.<sup>10</sup> In severe instances, the dilution may be so great that measurable clotting does not occur. In this case, the patient may experience postoperative bleeding due to a lack of plasma clotting factors. The adequacy of plasma clotting factors cannot be determined by the patient's hematocrit. Frequently, washed cells are administered to the patient without compensating for lost plasma. If the HR-ACT fails to clot following reversal of heparinization and a protamine titration indicates no heparin is present (Heparin Assay is not affected to the same degree by dilution as is an activated clotting time), a dilutional coagulopathy is a strong possibility. If a presurgical coagulation panel is available for the patient, either an activated partial thromboplastin time (aPTT) or a prothrombin time (PT), a comparison of the pre- and post-surgical values should indicate an extension. The degree of extension is patient-related and cannot be judged by reference to normal control clotting time ranges. An aPTT which changes from the low normal range to the high normal range is indicative of a significant change in patient coagulation characteristics.

**Warning:** Due to the formulation of the HR-ACT cartridge, it should not be used to monitor regional anticoagulation achieved by infusing citrate.

For instrument precautions and limitations, consult the operator's manual.

#### 12 EXPECTED VALUES

Clotting times using the Medtronic Low Range Activated Clotting Time Cartridge (LR-ACT), the Medtronic Recalcified Activated Clotting Time Cartridge (RACT) and the Medtronic High Range Activated Clotting Time Cartridge (HR-ACT) were determined on 3 normal volunteers. Table 3 is a sample of an individual response to heparin using 3 donors. The HR-ACT and the LR-ACT cartridges were run with fresh drawn whole blood, and the RACT cartridges were run with citrated whole blood (3.2% sodium citrate). Response to an increased concentration of heparin may vary between individuals. Each laboratory/facility should establish their own range based on their local demographics.

### SPECIFIC PERFORMANCE CHARACTERISTICS

Table 3.

Low Range ACT Cartridge		Recalcified ACT Cartridge		High Range ACT Cartridge	
Heparin Level <sup>a</sup>	Avg. Clot Time	Heparin Level	Avg. Clot Time	Heparin Level	Avg. Clot Time
Donor 1					
0.0	118	0.0	130	0.0	132
1.0	357	1.0	294	1.0	183
2.0	591	2.0	472	2.0	286
3.0	712	3.0	615	3.0	414
4.0	926	4.0	842	4.0	474
5.0	>999	5.0	>999	5.0	622
6.0	>999	6.0	>999	6.0	850
Donor 2					
0.0	106	0.0	105	0.0	109
1.0	258	1.0	221	1.0	145
2.0	387	2.0	310	2.0	214
3.0	652	3.0	404	3.0	282
4.0	790	4.0	509	4.0	334
5.0	879	5.0	624	5.0	389
6.0	>999	6.0	850	6.0	477
Donor 3					
0.0	106	0.0	121	0.0	116
1.0	254	1.0	308	1.0	164
2.0	426	2.0	440	2.0	264
3.0	516	3.0	572	3.0	323
4.0	648	4.0	804	4.0	427
5.0	734	5.0	>999	5.0	497
6.0	>999	6.0	>999	6.0	679

<sup>a</sup> Heparin level in units of heparin per ml of whole blood. Heparin was added to each blood sample just prior to running the cartridge.

## 13 SPECIFIC PERFORMANCE CHARACTERISTICS

The Medtronic Activated Clotting Cartridges in the clotting time range of 0–600 seconds, typically do not exceed a variation of 12% between duplicate samples for any one individual. In subjects with extreme sensitivity to heparin, an activated clotting time in excess of 600 seconds is generally not considered to be adequately reliable enough to predict heparin effect by coagulation times.

The LR-ACT cartridge (Catalog #402-01) is intended for therapeutic heparin monitoring via the activated clotting time of fresh drawn whole blood. Its increased sensitivity to heparin allows better control of heparin anticoagulation at therapeutic heparin levels.

The RACT cartridge (Catalog #402-02) is a general-purpose recalcified activated clotting time cartridge.

The instrument's operable range is up to 999 seconds.

## 14 QUALITY CONTROL

Two levels of controls should be run when the instrument is in operation to ensure adequate quality control of the instrument and cartridges.

Results should fall within the specified ranges listed by the manufacturer or ranges determined by the user.

Details concerning Quality Assurance can be found in the Medtronic QA Manual, (catalog # 300-29).

CLOTTrac® CWB Coagulation Controls (Catalog #550-01) are available to verify the correct functioning of the RACT and LR-ACT cartridges, and the Automated Coagulation Timers. LR-ACT cartridges require the addition of 0.1M Calcium Chloride (Catalog #550-11) (see calcium chloride package insert for directions).

## QUALITY CONTROL

CLOTtrac® LR Abnormal Controls (Catalog #550-09) are available to verify the correct functioning of the LR-ACT cartridges and the Automated Coagulation Timers in the extended clotting range. LR-ACT cartridges require the addition of one drop (40 microliters) of 0.1 M Calcium Chloride, Catalog #550-11, (see calcium chloride package insert for directions).

CLOTtrac® RACT Abnormal Controls (Catalog #550-10) are available to verify the correct functioning of the RACT cartridges and the Automated Coagulation Timers in the extended clotting range.

CLOTtrac® HR Coagulation Controls (Catalog #550-07) are available to verify the correct functioning of the HR-ACT cartridge and the Automated Coagulation Timers.

CLOTtrac® HR Abnormal Coagulation Controls (Catalog #550-08) are available to verify the performance of the HR-ACT cartridge and the Automated Coagulation Timers in the extended clotting range.

## DISCLAIMERS OF WARRANTY

### 15 DISCLAIMERS OF WARRANTY

#### THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO UNITED STATES CUSTOMERS ONLY:

- Le déni de garantie ci-après vaut uniquement pour les clients des Etats-Unis : /
- Der nachfolgende Haftungsausschluß gilt nur für Kunden in den Vereinigten Staaten: /
- La siguiente renuncia de responsabilidad se aplica sólo a los clientes de los Estados Unidos: /
- De volgende uitsluiting van garantie geldt uitsluitend voor klanten in de Verenigde Staten: /
- La seguente esclusione dalla garanzia è valida solo per i clienti negli Stati Uniti: /
- Seuraava takuu vastuuvaapauslauseke koskee vain Yhdysvaltojen asiakkaita: /
- Nedanstående friskrivning från ansvar är endast tillämplig på kunder i USA: /
- Η ακόλουθη δήλωση αποποίησης εγγύησης ισχύει μόνον για πελάτες από τις Ηνωμένες Πολιτείες: /
- Nedenstående ansvarsfriskrivelse er kun gældende for kunder i USA:
- A seguinte exoneração de responsabilidade aplica-se apenas aos clientes nos Estados Unidos:

#### DISCLAIMER OF WARRANTY

ALTHOUGH THE ACTIVATED CLOTTING TIME CARTRIDGES, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

#### THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

#### DISCLAIMER OF WARRANTY

ALTHOUGH THE ACTIVATED CLOTTING TIME CARTRIDGES, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Tabela 3.

Cartucho ACT de valores baixos		Cartucho ACT recalcificado		Cartucho ACT de valores altos	
Nível de heparina <sup>a</sup>	Tempo de coagulação médio	Nível de heparina	Tempo de coagulação médio	Nível de heparina	Tempo de coagulação médio
Dador 1					
0,0	118	0,0	130	0,0	132
1,0	357	1,0	294	1,0	183
2,0	591	2,0	472	2,0	286
3,0	712	3,0	615	3,0	414
4,0	926	4,0	842	4,0	474
5,0	>999	5,0	>999	5,0	622
6,0	>999	6,0	>999	6,0	850
Dador 2					
0,0	106	0,0	105	0,0	109
1,0	258	1,0	221	1,0	145
2,0	387	2,0	310	2,0	214
3,0	652	3,0	404	3,0	282
4,0	790	4,0	509	4,0	334
5,0	879	5,0	624	5,0	389
6,0	>999	6,0	850	6,0	477
Dador 3					
0,0	106	0,0	121	0,0	116

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Cartucho ACT de valores baixos		Cartucho ACT recalcificado		Cartucho ACT de valores altos	
Nível de heparina <sup>a</sup>	Tempo de coagulação médio	Nível de heparina	Tempo de coagulação médio	Nível de heparina	Tempo de coagulação médio
1,0	254	1,0	308	1,0	164
2,0	426	2,0	440	2,0	264
3,0	516	3,0	572	3,0	323
4,0	648	4,0	804	4,0	427
5,0	734	5,0	>999	5,0	497
6,0	>999	6,0	>999	6,0	679

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1. Bull BS, Korpman RA, Huse WM, Briggs BD. Heparin therapy during extra-corporeal circulation. I. Problems inherent in existing heparin protocol. *J Thorac Cardiovasc Surg.* 1975;69:674-684.
2. Bull BS, Huse MH, Brauer FS, Korpman RA. Heparin therapy during extra-corporeal circulation. II. The use of a dose-response curve to individualize heparin and protamine dosage. *J Thorac Cardiovasc Surg.* 1975;69:685-689.
3. Dauchot PJ, Berzina-Moettus L, Rabinovitch A, Ankeney JL. Activated coagulation and activated partial thromboplastin times in assessment and reversal of heparin-induced anticoagulation for cardiopulmonary bypass. *Anesth Analg.* 1983;62:710-719.
4. Effeney DJ, Goldstone J, Chin D, Krupski WC, Ellis RJ. Interoperative anticoagulation in cardiovascular surgery. *Surgery.* 1981;90:1068-1074.
5. Esposito RA, Culliford AT, Colvin SB, Thomas SJ, Lackner H, Spencer FC. The role of the activated clotting time in heparin administration and neutralization for cardiopulmonary bypass. *J Thorac Cardiovasc Surg.* 1983;85:174-185.
6. Hattersley PG. Activated coagulation time of whole blood. *JAMA.* 1966;196:436.
7. Hattersley PG, Mitsuoka JC, King JH. The use of the activated coagulation time (ACT) as a monitor of heparin therapy. In: Lundblad RL, Brown WV, Mann KG, Roberts HR, eds. *Chemistry and Biology of Heparin.* North Holland: Elsevier; 1981;659.
8. Henry JB. *Clinical Diagnosis and Management by Laboratory Methods.* Vol 1. 16<sup>th</sup> ed. Philadelphia, PA: WB Saunders Company; 1979.
9. Moorehead MT, Westengard JC, Bull BS. Platelet involvement in the activated coagulation time of heparinized blood. *Anesth Analg.* 1984;63:394-398.
10. Niinikoski J, Laato M, Laaksonen V, Jalonen J, Ingerg MV. Use of activated clotting time to monitor anticoagulation during cardiac surgery. *Scand J Thor Cardiovasc Surg.* 1984;18:57-61.
11. Reno WJ, Rotman M, Grumbine FC, Dennis LH, Mohler ER. Evaluation of the BART test (a modification of the whole-blood activated recalcification time test) as a means of monitoring heparin therapy. *Am J Clin Path.* 1974;61:78-84.
12. Roth JA, Cukingnan RA, Scott CR. Use of activated coagulation time to monitor heparin during cardiac surgery. *Ann Thorac Surg.* 1978;28:69-72.
13. Sirridge MS. *Laboratory Evaluation of Hemostasis.* 2<sup>nd</sup> ed. Philadelphia, PA: Lea and Febiger; 1974.
14. Wilhelmsson S, Lins L-E. Whole-blood activated coagulation time for evaluation of heparin activity during hemodialysis: A comparison of administration by single-dose and by infusion. *Clin Nephrol.* 1983;19:82-86.
15. Young DS, Pastaner LC, Gibberman V. Effects of drugs on clinical laboratory tests. *Clin Chem.* 1975;21:1D.

231036-00

Blank For Position Only. Do Not Print!

#### Europe

##### Europe/Africa/Middle East Headquarters

Medtronic Europe Sàrl  
Route du Molliau 31  
Case Postale  
CH-1131 Tolochenaz  
Switzerland  
Internet: [www.medtronic.co.uk](http://www.medtronic.co.uk)  
Tel. 41-21-802-7000  
Fax 41-21-802-7900

##### Medtronic E.C. Authorized Representative/Distributed by

Medtronic B.V.  
Earl Bakkenstraat 10  
6422 PJ Heerlen  
The Netherlands  
Tel. 31-45-566-8000  
Fax 31-45-566-8668

#### Asia-Pacific

##### Japan

Medtronic Japan  
Solid Square West Tower 6F  
580 Horikawa-cho, Sawai-ku  
Kawasaki, Kanagawa 210-0913  
Japan  
Tel. 81-44-540-6112  
Fax 81-44-540-6200

##### Australia

Medtronic Australasia Pty. Ltd.  
Unit 4/446 Victoria Road  
Gladesville NSW 2111  
Australia  
Tel. 61-2-9879-5999  
Fax 61-2-9879-5100

##### Asia

Medtronic International Ltd.  
Suite 1602 16/F, Manulife Plaza  
The Lee Gardens, 33 Hysan Avenue  
Causeway Bay  
Hong Kong  
Tel. 852-2891-4068  
Fax 852-2591-0313

#### Latin America

Medtronic Inc.  
710 Medtronic Parkway  
Minneapolis, MN 55432-5604  
USA  
Tel. 763-514-4000  
Fax 763-514-4879

#### North America

Canada  
Medtronic of Canada Ltd.  
6733 Kitimat Road  
Mississauga, Ontario L5N 1W3  
Canada  
Tel. 905-826-6020  
Fax 905-826-6620  
Toll-free in Canada:  
1-800-268-5346

#### United States

Manufacturer  
Medtronic  
710 Medtronic Parkway  
Minneapolis, MN 55432-5604  
USA  
Internet: [www.medtronic.com](http://www.medtronic.com)  
Tel. 763-514-4000  
Toll-free in the USA:  
1-800-328-2518

##### Medtronic Perfusion Systems

7611 Northland Drive  
Minneapolis, MN 55428-9947  
Tel. 763-391-9000  
Fax 763-391-9100  
Information Toll-free:  
1-800-328-3320  
Customer Service and Product  
Orders Toll-free: 1-800-854-3570  
Visit Medtronic at:  
[www.perfusionsystems.com](http://www.perfusionsystems.com)

\* 2 3 1 0 3 6 \*

© Medtronic, Inc. 2004  
All Rights Reserved

UC199600212b ML 231036-00  
May 2004