# Medtronic

## autoLog IQ™

Autotransfusion System

Autotransfusionssystem

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User Manual • Benutzerhandbuch • Manual del usuario • Manuel d'utilisation • Manuale per l'utente

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.

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autoLog IQ, Medtronic

## Contents

<b>1 Explanation of packaging, labeling, and device symbols</b>	<b>5</b>
<b>2 Introduction</b>	<b>7</b>
2.1 Advantages of autotransfusion over allogeneic transfusion	8
2.2 Principles of operation	8
<b>3 Intended use</b>	<b>9</b>
<b>4 Contraindications</b>	<b>9</b>
<b>5 Warnings and precautions</b>	<b>9</b>
<b>6 Adverse effects</b>	<b>11</b>
<b>7 Device description</b>	<b>12</b>
7.1 Device components	12
7.2 Items included with the device	18
7.3 Optional items not included with the device	18
7.4 Audible alert signals	18
7.5 Visual alert signals	18
<b>8 In vitro performance</b>	<b>19</b>
<b>9 Transportation and storage</b>	<b>19</b>
9.1 Transporting the device	19
9.2 Transporting the device between clinical sites	19
9.3 Storing the device	19
<b>10 Setting up the optional barcode scanner</b>	<b>20</b>
<b>11 Turning on the device power</b>	<b>20</b>
11.1 Connecting the device power cord	20
11.2 Turning on the device power	20
<b>12 Managing device settings</b>	<b>21</b>
12.1 Settings screen	21
12.2 Setting the time	21
12.3 Setting the alert theme	22
12.4 Initial password setup	22
12.5 Updating or resetting the password	23
12.6 Main screen overview	23
12.7 Function capabilities	24
<b>13 Operating the device</b>	<b>25</b>
13.1 Supplies required	25
13.2 Entering case data	25
13.3 Setting up the device for blood collection	26
13.4 Setting up the wash kit	28
13.5 Initiating blood processing	30
13.6 Resuming a case	31
<b>14 Managing case data</b>	<b>31</b>
14.1 Viewing case data	31
14.2 Searching and filtering case data	32
14.3 Exporting case data to a USB flash drive	33
14.4 Exporting case data to an external device	34
14.5 Deleting case data	35
<b>15 Ending a case</b>	<b>35</b>
<b>16 Cleaning the device</b>	<b>35</b>
<b>17 Special blood processing conditions</b>	<b>37</b>
17.1 Emergency setup using multiple suction/anticoagulant lines on 1 reservoir	37
17.2 Using 2 reservoirs in series	38
17.3 Using 2 reservoirs in parallel	40
17.4 Draining cardiopulmonary bypass circuit (CPB) contents to the device	41
17.5 Postoperative wound drainage	41
<b>18 Troubleshooting guide</b>	<b>41</b>
18.1 Power failure	41
18.2 Unusual behavior of the device	42

18.3 Normal startup prompts, audible and visual alert signals . . . . . 42

18.4 Normal operating prompts, audible and visual alert signals . . . . . 43

18.5 Error prompts, audible and visual alert signals . . . . . 47

18.6 Irrecoverable errors, audible and visual alert signals . . . . . 47

**19 Service and maintenance . . . . . 48**

19.1 Service interval . . . . . 48

19.2 Pump head . . . . . 48

19.3 Service and maintenance screens . . . . . 48

19.4 Notifications . . . . . 49

19.5 System information . . . . . 50

19.6 Error log . . . . . 50

19.7 Calibrate touchscreen . . . . . 50

19.8 Clean centrifuge . . . . . 50

19.9 Update software . . . . . 50

19.10 Export service information . . . . . 50

**20 Device disposal and ownership transfer . . . . . 51**

20.1 Deleting case data . . . . . 51

20.2 End of life disposition . . . . . 51

**21 Specifications . . . . . 51**

21.1 Technical specifications . . . . . 51


21.2 Electromagnetic compatibility declaration . . . . . 51

**22 Warranties . . . . . 53**

22.1 Equipment limited warranty for US customers . . . . . 53

22.2 Equipment limited warranty for customers outside the US . . . . . 54

### 1 Explanation of packaging, labeling, and device symbols

	Caution
	Fuse
	Equipotentiality
	Ground
	Power
	Stop
	Vacuum pump power
	Alternating current
	Humidity limitation
	Temperature limit
	Consult instructions for use
	Consult accompanying documents
	Do not resterilize
	Do not reuse
	Do not use if package is damaged
	Nonsterile
	Catalog number

	Lot number
	Serial number
	Manufacturer
	Date of manufacture
	Use-by date
	Sterilized using ethylene oxide
	Nonpyrogenic fluid path
	Quantity
	For US audiences only
	Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union Acts.
	Medical equipment with respect to electric shock, fire, and mechanical hazards only in accordance with ANSI/AAMI ES60601-1 (2012), CAN/CSA C22.2 No. 60601-1 (2014) E352357.
	Do not dispose of this product in the unsorted municipal waste stream. Dispose of this product according to local regulations. See <a href="http://recycling.Medtronic.com">http://recycling.Medtronic.com</a> for instructions on proper disposal of this product.
	This way up
	Keep dry
	Fragile, handle with care
	Manufactured in

	Corrugated recycles
	Protective earth ground
	Atmospheric pressure limitation
	China RoHS standard (SJ/T11364-2006)
	Contains di(2-ethylhexyl)phthalate (DEHP)
	Authorized representative in the European Community
	Heavy
	Do not push by the IV pole
	Storage temperature limit
	Transit temperature limit
	Storage humidity limitation
	Transit humidity limitation
	Storage atmospheric pressure limitation
	Transit atmospheric pressure limitation
	Type BF applied part

## 2 Introduction

Autotransfusion is a procedure in which the blood lost by, or removed from, a patient (autologous blood) is subsequently returned to the patient.

## 2.1 Advantages of autotransfusion over allogeneic transfusion

There are several benefits to autotransfusion:

- The risk of blood-transmitted diseases is eliminated.
- The related costs of transfusion are minimized.
- The risk of clerical errors is minimized.
- The use of autologous blood provides additional assurance during surgery on patients with multiple red blood cell antibodies or rare blood phenotypes.
- Valuable allogeneic blood is conserved.
- This process is generally accepted by patients with religious objections to allogeneic transfusion.

## 2.2 Principles of operation

Blood contains a suspension of formed elements of different densities. The autoLog IQ autotransfusion system (the device) operates by separating whole blood into its individual components by centrifugation.

**During operation, the following materials are recovered:**

- Washed, packed red blood cells
- A small volume of saline

**The following materials are less dense than red blood cells; they are removed during operation:**

- Lipids and fats
- Plasma-free hemoglobin
- Pharmacologic agents
- Irrigation solutions
- Anticoagulants
- Potassium
- Plasma
- Saline


As blood continues to enter the spinning bowl, the amassing red blood cell pack begins to occupy more of the bowl volume, and the excess plasma is pushed ahead of the red blood cells. When the total liquid volume of the bowl has been exceeded, the excess plasma exits the bowl through the effluent outlet to the waste bag via the connecting tubing.

To rid the red blood cell pack of contaminants, it is washed with isotonic saline (0.9% sodium chloride solution).

**Note:** Trace amounts of red blood cell contamination may cause the effluent to have a slightly reddish tint.

After washing, the clean, packed red blood cells are pumped from the bottom of the bowl to the holding bag, via the attached tubing. The blood should be transferred to a blood transfer bag and then to the patient.

**Figure 1.** Whole blood separation

1 Whole blood inlet  
2 Effluent outlet

3 Plasma  
4 Red blood cells

### 3 Intended use

The autoLog IQ autotransfusion system and Medtronic wash kit are intended for use in the collection, concentration, washing, and reinfusion of autologous blood. Areas of application may include, but are not limited to, the following:

- Surgeries including general, cardiovascular, orthopedic, vascular, plastic/reconstructive, obstetric/gynecologic, and neurological
- Postoperative treatment areas

### 4 Contraindications

The process of blood recovery is associated with few complications. The surgical team must consider the risks and relative contraindications of autotransfusion in any surgical procedure before proceeding with autotransfusion.

### 5 Warnings and precautions

- Read this user manual before using the device.
- Do not use the device for direct patient reinfusion (from the device directly into the patient). Adequate safeguards do not exist to protect the patient in this situation.
- Reinfusion of washed red blood cells can be carried out by gravity or pressure infusion after transferring processed blood to blood transfer bags. Do not directly reinfuse processed blood from the holding bag to the patient. Directly reinfusing the blood from the holding bag exposes the patient to the risk of air embolism. If a direct patient connection is required, follow the guidance for detecting and preventing air embolism in accordance with the current version of *AABB Standards for Perioperative Autologous Blood Collection and Administration*.
- Actual performance results may vary depending on many in-use variables. It is important to read and understand this user manual and understand the principles of cell washing before using the device for clinical operation. The responsibility for the use of this device in all cases belongs solely to the physician ordering its use.
- Failure to follow instructions may result in blood spillage, biohazard risk, slip hazard, low blood quality, infection, and hemolysis.
- Treat all blood and fluids in accordance with universal bloodborne pathogen precautions.
- *AABB Standards for Perioperative Autologous Blood Collection and Administration* recommends the expiration period for blood recovered intraoperatively, processed, and stored at room temperature, be no longer than 8 hours from the completion of processing. Intraoperative blood collected with processing can be stored for 24 hours between 1°C and 6°C (34°F and 43°F),

if the storage begins within 8 hours of completing processing. The transfusion of blood collected and processed under postoperative or posttraumatic conditions must begin within 8 hours of initiating the collection. Refer to the current version of *AABB Standards for Perioperative Autologous Blood Collection and Administration* for more information.

- Medtronic recommends the use of a blood transfusion filter, designed to retain particles that are potentially harmful to the patient, between the reinfusion container and the patient, in accordance with the current version of *AABB Standards for Perioperative Autologous Blood Collection and Administration*.
- Tandem collection reservoirs and leukocyte filtration can be used in relatively contraindicated cases; refer to the current version of *AABB Guidelines for Blood Recovery and Reinfusion in Surgery and Trauma* for more information.
- For orthopedic cases, always agitate the reservoir before processing the first cycle.
- Do not leave the device unmonitored, as it could present a biohazard risk if the waste or holding bag overflows, or could result in the loss of salvageable blood, among other concerns. It is the responsibility of the hospital to ensure that the individuals assigned to operate this device are well trained in the operation of the device and alert to potential problems.
- This device is intended for autotransfusion use in clinical patient care areas such as operating rooms, intensive care, and recovery rooms.
- Failure to use aseptic technique increases the risk of contaminating the disposables and harming the patient.
- Carefully monitor the patient's anticoagulation status in order to prevent complications. The basic concept of cell washing involves the removal of contaminated plasma and debris, while leaving red blood cells suspended in saline. Removal of large amounts of plasma during autotransfusion can cause patient hypovolemia. Since platelets and coagulation factors are contained in the plasma, this plasma removal may also reduce coagulation factors or platelet levels below normal levels. It is also possible that inadequate washing of salvaged blood may result in insufficient removal of anticoagulant or the development of coagulopathies upon the return of that blood to the patient.
- All autologous collected blood should be washed prior to reinfusion.
- Blood may be salvaged from body cavities, joint spaces, and other operative or trauma sites. If there is clinical evidence of sepsis, malignancy, or wound contamination, use tandem blood collection reservoirs and leukocyte filtration in accordance with the current version of *AABB Guidelines for Blood Recovery and Reinfusion in Surgery and Trauma*.
- Never transfuse blood that is suspected of having high hemolysis.
- Only use isotonic solutions that are labeled for intravenous use, such as 0.9% saline.
- To avoid hemolytic and subhemolytic stress and subsequent red blood cell lysis, vacuum pressure should be regulated at –80 mmHg to –120 mmHg, which is adequate for most surgical procedures. The vacuum setting may be temporarily increased to clear the field in the event of massive blood loss, then reduced to a lower setting; refer to the current version of *AABB Guidelines for Blood Recovery and Reinfusion in Surgery and Trauma* for more information.
- Applicable to products with phthalates identified on the product label: This product contains phthalates.
- Only use Medtronic wash kits and blood collection reservoirs with the device.
- The disposables are sterile and nonpyrogenic as long as package integrity has not been violated. Do not use the disposables if the package is damaged or open. Store all disposables in a dry place away from extreme temperatures.
- The disposable components utilized with this device are single use only. Do not reuse, reprocess, or resterilize the disposable components. Reuse, reprocessing, or resterilization may compromise the structural integrity of the disposables or create a risk of contamination of the disposables, which could result in patient injury, illness, or death. Only Medtronic sterilized disposable kits are approved for patient use with the device. Maintain a sterile field at the collection site. Use aseptic technique to minimize the possibility of contamination of the disposable components or the patient.
- The disposables must be used immediately after the removal of the protective packaging. Visually inspect the disposables. If any evidence of damage to the disposables is found during inspection or setup, do not use the disposable; return it to Medtronic for replacement. Do not use silicone, oils, or greases near the disposables.
- Dispose of the disposables in accordance with local hospital, administrative, and other government policies.
- Inspect the device and included items for shipping damage upon unboxing them. Confirm that all items are included in the box (Section 7.2). Contact your Medtronic service representative if damage is observed, or if an item is missing.
- Failure to fully secure connectors can lead to blood leakage.
- Close all clamps and use caution if transferring the disposables to a new device.
- The device must not be used in the presence of flammable anesthetic agents.
- Do not touch the bowl while it is spinning. Touching the moving bowl can lead to injury.
- Do not place items on the reservoir holder or collection reservoir.
- If using the rear reservoir holder, self-start will not be enabled, which could result in fluid overflowing the reservoir. Reservoir overflow could cause blood spillage, biohazard risk, and slip hazard.
- Do not use any solutions hotter than 42°C (108°F), because high heat can result in crenation and can destroy red blood cells.
- To avoid overheating the centrifuge, which could cause hemolysis and crenation, do not use the device at temperatures higher than 30°C (86°F).
- Do not use the device at temperatures lower than 15°C (59°F).
- Plastic materials used in the device and its disposable components may be sensitive to chemicals and certain detergents. Under certain adverse conditions, exposure to these chemicals (including vapors) can cause the plastics to fail or malfunction.
- Power loss or other failure during the **Fill** and **Wash** cycles may result in lower than normal blood quality parameters. Test the blood for hematocrit and contaminant removal to confirm suitability for patient infusion.
- Setting the vacuum at higher than desired settings, or unintentionally selecting **Max vacuum**, could result in hemolysis and tissue damage.
- Do not restrict the flow in any tubing line. If a tubing line is inadvertently clamped or kinked during operation, pressure may build up in the centrifuge bowl causing failure or leakage. Always examine the entire disposable kit to confirm that all tubing is free of

any kinks, twists, or flat areas. Examine the pump head and wash kit to confirm that all components are placed in the proper flow direction.

- The standard waste bag for the device holds approximately 10 L. Check the waste bag volume periodically and empty as required. The waste bag may be emptied at any time; however, 100 mL to 200 mL of fluid should be left in the bag to provide for proper expansion during filling and emptying. A full waste bag will cause back pressure, which may cause bowl leaks to occur.
- The device includes a centrifuge that rotates at 10,000 rpm. Parts that turn at high speeds may be dangerous. Follow safety rules related to the use of centrifuges. Do not open the centrifuge or remove the bowl before the centrifuge has completely stopped.
- Failure to respond to unusual behavior of the device could result in blood contamination or spillage.
- Choosing to wash a partially filled bowl may result in low blood quality and low red cell volume.
- Do not invert the bowl when examining for clots, as this could lead to fluid leakage.
- Do not service the device or update the device's software during use.
- Do not move the device by holding the IV pole, as this action could damage the device or cause tipping.
- Use caution when transporting the device by vehicle. Do not expose the device to extremes of temperature, shock, or vibration.
- Do not modify the device.
- Contact your Medtronic service representative, or a biomedical equipment technician trained and certified by Medtronic, to perform preventative maintenance and calibration of the device. Preventative maintenance and calibration should be performed on an annual basis. Blood spillage, biohazard risk, low blood quality, slip hazard, and other concerns, could result if preventative maintenance and calibration are not performed.
- Always turn off and unplug the device prior to working inside the cabinet or changing any fuses. Inside the cabinet are various electrical components and wiring. Physical contact with any of these components while the device is plugged in could result in severe electrical shock. For continued protection against risk of fire, replace fuses only with the same type and rating. Fuses should be replaced by a Medtronic service representative or Medtronic-trained and certified biomedical equipment technician. Internal grounding is provided for safety.
- Enclosure current leakage is a primary indicator of electrical shock hazard to personnel making contact with any exposed portion of the equipment. Each device is tested during the final quality inspection to verify enclosure current leakage is less than 100  $\mu$ A (under normal operating conditions) and less than 300  $\mu$ A (under single-fault conditions). Ensure that the current leakage is tested at least once a year, or more frequently as required, by the operating facility's biomedical engineering department or other qualified service technician. In addition, always test the enclosure current leakage and insulation after an event such as a fluid spill, after a major voltage surge in the power source has occurred, or after any device repair.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with a protective earth ground.
- The power cord may be used to disconnect the device from electrical power. Position the device and power cord so that the power cord can easily be accessed and disconnected.
- Although this device has been verified for EMC compliance, the potential exists that, in some situations, this device and other devices might electromagnetically interfere with each other. Take steps to minimize this possibility.
- The device should only be used in a hospital environment away from sources of high intensity EM disturbances, such as MRI devices.
- The emissions characteristics of this device make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this device might not offer adequate protection to radio frequency communication services. The user might need to take mitigation measures, such as relocating or reorienting the device.
- Do not stack or attach additional instrumentation to the device to minimize the potential for RF interference.
- Use of this device with accessories, transducers, and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this device could result.
- Power cord and USB cable length must not exceed 3 m (118 in).
- Do not touch a USB port and the patient at the same time.
- Use only a USB flash drive provided by Medtronic.
- Use appropriate personal protective equipment such as eye protection, mask, and gloves when cleaning the device and discarding used fluids. Failure to do so could result in exposure to bloodborne pathogens, blood, and cleaning fluids.
- While cleaning the centrifuge, failure to correctly install the tray for fluid collection will lead to spillage.
- While cleaning the centrifuge, failure to add sufficient volume of a 10% bleach solution or equivalent disinfectant solution for a sufficient duration can lead to ineffective cleaning.
- Do not use bleach to clean the touchscreen, as this could discolor the screen and result in failure.
- If it is suspected that fluid has entered the device, it should be immediately examined by a trained technician. Fluid ingress could result in electric shock to the user.
- Clean and maintain the device in accordance with this manual and hospital procedures.

## 6 Adverse effects

Proper or improper use of a cell salvage device can lead to the following risks to the patient and user:

Patient:

- Death
- Excessive bleeding

- Hemolysis
- Impaired liver function
- Loss of salvageable blood
- Infection
- Renal insufficiency and failure
- Low hematocrit
- Elevated serum lactate dehydrogenase levels
- Elevated total serum bilirubin levels
- Air emboli
- Febrile reaction
- Hypotension
- Increase in pulmonary and systemic vascular resistance
- Dyspnea
- Increased clotting
- Thromboembolism
- Tissue damage
- Arrhythmia
- Micro-shock

User:


- Exposure to biological fluids
- Biohazard
- Chemical exposure
- Burn injury
- Mechanical injury
- Electrical shock
- Slipping

## 7 Device description

### 7.1 Device components


The following figures detail the components of the device:

**Figure 2.** Device controls


- 1 Power: Turns the device power on and off
- 2 Power LED: Illuminates when the device power is on
- 3 Vacuum: Turns the device's internal vacuum on and off
- 4 Vacuum LED: Illuminates when the device's internal vacuum is on
- 5 Stop: Stops the blood processing sequence and returns blood to the reservoir
- 6 Alert light: Illuminates when a prompt is displayed on the touchscreen (for more information regarding visual alert signals, see Section 7.5)

**Figure 3.** Top panel

- |   |                            |
|---|----------------------------|
| 1 Centrifuge chamber (level sensor not shown) | 7 Pump lever               |
| 2 IV pole collar                              | 8 Pump outlet tubing guide |
| 3 IV pole knob                                | 9 Pump head                |
| 4 Centrifuge notches                          | 10 Pump tubing channel     |
| 5 Touchscreen                                 | 11 Manifold cover latch    |
| 6 Removable storage tray                      | 12 Manifold cover          |

**Figure 4.** Top panel with manifold cover open


- |                  |                                 |
|------------------|---------------------------------|
| 1 Air detector   | 4 Manifold cover locking sensor |
| 2 Tubing guides  | 5 Manifold sensor               |
| 3 Manifold valve | 6 Recessed cavity               |

**Figure 5.** Front panel


- |   |                                |
|---|--------------------------------|
| 1 Reservoir weight sensor (enables self-start mode) | 5 Speaker (under top panel)    |
| 2 Reservoir holder                                  | 6 Internal storage compartment |
| 3 IV pole   | 7 Device base                  |
| 4 Wash kit cradle arm                               |                                |

**Note:** The internal storage compartment can be used to store additional disposables and solutions as needed. Use caution when opening and closing the compartment, and while moving items to and from the compartment.

**Figure 6.** Back panel

- |   |  |
|---|--|
| 1 Manifold valve cavity fluid drain     | 7 Equipotential connector              |
| 2 Vacuum overflow trap quick-disconnect | 8 Power inlet module                   |
| 3 Vacuum overflow trap                  | 9 Power cord retainer clip             |
| 4 Wash kit storage brackets             | 10 Locking wheels (4)                  |
| 5 Power cord wrap                       | 11 Rear reservoir holder mounting base |
| 6 Wheel lock (4)                        | 12 Handle                              |

**Figure 7.** Handle positions

- 1 Upper storage position
- 2 Push position

- 3 Handle release button
- 4 Lower storage position

**Note:** Push the handle release button to move the handle into each position.


**Figure 8.** Left panel

- 1 Accessory USB ports (type A)
- 2 Barcode scanner holster (optional)


- 3 Barcode scanner cord retainer clip
- 4 Barcode scanner (optional)

**Figure 9.** Right panel

- |                            |                               |
|----------------------------|-------------------------------|
| 1 Waste bag pins           | 3 USB port (type B)           |
| 2 Wash kit storage holders | 4 Accessory USB port (type A) |

**Figure 10.** Bottom panel

- |                                 |                      |
|---------------------------------|----------------------|
| 1 Wash kit tray holder          | 3 Fan                |
| 2 Centrifuge drain (1/4-18 NPT) | 4 Vacuum pump outlet |

## 7.2 Items included with the device

The following items are included with the device:

- IV pole
- Reservoir holder
- Quick setup guide
- Quick setup video
- User manual
- Power cord
- USB flash drive

## 7.3 Optional items not included with the device

The following items are sold separately:

- Barcode scanner (used to simplify data entry; Chapter 10)
- Reservoir holder (attaches to the device rear reservoir holder mounting base; Figure 6)
- Medtronic cardiotomy reservoir holder E302 (attaches to an external IV pole)

**Warning:** Do not attach a reservoir to the device IV pole using a reservoir holder such as the Medtronic E302. The weight of a filled reservoir attached to the device IV pole may cause tipping.

## 7.4 Audible alert signals

The device emits an audible alert signal when a prompt appears on the touchscreen. Audible alert signals are divided into 3 categories: notify alert, low alert, and medium alert; all require a response. See Section 13.4, Chapter 18, and Section 19.4 for the full list of prompts and corresponding audible alert signals.

- Notify alert: 2 beeps
- Low alert: 2 beeps, repeating
- Medium alert: 3 beeps, repeating

To adjust and test the audible alert signal volume level, go to the **Settings** screen, and press **+** or **-** in the **Volume** field.


## 7.5 Visual alert signals

The alert light (Figure 2) on the device illuminates to indicate that there are prompts that may require a response. Visual alert signals are divided into 2 categories: low alert and medium alert. The visual alert signal categories align with the audible signal categories. See Section 13.4, Chapter 18, and Section 19.4 for the full list of prompts and corresponding visual alert signals.

- Low alert: alert light is illuminated, and not flashing
- Medium alert: alert light is flashing

## 8 In vitro performance

Figure 11. Typical hematocrit

1 Inlet hematocrit (Hct %)

2 Achievable hematocrit, mean (Hct %)

The data in Figure 11 are from in vitro testing using human blood and standard wash protocol. Clinical results may vary.

## 9 Transportation and storage

### 9.1 Transporting the device

**Warning:** Failure to adhere to the following loading conditions when transporting the device may lead to tipping:

- When using the bottom IV pole hooks, do not attach more than 5 L of liquid, and do not attach any bag higher than 48.3 cm (19 in) above the device.
- When using the top IV pole hooks, do not attach more than 3 L of liquid, and do not attach any bag higher than 86.4 cm (34 in) above the device.
- Do not attach a collection reservoir.
- Do not attach a waste bag.

**Caution:** Do not move the device by holding the IV pole, as this action could damage the device or cause tipping.

- Unplug the device power cord before transporting the device.
- Do not drop the device.
- Do not push, tip, or lean on the IV pole.
- If the device is tipped over or dropped, ensure that a trained service technician examines the device before use.
- Unlock the wheels if necessary.
- Only move the device while holding the handle. Rotate the handle to the push position (Figure 7).
- Use caution when transporting the device when fluids are present.
- See Chapter 21 for transit and storage environmental limits.

### 9.2 Transporting the device between clinical sites

- The device stores data after a password is set. Refer to institutional policies regarding the handling of protected health information and data security before transporting the device to another clinical site.
- Before moving the device, retract the IV pole and secure it by tightening the IV pole knobs.
- Lock the wheels before lifting the device into a vehicle.
- If it is necessary to lift the device, it should be done by a minimum of 2 people. Lift the device by holding it at the device base (Figure 5) or the power cord wrap (Figure 6). Do not lift the device by holding any other points. Use caution when lifting the device.
- The device may be transported without its original packaging, either standing upright or lying on its back.

**Warning:** Use caution when transporting the device by vehicle. Do not expose the device to extremes of temperature, shock, or vibration.

### 9.3 Storing the device

- Store the device according to the environmental conditions specified in Chapter 21.
- When not in use, turn off the device power and disconnect the power cord from the supply mains. Store the power cord by wrapping it around the power cord wrap on the back of the device (Figure 6).
- Store the device in a cool, dry, and well-ventilated area.

- Refer to institutional policies regarding the handling of protected health information and data security before storing the device.
- Store the USB flash drive in a secure location.

## 10 Setting up the optional barcode scanner

The barcode scanner (sold separately) can be used to enter the following information into the device:

- Patient and operator IDs
- Blood collection and processing disposables lot data
- Heparin lot data
- ACD-A lot data

To set up the barcode scanner, perform the following steps:


1. To install the holster, use the tool included in the barcode scanner kit to remove the screw from the left side of the device (Figure 8).
2. Attach the holster to the left side of the device, and then reinstall the screw using the tool.
3. Place the scanner in the holster, in an upright position.
4. Connect the scanner cord to a USB port (type A). See Figure 8.
5. Secure the cord using the barcode scanner cord retainer clip.
6. To scan data into the device, see Section 13.2.

## 11 Turning on the device power

### 11.1 Connecting the device power cord

1. Remove the power cord retainer clip (Figure 6) using the thumbscrew (Figure 12).
2. Insert the connector end of the power cord into the device power inlet module, and then insert the cord into the retainer clip.
3. Reinstall the power cord retainer clip using the thumbscrew.

**Figure 12.** Installing the power cord retainer clip and power cord

**Warning:** The power cord may be used to disconnect the device from electrical power. Position the device and power cord so that the power cord can easily be accessed and disconnected.


4. Insert the plug end of the power cord into a grounded, hospital-grade supply mains.

### 11.2 Turning on the device power

1. Press  $\text{\textcircled{P}}$  to turn on the device. The power LED illuminates, the alert light flashes, and the speaker emits 2 beeps.

**Note:** Press **VAC** to activate the vacuum at the default setting, during startup. The default vacuum setting can be changed using the **Settings** screen (Figure 14).

**Figure 13.** Turning on the power to the device

2. If turning on power for the first time, select the desired language.
3. Set the password if desired.


**Note:** A password must be set in order to store and access case data.

## 12 Managing device settings

### 12.1 Settings screen

Press **Settings** to access device setup, security, and user interface options.

**Figure 14.** Settings screen


- 1 Choose the device interface language
- 2 Set and update the password (required for storing and accessing case data)
- 3 Set the device name
- 4 Enable an audible click every time the screen is touched (a check mark appears when selected)
- 5 Enable a visual indicator every time the screen is touched (a check mark appears when selected)
- 6 Set the current time and the local time zone
- 7 Set the current date
- 8 Set the audible alert pattern
- 9 Adjust the audible alert volume
- 10 Adjust the touchscreen brightness
- 11 Adjust the default vacuum setting
- 12 Set the default maximum vacuum setting

**Note:** The factory default vacuum setting is -120 mmHg.

### 12.2 Setting the time

1. Select **Time** to display the time menu.
2. Select your **Time zone**.
3. Set the **Hours** and **Minutes**.
4. Select the desired **Time format**.
5. Select **Daylight savings** if applicable. A check mark appears when selected.

**Figure 15.** Setting the time screen

### 12.3 Setting the alert theme


1. Select **Alert theme** to display the alert theme menu.
2. Select **Theme One**, **Theme Two**, or **Theme Three**.
3. To sample each alert, press  next to the desired alert. Press  to stop the alert.
4. Adjust the audible signal volume as needed using the **Volume** field.

Figure 16. Setting the alert theme screen

### 12.4 Initial password setup

A password must be set in order to store and access case data.

The device stores data for up to 100 cases. For more information, see Chapter 14.

Refer to institutional policies regarding the handling of protected health information and data security.

1. Press **Set** in the **Password** field.
2. Correctly enter the password twice in the designated text fields, using the keyboard. Passwords must contain at least 6 characters.
3. Press .


Figure 17. Password entry keyboard screen

- |  |  |
|--|--|
| 1 Capitalization lock                                  | 5 Symbol, punctuation mark, and other character access       |
| 2 Dismiss keyboard                                     | 6 Enter; changes from gray to green when password is entered |
| 3 Number, punctuation mark, and other character access |  |
| 4 Delete   |  |

4. Press **Set**.


**Figure 18.** New password entry screen

**12.5 Updating or resetting the password**

1. Press **Update** in the **Password** field.
2. If the old password is known, enter it in the **Old password** text field.
3. Correctly enter the password twice in the designated text fields, using the keyboard. Passwords must contain at least 6 characters.
4. Press .
5. Press **Update**.


**Figure 19.** Updating or deleting the password screen

**12.6 Main screen overview**

The **Main** screen is divided into 4 areas: current case data, screen selection, user controls, and prompts.

**Figure 20.** Main screen

1 Current case data area	3 User controls area
2 Screen selection area	4 Prompts area

### Current case data area

- **Fill:** The volume of blood transferred from the collection reservoir to the centrifuge bowl in the current processing cycle
- **Wash:** The volume of saline used for the current processing cycle
- **Empty:** The volume of washed red blood cells transferred to the holding bag in the current processing cycle
- **Bowl (X):** The number of completed **Fill**, **Wash**, and **Empty** processing cycles
- **Total:** The total **Fill**, **Wash**, and **Empty** volumes for a case

### Screen selection area

- **Main:** Contains current case data, user controls, and prompts
- **Data:** Allows the user access to data stored in the device (Chapter 14)
- **Settings:** Allows the user to adjust the device settings (Chapter 12)
- **|||||:** Allows the user to enter patient ID, operator ID, and disposable lot, heparin lot, and ACD-A lot data (Section 13.2)
- **🔍:** Allows the user to access the **Service and maintenance** screen (Section 19.3)

**Note:** When **🔍** is highlighted in yellow, it indicates the presence of a notification. Click **🔍** to view the notification. See Section 19.4.

### User controls area

- **Vacuum:** Allows the user to set the vacuum setting. The vacuum setting can also be adjusted during a case. The default and maximum vacuum settings can be preset using the **Settings** screen (Section 12.1).
- **Wash volume:** The default wash volume is 250 mL. The wash volume can also be adjusted during a case; adjustments are effective immediately. If another volume is desired, select **Fast** or **Emergency** from the **Wash volume** menu on the **Main** screen. **Fast** uses 100 mL of wash solution. **Emergency** does not use any wash solution.

**Note:** Selecting **Emergency** causes a notification to appear every 10 min. If the notification is postponed, the notification will reappear at the interval selected in the **Reset alert** field. See Section 19.4.

- **🟢:** Allows the user to begin or resume automated blood processing, and to control the manual **FCT** functions. **🟢** changes to **🟡** during operation; press **🟡** to pause the applicable operation.

**Note:** If the centrifuge is operating when **🟡** is pressed, the paused state persists for a maximum of 60 s. After 60 s, the device automatically returns the bowl contents to the reservoir.

- **Function:** Allows the user to manually transfer volumes or remove air from the holding bag (Figure 21).
- **🔇:** Allows the user to mute and unmute the audible signals. When muted, the button icon changes, providing a visual indication of the selection.

**Note:** The device automatically reverts to the unmuted state after 60 s.

- **Date & Time:** Displays the local date and time. The flashing **:** indicates that the user interface is operating normally.

**Note:** The date and time can be adjusted on the **Settings** screen (Chapter 12).

### Prompts area

The Prompts area presents prompts and notifications. Prompts are messages that require user action. Notifications provide information about device conditions that are important but do not stop device operation. See Section 18.4 and Section 19.4 for prompts and notifications that may appear during operation.

## 12.7 Function capabilities


**Note:** **Function** capabilities are only enabled when the device is in Stop mode.

1. Press **🔴** while the pump is running, to empty the bowl to the reservoir. Wait for the bowl to empty completely, or press **🔴** again to immediately enable all **Function** capabilities.

- Press **FCT** to allow for the manual transfer of volumes or to remove air from the holding bag. Choose from the 3 options shown in Figure 21.

**Figure 21. FCT capabilities screen**

- FCT 1: Empty the bowl to the reservoir. Allows for the manual transfer of blood to the reservoir.
- FCT 2: Empty the bowl to the holding bag. Allows for the manual transfer of blood to the holding bag.
- FCT 3: Invert the holding bag. Allows for the manual removal of air from the holding bag.

- Press and hold to perform the selected **Function**.

**Warning:** Failure to press and hold when inverting the holding bag can result in air not being removed from the bag.

- Press at any time to return to the **Main** screen.

## 13 Operating the device

### 13.1 Supplies required

**Warning:** Only use Medtronic wash kits and blood collection reservoirs with the device.

**Note:** The Medtronic disposables are sterilized using ethylene oxide.

The following supplies are required:

- Medtronic wash kit that includes the following:
  - Centrifuge bowl (135 mL capacity) with tubing
  - Holding bag (1000 mL capacity)
  - Waste bag (10 L capacity)
- Medtronic blood collection reservoir (4 L capacity)
- Adequate quantity of 0.9% isotonic saline for washing
- Adequate quantity of anticoagulant solution
- Blood transfer bags
- Suction/anticoagulant line

### 13.2 Entering case data

- Press .

**Figure 22. Case data entry screen**

The screenshot shows the 'Case data entry' screen. At the top, there are three tabs: 'Main', 'Data', and 'Settings'. Below the tabs, there are several input fields for case data: 'Patient ID', 'Operator ID', 'Reservoir lot #', 'Washkit lot #', 'Suction line lot #', 'Notes', 'Heparin lot #', and 'ACD-A lot #'. The 'Notes' field is a larger text area. At the bottom right, there is a speaker icon and the date and time '20JAN16 12 24 PM'.

- Press the text field corresponding to the type of data being entered.
- To manually enter data, use the keyboard to enter the data into the selected text field.

**Note:** If using multiple disposables for a case, record the additional data in the **Notes** text field using the barcode scanner or keyboard.

4. Press  on the keyboard to enter the data in the device.
5. To scan data into the device using the optional barcode scanner (Chapter 10), press  on the screen.
6. Press the text field corresponding to the type of data being scanned.
7. To use the scanner, aim it at a barcode. Ensure that the red scanner line fully covers the barcode, and then press the trigger.


**Note:** The scanner can read the following barcode formats, at a minimum: aztec, code 128, code 39, data matrix ECC 200, and PDF 417.

**Note:** The scanner cannot read custom barcode formats.

8. After a successful scan, the scanner beeps, and the green LED illuminates. If unsuccessful, the red LED illuminates.
  9. Press  on the screen to record the scanned data in the device.
- Note:** For additional barcode scanner information, refer to the literature included with the scanner.
10. If desired, press **Main** to return to the **Main** screen.

### 13.3 Setting up the device for blood collection

**Figure 23.** Device components

- |  |   |
|--|---|
| 1 Anticoagulant solution               | 8 Step-down connector   |
| 2 Saline bag (1 or more)               | 9 Suction/anticoagulant line to patient   |
| 3 Holding bag                          | 10 Blood collection reservoir   |
| 4 Manifold (3-way tubing assembly)     | 11 Vacuum line to the vacuum overflow trap of the device or to external regulated vacuum source |
| 5 Line to holding bag                  | 12 Waste bag  |
| 6 Line from saline bag                 |   |
| 7 Line from blood collection reservoir |   |

**Warning:** Failure to adhere to the following loading conditions when using the device may lead to tipping:

- When using the bottom IV pole hooks, do not attach more than 6 L of fluid, and do not attach any bag higher than 48.3 cm (19 in) above the device.
- When using the top IV pole hooks, do not attach more than 3 L of fluid, and do not attach any bag higher than 86.4 cm (34 in) above the device.
- 2 collection reservoirs may be filled with no more than 4 L of fluid in each.
- 1 waste bag may be filled with no more than 7 L of fluid.

**Note:** If the device must be moved, different loading conditions apply (Section 9.1).

**Warning:** Failure to use aseptic technique increases the risk of contaminating the disposables and harming the patient.

**Warning:** Do not touch a USB port and the patient at the same time.

**Warning:** Failure to follow instructions may result in blood spillage, biohazard risk, slip hazard, low blood quality, infection, and hemolysis.

**Warning:** Do not use any solutions hotter than 42°C (108°F), because high heat can result in crenation and can destroy red blood cells.

**Warning:** Failure to fully secure connectors can lead to blood leakage.

**Warning:** Treat all blood and fluids in accordance with universal bloodborne pathogen precautions.


**Caution:** Do not touch the bowl while it is spinning. Touching the moving bowl can lead to injury.

**Note:** If using tandem or parallel collection reservoirs, see Chapter 17.

**Note:** It is not necessary to turn on the device power in order to set up the disposables.

1. Position the device on a level surface.
2. Examine the device for damage and cleanliness, and clean it if necessary (Chapter 16).
3. Lock the wheels.
4. Position the IV pole, and then tighten the IV pole knobs. Ensure that the maximum height indicator on the lower half of the IV pole is aligned with the top of the IV pole collar. See Figure 24. Do not extend the IV pole beyond the maximum height indicator.

**Figure 24.** IV pole maximum height indicator

5. Attach the reservoir holder to the reservoir weight sensor or rear reservoir holder mounting base.
 

**Warning:** If using the rear reservoir holder, self-start will not be enabled, which could result in fluid overflowing the reservoir. Reservoir overflow could cause blood spillage, biohazard risk, and slip hazard.

**Caution:** Do not place items on the reservoir holder or collection reservoir.
6. Remove the step-down connector from the collection reservoir.
7. Attach the step-down connector to the bottom of the collection reservoir. Ensure that the connector is fully secured.
8. Close the clamp on the step-down connector before the kit setup is complete. Ensure that the clamp is fully closed.
9. Place the collection reservoir into the reservoir holder.
10. Attach a vacuum line to the yellow-capped port on the collection reservoir.
11. Attach the other end of the vacuum line to the vacuum overflow trap of the device.
 

**Note:** The vacuum line also can be connected to an external regulated vacuum source, if necessary.
12. Aseptically pass the suction/anticoagulant line into the sterile field.
13. Attach the blue-capped suction/anticoagulant line connector to either the blue-capped or white-capped inlet port on the top of the collection reservoir.
14. Close the roller clamp on the suction/anticoagulant line. Spike the anticoagulant container. If the container is nonvented, open the vent cover on the drip chamber.
15. If it has not already been done, turn on the device power by pressing  $\downarrow$ .
 

**Note:** The vacuum can be enabled immediately after turning on the device power, and can be adjusted on the **Main** screen. If the device power is turned on within 6 hours of previous shutdown or power loss, the message **Resume previous case? If Yes, bowl contents will be pumped to reservoir. Make sure the reservoir is connected.** appears.

  - To continue the previous case, select **Yes**. The bowl contents will be pumped to reservoir. Make sure the reservoir is connected.
  - To start a new case, select **No**.

**Note:** The device defaults to a new case if the device power is turned on 6 or more hours after shutdown or power loss. See Section 13.6 for more information.
16. Turn on the vacuum by pressing **VAC**.

17. If necessary, adjust the vacuum setting using the **Vacuum** drop-down menu on the **Main** screen. The vacuum setting can be increased or decreased manually using **−** and **+**, or it can be set to the preprogrammed **Default** or **Max vacuum** settings. See Section 12.1 to set the **Default** and **Max vacuum** settings.


**Warning:** Setting the vacuum at higher than desired settings, or unintentionally selecting **Max vacuum**, could result in hemolysis and tissue damage.

**Note:** The factory default vacuum setting is  $-120$  mmHg.

**Note:** After selecting **Max vacuum**, the message **Confirm max vacuum?** appears. Select **Yes** to confirm or **No** to decline.

**Note:** Selecting **Max vacuum** causes a notification to appear every 10 min. If the notification is postponed, the notification will reappear at the interval selected in the **Reset alert** field.

**Figure 25.** Vacuum settings screen

18. Before blood collection, open the roller clamp and prime the collection reservoir with a minimum of 200 mL of anticoagulant solution.
19. Reduce the anticoagulant flow to a ratio of approximately 15 mL of solution to 100 mL of blood.
20. Monitor the collection reservoir periodically for appropriate anticoagulation.
21. Prime the collection reservoir with anticoagulant solution each time the collection reservoir is emptied.

### 13.4 Setting up the wash kit

See Table 1 for prompts associated with the setup of the disposables.

1. Attach the saline wash solution containers to a lower IV pole hook.
2. Rotate the blue wash kit cradle arm as shown in Figure 26.
3. Place the wash kit in the cradle arm and remove the wash kit cover.

**Figure 26.** Wash kit cradle arm

4. To enter the patient and operator IDs, and the disposable lot, heparin lot, and ACD-A lot data into the device, see Section 13.2.
5. Remove the waste bag from the wash kit.
6. Confirm that the drain valve is closed on the waste bag.
 

**Warning:** Failure to close the drain valve can lead to blood leakage.
7. Install the waste bag onto the waste bag pins on the right-facing side of the device. Ensure that the waste bag volume label faces away from the device.
8. Remove the holding bag from the wash kit and hang it on an upper IV pole hook. Close the clamps on the 2 red-capped outlet lines that each include a membrane port.
9. Place the bowl in the centrifuge chamber with the yellow-capped waste line outlet tubing facing the waste bag. Refer to the markings on the top panel.
10. Align the centrifuge markings with the vertical assemblies on the bowl.
11. To lock the bowl, simultaneously press down on the bowl and turn it clockwise until a click is heard.
12. Insert the yellow-capped connector into the yellow-capped inlet port on the waste bag.

13. Confirm that there are no kinks in the waste line.  
**Note:** If desired, the wash kit tray may be removed and the cradle arm rotated to its original position.
14. Open the manifold cover.  
**Note:** See the top panel components in Figure 3 and Figure 4.
15. Place the manifold into the recessed cavity.
16. Place the 3 manifold tubing lines into the tubing guides, and ensure that each tube is fully seated in each guide.
17. Place the single tubing line of the manifold into the air detector, and ensure that it is fully seated.
18. Close the manifold cover.  
**Note:** There is a feature on the manifold cover that pushes the tubing into the air detector when closing the cover.
19. Rotate the manifold cover latch until the lock indicators are aligned.
20. Pull the pump lever and hold it open.
21. Stretch the pump tubing around the pump head and over the pump outlet tubing guide. Release the pump lever.
22. Close the saline line clamps.
23. Spike the middle manifold tubing into the saline bag or bags.  
**Note:** If using only 1 saline bag, confirm that the clamp on the unused saline line is closed.
24. Open the saline clamp or clamps.
25. Connect the left-hand tube of the manifold to the blue-capped step-down connector at the bottom of the reservoir, and then open the clamp on the step-down connector.
26. Respond to prompts on the touchscreen.

**Table 1.** Setup of disposables prompts, audible alert signals, and visual alert signals

Prompt	Response	Audible alert signal	Visual alert signal
Kit is not installed. Lock the bowl cover.	The sensor has detected that the bowl cover is not locked.  <ol style="list-style-type: none"> <li>1. Place the bowl in the centrifuge chamber with the yellow-capped waste line outlet tubing facing the waste bag. Refer to the markings on the top panel.</li> <li>2. Align the centrifuge markings with the vertical assemblies on the bowl.</li> <li>3. To lock the bowl, simultaneously press down on the bowl and turn it clockwise until a click is heard.</li> </ol>	None	None
Insert the kit in the valve and the pump. Press  to resume. <sup>a</sup>	The kit is not properly installed in either the manifold valve or the pump.  <ol style="list-style-type: none"> <li>1. Open the manifold cover.</li> <li>2. Place the manifold into the recessed cavity.</li> <li>3. Place the 3 manifold tubing lines into the tubing guides, and ensure that each tube is fully seated in each guide.</li> <li>4. Place the single tubing line of the manifold into the air detector, and ensure that it is fully seated.</li> <li>5. Close the manifold cover. <b>Note:</b> There is a feature on the manifold cover that pushes the tubing into the air detector when closing the cover.</li> <li>6. Rotate the manifold cover latch until the lock indicators are aligned.</li> <li>7. Pull the pump lever and hold it open.</li> <li>8. Stretch the pump tubing around the pump head and over the pump outlet tubing guide.</li> <li>9. Release the pump lever.</li> </ol>	None	None
Kit is not installed. Lock the manifold cover. Press  to resume. <sup>a</sup>	A sensor has detected that the manifold cover is not locked.  <ol style="list-style-type: none"> <li>1. Open the manifold cover.</li> <li>2. Place the manifold into the recessed cavity.</li> <li>3. Place the 3 manifold tubing lines into the tubing guides, and ensure that each tube is fully seated in each guide.</li> <li>4. Place the single tubing line of the manifold into the air detector, and ensure that it is fully seated.</li> <li>5. Close the manifold cover.</li> </ol>	None	None

**Table 1.** Setup of disposables prompts, audible alert signals, and visual alert signals (continued)

Prompt	Response	Audible alert signal	Visual alert signal
	<p><b>Note:</b> There is a feature on the manifold cover that pushes the tubing into the air detector when closing the cover.</p> <p>6. Rotate the manifold cover latch until the lock indicators are aligned.</p>		
Connect the saline line and reservoir. Press  to resume.	<p>1. Spike the middle manifold tubing into the saline bag or bags.</p> <p><b>Note:</b> If using only 1 saline bag, confirm that the clamp on the unused line is closed.</p> <p>2. Open the saline clamp or clamps.</p> <p>3. Connect the left-hand tube of the manifold to the blue-capped step-down connector at the bottom of the reservoir, and then open the clamp on the step-down connector.</p> <p>4. Press .</p>	None	None
Air in saline line. Press  to resume.	<p>If the device detects air during the prime or wash phases, this message appears.</p> <p>1. Confirm that the saline bags are mounted and spiked. Mount and spike the saline bags if necessary.</p> <p>2. If the previous step does not resolve the message, reset the alert multiple times.</p> <p>3. If previous steps do not resolve the message, ensure that the tubing is correctly inserted into the air detector.</p> <p>4. If the previous steps do not resolve the message, clean the air detector and tubing with water, and then dry.</p> <p>5. If the previous steps do not resolve the message, use a different device if available. Contact your Medtronic service representative.</p> <p><b>Warning:</b> Close all clamps and use caution if transferring the disposables to a new device.</p>	Low	Low
Failed to move the valve. Check the kit installation. Press  to resume.	<p>1. Open the manifold cover.</p> <p>2. Place the manifold into the recessed cavity.</p> <p>3. Place the 3 manifold tubing lines into the tubing guides, and ensure that each tube is fully seated in each guide.</p> <p>4. Place the single tubing line of the manifold into the air detector, and ensure that it is fully seated.</p> <p><b>Note:</b> There is a feature on the manifold cover that pushes the tubing into the air detector when closing the cover.</p> <p>5. Close the manifold cover.</p> <p>6. Rotate the manifold cover latch until the lock indicators are aligned.</p> <p>7. Press .</p>	Medium	Medium

<sup>a</sup> Prompt will have different audible and visual alert signals when issued during operation.

## 13.5 Initiating blood processing

### 13.5.1 Blood processing

**Warning:** Do not apply pressure to the holding bag, as this can lead to air embolism.

**Note:** Respond to prompts that appear during operation. See Section 18.4 and Section 19.4 for a list of prompts that may appear during operation.

1. The message **Machine stopped. Press  to enable self-start.** appears.
2. Press  to activate the self-start mode. The message **Self-start enabled. Press  to manually start.** appears.
 


**Note:** The device automatically starts when it detects approximately 800 mL of blood in the reservoir for at least 5 s.
3. If it is desired to begin the fill before 800 mL is collected, press .
4. The default wash volume is 250 mL. The wash volume level can be adjusted during a case; adjustments are effective immediately. If a different wash volume is desired, select one from the **Wash volume** menu on the **Main** screen.

Figure 27. Wash volume screen

### 13.5.2 Pausing and stopping blood processing

If necessary to pause a **Fill**, **Wash**, or **Empty** phase, press **⏸**. During the **Fill** and **Wash** phases, the centrifuge continues to spin, and the pump pauses for 60 s. Press **▶** to resume. The device times-out and returns the bowl contents to the reservoir if **⏸** is not selected within 60 s.

If necessary to stop a **Fill**, **Wash**, or **Empty** phase, press **⏹**. This automatically returns the bowl contents to the reservoir.

**Note:** **Function** capabilities are only enabled when the device is in Stop mode.

### 13.5.3 Blood transfer

1. If necessary, empty the holding bag line by pressing FCT 2. The following prompt appears: **FCT 2: Empty the bowl to the holding bag. Press and hold **▶** to run the pump. Press **⏹** when complete.** (Section 12.7).

2. Transfer blood from the holding bag to a transfer bag, transfer blood in the holding bag, or transfer blood by direct patient connection.

**Warning:** If a direct patient connection to the processing device is required, take additional measures to detect and prevent air embolism (refer to the current version of *AABB Standards for Perioperative Autologous Blood Collection and Administration*).

**Note:** If using a transfer bag, remove all air from the transfer bag line by squeezing the transfer bag, before transferring it to the anesthesiologist.

**Note:** If using the holding bag to transfer blood, use FCT 3 to remove air from the bag, before transferring it to the anesthesiologist.

3. Place a new holding bag on the IV pole and connect it to the device if applicable.
4. Press **▶** to continue.

### 13.6 Resuming a case

If the device power has been turned off during a case, resume the case by restarting the device within 6 hours. To resume a case, perform the following steps:

1. Press **⏻** to turn on the device. The power LED illuminates, the alert light flashes, and the speaker emits 2 beeps.
2. The following message appears: **Resume previous case? If Yes, bowl contents will be pumped to reservoir. Make sure the reservoir is connected.** To resume the previous case, select **Yes**.

**Note:** The device defaults to a new case if the device power is turned on 6 or more hours after shutdown or power loss.

## 14 Managing case data

A password must be set in order to store and access case data.

The device stores data for up to 100 cases. A notice appears when 90 cases have been stored, to inform the user that the storage is nearly full. If additional case data is stored after the limit of 100 cases has been met, 1 old case will be deleted for each new case added.

Refer to institutional policies regarding the handling of protected health information and data security.

### 14.1 Viewing case data

1. To view current case data, press **Data**.

**Note:** To view all case data, press **Login**, enter the password, and press **✓**. If a password has not been set, see Section 12.4.

**Note:** After the password has been entered, cases appear chronologically, with the most recent case on top. A scroll bar appears if more than 6 cases are displayed.

2. To filter the data, see Section 14.2.
3. Select the desired cases by pressing anywhere in the applicable rows.

**Note:** Press **Select all** to see data for all stored cases. Press **Select none** to undo this.

4. Press **Details**.

**Figure 28.** Case list screen (password entered)

Case date	Patient ID	Bowls	Fill	Wash	Empty
20JAN16	123456789	3	1607	750	405
20JAN16	234567899	2	1002	500	270
19JAN16	172310894	1	500	250	135
19JAN16	523125584	2	1105	500	270
18JAN16	937585721	3	1475	750	405
18JAN16	321654987	1	472	250	135

5. Case date, bowls, and data entered by the user appear. Press **Details** to view bowl data for the case.

**Figure 29.** Case date, bowls, and data entered by user screen

Case date	Patient ID	Operator ID
20JAN16	123456789	223456789
Washkit lot #	Suction line lot #	Reservoir lot #
75881868	10233729	116909914
Heparin lot #	ACD-A lot #	Bowls
4456789	5567891	3
Notes		

6. The volumes for each bowl appear. Press **Start times** to display the **Fill**, **Wash**, and **Empty** start times of each bowl. Press **Volumes** to return to the bowl volume data screen.
7. Press **Details** to return to the case date, bowls, and data entered by user screen.
8. Press **Previous** to access the previous case record.
9. Press **Next** to access the next case record.
10. Press **Back** or **Data** to return to the case list screen.

**Figure 30.** Bowl volume data screen

Bowl	Time	Fill	Wash	Empty	Incomplete fill?
1	09:00 AM	520	250	90	No
2	11:00 AM	535	250	135	No
3	11:15 AM	552	250	100	No

Total volume mL:      1607      750      325

## 14.2 Searching and filtering case data

1. From the case list screen, press **Filter** to narrow results.

- To filter by date, press anywhere in **Oldest case** or **Newest case**.

**Figure 31.** Data filter screen

- Change the date as desired. Press  to close the menu.
- For all other text fields, enter text by pressing the text field and using the keyboard. Clear the text fields by pressing .
- Press **Back** to apply the filter.
- To clear the filter and return to the case list screen, press .

### 14.3 Exporting case data to a USB flash drive

#### 14.3.1 Exporting case data

**Caution:** Use only a USB flash drive provided by Medtronic.

- Plug a USB flash drive into one of the accessory USB ports (type A). See Figure 8 and Figure 9 for the location of the USB ports.
- From the case list screen, select the desired cases and press **Download to USB**.

**Note:** If the USB flash drive was not installed, or was improperly installed, the message **Download failed** appears. Properly install a USB drive and press **Download to USB**.

- After download, select **Yes** to delete the downloaded cases from the device. Select **No** to keep the downloaded cases stored in the device.
- Remove the USB flash drive.

**Figure 32.** Download is complete screen

Case date	Patient ID	Bowls	Fill	Wash	Empty
20JAN16	12345678				405
20JAN16	23456789				270
19JAN16	17231085				135
19JAN16	52312556				270
18JAN16	93758572				405
18JAN16	321654987	1	472	250	135

#### 14.3.2 Viewing exported case data

Exported case data is saved on the USB flash drive in a standard comma-separated values (CSV) text file, which can be opened and viewed in most spreadsheet software applications. To view exported case data, browse to the file location on the USB flash drive containing the file, and then open the file using a spreadsheet software application. The case data will appear in a tabular format (Figure 33), with the following headings:

- Print Date
- Device SN
- System name
- Case Date
- Patient ID
- Operator ID
- Reservoir Lot #
- Wash Kit Lot #

- Suction Line Lot #
- Heparin Lot #
- ACD-A Lot #
- Notes
- Bowl #
- Fill (mL)
- Fill start time
- Wash (mL)
- Wash start time
- Empty (mL)
- Empty start time
- Incomplete Fill?
- Total
- End of data

**Figure 33.** Case data file (example only)

Print Date	Device SN	System name						
13-Oct-16	ATLGIQ							
Case Date	Patient ID	Operator ID	Reservoir Lot #	Wash Kit Lot #	Suction Line Lot #	Heparin Lot #	ACD-A Lot #	Notes
20-Jan-16	123456789	223456789	116909914	75881868	10233729	4456789	5567891	Ortho case.
Bowl #	Fill (mL.)	Fill start time	Wash (mL.)	Wash start time	Empty (mL.)	Empty start time	Incomplete Fill?	
1	520	9:00	250	9:01	90	9:03	No	
2	535	11:00	250	11:01	135	11:03	No	
3	552	11:15	250	11:16	100	11:18	No	
Total	1607		750		325			
End of data								

#### 14.4 Exporting case data to an external device

Case data can be exported to an external device that has been configured with software that supports device case data exports. Contact your Medtronic service representative for support.

**Note:** When the device is connected to an external device, it is more susceptible to electrostatic discharge, which may result in the interruption of data export. Data export may be restarted if interrupted.

1. Connect a USB cable (type A male to type B male) from the type B USB port on the device (Figure 9) to a type A port on the external device.

**Warning:** USB cable length must not exceed 3 m (118 in).


**Note:** Ensure that any PC, external device, or external information technology equipment connected to the type B USB port on the device is certified equipment that complies with IEC 60950-1.

2. From the case list screen, select the desired cases and press .
3. The message **Selected cases received?** appears.
  - Select **Yes** if the data was successfully received.
  - Select **No** if the export was not received completely. The message **Reset the connection and try again.** appears. Select **OK.** The message **Download failed** appears. Select **OK.** To retry the export, check the USB cable connection, and then press .
4. After successful export, The message **Download is complete. Delete selected cases?** appears.
  - Select **Yes** to delete the exported cases from the device.
  - Select **No** to keep the exported cases stored in the device.
5. Remove the USB cable.

## 14.5 Deleting case data

1. To delete cases, select desired cases, and then press **Delete** on the case list screen.
2. Press **Yes** to accept the deletion. Press **No** to reject the deletion.

**Figure 34.** Delete selected cases screen

## 15 Ending a case

1. If necessary, turn off the vacuum by pressing **VAC**.
2. Document the case results.
3. Turn off the device by pressing **⏻** for 1.5 s. Unplug the power cord from the supply mains.
4. Cover all exposed ports and connectors of the disposables.
5. Close all tubing clamps.
6. Unlock and open the manifold cover.
7. Pull the pump lever and remove the tubing from the pump.
8. Remove the reservoir, centrifuge bowl, tubing harness, waste bag, vacuum line, and all solution containers from the device.
  - Warning:** Only clamp the waste line after the disposables have been removed.
9. Follow hospital procedures when disposing of biohazardous waste.
10. Clean the device in accordance with hospital protocols and as described in Chapter 16.

## 16 Cleaning the device

**Warning:** Use appropriate personal protective equipment such as eye protection, mask, and gloves when cleaning the device and discarding used fluids. Failure to do so could result in exposure to bloodborne pathogens, blood, and cleaning fluids.

**Warning:** If it is suspected that fluid has entered the device, it should be immediately examined by a trained technician. Fluid ingress could result in electric shock to the user.

**Warning:** Clean and maintain the device in accordance with this manual and hospital procedures.

Component/area	Cleaning
External	<ol style="list-style-type: none"> <li>1. Any time the external portions of the device become dirty (for example, blood spills), they should be cleaned according to approved hospital protocols with a 10% bleach solution or equivalent disinfectant solution.</li> <li>2. After cleaning, wipe the device with a cloth and water to remove any disinfectant solution residue. Do not use alcohol.</li> <li>3. Allow the device to dry.</li> </ol>
Level sensor	<p>If blood or other solutions are spilled in the centrifuge chamber, or dust builds up on the level sensor during use, the sensor can become inoperative.</p> <ol style="list-style-type: none"> <li>1. To clean the level sensor, use a foam-tipped swab dipped in a mild detergent and water.           <ul style="list-style-type: none"> <li><b>Caution:</b> Do not use liquids that may leave a film, such as alcohol. Do not use an abrasive cloth.</li> <li><b>Note:</b> To clean blood spills, use a 10% bleach solution or equivalent disinfectant solution.</li> </ul> </li> <li>2. Carefully wipe across the sensor to remove any blood, dust, or debris.</li> <li>3. Carefully wipe across the sensor with another foam-tipped swab to remove any remaining moisture.</li> </ol>
Air detector and valve assembly	<ol style="list-style-type: none"> <li>1. To clean the air detector or valve assembly, use a foam-tipped swab dipped in a mild detergent and water.           <ul style="list-style-type: none"> <li><b>Caution:</b> Do not use liquids that may leave a film, such as alcohol. Do not use an abrasive cloth.</li> </ul> </li> </ol>

Component/area	Cleaning
	<p><b>Note:</b> To clean blood spills, use a 10% bleach solution or equivalent disinfectant solution.</p> <ol style="list-style-type: none"> <li>Carefully wipe across the air detector or valve assembly to remove any blood, dust, or debris.</li> <li>Carefully wipe across the air detector or valve assembly with another foam-tipped swab to remove any remaining moisture.</li> </ol>
Centrifuge	<ol style="list-style-type: none"> <li>If blood or fluids spill in the centrifuge chamber, turn off the device.</li> <li>Remove the wash kit.</li> <li>To collect disinfectant solutions during the cleaning process, slide the empty tray from the wash kit into the wash kit tray holder under the device, from the right side. Tilt the tray slightly upward in order to insert the edges of the tray into the wash kit tray holder.</li> <li>Confirm that the wash kit tray is correctly installed. <ul style="list-style-type: none"> <li><b>Warning:</b> While cleaning the centrifuge, failure to correctly install the tray for fluid collection will lead to spillage.</li> <li><b>Note:</b> The centrifuge drain is equipped with female pipe threading (1/4-18 NPT). If desired, a customized fluid container can be used to collect fluids, instead of using the wash kit tray.</li> </ul> </li> <li>Turn on the device.</li> <li>Press  on the touchscreen to access the <b>Service and maintenance</b> screen.</li> <li>Press <b>Clean centrifuge</b>.</li> <li>Press  to start the centrifuge rotating. <ul style="list-style-type: none"> <li><b>Warning:</b> Failure to press  will result in the centrifuge not rotating and a failure to effectively clean the centrifuge.</li> </ul> </li> <li>Over a 1 min to 2 min period, slowly pour 100 mL of a 10% bleach solution or equivalent disinfectant solution into the centrifuge chamber while it is rotating until it starts to overflow. During rotation, fluid will spin out of the centrifuge and into the tray. Continue adding solution until the centrifuge is clean. <ul style="list-style-type: none"> <li><b>Warning:</b> While cleaning the centrifuge, failure to add sufficient volume of a 10% bleach solution or equivalent disinfectant solution for a sufficient duration can lead to ineffective cleaning.</li> </ul> </li> <li>Over a 1 min to 2 min period, rinse the centrifuge by slowly pouring 100 mL of water into the centrifuge chamber.</li> <li>Press . <ul style="list-style-type: none"> <li><b>Caution:</b> Failure to press  can lead to minor injury if the chamber is contacted while it is spinning.</li> </ul> </li> <li>Dispose of the used cleaning solutions in accordance with hospital procedures for the disposal of biohazardous waste.</li> <li>Dry the chamber with a soft cloth. Do not use alcohol to clean the centrifuge chamber. Remove and empty the tray.</li> </ol>
Pump head	<p>Clean the exterior and interior of the pump head after every use. To clean the pump head, perform the following procedure:</p> <ol style="list-style-type: none"> <li>Turn off the device and unplug the power cord from the supply mains.</li> <li>Pull back the pump lever with one hand. <ul style="list-style-type: none"> <li><b>Caution:</b> Be careful not to injure fingers in the pump lever.</li> </ul> </li> <li>With the other hand, lift the pump head up and off of the shaft.</li> <li>Release the pump lever.</li> <li>Clean the spilled fluids from the pump head as much as possible, using towels and foam-tipped swabs dampened with a 10% bleach solution or equivalent disinfectant solution. Dry the area thoroughly. <ul style="list-style-type: none"> <li><b>Caution:</b> Failure to thoroughly dry the area can result in corrosion.</li> </ul> </li> <li>Ensure that the pump head rollers move freely. Clean the pump head rollers with soap and water or with water alone. Dry the area thoroughly.</li> <li>To reinstall the pump head onto the shaft, pull back the pump lever with one hand.</li> <li>With the other hand, rotate the pump head while applying pressure until the keyed slot in the pump head aligns with the shaft and the pump head begins to slide onto the shaft. <ul style="list-style-type: none"> <li><b>Note:</b> The dimple on the pump head aligns with the shaft when the pump head slides onto the shaft.</li> </ul> </li> <li>Release the pump lever.</li> </ol>
Vacuum overflow trap	<ol style="list-style-type: none"> <li>After every use, examine the vacuum overflow trap for fluids or contamination.</li> <li>If fluids are present or if the vacuum overflow trap is dirty, remove, disassemble, and clean the vacuum overflow trap according to approved hospital protocols with a 10% bleach solution or equivalent disinfectant solution.</li> <li>Rinse with water and allow the components to dry prior to reassembly and reconnection to the device.</li> </ol>


Component/area	Cleaning
Storage tray	<ol style="list-style-type: none"> <li>1. Any time the storage tray becomes dirty (for example, blood spills), it should be removed and cleaned according to approved hospital protocols with a 10% bleach solution or equivalent disinfectant solution.</li> <li>2. After cleaning, wipe the tray with a cloth and water to remove any cleaning solution residue. Do not use alcohol.</li> <li>3. Allow the storage tray to dry before reinstalling.</li> </ol>
Internal storage compartment	<ol style="list-style-type: none"> <li>1. Whenever the internal storage compartment becomes dirty (for example, blood spills), it should be removed and cleaned.</li> <li>2. To remove the storage compartment, open it.</li> <li>3. Retract the spring plunger in the upper left corner on the inside of the compartment, and open the compartment completely.</li> <li>4. Lift the compartment off of the 2 pins in the device front panel.</li> <li>5. Clean the compartment and the area inside the device from which the compartment was removed according to approved hospital protocols using a 10% bleach solution or equivalent disinfectant solution.</li> <li>6. After cleaning, wipe with a cloth and water to remove any cleaning solution residue. Do not use alcohol.</li> <li>7. Allow the storage compartment and device to dry before reinstalling.</li> <li>8. To reinstall the compartment, align the sockets on the bottom of compartment with the 2 pivot pins in the device front panel.</li> <li>9. While closing the compartment, retract the spring plunger in the upper left corner of the compartment. Release the spring plunger once it has passed the compartment stop plate.</li> </ol>
Touchscreen	<p>Clean the touchscreen with a soft, lint-free cloth.</p> <p><b>Warning:</b> Do not use bleach to clean the touchscreen, as this could discolor the screen and result in failure.</p>
Barcode scanner	<p>For barcode scanner cleaning instructions, refer to the literature included with the scanner.</p>

## 17 Special blood processing conditions

**Note:** The following steps should be followed in conjunction with normal setup steps in Section 13.3.

### 17.1 Emergency setup using multiple suction/anticoagulant lines on 1 reservoir

The setup for this configuration is similar to that used for a normal case.


**Figure 35.** Emergency setup

- |  |  |
|--|--|
| 1 Anticoagulant solution containers      | 4 Step-down connector to device          |
| 2 Suction/anticoagulant lines to patient | 5 Vacuum line to regulated vacuum source |
| 3 Blood collection reservoir             |  |

1. Collect additional supplies necessary for this configuration:
  - Suction/anticoagulant line
  - Adequate quantity of anticoagulant solution
2. Attach the vacuum line to the yellow-capped port on the collection reservoir.
3. Attach the other end of the vacuum line to a regulated vacuum source.
4. Attach the blue-capped connector of the additional suction/anticoagulant line to one of the available inlet ports on the collection reservoir.
5. Close the roller clamp on the suction/anticoagulant line. If the anticoagulant container is nonvented, spike the anticoagulant container and open the vent cover on the drip chamber.
6. Spike the additional anticoagulant container for the second suction/anticoagulant line.

### 17.2 Using 2 reservoirs in series

**Note:** If necessary to transition to 2 collection reservoirs midcase, follow the steps in this section.

**Figure 36.** Using 2 collection reservoirs in series

- |   |  |
|---|--|
| 1 Anticoagulant solution                | 6 Step-down connector from secondary collection reservoir to primary collection reservoir inlet port |
| 2 Suction/anticoagulant line to patient | 7 Vacuum line to regulated vacuum source   |
| 3 Primary blood collection reservoir    | 8 Vacuum line to device  |
| 4 Step-down connector to device         |  |
| 5 Secondary blood collection reservoir  |  |

1. Collect additional supplies necessary for this configuration:

- Blood collection reservoirs
- Suction/anticoagulant line
- Adequate quantity of anticoagulant solution
- Vacuum source, if necessary
- Vacuum line, if necessary
- Medtronic cardiotomy reservoir holder E302, which attaches to an external IV pole

**Warning:** Do not attach a reservoir to the device IV pole using a reservoir holder such as the Medtronic E302. The weight of a filled reservoir attached to the device IV pole may cause tipping.


2. Attach a second reservoir holder to an external IV pole, above the primary collection reservoir.
3. Attach a step-down connector to the bottom of the secondary collection reservoir. Ensure that the connector is fully secured.
4. Place the secondary collection reservoir into the holder above the primary reservoir holder.
5. Connect the step-down connector from the secondary collection reservoir to the secondary collection reservoir inlet port on top of the primary collection reservoir, and close the clamp on the step-down connector.
6. If using 1 regulated vacuum source, connect tubing from a dual vacuum assembly to the yellow-capped ports on each collection reservoir. Connect the single tubing from the dual vacuum assembly to the vacuum source.
 

**Note:** The dual vacuum assembly is a "T" shaped tubing component, included with the tandem reservoir Y assembly, which connects 1 vacuum pump to the vacuum ports of the 2 reservoirs.
7. If using a secondary regulated vacuum source, attach the second vacuum line to the yellow-capped port on the secondary collection reservoir. Attach the other end of the second vacuum line to the secondary regulated vacuum source.
8. Attach the blue-capped connector of an additional suction/anticoagulant line to one of the available inlet ports of the secondary collection reservoir.
9. Close the roller clamp on the additional suction/anticoagulant line. If the anticoagulant container is nonvented, spike the anticoagulant container and open the vent cover on the drip chamber.
10. Spike an additional anticoagulant container for the second suction/anticoagulant line.

11. Before blood collection, open the roller clamp and prime the secondary collection reservoir with a minimum of 200 mL of anticoagulant solution.
12. Reduce the anticoagulant flow to each collection reservoir to a ratio of approximately 15 mL of solution to 100 mL of blood.
13. Monitor both collection reservoirs periodically for appropriate anticoagulation.
14. To move blood from the secondary collection reservoir to the primary collection reservoir, open the clamp on the connector to the primary collection reservoir, and momentarily clamp the suction/anticoagulant line on the secondary collection reservoir, using a hemostat.
15. To pause the movement of blood from the secondary collection reservoir to the primary collection reservoir, open the hemostat on the suction/anticoagulant line to the secondary collection reservoir, and close the clamp on the connector to the primary collection reservoir.
16. To maintain the desired vacuum setting during processing using 1 vacuum source, either clamp one of the suction/anticoagulant lines not in use, or increase the vacuum setting to maximum.

### 17.3 Using 2 reservoirs in parallel

Figure 37. Using 2 collection reservoirs in parallel

- |   |   |
|---|---|
| 1 Anticoagulant solution                | 5 Secondary blood collection reservoir                  |
| 2 Suction/anticoagulant line to patient | 6 Vacuum line to the vacuum overflow trap of the device |
| 3 Primary blood collection reservoir    | 7 Vacuum line to external regulated vacuum source       |
| 4 Tandem reservoir Y assembly           |   |

1. Collect additional supplies necessary for this configuration:
  - Blood collection reservoirs
  - Suction/anticoagulant line
  - Adequate quantity of anticoagulant solution
  - Tandem reservoir Y assembly
  - Vacuum source, if necessary
  - Vacuum line, if necessary
  - 1 reservoir holder: either the type that attaches to the device rear reservoir holder mounting base (Figure 6) or the Medtronic cardiotomy reservoir holder E302, which attaches to an external IV pole.
- Warning:** Do not attach a reservoir to the device IV pole using a reservoir holder such as the Medtronic E302. The weight of a filled reservoir attached to the device IV pole may cause tipping.
2. Attach a second reservoir holder to the device rear reservoir holder mounting base or to an external IV pole.
3. Place the secondary collection reservoir into the secondary reservoir holder.
4. Do not attach the step-down connector to either collection reservoir. Instead, use a tandem reservoir Y assembly to connect the primary and secondary collection reservoirs. Connect the single tubing to the secondary collection reservoir and the Y assembly to the primary collection reservoir. Ensure that the connectors are fully secured.

5. Close the clamp on the bottom of the secondary collection reservoir.
6. If using 1 regulated vacuum source, connect tubing from a dual vacuum assembly to the yellow-capped ports on each collection reservoir. Connect the remaining tubing from the dual vacuum assembly to the vacuum source.  
**Note:** The dual vacuum assembly is a "T" shaped tubing component, included with the tandem reservoir Y assembly, which connects 1 vacuum pump to the vacuum ports of the 2 reservoirs.
7. If using a secondary regulated vacuum source, attach the second vacuum line to the yellow-capped port on the secondary collection reservoir. Attach the other end of the second vacuum line to the secondary regulated vacuum source.
8. Attach the blue-capped connector of an additional suction/anticoagulant line to one of the available inlet ports of the secondary collection reservoir.
9. Close the roller clamp on the additional suction/anticoagulant line. If the anticoagulant container is nonvented, spike the anticoagulant container and open the vent cover on the drip chamber.
10. Spike an additional anticoagulant bag for the second suction/anticoagulant line.
11. Before blood collection, open the roller clamp and prime the secondary collection reservoir with a minimum of 200 mL of anticoagulant solution.
12. Reduce the anticoagulant flow to each collection reservoir to a ratio of approximately 15 mL of solution to 100 mL of blood.
13. Monitor both collection reservoirs periodically for appropriate anticoagulation.
14. To move blood from the secondary collection reservoir to the primary collection reservoir, open the clamp on the connector to the primary collection reservoir, and momentarily clamp the suction/anticoagulant line on the secondary collection reservoir, using a hemostat.
15. To pause the movement of blood from the secondary collection reservoir to the primary collection reservoir, open the hemostat on the suction/anticoagulant line to the secondary collection reservoir, and close the clamp on the connector to the primary collection reservoir.
16. To maintain the desired vacuum setting during processing using 1 vacuum source, either clamp one of the suction/anticoagulant lines not in use, or increase the vacuum setting to maximum.

#### 17.4 Draining cardiopulmonary bypass circuit (CPB) contents to the device

1. Connect the CPB circuit to an inlet port on a collection reservoir on the device. Ensure that the connector is fully secured.
2. Drain the contents of the CPB circuit to the collection reservoir.
3. Process the blood. See Section 13.5.

#### 17.5 Postoperative wound drainage

The following steps are only necessary if they were not completed previously for intraoperative blood processing.

The bowl and tubing setup and the operating instructions for this case are the same as those for a standard case (Section 13.4).

1. Open the sterile collection reservoir and remove the step-down connector.
2. Attach the step-down connector to the bottom of the collection reservoir. Ensure that the connector is fully secured.
3. Place the collection reservoir into a reservoir holder. Attach a vacuum line to the yellow-capped port on the collection reservoir.
4. Attach the other end of the vacuum line to a regulated vacuum source. Use an appropriate, and preferably intermittent, vacuum source adjusted to -80 mmHg or less.
5. Aseptically pass the suction/anticoagulant line into the sterile field.
6. Attach the blue-capped connector of the suction/anticoagulant line to either the blue-capped or white-capped inlet port on the top of the collection reservoir.
7. Close the roller clamp on the suction/anticoagulant line. Spike the anticoagulant container. If the anticoagulant container is nonvented, open the vent cover on the drip chamber.
8. Before blood collection, open the roller clamp and prime the collection reservoir with a minimum of 200 mL of anticoagulant solution.
9. Reduce the anticoagulant flow to a ratio of approximately 15 mL of solution to 100 mL of blood.
10. Connect the universal Y connector to the suction/anticoagulant line. When the device is turned off, stop the anticoagulant flow. Start the anticoagulant flow when the vacuum is turned on.
11. Begin blood collection. After adequate blood volume has been collected, it can be processed by the device.
12. If simultaneous blood collection is needed for 2 drain lines, a second collection reservoir, suction/anticoagulant line with anticoagulant container, and vacuum source may be used (with reservoirs in either series or parallel).

### 18 Troubleshooting guide

**Note:** This guide is intended to provide a simple description of common problems that can occur during use of the device, and of the messages that appear during its operation. A brief description of the cause of the problem or message is also included.

#### 18.1 Power failure

Normally, each time the device power is turned on, parameters are initialized to prepare it for use with a new patient. However, after a power failure or interruption (and if the device is restarted within 6 hours), the following message appears: **Resume previous case? If Yes, bowl contents will be pumped to reservoir. Make sure the reservoir is connected.**

1. To continue the previous case, select **Yes**.
2. To start a new case, select **No**.

**Note:** The device defaults to a new case if the device power is turned on 6 or more hours after shutdown or power loss.

## 18.2 Unusual behavior of the device

**Warning:** Failure to respond to unusual behavior of the device could result in blood contamination or spillage.

Monitor the device during operation, and if unusual behavior is observed, respond as described in the following table.

Unusual behavior	Response
The centrifuge bowl fills with blood and then blood returns to the collection reservoir, repeatedly.	<ol style="list-style-type: none"> <li>1. Press  to place the device in Stop mode.</li> <li>2. Determine if the input volume is diluted.</li> <li>3. If the input volume is diluted, transition to manual start for blood processing.</li> <li>4. Examine the inner and outer surfaces of the reservoir filter for clots.</li> <li>5. If clots are detected on the outside of the reservoir filter, replace the reservoir.</li> <li>6. If clots are detected on the inside of the reservoir filter, prime the reservoir with at least 200 mL of anticoagulant solution and agitate the reservoir.</li> <li>7. Restart blood processing.</li> <li>8. Frequently agitate the reservoir, increasing anticoagulation, for the remaining duration of the case.</li> </ol>
The bowl is vibrating or making excessive noise.	<ol style="list-style-type: none"> <li>1. Press  to place the device in Stop mode.</li> <li>2. Remove the centrifuge bowl.</li> <li>3. Examine the bowl for damage, clots, or blood contaminated with particulates.</li> <li>4. Discard and replace the bowl if it is damaged or contaminated with particulates.</li> <li>5. Reinstall the bowl.</li> <li>6. Examine the inner and outer surfaces of the reservoir filter for clots.</li> <li>7. If clots are detected in the reservoir filter, replace the reservoir, centrifuge bowl, and tubing harness.</li> <li>8. If no clots are detected in the reservoir filter, replace the centrifuge bowl and tubing harness.</li> <li>9. Prime the new disposables with at least 200 mL of anticoagulant solution.</li> <li>10. Restart blood processing.</li> <li>11. Frequently agitate the reservoir, increasing anticoagulation, for the remaining duration of the case.</li> <li>12. If the previous steps do not resolve the message, use another device if available. Contact your Medtronic service representative.</li> </ol> <p><b>Warning:</b> Close all clamps and use caution if transferring the disposables to a new device.</p>
The reservoir is making a cracking noise as a result of negative pressure.	<ol style="list-style-type: none"> <li>1. Turn off the vacuum or reduce the vacuum setting.</li> <li>2. Examine the suction/anticoagulant line for obstructions.</li> <li>3. Vent the reservoir.</li> </ol>
The barcode scanner does not work.	<ol style="list-style-type: none"> <li>1. Plug the barcode scanner cord into a different USB port (type A).</li> <li>2. Confirm that the barcode type can be read by the scanner. See Section 13.2.</li> <li>3. Contact your Medtronic service representative.</li> </ol>

## 18.3 Normal startup prompts, audible and visual alert signals



Normal prompts may appear during device operation or startup, and they typically can be resolved by making a quick adjustment of the device or the disposable kit.

Prompt	Response	Audible alert signal	Visual alert signal
Resume previous case? If Yes, bowl contents will be pumped to reservoir. Make sure the reservoir is connected.	<ol style="list-style-type: none"> <li>1. To continue the previous case, select <b>Yes</b>.</li> <li>2. To start a new case, select <b>No</b>.</li> </ol> <p><b>Note:</b> If <b>No</b> is selected in error, case data can be retrieved by following steps in Section 14.1.</p> <p><b>Note:</b> The device defaults to a new case if the device power is turned on 6 or more hours after shutdown or power loss.</p>	None	None
Connect the saline line and reservoir. Press  to resume.	<ol style="list-style-type: none"> <li>1. Spike the middle manifold tubing into the saline bag or bags.</li> </ol> <p><b>Note:</b> If using only 1 saline bag, confirm that the clamp on the unused line is closed.</p> <ol style="list-style-type: none"> <li>2. Open the saline clamp or clamps.</li> <li>3. Connect the left-hand tube of the manifold to the blue-capped step-down connector at the bottom of the reservoir.</li> <li>4. Press .</li> </ol>	None	None

Prompt	Response	Audible alert signal	Visual alert signal
Machine stopped. Press  to enable self-start.	This message appears after the disposable kit is installed. Enable self-start by pressing  . The device automatically starts when it detects approximately 800 mL of blood in the reservoir for at least 5 s.	None	None
Self-start enabled. Press  to manually start.	Press  to cause a manual start.	None	None
Vacuum pressure detected at startup. Make sure the suction line is not clamped. Press  to resume.	This message appears if existing vacuum pressure is detected in the reservoir when vacuum is turned on at device startup. <ol style="list-style-type: none"> <li>1. Confirm that vacuum and suction lines are not clamped or immersed in liquid.</li> <li>2. To continue, press .</li> <li>3. Press <b>VAC</b> to restore vacuum.</li> </ol> <p><b>Note:</b> If prompt returns, the device can still be used if an external regulated vacuum source is connected to the reservoir. Contact your Medtronic service representative.</p>	Medium	Medium

#### 18.4 Normal operating prompts, audible and visual alert signals

Normal prompts may appear during device operation or startup, and they typically can be resolved by making a quick adjustment of the device or the disposable kit.

Prompt	Response	Audible alert signal	Visual alert signal
Air in reservoir line. Press  to try again or select Return, Concentrate, or Wash.	The <b>Fill</b> phase emptied all the contents of the reservoir. <ol style="list-style-type: none"> <li>1. If possible, remove air from the reservoir line. Press  to continue filling.</li> <li>2. Depending on the clinical use condition, select <b>Return</b>, <b>Concentrate</b>, or <b>Wash</b>.</li> </ol> <p><b>Warning:</b> Choosing to wash a partially filled bowl may result in low blood quality and low red cell volume.</p>	Low	Low
Air in holding bag line. Bowl is not full. Press  to try again or select Return or Wash.	This message appears after selecting <b>Concentrate</b> , and after the red blood cells from the holding bag are consumed before the bowl is full. <ol style="list-style-type: none"> <li>1. Unclamp the holding bag line if necessary.</li> <li>2. Press and hold  to advance blood to the bowl.</li> <li>3. If the previous steps do not resolve the message, select <b>Return</b> to return the blood to the reservoir or <b>Wash</b> to wash the contents of the bowl.</li> <li>4. The centrifuge continues to spin for 1 min, and blood is returned to the reservoir if no selection is made.</li> </ol> <p><b>Warning:</b> Choosing to wash a partially filled bowl may result in low blood quality and low red cell volume.</p>	Low	Low
Kit is not installed. Lock the bowl cover.	A sensor has detected that the bowl cover is not locked. <ol style="list-style-type: none"> <li>1. Place the bowl in the centrifuge chamber with the yellow-capped waste line outlet tubing facing the waste bag. Refer to the markings on the top panel.</li> <li>2. Align the centrifuge markings with the vertical assemblies on the bowl.</li> <li>3. To lock the bowl, simultaneously press down on it and turn it clockwise until a click is heard.</li> </ol>	None	None
Kit is not installed. Lock the manifold cover. Press  to resume. <sup>a</sup>	A sensor has detected that the manifold cover is not locked. <ol style="list-style-type: none"> <li>1. Open the manifold cover.</li> <li>2. Place the manifold into the recessed cavity.</li> </ol>	Medium	Medium

Prompt	Response	Audible alert signal	Visual alert signal
	<ol style="list-style-type: none"> <li>3. Place the 3 manifold tubing lines into the tubing guides, and ensure that each tube is fully seated in each guide.</li> <li>4. Place the single tubing line of the manifold into the air detector, and ensure that it is fully seated.</li> </ol> <p><b>Note:</b> There is a feature on the manifold cover that pushes the tubing into the air detector when closing the cover.</p> <ol style="list-style-type: none"> <li>5. Rotate the manifold cover latch until the lock indicators are aligned.</li> </ol>		
Insert the kit in the valve and the pump. Press  to resume. <sup>a</sup>	<p>The kit is not properly installed in either the manifold valve or the pump.</p> <ol style="list-style-type: none"> <li>1. Open the manifold cover.</li> <li>2. Place the manifold into the recessed cavity until the connector is fully seated in the cavity.</li> <li>3. Place the 3 manifold tubing lines into the tubing guides, and ensure that each tube is fully seated in each guide.</li> <li>4. Place the single tubing line of the manifold into the air detector, and ensure that it is fully seated.</li> </ol> <p><b>Note:</b> There is a feature on the manifold cover that pushes the tubing into the air detector when closing the cover.</p> <ol style="list-style-type: none"> <li>5. Rotate the manifold cover latch until the lock indicators are aligned.</li> <li>6. Pull the pump lever and hold it open.</li> <li>7. Stretch the pump tubing around the pump head and over the pump outlet tubing guide.</li> <li>8. Release the pump lever.</li> </ol>	Medium	Medium
Bowl may not be full. Wash it?	<p>This message occurs when selecting <b>Wash</b>, after the prompts <b>Air in reservoir line. Press  to try again or select Return, Concentrate, or Wash. or Air in holding bag line. Bowl is not full. Press  to try again or select Return or Wash.</b></p> <ol style="list-style-type: none"> <li>1. Select <b>Yes</b> to wash the bowl.</li> </ol> <p><b>Warning:</b> Choosing to wash a partially filled bowl may result in low blood quality and low red cell volume.</p> <ol style="list-style-type: none"> <li>2. Selecting <b>No</b> stops the bowl and returns the blood to the reservoir.</li> </ol>	Low	Low
Waste bag may be full. Empty if necessary. Press  to resume.	<p>The device estimates the volume in the waste bag and a notice appears when the bag is full.</p> <ol style="list-style-type: none"> <li>1. Empty the waste bag.</li> </ol> <p><b>Note:</b> Leave 100 mL to 200 mL of fluid in the bag to provide for proper expansion during filling and emptying.</p> <ol style="list-style-type: none"> <li>2. Press  to resume.</li> </ol>	Medium	Medium
Unlock then relock the bowl cover. Press  to resume.	<p>The bowl moved during rotation and is unlocked.</p> <ol style="list-style-type: none"> <li>1. Unlock and remove the bowl.</li> <li>2. Examine the centrifuge chamber for debris and liquid.</li> <li>3. Remove the debris or liquid if necessary.</li> <li>4. Examine the bowl for damage.</li> <li>5. Discard and replace the bowl if it is damaged.</li> <li>6. If clots are observed in the bowl, discard and replace the wash kit.</li> <li>7. Confirm that the bowl sensor and sensor spring are moving freely and not damaged or dirty.</li> <li>8. Reinstall the bowl.</li> </ol>	Medium	Medium

Prompt	Response	Audible alert signal	Visual alert signal
	<p>9. If the message is not resolved after the previous steps, and if the bowl was not replaced, then replace the bowl.</p> <p>10. If the message is still not resolved, use another device if available. Contact your Medtronic service representative.</p> <p><b>Warning:</b> Close all clamps and use caution if transferring the disposables to a new device.</p>		
Optics obstructed. Empty the bowl or clean the optics sensor. Press  to resume.	<p>The level sensor detects blood at the top of the bowl when there should be no blood in the bowl.</p> <ol style="list-style-type: none"> <li>1. Examine the bowl for blood.</li> <li>2. If there is blood in the bowl, manually return the blood to the reservoir using the <b>FCT</b> key (Figure 21).</li> <li>3. If the previous step does not resolve the message, unlock the bowl.</li> <li>4. If clots are observed in the bowl, discard and replace the wash kit.</li> <li>5. If clots are observed inside the reservoir filter, prime the reservoir with anticoagulant.</li> <li>6. Agitate the reservoir vigorously.</li> <li>7. If clots are observed outside the reservoir filter, replace the reservoir.</li> <li>8. If the previous steps do not resolve the message, document all fluid volumes.</li> <li>9. Clean and dry the optic sensors using an approved cleaning solution (Chapter 16).</li> <li>10. If the previous steps do not resolve the message, use another device if available. Contact your Medtronic service representative.</li> </ol> <p><b>Warning:</b> Close all clamps and use caution if transferring the disposables to a new device.</p>	Low	Low
Air in saline line. Press  to resume.	<p>This message appears if the device detects air during the prime or wash phases.</p> <ol style="list-style-type: none"> <li>1. Confirm that the saline bags are mounted and spiked. Mount and spike the saline bags if necessary.</li> <li>2. If the previous step does not resolve the message, reset the alert as needed.</li> <li>3. If the previous steps do not resolve the message, ensure that the tubing is correctly inserted into the air detector.</li> <li>4. Clean the air detector and tubing with water, and then dry.</li> <li>5. If the previous steps do not resolve the message, use another device if available. Contact your Medtronic service representative.</li> </ol> <p><b>Warning:</b> Close all clamps and use caution if transferring the disposables to a new device.</p>	Low	Low
Emptying time was too short. Check the holding bag line for air. Press  to resume.	<p>Examine the disposables for air leaks.</p> <ol style="list-style-type: none"> <li>1. Confirm that there is flow to the holding bag.</li> <li>2. If no flow is observed, remove the centrifuge bowl and examine it for clots.</li> </ol> <p><b>Warning:</b> Do not invert the bowl when examining for clots, as this could lead to fluid leakage.</p> <ol style="list-style-type: none"> <li>3. If clots are observed in the bowl, discard and replace the wash kit.</li> <li>4. If no clots are observed, clean the air detector and tubing with water, and then dry.</li> </ol>	Low	Low
Emptying time was too long. Check the holding bag line for occlusions. Press  to resume.	<p>This message appears if the holding bag line is clamped and the device is unable to move blood volume from the bowl, or if the air detector has failed.</p>	Low	Low

Prompt	Response	Audible alert signal	Visual alert signal
	<ol style="list-style-type: none"> <li>1. Examine the holding bag line for closed clamps and other occlusions.</li> <li>2. Examine the holding bag. Empty or replace the holding bag if it is full.</li> <li>3. Clean the air detector and tubing with water, and then dry.</li> <li>4. If the message is not resolved, use another device if available. Contact your Medtronic service representative.</li> </ol> <p><b>Warning:</b> Close all clamps and use caution if transferring the disposables to a new device.</p>		
Holding bag may be full. Empty or change if necessary. Press  to resume.	<p>The device estimates the volume pumped to the holding bag and a notice appears when the bag is full.</p> <ol style="list-style-type: none"> <li>1. Empty or replace the holding bag.</li> <li>2. Press .</li> </ol>	Medium	Medium
Machine paused. Press  to resume.	<p>The pump has been manually paused using .</p> <ol style="list-style-type: none"> <li>1. Resolve the issue.</li> <li>2. Press  to resume.</li> </ol>	None	Low
Machine stopped. Press  to enable self-start.	<p>This message appears after the disposable kit is installed. Enable self-start by pressing . The device automatically starts when it detects approximately 800 mL of blood in the reservoir for at least 5 s.</p>	None	None
Self-start enabled. Press  to manually start.	<p>Press  to cause a manual start.</p>	None	None
The pump is running too slow. Check the kit installation. Press  to resume.	<p>Tubing is binding on the pump.</p> <ol style="list-style-type: none"> <li>1. Confirm that the pump is moving. If it is moving, proceed to the next step. If not, contact your Medtronic service representative.</li> <li>2. Ensure that the tubing is correctly inserted in the pump tubing channel.</li> <li>3. If the previous steps do not resolve the message, remove the tubing from the pump tubing channel.</li> <li>4. Remove the pump head.</li> <li>5. Confirm that the pump head rollers move freely, and are free of burrs and other debris.</li> <li>6. Clean the pump head and rollers if necessary using an approved cleaning solution (Chapter 16).</li> <li>7. Reinstall the pump head.</li> <li>8. If the previous steps do not resolve the message, use another device if available. Contact your Medtronic service representative.</li> </ol> <p><b>Warning:</b> Close all clamps and use caution if transferring the disposables to a new device.</p>	Medium	Medium
Centrifuge is running too slow. Check the bowl installation. Press  to resume.	<p>The bowl may be binding in the centrifuge.</p> <ol style="list-style-type: none"> <li>1. Unlock then relock the centrifuge bowl to ensure that it is correctly installed.</li> <li>2. Remove the centrifuge bowl.</li> <li>3. Examine the bowl for damage.</li> <li>4. Discard and replace the bowl if it is damaged.</li> <li>5. Examine the centrifuge chamber for debris and liquid.</li> <li>6. Examine the centrifuge chamber for binding by moving the chamber by hand. If it is binding, contact your Medtronic service representative.</li> <li>7. Remove the debris or liquid if necessary.</li> <li>8. Reinstall the bowl and confirm normal operation.</li> <li>9. If the previous steps do not resolve the message, document all fluid volumes.</li> </ol>	Medium	Medium

Prompt	Response	Audible alert signal	Visual alert signal
	10. Turn off the device, then turn it on. 11. If the previous step does not resolve the message, use another device if available. Contact your Medtronic service representative. <b>Warning:</b> Close all clamps and use caution if transferring the disposables to a new device.		
Reservoir weight not detected. Self-start not available. Press  to manually start.	The weight sensor does not detect the reservoir. 1. Mount the reservoir in the reservoir holder (enabling self-start) if self-start mode is desired. 2. Or, press  to manually start.	Notify	None
FCT 1: Empty the bowl to the reservoir. Press and hold  to run the pump. Press  when complete.	FCT 1 was selected while the device was in Stop mode. These functions enable the pump to be manually engaged.	None	None
FCT 2: Empty the bowl to the holding bag. Press and hold  to run the pump. Press  when complete.	FCT 2 was selected while the device was in Stop mode. These functions enable the pump to be manually engaged.	None	None
FCT 3: Invert the holding bag. Press and hold  to remove air. Press  when complete.	FCT 3 was selected while the device was in Stop mode. These functions enable the pump to be manually engaged.	None	None
Vacuum pressure exceeded the setting. Make sure the suction line is not clamped. Press  to resume.	This message appears if the device vacuum system could not regulate the vacuum pressure, possibly because of a clamp on the vacuum line between the fluid trap and the reservoir. 1. Confirm that vacuum and suction lines are not clamped or immersed in liquid. 2. To continue, press  . 3. Press <b>VAC</b> to restore vacuum. 4. If the previous steps do not resolve the issue, turn off the power, then turn it on, and retry.  <b>Note:</b> If prompt returns, the device can still be used if an external regulated vacuum source is connected to the reservoir. Contact your Medtronic service representative.	Medium	Medium

<sup>a</sup> Prompt will have different audible and visual alert signals when issued during setup.

### 18.5 Error prompts, audible and visual alert signals

Error prompts may appear during device operation or startup.

Error	Response	Audible alert signal	Visual alert signal
The vacuum system failed. Vacuum is not available. Press  to resume.	This message appears if the device vacuum system has failed. 1. Confirm that vacuum and suction lines are not clamped or immersed in liquid. 2. Turn off the power, then turn it on, and retry. 3. If the previous steps do not resolve the issue, contact your Medtronic service representative. <b>Note:</b> The device can still be used if an external regulated vacuum source is connected to the reservoir.	Medium	Medium

### 18.6 Irrecoverable errors, audible and visual alert signals

These errors cannot be recovered by the device or user. The device should be taken out of service until it is examined by a Medtronic service representative or biomedical equipment technician trained and certified by Medtronic.

Error	Response	Audible alert signal	Visual alert signal
Failure of any of the following components will produce a message containing the component name, the error code, and a statement to contact Medtronic for service: <ul style="list-style-type: none"> <li>• Air detector</li> <li>• Centrifuge</li> </ul>	Contact your Medtronic service representative, and provide the error code displayed on the screen.  <b>Warning:</b> Close all clamps and use caution if transferring the disposables to a new device.	Medium	Medium

Error	Response	Audible alert signal	Visual alert signal
<ul style="list-style-type: none"> <li>• Diagnostic hardware</li> <li>• External A/D convertor</li> <li>• Level sensor</li> <li>• Weight sensor</li> <li>• Device microprocessor</li> <li>• Roller pump controller</li> <li>• Power supply</li> <li>• Valve</li> <li>• GUI communications</li> <li>• Keypad</li> <li>• Temperature sensor</li> <li>• SD card</li> <li>• NVRAM</li> <li>• Vacuum system</li> <li>• System controller<sup>a</sup></li> <li>• System error<sup>b</sup></li> </ul>			

<sup>a</sup> No visual alert signal is issued.

<sup>b</sup> No audible alert signal is issued; message to contact Medtronic does not appear.

## 19 Service and maintenance

### 19.1 Service interval

A notification appears when maintenance is due.

**Warnings** – Read and follow all of the following warnings when servicing or maintaining the device.

- Contact your Medtronic service representative, or a biomedical equipment technician trained and certified by Medtronic, to perform preventative maintenance and calibration of the device. Preventative maintenance and calibration should be performed on an annual basis. Blood spillage, biohazard risk, low blood quality, slip hazard, and other concerns, could result if preventative maintenance and calibration are not performed.
- Always turn off and unplug the device prior to working inside the cabinet or changing any fuses. Inside the cabinet are various electrical components and wiring. Physical contact with any of these components while the device is plugged in could result in severe electrical shock. For continued protection against risk of fire, replace fuses only with the same type and rating. Fuses should be replaced by a Medtronic service representative or Medtronic-trained and certified biomedical equipment technician. Internal grounding is provided for safety.
- Enclosure current leakage is a primary indicator of electrical shock hazard to personnel making contact with any exposed portion of the equipment. Each device is tested during the final quality inspection to verify enclosure current leakage is less than 100  $\mu\text{A}$  (under normal operating conditions) and less than 300  $\mu\text{A}$  (under single-fault conditions). Ensure that the current leakage is tested at least once a year, or more frequently as required, by the operating facility's biomedical engineering department or other qualified service technician. In addition, always test the enclosure current leakage and insulation after an event such as a fluid spill, after a major voltage surge in the power source has occurred, or after any device repair.

**Note:** Contact your Medtronic service representative for replacement packaging prior to returning the device to Medtronic.


### 19.2 Pump head

- The tubing guides and pump head should be examined periodically for burrs and sharp edges that could lead to tubing damage and failure. Contact your Medtronic service representative if any defective or damaged parts are observed.
- Ensure that the pump head rollers move freely.

### 19.3 Service and maintenance screens

1. To view service and maintenance screens, press .
2. Select the desired field.


**Figure 38. Service and maintenance screen**

- 1 Current device notifications
- 2 Current device information such as software operating system, graphical user interface, and system controller versions; serial number; and device, centrifuge, roller pump, and vacuum pump run times
- 3 Logged device errors
- 4 Allows the user to calibrate the touchscreen
- 5 Allows the user to perform the centrifuge cleaning operation
- 6 Allows the service technician to perform authorized software updates
- 7 **Export service information** (not shown): Allows the user to export notifications, device information, and error logs

## 19.4 Notifications

Notifications provide information about device conditions that are important but do not stop device operation. The presence of a notification is indicated by a yellow outline on the  and the **Reset alert** field (Figure 39), and on the **Notifications** field of the **Service and maintenance** screen (Figure 38). The following information describes the appearance and behavior of device notifications:

- A new notification, or a postponed notification that has reappeared, causes yellow outlines to appear in the screen areas described previously.
- Pressing the **Reset alert** field allows the user to postpone a yellow-outlined notification. Postponing the notification causes the yellow outlines to disappear from the screen areas described previously. The following reset options are available:
  - **Minutes:** Select this option, and then use  and  to set the number of minutes before a notification reappears.
  - **Cases:** Select this option, and then use  and  to set the number of cases before a notification reappears.
  - **Forever:** For applicable notifications, select this option and the notification will not reappear.
- Notifications disappear entirely once the underlying condition is corrected. For example, the notification **Max vacuum is still set. Adjust if necessary.** disappears after the vacuum setting is changed.
- Active notifications, whether yellow-outlined or not, can be accessed by pressing **Notifications** from the **Service and maintenance** screen.
- If there are active notifications, and if there are no higher-priority prompts also active, the notifications will continuously appear (scroll) in the Prompts area of the **Main** screen.


The following table contains a list of device notifications:

**Table 2.** Notifications, audible and visual alert signals

Notification	Description	Audible alert signal	Visual alert signal
Maintenance due in __ days.	The device begins counting down when maintenance is due in 30 days, and updates the counter daily. Contact your Medtronic service representative to perform maintenance on the device.	None	Low
The real-time clock failed. See error log for details. Contact Medtronic for service.	Case records will not contain an accurate time stamp. Contact your Medtronic service representative.	None	Low
Max vacuum is still set. Adjust if necessary.	<b>Max vacuum</b> has been set for 10 min or more. If <b>Max vacuum</b> is not desired, change the vacuum setting as shown in Section 13.3.	None	Low
Emergency wash is still set. Adjust if necessary.	<b>Emergency wash</b> has been set for 10 min or more. If <b>Emergency wash</b> is not desired, change the <b>Wash volume</b> setting as shown in Section 13.5.1.	None	Low
Case records will not be stored unless a password is set in Settings.	Set the password, in order to store case records as shown in Section 12.4.	None	Low

**Table 2.** Notifications, audible and visual alert signals (continued)

Notification	Description	Audible alert signal	Visual alert signal
The SD card failed. See error log for details. Contact Medtronic for service.	The service log is not functional. Contact your Medtronic service representative.	None	Low
The vacuum system failed. See error log for details. Contact Medtronic for service.	The vacuum system is not functional. Contact your Medtronic service representative.	None	Low
Storage available for only ___ more case records. Export and delete records before they are overwritten.	Device internal data storage contains 90 or more case records. Delete records as described. See Section 14.5.	None	Low
The GUI communications failed. See error log for details. Contact Medtronic for service.	There is a problem with the internal electronics. Contact your Medtronic service representative.	None	Low
The keypad failed. See error log for details. Contact Medtronic for service.	The VAC or  control is not functioning properly. Contact your Medtronic service representative.	None	Low

**Figure 39.** Notifications screen

### 19.5 System information

To view system information (software versions and runtime statistics) press **System information** from the **Service and maintenance** screen.

### 19.6 Error log

To view device errors, press **Error log** from the **Service and maintenance** screen.

### 19.7 Calibrate touchscreen

**Note:** This function is not available during blood processing.

**Note:** If the touchscreen is inoperable, press and hold  for 10 s to calibrate the touchscreen. Follow the instructions that appear on the screen after 10 s.

1. Press  to stop processing if necessary.
2. Wait until the device is fully stopped. The message **Machine stopped. Press  to enable self-start.** appears.
3. Press **Calibrate touchscreen** from the **Service and maintenance** screen.
4. Follow the instructions on the screen.

### 19.8 Clean centrifuge

For instructions for cleaning the centrifuge, press **Clean centrifuge** from the **Service and maintenance** screen. Follow the instructions on the screen. See Chapter 16 for full cleaning instructions.

### 19.9 Update software

Contact your Medtronic service representative to perform software updates.

### 19.10 Export service information

1. To export information for service, press **Export service information** from the **Service and maintenance** screen.
2. Insert a USB flash drive into a USB port (type A).  
**Caution:** Use only a USB flash drive provided by Medtronic.
3. Press **Download to USB**.
4. The message **Download is complete.** appears. Press **OK**.

**Note:** If the download is unsuccessful, the message **Download failed** appears.

## 20 Device disposal and ownership transfer

### 20.1 Deleting case data

Delete all case data before disposing or transferring ownership of the device (Section 14.5).

### 20.2 End of life disposition

Contact your Medtronic service representative before disposing of the device.

## 21 Specifications

### 21.1 Technical specifications

<b>Electrical classification</b>	Class I, Type BF (suction/anticoagulant line), Ordinary, Continuous operation
<b>Power</b>	
Voltage	100 V~ to 240 V~
Frequency	50 Hz / 60 Hz
Phase	Single
Current	10 VA to 425 VA
Fuses	7 A / 250 V slow blow, 3AG, 200 A breaking capacity (Littelfuse 0313007.MXP or equivalent)
Power cord	3 prong hospital grade connector (varies by geography); length must not exceed 3 m (118 in)
<b>Speed, flow rate, and pressure</b>	
Centrifuge	0 rpm to 10 000 rpm ( $\pm 5\%$ )
Pump	0 mL/min to 1 000 mL/min ( $\pm 7\%$ )
Vacuum	-10 mmHg to -370 mmHg $\pm ( 5\%  + 8 \text{ mmHg})$
<b>Weight sensor</b>	Self-start: 800 mL $\pm 400$ mL
<b>Dimensions</b>	
Width	69 cm (27 in)
Height (without IV pole)	80.5 cm (31.7 in)
Depth	42.5 cm (16.7 in)
<b>Weight (device including IV pole)</b>	50 kg (110 lb)
<b>IP rating</b>	IPX1
<b>Temperature limit</b>	
Operational	15°C to 30°C (59°F to 86°F)
Storage (clinic)	15°C to 30°C (59°F to 86°F)
Storage (warehouse)	15°C to 30°C (59°F to 86°F)
Transit	-35°C to 60°C (-31°F to 140°F)
<b>Humidity range</b>	
Operational	25% to 70% noncondensing
Storage (clinic)	25% to 70%
Storage (warehouse)	10% to 90%
Transit	10% to 90%
<b>Pressure range</b>	
Operational	80 kPa to 101 kPa (11.6 psi to 14.6 psi)
Storage (clinic)	80 kPa to 101 kPa (11.6 psi to 14.6 psi)
Storage (warehouse)	80 kPa to 101 kPa (11.6 psi to 14.6 psi)
Transit	59.5 kPa to 101 kPa (8.6 psi to 14.6 psi)

### 21.2 Electromagnetic compatibility declaration

- IEC 60601-1-2, Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests.
- Class A, Group 1

**Warning:** The emissions characteristics of this device make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this device might not offer adequate protection to radio frequency communication services. The user might need to take mitigation measures, such as relocating or reorienting the device.

**Warning:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this device could result.

**Warning:** Do not stack or attach additional instrumentation to the device, in order to minimize the potential for RF interference.

**Warning:** Power cord and USB cable length must not exceed 3 m (118 in).

IEC 60601-1-2 table 4 – enclosure port		
Phenomenon	Basic EMC standard or test method	Immunity test levels
		Professional healthcare facility environment
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Radiated RF EM fields	IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See IEC 60601-1-2 table 9
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz

IEC 60601-1-2 table 5 – input AC power port (1 of 2)		
Phenomenon	Basic EMC standard	Immunity test levels
		Professional healthcare facility environment
Electrical fast transients / bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency
Surges Line-to-line	IEC 61000-4-5	±0.5 kV, ±1 kV
Surges Line-to-ground	IEC 61000-4-5	±0.5 kV, ±1 kV, ±2 kV
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage dips	IEC 61000-4-11	0% $U_T$ ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% $U_T$ ; 1 cycle and 70% $U_T$ ; 25/30 cycles <sup>a</sup> Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0% $U_T$ ; 250/300 cycle <sup>a</sup>

<sup>a</sup> Significant power loss may cause the device to shut off.

IEC 60601-1-2 table 8 – signal input/output parts port		
Phenomenon	Basic EMC standard	Immunity test levels
		Professional healthcare facility environment
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transients / bursts	IEC 61000-4-4	±1 kV 100 kHz repetition frequency
Surges Line-to-ground	IEC 61000-4-5	±2 kV
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

IEC 60601-1-2 table 9 – test specifications for enclosure port immunity to RF wireless communications equipment						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 – 470	GMRS 460, FRS 460	FM ±5 kHz deviation 1 kHz sine	2	0.3	28
710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						

IEC 60601-1-2 table 9 – test specifications for enclosure port immunity to RF wireless communications equipment						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1 720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1 845						
1 970						
2 450	2 400 – 2 570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5 240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5 500						
5 785						

## 22 Warranties

### 22.1 Equipment limited warranty for US customers<sup>1</sup>

(US customers only)

- A. This LIMITED WARRANTY provides the following assurance to the purchaser of the autoLog IQ autotransfusion system, hereafter referred to as the "Equipment":
- (1) Should the Equipment fail to function within normal tolerances due to a defect in materials or workmanship within a period of 1 year, commencing with the delivery of the Equipment to the purchaser, Medtronic will at its option: (a) repair or replace any defective part or parts of the Equipment; (b) issue a credit to the purchaser equal to the Purchase Price, as defined in Subsection A(2), against the purchase of the replacement Equipment; or (c) provide a functionally comparable replacement Equipment at no charge.
  - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement, Equipment.
- B. To qualify for this repair, replacement, or credit set forth in Section A, the following conditions must be met:
- (1) The Equipment must be returned to Medtronic within 60 days after discovery of the defect (Medtronic may, at its option, repair the Equipment on site).
  - (2) The Equipment must not have been repaired or altered outside of Medtronic's factory in any way which, in the judgment of Medtronic, affects its stability and reliability. The Equipment must not have been subjected to misuse, abuse, or accident.
- C. This LIMITED WARRANTY is limited to its express terms. In particular:
- (1) Except as expressly provided by this LIMITED WARRANTY, MEDTRONIC IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE, OR MALFUNCTION OF THE EQUIPMENT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.
  - (2) This LIMITED WARRANTY is made only to the purchaser of the Equipment. AS TO THE PURCHASER OF THE EQUIPMENT AND ALL OTHERS, MEDTRONIC MAKES NO OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOM, OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PURCHASER SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
  - (3) The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid. This LIMITED WARRANTY gives the purchaser specific legal rights. The purchaser may also have other rights which vary from state to state.
  - (4) No person has any authority to bind Medtronic to any representation, condition, or warranty except this LIMITED WARRANTY.

<sup>1</sup> This LIMITED WARRANTY is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States. Areas outside the United States should contact a local Medtronic representative for exact terms of the LIMITED WARRANTY.

## 22.2 Equipment limited warranty for customers outside the US<sup>2</sup>

(For countries outside the US)

- A. This LIMITED WARRANTY provides the following assurance to the purchaser of the autoLog IQ autotransfusion system, hereafter referred to as the "Equipment", that should the Equipment fail to function within normal tolerances due to a defect in materials or workmanship within a period of 1 year, commencing with the delivery of the Equipment to the purchaser, Medtronic will at its option: (a) repair or replace any defective part or parts of the Equipment; (b) issue a credit equal to the original Equipment purchase price (but not to exceed the value of the replacement Equipment), against the purchase of replacement Equipment; or (c) provide functionally comparable replacement Equipment at no charge.
- B. To qualify for the repair, replacement, or credit set forth in Section A, the following conditions must be met:
  - (1) The Equipment must be returned to Medtronic within 60 days after discovery of the defect (Medtronic may, at its option, repair the Equipment on site).
  - (2) The Equipment must not have been repaired or altered by someone other than Medtronic in any way which, in the judgment of Medtronic, affects its stability and reliability.
  - (3) The Equipment must not have been subjected to misuse, abuse, or accident.
- C. This LIMITED WARRANTY is limited to its express terms. In particular, Medtronic is not responsible for any incidental or consequential damages based on warranty, contract, tort, or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

---


<sup>2</sup> This LIMITED WARRANTY is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. Areas outside the United States should contact a local Medtronic representative for exact terms of the LIMITED WARRANTY.

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© 2018 Medtronic
Medtronic B.V.
Medtronic of Canada Ltd