

INSTRUCTIONS FOR USE

VO 100 Venous Oximetry



ENGLISH

Mespere
LIFE SCIENCES

Table of Contents

Section 1: Introduction	4
About this Manual	4
Paper Copy	4
Download	4
Previous Versions	4
Intended Use	4
Description	4
Indications for Use	4
Contraindications	4
Section 2: Safety Information	5
Warnings	5
Cautions	6
Section 3: Setup	7
Installation	7
Display Software Setup	7
Turning the Computer On and Off:	7
Monitor Display	8
VO 100 Venous Oximetry Sensor and Isolator	10
Sensor Adhesive	10
Section 4: System Operation	12
Setup for Non-Invasive and Continuous Monitoring of Jugular Venous Oximetry	12
Obtaining Optimal Measurement (i.e. adjusting patient incline angle)	12
Computer Functions	14
System Messages	17
Patient Environment	18
System Trend Display	18
Reading the Trend Data Display	18
Clearing Trend Information	18
Essential Performance	18
Sensor Performance Considerations	18
Factory Default Settings	18

Display Software	18
Section 5: Cleaning, Disinfection, Maintenance, and Disposal	19
Cleaning and Disinfection	19
Maintenance	19
Inspecting the Equipment and Accessories	19
Returning the System for Repairs and/or Service	19
Disposal	19
Section 6: Troubleshooting	19
For Technical Assistance	20
Section 7: Specifications	21
Performance	21
Electrical	21
Environmental	22
Physical Characteristics	22
Compliance	23
Warranty	23
Lifetime	23
Manufacturer’s Declaration	23
Section 8: Appendices	28
Appendix 1: List of Abbreviations	28
Appendix 2: List of Symbols	28
Appendix 3: List of VO 100 Venous Oximetry Components	31

Section 1: Introduction

About this Manual

Paper Copy

If a paper copy of this manual is needed, please contact your distributor or Mespere LifeSciences directly.

Download

The electronic version of this manual is available at <http://www.mespere.com/resources.html>. Please contact your distributor or Mespere LifeSciences directly to obtain the login information to download this manual.

Previous Versions

Please contact Mespere LifeSciences directly to obtain previous versions of this manual.

Intended Use

The VO 100 Venous Oximetry is intended to be used by healthcare professionals for the assessment of venous blood oxygenation through jugular vein (SjvO₂) of adult individuals.

Description

Mespere LifeSciences Inc. model VO 100 Venous Oximetry System is intended to be used by healthcare professionals for the assessment of hemodynamic parameters such as venous blood oxygenation from the jugular vein (SjvO₂) value (in percentage, %), Blood Volume Index (BVI) (in milliliter of blood per liter of tissue, %), a calculated Jugular Vein Distention Index based upon BVI (JVD^{index}) (in percentage, %), and corresponding plethysmographic waveform using non-invasive near infrared and proprietary signal-processing technologies. The SjvO₂, BVI, and JVD^{index} are displayed at a refresh rate of approximately one (1) second.

The saturation (SjvO₂) is calculated through the ratio of the amplitude attenuation of two different near infrared wavelengths. Blood Volume Index (BVI) is calculated through the sum of the amplitude attenuation of the two wavelengths. Jugular Vein Distention Index is calculated as the relative change in BVI using the following formula:

$$JVD^{index} = \frac{BVI_{max} - BVI_{min}}{BVI_{min}}$$

Indications for Use

Health Canada/CE Marking:

The Mespere VO 100 can measure the venous blood oxygenation (SjvO₂), static blood volume (BVI), and the relative change in blood volume (JVD^{index}) through the jugular vein. The VO 100 Venous Oximetry should not be used as the sole basis for diagnosis or therapy.

Contraindications

- Do NOT use the VO 100 Venous Oximetry on patients:
 - Undergoing photodynamic therapy (PDT).
 - Known to have allergic reactions to medical-grade adhesive tape.
 - With central vein stenosis.

Section 2: Safety Information

Warnings

WARNING



- Do NOT use the VO 100 Venous Oximetry unless on the order of a physician.
 - The VO 100 Venous Oximetry is intended for use by healthcare professionals only.
 - As per Standard IEC 60601-1 Section 11, it has been disclosed that the maximum temperature of the applied parts for the VO 100 Venous Oximetry is 43.2 degrees Celsius. It is recommended that the VO 100 Venous Oximetry is not left running on the patient when it is not in use so that the applied parts will not reach the maximum temperature. Caution must be taken for patients who are taking medication that would make their skin more susceptible to damage from sun or light.
 - The VO 100 Venous Oximetry may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the VO 100 Venous Oximetry or shielding the location.
 - To avoid the risk of electric shock, the DM-10-S Power Supply (PS-03) must be connected to an AC outlet with protective earth.
 - Do NOT re-use sensor adhesives. Doing so may increase the risk of infection or cross-contamination.
 - Always unplug the system before cleaning. Refer to **Section 5: Cleaning, Disinfection, Maintenance, and Disposal** for more information.
 - Any changes or modification to this equipment not expressly approved by Mespere LifeSciences Inc. may cause bodily harm and void your authority to use this equipment.
 - The use of accessories and cables other than those specified, with the exception of accessories and cables qualified and sold by Mespere LifeSciences Inc. may result in increased emissions or decreased immunity of the equipment and may cause the VO 100 Venous Oximetry to be non-compliant with the requirements of IEC 60601-1-2:2014.
 - Pins of connectors should not be touched and that connections should not be made to these connectors unless ESD precautionary procedures are used.
 - The operator should NOT touch the computer and the patient at the same time.
 - The VO 100 Venous Oximetry should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment should be observed to verify normal operation in the configuration in which it will be used.
-

Cautions

CAUTION



- Read and understand this entire manual carefully before using the VO 100 Venous Oximetry.
 - The VO 100 Venous Oximetry needs precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in **Section 7: Specifications**.
 - Portable and mobile RF communications equipment may affect the VO 100 Venous Oximetry.
 - Do NOT sterilize any part of the VO 100 Venous Oximetry using any process.
No part of the VO 100 Venous Oximetry is to be sterilized. Cleaning must be performed as described in **Section 5: Cleaning, Disinfection, Maintenance, and Disposal**.
 - A cable connects the VO 100 Venous Oximetry sensor and isolator to the computer.
Do NOT twist or pull on this connection.
Do NOT allow this cable to become tangled in other equipment such as bed frames.
-

Section 3: Setup

Installation

Before the VO 100 Venous Oximetry can be used in a clinical setting, it should be inspected and properly setup.

Display Software Setup

The display software can be pre-loaded by distributor or Mespere LifeSciences to a computer. If the display software is not preloaded, please refer to the Display Software Installation document.

Procedure

1. Download VO 100 Venous Oximetry Installer file provided to you via Dropbox or by a Mespere LifeSciences representative.
2. Plug an external hard drive, thumb drive, or USB drive into your computer.
3. Place the VO 100 Venous Oximetry Installer file onto the external hard drive, thumb drive, or USB stick.
4. Plug external hard drive, thumb drive, or USB drive into the Computer. A pop-up window should appear showing the contents of this drive.
5. Open the VO 100 Venous Oximetry Installer file.
6. Allow the VO 100 Venous Oximetry Installer file to make changes to your computer by clicking "Yes".
7. Click "Install" on the VO 100 Venous Oximetry Installer Setup Wizard.
8. Once the VO 100 Venous Oximetry Installation is complete press "Finish". VO 100 Venous Oximetry Display Software is now installed on the computer.

Turning the Computer On and Off:

Once the computer is turned on the Mespere home screen will appear. If the display software is not pre-loaded please refer to the Display Software Installation document.

Procedure

1. Turn the computer on by pushing the power button.
2. Once the Mespere home screen appears, the system is ready for use.
3. Turn the monitor off by pushing the virtual power button on the home screen

Monitor Display

The VO 100 Venous Oximetry displays the plethysmographic waveform, S_{jv}O₂, BVI, and JVD^{index} values in real time.




Control/Indicator	Description
Jugular Venous Oxygen Saturation (S _{jv} O ₂) Number	Shows the numerical value of venous blood oxygenation from the jugular vein in percentage (%)
Pulse Strength Index (PI)	Indicates the strength of the pulse signal. When the signal is less than 1, indicating a weak signal, the number will be displayed in red
Blood Volume Index (BVI)	Shows the index of the static blood volume in the jugular vein beneath the sensor in mL of blood per L of tissue (%)
Jugular Vein Distention Index (JVD ^{index})	<p>Respiratory Variation: relative change in jugular vein cross section/blood volume index throughout a respiratory cycle.</p> <p>Fluid Challenge: relative change in jugular vein cross section/blood volume index during a fluid challenge protocol.</p> $JVD^{index} = \frac{BVI_{max} - BVI_{min}}{BVI_{min}}$
Sensor Selection	Choose which sensor (R or L) will be used for the displayed JVD ^{index} calculations (R is suggested). If only one VO100 sensor is connected it will automatically be chosen. This can be switched throughout the monitoring session.
Trend Plot	The S _{jv} O ₂ , BVI, and JVD ^{index} data trends
Waveform	The plethysmographic waveform detected by the system
ID – Patient Folder	Once a patient ID is entered, a patient folder will be created using this ID and all patient data will be saved in this folder. If a patient ID is not entered the files will be saved in the temp folder using the current time stamp as the file name.
New Measurement	Pressing this button will clear the stopped data on the screen
Start	Pressing this button will start a new monitoring session
Stop – Auto Save	Press this button to stop the current monitoring session. Once the monitoring session is stopped an excel file, screenshot, marked events summary, and summary pdf will automatically be saved to the patient's folder.

Control/Indicator	Description
Event Marker	Mark an event by displaying a vertical line on the trend window. Enter event details and/or press OK once complete. Marked events can be reviewed at any time under the Review – Trend window.
Review	Trend: view the complete StO ₂ , BVI, and JVD ^{index} trend summary and marked events from the entire monitoring session. Press Save Review on the top toolbar to save a pdf summary and screenshot to the patient's folder Compare Snap Shot: Select two Snap Shot files to compare the data and create a comparison report. Press Save Comparison on the top toolbar to save a pdf summary to the patient's folder.
Snap Shot	Create a report based on the average of the last 30 seconds of data. An excel file and pdf summary will automatically be saved to the patient's folder.
Trend Toggle	Change the displayed data trend (SjvO ₂ and BVI)
Export Data	Export previously saved patient files such as Snap Shots, Review Summaries, Excel Files, and Screenshots to a removable drive. Patient data is saved by patient ID. Locate the patient's folder and copy desired files to removable drive. If a patient ID was not entered, the patient's data files will be located in the temp folder and can be identified using the time stamps as the file name.
Help	Help Center
Trend lines	The SjvO ₂ , BVI, and JVD ^{index} trend.
Help	Help Center
Data files	Previously saved data files.
Sensor Life	Provides the number of remaining hours on the sensor.
Demo Mode	Opens the Demo Application.
Save	Saves the screenshot, including SjvO ₂ , BVI, and JVD ^{index} values and plethysmographic waveform, to a user-specified file.
Start/Stop	Clicking this button controls the start or stop measurement. Depending on the program state, its function will either start or stop the measurement.
Exit	Returns back to the home screen.
Battery Status	Provides the status of the battery.


VO 100 Venous Oximetry Sensor and Isolator

The VO 100 Venous Oximetry Sensor consists of a sensor cable and sensor head. The sensor head contains the LEDs and PDs. The isolator consists of a cable connector, an isolator, flexible cable, and USB connector. The isolator connects the VO 100 Venous Oximetry sensor to the computer or USB hub.



Procedure

1. Open the clear case on the isolator cable connector.
2. Plug the cable connector from the VO 100 Venous Oximetry sensor into the cable connector of the isolator.
3. Close the clear case on the cable connector.



Sensor Adhesive

The VO 100 Venous Oximetry sensor adhesives are single-use, and use a medical-grade adhesive and foam.

The sensor adhesive clips onto the VO 100 Venous Oximetry sensor head and is used to adhere the VO 100 Venous Oximetry sensor to the patient.



WARNING



- The VO 100 Venous Oximetry sensor adhesives are intended for single-use only.
 - Do NOT re-use the VO 100 Venous Oximetry sensor adhesives. Doing so may increase the risk of infection or cross-contamination.
-

Section 4: System Operation

Setup for Non-Invasive and Continuous Monitoring of Jugular Venous Oximetry

Monitor Setup

1. Choose the S_{ijv}O₂ application from the home screen
2. Connect the VO 100 Venous Oximetry sensor to the isolator via the cable connector.
3. Plug the VO 100 Venous Oximetry sensor that is connected to the isolator into the monitor.

Sensor Setup

4. Snap a disposable sensor adhesive onto the sensor head with the yellow dot adjacent to the logo.

Sensor Placement

1. Identify the patient's right internal or external jugular vein (IJV or EJV).
2. Position the right sensor (R) so that the arrow points towards the patient's head and the midline of the sensor is directly over the IJV or EJV.
3. Repeat with the left sensor (L) if conducting bilateral readings.



Patient Setup

4. Place the patient at a comfortable inclination (between 0 and 30 degrees is suggested).
5. Slightly tilt the patient's head to the left, and have the patient remain quiet and keep their neck relaxed.

Measuring Jugular Venous Oximetry (S_{ijv}O₂)

6. Press the Start button on the monitor.
7. Wait for instructions to appear in the message center on the screen.

Measuring Jugular Vein Distention Index (JVD^{index})

8. Choose which sensor (R or L) will be used for the displayed JVD^{index} calculations (R is suggested) by pressing the radio button. A black dot will be in the radio button indicating the selected sensor. If only one VO100 sensor is connected it will automatically be chosen.
9. Select which method will be used for the JVD^{index} calculation, Respiratory Variation or Fluid Challenge, by clicking on the box to the right. A check mark will appear in the box indicating the selected method.
10. Sensor and JVD^{index} method can be switched throughout the monitoring session using the radio and check box buttons.

Obtaining Optimal Measurement (i.e. adjusting patient incline angle)

Follow the instructions on the screen, recline or incline the patient until a stable S_{ijv}O₂ reading and clear waveform appear on the screen. The VO 100 sensor uses Near Infrared Spectroscopy (NIRS) technology to measure jugular venous oxygenation. The VO 100 sensor shines near infrared photons into the tissue containing jugular venous blood vessels via a sensor array. The sensor array analyzes the diffusely reflected photons and is then able to determine jugular venous oxygenation. In order to get an exact reading the sensor must detect the height of the venous blood column beneath the sensor. Changing the

degree of inclination for the patient will help position the jugular venous blood column beneath the VO 100 sensor in order to achieve a stable reading.

If PI is less than 1 and in red:

1. Ensure that the sensor is placed directly over the right or left jugular vein.
2. The reading may not be reliable and must be used with the physician's judgement of the waveform.

Computer Functions

ID – Enter Patient ID

1. Enter a patient ID into the white input box anytime throughout the monitoring session.
2. A patient folder will be created using this ID and all patient data will be saved in this folder. If a patient ID is not entered the file will be saved in the temp folder using the current time stamp as the file name.

New Measurement (New Meas.)

1. Press New Meas. to clear the previously stopped data showing on the screen.

Start

1. Press Start to begin a new monitoring session.

Stop – Auto Save

1. Press Stop to stop the current monitoring session.
2. Once the monitoring session is stopped an excel file, screenshot, and summary pdf will automatically be saved to the patient's folder.

Trend

1. Press the Trend button.
2. A dropdown menu will appear to change the displayed trend in the trend plot; SjvO2 or BVI.

Event Marker

1. Press the Event Marker button.
2. Enter event details and press OK.
3. A vertical line will appear in the trend window marking an event.

Change Event Description

1. Press the Review button and select Trend.
2. The Review – Trend screen will appear. Press View Events from the top toolbar.
3. A pop-up window will appear with the marked events.
4. Click the event box and enter the appropriate description.

Review - Trend

1. Press the Review button on the monitor. A dropdown menu will appear and select Trend.

2. A window will appear with a complete SjvO₂, BVI, JVD^{index} trend summary and marked events from the entire monitoring session.
3. Press Save Review on the top toolbar to save a pdf summary and screenshot to the patient's folder.

Review – Compare Snap Shot

1. Press the Review button on the monitor. A dropdown menu will appear and select Compare Snap Shot.
2. A pop-up window will appear. Select two Snap Shot files to compare the data and create a comparison report.
3. Press Save Comparison on the top toolbar to save a pdf summary to the patient's folder.

Stop – Auto Save

3. Press the Stop button located on the top toolbar once the monitoring session is complete.
4. A pop-up window will appear showing that an excel file, screenshot, marked events, and summary pdf will automatically be saved to the patient's folder.

Snap Shot

1. Press the Snap Shot button located on the top toolbar.
2. A report will be created based on the average of the last 30 seconds of data. An excel file and pdf summary will automatically be saved to the patient's folder.

Export Data – Snap Shots, Review Summaries, Excel Files, and Screenshots

1. Press the Export Data button on the top toolbar.
2. Locate the patient's folder. If a patient ID was entered, a folder was created using this ID and all patient data will be saved in this folder. If a patient ID is not entered the file will be saved in the temp folder using the current time stamp as the file name.
3. Plug in an external memory source to the USB port on the right side of the computer.
4. Click the file until a menu appears. Select Copy
5. Open the external memory source and hold the screen until a menu appears. Select Paste. The copied files will be transferred to the connected removable drive.

Deleting saved files

1. Press Export Data button on the top toolbar.
2. Locate the patient folder.

3. Click the file until a menu appears.
4. Select delete from the menu.

Checking Sensor Life

1. Plug the sensor into the USB port on the right side of the computer.
2. Press the Help button.
3. Press the sensor life button.
4. A pop-up window will appear indicating the number of hours remaining for that specific sensor.

System Messages

The guidance messages on the Display are shown in Table 1:

Table 1: Informative system messages

Message	Cause and Resolution
"VO100 Sensor is disconnected"	<p>The VO 100 Venous Oximetry is detecting that the sensor is not connected to the computer.</p> <ul style="list-style-type: none"> - Ensure that the VO 100 sensor is correctly connected to the isolator. - Ensure that the isolator is correctly connected to the computer.
"VO100 Sensor Expired. Please Replace!"	<p>The VO 100 Venous Oximetry is detecting that the attached VO100 sensor has expired.</p> <ul style="list-style-type: none"> - Replace the connected VO100 sensor with a new one or one that hasn't expired
"Check VO100 Sensor Patch Adhesion"	<p>This may if the sensor has loose adhesion to the patient's neck.</p> <ul style="list-style-type: none"> - Ensure that the sensor is adhered to the patient's neck.
"Weak Pulse"	<p>The VO 100 Venous Oximetry is detecting that the PI is less than 1.0 therefore a defined waveform may not appear.</p> <ul style="list-style-type: none"> - Adjust the patient angle in order to position the taper of the jugular vein beneath the LEDs and PDs for a defined waveform.
"Demo Mode"	<p>The VO 100 Venous Oximetry is running in demo mode.</p> <ul style="list-style-type: none"> - Press "Exit Demo Mode" or connect a sensor to exit the demo application.
"Please wait"	<p>The VO 100 Venous Oximetry is loading. Please wait until this message has closed, before continuing.</p>
"__ hours of usage has elapsed"	<p>The attached VO 100 sensor has been used for this many hours.</p>

Patient Environment

The sensor adhesives are the applied parts of the VO 100 Venous Oximetry that will come in contact with the patient during normal use.

The sensor and sensor adhesive are components of the VO 100 Venous Oximetry that will be in the patient environment.

NOTE

The patient environment is defined as the area that is 1.5m from the patient in their bed or chair.

CAUTION

The isolator chip and the computer should be kept away from the patient environment.

System Trend Display

The history of the S_{jv}O₂, BVI, and JVD^{index} trend is displayed on the computer.

Reading the Trend Data Display

The S_{jv}O₂, BVI, and JVD^{index} trend data is displayed on the computer, with the x axis representing number of data points and the y axis displaying the history S_{jv}O₂, BVI, and JVD^{index} data. The current version of the VO 100 Venous Oximetry does not allow the user to display the S_{jv}O₂, BVI, and JVD^{index} trend data in another unit.

Clearing Trend Information

All history VO 100 Venous Oximetry data is erased when a new measurement is started by pressing the Start button.

Essential Performance

The VO 100 Venous Oximetry has no essential performance characteristics.

Sensor Performance Considerations

CAUTION

Patients must not move excessively during measurements. Excessive talking and motions such as yawning, deep-breathing and the like should be discouraged while measurements are being taken.

Factory Default Settings

Display Software

The computer will include the Mespere Display Software (OS-VO100) with software version VO100.EN.01.00.046

Section 5: Cleaning, Disinfection, Maintenance, and Disposal

Cleaning and Disinfection

The VO 100 Venous Oximetry requires no preventive maintenance. Only regular cleaning and disinfection based on the institutional or jurisdictional guidelines is required. In the absence of guidelines, it is recommended that cleaning is performed after each use of the VO 100 Venous Oximetry.

WARNING



- Always unplug the VO 100 Venous Oximetry before cleaning.
- Allow the VO 100 Venous Oximetry to dry completely before re-connecting power.

CAUTION



- Do NOT immerse any part of the VO 100 Venous Oximetry in any liquid for any reason.
- Cleaning should be performed using institutional or jurisdictional guidelines.
- Do NOT subject any component of the system to sterilization (e.g. ETO, steam). Doing so may damage the components of the VO 100 Venous Oximetry.

Maintenance

The VO 100 Venous Oximetry requires no preventive maintenance.

Inspecting the Equipment and Accessories

Prior to each use of the VO 100 Venous Oximetry, it is suggested that the user checks that all connections to the computer, docking station, isolator, and sensor are correct and secure. For instructions on how to make these connections, please refer to **Section 3: Setup**.

Returning the System for Repairs and/or Service

Contact your distributor for information on repairs and/or services for the VO 100 Venous Oximetry.

Disposal

Comply with local laws for disposal of VO 100 Venous Oximetry.

CAUTION



Product is to be taken to separate collection at end of product life. Do NOT dispose of the product as unsorted municipal waste.

Section 6: Troubleshooting

The following table provides basic guidance in the event of a problem with the VO 100 Venous Oximetry.

Table 2: List of Problems for Troubleshooting

Problem	Possible Causes	Possible Remedies
Computer will not power on	<ul style="list-style-type: none"> - No power 	<ul style="list-style-type: none"> - Check connection of the power supply to the computer. - Check connection of Power Supply to AC power (wall outlet). - Try alternate outlet.
Computer cannot communicate sensor	<ul style="list-style-type: none"> - Connection of sensor 	<ul style="list-style-type: none"> - Ensure sensor is connected correctly to the isolator and computer. - Try re-starting software by powering down the computer, then turn the computer on.
VO 100 Venous Oximetry does not display S _{ijv} O ₂ value, or plethysmographic waveform	<ul style="list-style-type: none"> - Patient angle - Patient status - sensor placement 	<ul style="list-style-type: none"> - Be sure to follow the on-screen instructions for adjusting the patient angle. Refer to Section 3: Setup. - Check the sensor placement and orientation. - The VO 100 Venous Oximetry may not work with all patients due to anatomical or health status.
"Check VO100 Sensor Patch adhesion"	<ul style="list-style-type: none"> - Sensor has loose adhesion to the patient's neck. 	<ul style="list-style-type: none"> - Ensure that the sensor is adhered to the patient's neck.
" VO100 Sensor is disconnected"	<ul style="list-style-type: none"> - The VO 100 Venous Oximetry is detecting that the sensor is not connected to the isolator or the isolator is not connected to the computer. 	<ul style="list-style-type: none"> - Ensure that the sensor is correctly connected to the computer. - Ensure that the isolator is correctly connected to the computer.
"Please Wait"	<ul style="list-style-type: none"> - The VO 100 Venous Oximetry is loading 	<ul style="list-style-type: none"> - Before continuing wait until this message disappears

For Technical Assistance

In the event that your problem is not resolved with the help of Table or if you have additional problems to report, contact your distributor or Mespere LifeSciences directly.

Section 7: Specifications

Performance

Table 3: Performance specifications

Parameter	Specification	Notes
SjvO ₂ Measurement Range	0 – 100%	
Accuracy of SjvO ₂	± 2.7%	1,2
Resolution of SjvO ₂	1%	
Accuracy of BVI	± 4.3 mL of blood per L of tissue	3
Resolution of BVI	0.1 mL/L	

NOTES:

- 1 Measurement range is 45% - 91%
- 2 Accuracy is based on a validation study comparison with external jugular venous blood samples and a Co-Oximeter on healthy volunteers. In the SjvO₂ range of 45-91%, the accuracy was ± 2.7% in the study STP-9100001 at the University of California San Francisco, UCSF (San Francisco, California, USA).
- 3 Accuracy is based on a validation study comparison against ultrasound measurements during an inclination protocol on healthy volunteers. The accuracy was ± 4.3ml/L in the study STP-9200002 at the Wayne State University iBio Center (Detroit, Michigan, USA).

Electrical

Table 4: Electrical specifications

Parameter	Specification	Notes
Computer Specifications	CPU: Quad Core Intel® Atom™x7-Z8700 processor (2MB Cache, 1.6GHz with Intel Burst™ technology up to 2.4GHz) or faster RAM: 2G or more Storage: 64GB or more USB: 5Vdc, 500mA, USB 2.0 port OS: WIN 10	
Power Supply	AC Input: 100-240V, ~0.3A, 50-60Hz DC Output: 5.2V ± 2.5A	
Measurement NIR wavelength	905±5 and 660±3nm	
Base Sampling Modulation Rate	2 kHz	
NIR LED duty cycle	7.5 ms	
NIR power	Peak	14 mw
	Average	0.9 mw

NOTES:

- 3 Peak power indicates instantaneous electrical power dissipated by one NIR LED when 'on'. Average power reflects the time-averaged electrical power dissipated by one NIR LED when the system is running normally at a sampling rate of 2kHz and an NIR LED duty cycle of 9.1%.
- 4 Power figures are given for one NIR LED. Only one of the two NIR LEDs is active at any given time.

NOTE

The Computer and the device form a ME system.

Environmental

Table 5: Environmental specifications

Parameter		Specification
Temperature	Operating	10 °C – 38 °C
	Transportation/Storage	-20 °C – 50 °C
Humidity	Operating	30% – 70%
	Transportation/Storage	20% – 80%
Air Pressure	Operating	86 – 106 kPa
	Transportation/Storage	86 – 106 kPa

Physical Characteristics

Table 6: Physical characteristics specifications

Parameter		Specification
Dimensions of VO 100 Sensor	L × W × D	54mm × 44mm × 8.9mm
Dimensions of Isolator	L × W × D	33mm × 33mm × 9mm
Data Display		SjvO ₂ (%), BVI, JVD ^{index} and plethysmographic waveform and, SjvO ₂ , BVI, JVD ^{index} trend.
Display update rate		1 second

Compliance

Table 7: Compliance

Parameter	Specification
Safety Standard for Medical Equipment	IEC 60601-1 edition 3 Amendment 1
EMC Standard	IEC 60601-1-2 Class A
Degree of protection	Type BF-Applied part
Mode of Operation	Continuous

Warranty

Mespere LifeSciences Inc. warrants to the Purchaser for a period of one (1) year from the date of purchase for the monitor, and six (6) months or 100 hour operation hours for the sensor, whichever comes first. Mespere warrants that the Product delivered is free from defects in workmanship or materials and the Product will perform as labelled in the directions for use. Mespere LifeSciences will work with the end-users distributor, or directly with the end user to determine warranty coverage for issues, should they arise, and if the issue is not covered by warranty, will provide a quote with respect to the repair costs to rectify the issue.

EXCLUSIONS

The warranty does not apply in the following conditions:

- 1) The product is modified without authorization from Mespere LifeSciences Inc.
- 2) Used with devices not supplied or authorized by Mespere LifeSciences Inc.
- 3) Used with sensor not supplied by or approved by Mespere LifeSciences Inc.
- 4) Used in conditions other than mentioned in the Instructions For Use
- 5) Products used for testing or demonstration purposes
- 6) Mishandling of equipment
- 7) Accident, fire, water, vandalism, weather or any act of God

Lifetime

Product Code	Description	Lifetime
SEN-VO100	VO 100 Venous Oximetry Sensor	1 year
SA-VO100	VO 100 Adhesives (Single use)	1 year
DM-10	Mespere Computer	3 years
MM-01	Monitor Mount	3 years
IS-01	Mespere Isolator	3 years

Manufacturer's Declaration

WARNING



The use of accessories and cables other than those specified with the exception of accessories and cables qualified and sold by Mespere LifeSciences Inc. may result in increased emissions or decreased

immunity of the equipment and may cause the system to be non-compliant with the requirements of IEC 60601-1-2:2014.

Table 8: Electromagnetic Emissions

The VO 100 Venous Oximetry is intended for use in the electromagnetic environment specified below. The customer or the user of the VO 100 Venous Oximetry should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The equipment is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the equipment or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 9: Electromagnetic Immunity

The VO 100 Venous Oximetry is intended for use in the electromagnetic environment specified below. The customer or the user of the VO 100 Venous Oximetry should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±4 kV contact ±4 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines Input/output lines are exempted because none are more and 3 meter in length	Mains power quality should be that of a typical commercial/residential or hospital environment.

Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial/residential or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 U _T = 240 Vac and 120Vac	<5 % U _T (>95 % dip in U _T) for 0.5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles 70 % U _T (30 % dip in U _T) for 25 cycles <5 % U _T (>95 % dip in U _T) for 5 sec	<5 % U _T (>95 % dip in U _T) for 0.5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles 70 % U _T (30 % dip in U _T) for 25 cycles <5 % U _T (>95 % dip in U _T) for 5 sec	Mains power quality should be that of a typical commercial/residential or hospital environment. If the user of the equipment requires continued operation during power mains interruptions, it is recommended that the equipment be powered from an uninterruptible power supply or a battery.
Power frequency (50 Hz/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial/residential or hospital environment.

Table 10: Electronic Immunity for Portable and Mobile RF Equipment

The VO 100 Venous Oximetry is intended for use in the electromagnetic environment specified below. The customer or the user of the VO 100 Venous Oximetry should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the equipment including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \frac{3.5}{v_1} \sqrt{P}$ $d = \frac{3.5}{v_1} \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

			$d = \frac{3.5}{v_1} \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p>
--	--	--	--

NOTE



At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE



These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey¹ should be less than the compliance level in each frequency range².

Interference may occur in the vicinity of known RF transmitting devices and equipment marked with the following symbol:



Table 11: Recommended Separation Distances

The ME equipment is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ME equipment can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ME equipment as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Power of Transmitter W	Separation Distance According to Frequency of Transmitter m		
	150 kHz to 80 MHz $d = \frac{3.5}{v_1} \sqrt{P}$	80 MHz to 800 MHz $d = \frac{3.5}{E_1} \sqrt{P}$	800 MHz to 2.5 GHz $d = \frac{7}{E_1} \sqrt{P}$
0.01	0.12	0.12	0.24

¹ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.

² Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

0.1	0.37	0.37	0.74
1	1.17	1.17	2.34
10	3.69	3.69	7.38
100	11.7	11.7	23.4

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE



At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE



These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Table 12: Cable Compliance

The following cables comply with:

- RF emissions, CISPR 11, Class A/Group 1
- Electrostatic discharge (ESD), IEC 61000-4-2
- Radiated RF, IEC 61000-4-3
- Electric fast transient/burst, IEC 61000-4-4
- Conducted RF IEC 61000-4-6

Cables and Accessories	Maximum Length
Cable of the VO 100 Venous Oximetry sensor (SEN-VO100)	1.5m
Cable of isolator (IS-01)	1.5m

Section 8: Appendices



Appendix 1: List of Abbreviations




















Table 13: Glossary of Abbreviations


Acronym	Meaning
SjvO ₂	Venous Blood Oxygenation from jugular vein This is the oxygen saturation of the venous blood in the jugular vein.
BVI	Blood Volume Index The sum of the amplitude attenuation of the two wavelengths. It is the blood to tissue volume ratio of regional tissue. For example; if BVI is 10% that means that blood fills 10% of the tissue volume.
JVD ^{index}	Jugular Vein Distention Index The relative change in jugular vein cross section/BVI which indicates jugular vein distensibility.
EJV	External Jugular Vein This is the vein in the side of the neck in which the Neck Patch is placed.
IJV	Internal Jugular Vein As differentiated from the EJV used by the Mespere VO 100 System.
BVI	Blood Volume Index Index of the static blood volume in the jugular vein beneath the sensor in mL of blood per L of tissue (%).
NIR	Near-Infra-Red Light having a wavelength between 0.7- and 5-microns (700nm to 5000nm).

Appendix 2: List of Symbols

Table 14: Glossary of Symbols

Symbol	Meaning
	“Note” message used to convey tips, helpful information and instructive statements.
	“Caution” message used to call attention to installations, practices, and operations that, if not properly performed or adhered to, could result in damage to the equipment.

Symbol	Meaning
	“Warning” message used to call attention to installations, practices, and operations that, if not properly performed or adhered to, could result in personal injury to the patient and/or user.
	Mark of conformity to European Union Medical Device Directive.
	Federal law (U.S.) restricts this device to sale by or on the order of a physician.
	Direct current (DC)
	Power indicator, located onto or near power switches.
	ESD sensitivity. This symbol may appear adjacent ports on the equipment that have heightened sensitivity to static discharge.
	Type BF applied part. This symbol indicates that the applied part is suitable for connection to the patient and that it delivers Type BF isolation.
	Consult Manual
	Class II Symbol
	European Authorized Representative
	Lot Number
	Reference/Model Number
	Serial Number
	Part number (re-order number)
	Manufacturer
	Consult Instructions For Use
	Single Use Only
	Product is to be taken to separate collection at end of product life. Do NOT dispose of the product as unsorted municipal waste.
	Date of Manufacture

Symbol	Meaning
	Storage temperature range: -20°C to 50°C

NOTE



Not all symbols shown in Table may appear on the equipment or in this manual.

Appendix 3: List of VO 100 Venous Oximetry Components

Table 15: List of VO 100 Venous Oximetry Components

Description	Component Part Number
VO 100 Venous Oximetry Sensor	SEN-VO100
VO 100 Adhesives (Single use)	SA-VO100
VO 100 Operating System	OS-VO100

Contact:

Mespere LifeSciences Inc.
180 Frobisher Drive, Suite 1C
Waterloo, Ontario
N2V 2A2
Tel: (519) 884-7575
Fax: (519) 884-8810
Email: info@mespere.com
www.mespere.com

Authorized Representative:



Ms Lu Lu
Karl-Jatho-Weg 10,
71034 Boeblingen, Germany
Phone: +0049 17695706210

VO-100-IU_EN Rev S
Date of issue 2021-02
Printed in Canada

© 2021 Mespere LifeSciences Inc.

All rights reserved. No parts of this publication may be reproduced, in any form whatsoever, without the prior written permission of Mespere LifeSciences, Inc.

