Medtronic

MYOthem XP[®]

Cardioplegia Delivery System
Système d'infusion de cardioplégie
Kardioplegie-Infusionssystem
Sistema de administració para cardioplejía
Cardioplegietoedieningssysteem
Sistema di somministrazione della cardioplegia
Kardioplegisystem
Kardioplegiajärjestelmä
Kardioplegitillförselsystem
Σύστημα παροχής καρδιοπληγίας
Cardioplegi leveringssystem
Układ podawania kardiopleginy
Sistema de administração de cardioplegia
Kardiyopleji İletim Sistemi
Systém pro podávání kardioplegie
Kardioplégiát beadó rendszer
Kardioplegický privodný systém
Kardiopleegia tekitamise süsteem
Система за доставяне при сърдечен удар
Sustav za provedbu kardioplegije
Kardioplégijas šķīduma pievades sistēma
Kardiopleginė tiekimo sistema
Sistem de livrare cardioplegie
Система для кардиоплегии
Sistem za kardioplegiju
Sistem za dovajanje raztopine za kardioplegijo

Instructions For Use ▪ Mode d'emploi ▪ Gebrauchsanweisung ▪ Instrucciones de uso
Gebruiksaanwijzing ▪ Istruzioni per l'uso ▪ Bruksanvisning ▪ Käyttöohjeet ▪ Bruksanvisning
Οδηγίες χρήσης ▪ Brugsanvisning ▪ Instrukcja użytkowania ▪ Instruções de utilização
Kullanım Talimatları ▪ Návod k použití ▪ Használati utasítás ▪ Pokyny pre používanie
Kasutusjuhend ▪ Инструкции за употреба ▪ Upute za upotrebu ▪ Lietošanas pamācība
Naudojimo instrukcijos ▪ Instrucțiuni de utilizare ▪ Инструкция по эксплуатации
Uputstva za upotrebu ▪ Navodila za uporabo

USA Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Explanation of symbols on package labeling / Explicación des symboles des étiquettes sur l'emballage / Erklärung der Symbole auf dem Verpackungsetikett / Explicación de los símbolos del etiquetado del envase / Verklaring van de symbolen op de verpakkingslabels / Spiegazione dei simboli sulle etichette della confezione / Forklaring av symboler på emballasjen / Pakkauksen merkintöjen selitykset / Förklaring av symboler på förpackningsetiketten / Επεξήγηση των συμβόλων στην ετικέτα της συσκευασίας / Forklaring af symboler på emballagens mærkater / Objašniená symboli z najdujacych się na etykietach opakowania / Explicação dos símbolos nas etiquetas da embalagem / Ambalaj etiketindeki sembollerin açıklamaları / Vysvětlení symbolů uvedených na štítcích obalu / A csomagoláson látható szimbólumok jelentése / Vysvětlivky k symbolom na označení balenia / Pakendi siltidel olevate tingmärkide tähendus / Обяснение на символите от етикетите на опаковката / Objašnjenje simbola na etiketi na ambalaži / Iepakojuama etiķetē attēlotu simbolu skaidrojums / Simbolių paaiškinimai pakuočių etiketėje / Explicarea simbolurilor cuprinse în etichetele de pe ambalaj / Объяснение обозначений на наклейках продукта / Objašnjenje simbola na oznakama na pakovanju / Razlaga znakov na ovojnicini

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Cardioplegia Delivery System

Description

The MYOthern XP® systems are designed to mix arterial blood from an oxygenator with asanguineous cardioplegia solution in specific ratios, depending upon the system chosen. The MYOthern XP® Cardioplegia Delivery System, depending on the model selected, may be used with either a standard roller pump or the Medtronic® CSS™ Cardioplegia Safety System (Model 990 instrument). The Medtronic CSS system is a micro-computer controlled cardioplegia delivery instrument featuring dual, independently operating pumps. The dual-pump design allows desired blood and solution ratios to be delivered based on the relative speeds of each pump (for detailed information, refer to the CSS instrument Operator and Reference Manual). The specifications provided by each standard system and their corresponding product designations are presented below in Table 1. The patient delivery line is packaged separately to allow for ease of sterile transfer to the operative field prior to connection to the delivery outlet of the cardioplegia system. Depending on the system selected, a pressure monitoring line, three-way stopcock, and gauge protector, as well as a pressure relief line with pressure relief valve may be preconnected.

Table 1.

MYOthern XP® Model	Delivery Ratio ^a	Tubing ID Dimensions		Pump Type
		Blood Line	Cardioplegia Line	
41, 41-B	4:1	0.6 cm (1/4 in)	0.3 cm (1/8 in)	Standard Roller Pump
CSS Custom Only	Set on CSS Instrument	0.6 cm (1/4 in)	0.3 cm (1/8 in)	CSS 990 Instrument

^a **NOTE:** MYOthern XP®-B (Bridge) series provides the ability to administer 100% blood as prescribed by physician.

Specifications

Bellows:	Stainless Steel
Recommended Blood Flow Rate:	Up to 500 mL/min.
Filter Screen:	150 micron
Heat Exchanger Volume:	44 mL
Heat Exchanger Blood Contact Surface Area:	0.06 m ² (95 in ²)
Size:	Adult
Inlet Connector:	0.6 cm (1/4 in)
Outlet Connector:	0.5 cm (3/16 in)
Pressure Monitor/Pressure Relief Ports:	Standard Female Luer Lock
Waterports:	Hansen Series 3-ST Quick Disconnect [1.3 cm (1/2 in)]
Water phase pressure limit:	40 psi
Blood phase pressure limit (post prime):	750 mmHg
Maximum water temperature:	42°C

Indications for Use

The MYOthern XP® Cardioplegia Delivery System is a device intended for the mixing, warming/cooling and delivery of oxygenated blood/cardioplegia solution in a predetermined ratio.

The Medtronic® MYOthern XP® Cardioplegia Delivery System (when using the MYOthern XP® CSS Custom Models) is intended only for mixing and delivery of oxygenated blood and/or asanguineous cardioplegia solution with the CSS instrument, Model 990. The MYOthern XP® Heat Exchanger (when part of the MYOthern XP® CSS Custom Models) is used with the instrument to warm, cool, and deliver oxygenated blood and/or asanguineous cardioplegia solution.

Contraindications

This device used for any other purposes than for the indicated use is the responsibility of the user.

Warnings

Read all Warnings, Precautions, and Instructions for Use carefully prior to use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**

- Only persons thoroughly trained in cardiopulmonary bypass procedures should use this device. Operation of each device requires constant supervision by qualified personnel for patient safety.
- Each device has been carefully manufactured, tested and packaged; however, the state of the art has not been developed to the point that Medtronic® is able to ensure that the device will not leak, crack, or fail during use. Perfusion must be carefully and constantly monitored.
- Each device is intended for single use only. Do not reuse or resterilize. EO sterilized.

- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use if the package is opened or damaged or the protective caps are not in place.
- Tubing should be attached in such a manner as to prevent kinks or restrictions that may alter blood or water flow.
- Do not allow alcohol, alcohol-based fluids; anesthetic fluids (such as isoflurane), or corrosive solvents (such as acetone) to come into contact with the device as they may jeopardize the structural integrity.
- If air bubbles and/or leaks are observed during priming and/or operation, these conditions may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the device if these conditions are observed.
- All gas emboli must be cleared from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- This system is intended for use with the MYOthem XP® Cardioplegia holder to allow for secure mounting to a vertical pole during operation. Inspection and maintenance of this holder should be performed routinely.
- Never restrict the water outflow of the heat exchanger.
- Direction of water flow through heat exchanger should be countercurrent with blood flow for optimal thermal efficiency.
- It is recommended that appropriate clinical laboratory procedures prior to use of this system examine the presence of cold agglutinins in the patient's blood. The temperature of the blood/cardioplegia solution should not be less than the lowest temperature at which negative cold agglutinins were determined.
- Never exceed water inlet pressure of 40 psi.
- Do not exceed blood side pressures of 750 mmHg post prime.
- The heat exchanger must be in a 10° inclined position as provided by the standard MYOthem XP® holder position at all times to help avoid the introduction of air into the delivery line.
- Disinfectants must not be used in the heater/cooler system when the heat exchanger is in use. If disinfectants are used in the heater/cooler system prior to use, the system must be thoroughly flushed.
- The heat exchanger is extremely efficient. It is recommended that water temperature never be allowed to exceed 42°C as thermal blood damage may result.
- The luer cap on the monitoring port of this device is VENTED.
- The cardioplegia pump must be fully occlusive.
- A FULLY OCCLUSIVE ROLLER HEAD must be used while operating this system. Failure to fully occlude both tubing segments within the roller head may result in backflow of cardioplegia solution into the oxygenator and inaccurate mixing ratios of blood to asanguineous cardioplegia solution.
- Do not start cardioplegia pump unless arterial pump is on.
- Do not stop arterial pump unless cardioplegia pump is off.
- Arterial flow must always exceed cardioplegia flow.
- Ensure that the system outflow line (patient delivery line) is not clamped while the pump is on.
- When using with membrane oxygenators the arterial pump must be delivering flow through the oxygenator at a flow greater than the flow through the cardioplegia circuit. This will help prevent the introduction of air into the perfusion circuit.
- Use of a bottle for cardioplegia solution requires the presence of a functional standard air vent in the bottle, which extends at least 2.5 cm (1 in) above the top of the cardioplegia solution fluid level.
- The asanguineous cardioplegia solution should be maintained at an appropriate level within the bag or bottle at all times to avoid air entry into the system. Air entry into the system may result in air embolism to the patient.
- The concentrations of additives to the asanguineous cardioplegia solution should be adjusted to achieve the desired concentrations in the blood/cardioplegia solution being delivered to the patient.
- If a stopcock has been inserted between the patient end of the delivery line and the cardioplegia cannula, care should be taken to ensure that the stopcock is correctly oriented. If the stopcock becomes closed during delivery of cardioplegia solution, back pressure could build-up resulting in bursting of the heat exchanger housing or tubing connections. The stopcock should be tested to ensure that acceptable resistance is generated at desired flow rate when used with a specific cannula. Pressure monitoring of the system during operation should be maintained to determine system performance and pressure drop generated with the cannula used for the procedure.
- Possible side effects include, but are not limited to infections, mechanical failure, hemolysis, air embolism, blood loss, circulatory compromise and thromboembolic phenomena. These are potential side effects with all extracorporeal blood systems.
- Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

NOTE: Additional warnings and cautions applicable to specific procedures can be found at appropriate places in these Instructions for Use.

Precautions

- Refer to package labeling for storage temperature requirements.
- Use aseptic technique in all procedures.
- A strict anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.

Warnings and Precautions (Specific to the CSS™ Custom Instrument)

- MYOtherm XP® CSS Custom Models are to be operated only with the Medtronic® Model 990 CSS instrument. There is no safety or performance data known to Medtronic that establishes compatibility of any other manufacturer's devices or components to the Medtronic system.
- MYOtherm XP® CSS Custom Models should only be placed in a conventional (standard) roller pump in an emergency situation. Refer to the CSS instrument Operator and Reference Manual for emergency procedures.
- Ensure that the CSS set tubing is properly loaded in the roller heads of the CSS instrument (see the CSS instrument Operator and Reference Manual for loading instructions). Potential consequences of a misloaded set include: backflow of cardioplegia solution into the oxygenator, inaccurate mixing ratios of blood to asanguineous cardioplegia solution, inaccurate flow and pressure readings, increased hemolysis, or spallation of the tubing material.
- If the dual-lumen patient line is used (CSS Custom Models), the transducer protector provided (with fluid barrier) must be used on the external pressure line to protect against fluid reaching the CSS instrument's pressure transducer. Inaccurate pressure readings may result if the transducer protector is not used.
- If the dual-lumen patient line is used (CSS Custom Models), overtightening of the luer lock may cause the luer lock to crack and give inaccurate pressure readings.

Installation and Setup of MYOtherm XP®

WARNING: The setup and use of the MYOtherm XP® Cardioplegia Delivery System is the responsibility of the attending clinician.

System Setup (Figure 2)

1. Remove the MYOtherm XP® Cardioplegia Delivery System components from the packaging carefully to ensure sterile fluid pathway.

WARNING: Ensure aseptic technique is used during all stages of setup and use of this system.

WARNING: Before removing the MYOtherm XP® Cardioplegia Delivery System from the packaging, inspect the packaging and product for damage. If the package or product is damaged, do not use as sterility of the device may have been compromised and/or performance may be affected.

2. Securely mount the MYOtherm XP® holder to an appropriate mast near the pump head used. Verify that this attachment is secure and the mast is capable of safely supporting the system's weight. Orient the mounting bracket so the outlet chamber is in full view during operation.

WARNING: If the holder is mounted to anything other than a vertical mast, care must be taken to ensure the MYOtherm XP® is at the proper angle (10°) to maintain optimal air handling capability and avoid air embolism.

3. Slide the heat exchanger into the heat exchanger holder arms. The holder will securely seat the heat exchanger when inserted correctly.
4. Securely connect the Hansen Series 3-ST connectors from the water source to the inlet and outlet ports of the heat exchanger as labeled.
5. Perform a water leakage check by circulating water through the heat exchanger while examining the blood path chamber for evidence of water. This procedure must be performed prior to priming this system.

Should there be any question of system integrity, DO NOT USE.

WARNINGS:

A. Do not restrict the water outflow of the heat exchanger.

B. Never exceed water inlet pressure of 40 psi.

C. As this heat exchanger is extremely efficient, it is recommended that water temperature never be allowed to exceed 42°C as thermal blood damage may result.

D. Do not exceed blood side pressures of 750 mmHg, post prime.

6. Install tubing in the pump head (using the appropriate head inserts) with the asanguineous cardioplegia solution tubing on top of the blood line within the roller pump. Adjust the pump head occlusion to achieve full occlusion of both tubing segments. Refer to the CSS™ Instrument's Operator and Reference manual for setup instructions if using one of the MYOtherm XP® Model 990 systems.

WARNING: A fully occlusive roller head must be used while operating this system. Failure to fully occlude both tubing segments within the roller head may result in cardioplegia solution entering the oxygenator circuit and inaccurate mixing ratios of blood to cardioplegia solution.

7. Confirm secure connection of the male luer-lock fitting of the pressure monitoring line to the pressure monitoring site at the top of the heat exchanger. Rotate the three-way stopcock to measure pressure from the heat exchanger.

NOTE: If the pressure monitoring site on the heat exchanger is not used, place a nonvented cap on this port.
8. Confirm secure connection of the pressure relief line to any available luer port on a vented cardiotomy reservoir.

NOTE: MYOthern XP® CSS Custom Models (designed for use with the CSS instrument) do not contain the pressure relief line. The CSS instrument has built-in systems to guard against overpressurization. The line connected to the pressure monitoring port may be connected to a port on a cardiotomy reservoir for priming/recirculation purposes.

NOTE: The pressure relief and pressure monitoring lines are preattached to the heat exchanger.
9. Attach the oxygenated blood line to the extracorporeal circuit and clamp securely prior to priming the extracorporeal circuit.
10. Securely close both clamps and insert the first IV spike of the cardioplegia line into the appropriate bag/bottle containing the cardioplegia/priming solution. If only one spike is to be utilized, the second line clamp should be moved down as close as possible to the “Y” Connector and securely clamped. If both spiked lines are to be used, the second line clamp should be moved up and securely clamped next to the spike to facilitate the prime.

WARNING: Use of a bottle for cardioplegia solution requires the presence of a functional standard air vent in the bottle that extends at least 2.5 cm (1 in) above the top of the cardioplegia solution fluid level.
11. Attach the temperature probe to the temperature monitoring adaptor at the outlet port of the heat exchanger.



Figure 2. System Setup

MYOtherm XP® for CSS Custom Models


Figure 2. System Setup (cont.)

Priming

1. Prior to priming the cardioplegia system, the oxygenator and extracorporeal tubing circuit should be primed, recirculated, and debubbled.
2. Priming of the cardioplegia system may be completed either before; or immediately after the commencement of cardiopulmonary bypass. If FINAL priming of the cardioplegia circuit is to be deferred until after the commencement of bypass, the system MUST be primed to a point beyond the roller pump or CSS instrument. This avoids inadvertent introduction of air into the oxygenator through the blood line from the cardioplegia system.

NOTE: When utilizing both spiked lines of the MYOthem XP® Cardioplegia System, retro fill second line up to second clamp and then close this clamp; spike second cardioplegia solution bag/bottle; gently tap spike and tubing to remove any bubbles; unclamp second line only when ready to administer.

2a. On MYOthem XP®-B Series only:

Ensure that all “bridge” clamps are opened momentarily to allow tubing bridge to be properly primed.

NOTE: Complete priming of this system with asanguineous pump prime will result in the delivery of approximately 190 mL of asanguineous cardioplegia solution before the direct delivery of blood/ cardioplegia solution.

3. To prime the heat exchanger system:

3a. Ensure that the stopcock is off on the pressure monitoring line.

NOTE: The extracorporeal circuit must be actively recirculated during priming of the cardioplegia system.

3b. Release clamp from blood line; rotate pump slowly; release clamp from cardioplegia line. Fill system to a maximum flow of 100 mL/min.

NOTE: Ensure that all bubbles are removed from both lines proximal to heat exchanger prior to wetting heat exchanger.

3c. Ensure the system pressure monitoring line and the pressure relief line are securely connected.

3d. Stop the pump when the heat exchanger is appropriately filled. Gently tap the heat exchanger to dislodge any remaining bubbles and pass them to the outlet chamber for removal.

DO NOT USE clamps or other tools to tap on the device to dislodge bubbles to avoid damaging the seals and housing.

3e. Turn the cardioplegia pump on at a low (<100 mL) flow and momentarily clamp outlet line to activate pressure relief valve until all visible air is purged through the pressure relief line and then stop pump. Open stopcock to vent any remaining air in heat exchanger.

NOTE: This step does not apply to MYOthem XP® CSS Custom Models.

3f. Securely attach pressure gauge to the transducer protector.

NOTE: This step does not apply to MYOthem XP® CSS Custom Models. The CSS instrument automatically detects and reports system internal pressure.

3g. After priming the system and checking the system to ensure it is air free, open the pressure monitoring line to the manometer by turning the stopcock to the appropriate position. Ensure the monitoring line is filled with fluid to the stopcock.

If this monitoring line is not used, turn the stopcock off on the pressure monitoring line.

4. Remove the tubing clamp from the system outflow line.

CAUTION: If outlet line IS NOT unclamped, pressure buildup >750 mmHg can occur causing the pressure relief valve to activate.

5. Prime the extension/delivery line by slowly advancing the system pump head. Ensure the delivery line is connected securely to the outflow line of the MYOthem XP® and that this connection is bubble free.

WARNINGS:

A. This system must be free of air bubbles prior to administration of cardioplegia to the patient.

B. The heat exchanger must be in a 10° inclined position as provided by the Standard MYOthem XP® holder position at all times to help avoid introduction of air into the delivery line.

C. The various components and connections of the system should be checked for leaks during priming and operation. Leakage of any connections or components may result in air embolism and/or solution loss.

6. After priming the complete system, verify that the pump has been properly occluded and all clamps have been removed prior to operation.

System Operating Procedure

1. Cold water (2-4°C) recirculation through the heat exchanger should be initiated 1 to 2 minutes before the delivery of cold blood/cardioplegia solution to the patient.

WARNING: If cold agglutinins have been demonstrated in the patient's blood, the blood temperature MUST be maintained at a level greater than that at which agglutination was demonstrated.

WARNING: When using with membrane oxygenators, the arterial pump must be delivering flow through the oxygenator at a flow greater than the flow through the cardioplegia circuit. This will help prevent the introduction of air into the perfusion circuit.

2. To initiate blood/cardioplegia solution delivery to the patient, remove any tubing clamps from the delivery line, start the system pump and adjust to the desired rate of infusion.
3. Stop the system pump to terminate the infusion of solution after achieving desired volume administration.

NOTE: A clamp should always be placed on the delivery line distal to the outlet chamber at the conclusion of each cardioplegia administration.

4. Cold water recirculation through the heat exchanger should be discontinued when the system pump has been stopped to avoid excessive cooling of the residual contents of the heat exchanger.
5. If additional blood/cardioplegia solution delivery is required, repeat Steps 1 through 4.

WARNING: The asanguineous cardioplegia solution should be maintained at an appropriate level within the bag or bottle at all times to avoid air entry into the system. Air entry into the system may result in air embolism to the patient.

6. Should replacement of the asanguineous cardioplegia solution bag (or bottle) become necessary:

6a. Ensure that the system administration pump IS NOT RUNNING.

6b. Clamp the asanguineous cardioplegia solution line between the supply bag/bottle and the pump head.

6c. Replace the bag/bottle with a full bag/bottle.

6d. Remove the tubing clamp from the asanguineous cardioplegia solution supply line.

6e. Clear the solution supply line of any entrapped air by forcing it back into the bag/bottle PRIOR TO restarting system pump flow.

WARNING: If a stopcock has been inserted between the patient end of the delivery line and the cardioplegia cannula, care should be taken to ensure that the stopcock is correctly oriented. If the stopcock becomes closed during delivery of cardioplegia solution, back pressure could build-up resulting in bursting of the heat exchanger housing or tubing connections. The stopcock should be tested to ensure that acceptable resistance is generated at desired flow rate when used with a specific cannula. Pressure monitoring of the system during operation should be maintained to determine system performance and pressure drop generated with the cannula used for the procedure. To prevent device damage and subsequent leakage, DO NOT EXCEED a pressure of 750 mmHg in the fluid path.

IMPORTANT NOTICE – LIMITED WARRANTY (U.S. CUSTOMERS ONLY)

- A. This **LIMITED WARRANTY** provides the following assurance to the patient who receives a Medtronic® MYOthern XP® Cardioplegia Delivery System, hereafter referred to as the “Product.”
- (1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its “Use By” date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the Purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
 - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its “Use By” date.
 - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
 - (3) The Product must not have been altered or subjected to misuse, abuse or accident.
 - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) Except as expressly provided by this **LIMITED WARRANTY**, Medtronic is not responsible for any direct, incidental or consequential damages based on any defect, failure or malfunction of the product, whether the claim is based on warranty, contract, tort or otherwise.
 - (2) THIS **LIMITED WARRANTY** IS MADE ONLY TO THE PATIENT IN WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS **LIMITED WARRANTY** SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
 - (3) The exclusions and limitation set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid. This **LIMITED WARRANTY** gives the patient specific legal rights. The patient may also have other rights which vary from state to state.
 - (4) No person has any authority to bind Medtronic to any representation, condition or warranty except this **LIMITED WARRANTY**.
- NOTE:** This **LIMITED WARRANTY** is provided by Medtronic, Inc. 710 Medtronic Parkway, Minneapolis, MN 55432-5604. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for exact terms of the **LIMITED WARRANTY**.

IMPORTANT NOTICE – LIMITED WARRANTY (FOR COUNTRIES OUTSIDE THE U.S.)

- A. This **LIMITED WARRANTY** provides assurance for the purchaser who receives a Medtronic® MYOthern XP® Cardioplegia Delivery System, hereafter referred to as the “Product”, that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.
The Warnings contained in the product labeling are considered an integral part of this **LIMITED WARRANTY**. Contact your local Medtronic® representative to obtain information on how to process a claim under this **LIMITED WARRANTY**.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its “Use By” date.
 - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
 - (3) The Product may not have been used for any other patient.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
 - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.

- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

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Alleviating Pain · Restoring Health · Extending Life

Europe

Europe/Africa/Middle East Headquarters

Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale 84
CH - 1131 Tolochenaz
Switzerland
Internet: www.medtronic.co.uk
Tel. 41-21-802-7000
Fax 41-21-802-7900

Medtronic E.C. Authorized Representative/Distributed by Medtronic B.V.

Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. 31-45-566-8000
Fax 31-45-566-8668

Asia-Pacific

Japan

Medtronic Japan
Comodio Shiodome 5F
2-14-1 Higashi-Shimbashi, Minato-ku
Tokyo 105-0021
Japan
Tel. 81-3-6430-2011
Fax 81-3-6430-7140

Australia

Medtronic Australasia Pty. Ltd.
Business Address:
97 Waterloo Road
North Ryde NSW 2113
Mailing Address:
PO Box 945
North Ryde NSW 1670
Australia
Tel. 61-2-9857-9000
Fax 61-2-9878-5100

Asia

Medtronic International Ltd.
Suite 1602 16/F, Manulife Plaza
The Lee Gardens, 33 Hysan Avenue
Causeway Bay
Hong Kong
Tel. 852-2891-4068
Fax 852-2591-0313

Americas

Latin America

Medtronic Latin America
3750 NW 87th Avenue
Suite 700
Miami, FL 33178
Tel. 305-500-9328
Fax 786-709-4244

Canada

Medtronic of Canada Ltd.
6733 Kitimat Road, Mississauga
Ontario L5N 1W3
Canada
Tel. 905-826-6020
Fax 905-826-6620
Toll-free in Canada:
1-800-268-5346

United States

Manufacturer:

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Internet: www.medtronic.com
Tel. 763-514-4000
Fax 763-391-9100
Toll-free in the USA:
1-800-328-2518
(24-hour consultation service)

Manufacturing Facility:

7611 Northland Drive
Minneapolis, MN 55428-9947
Tel. 763-391-9000
Fax 763-391-9100
Customer Service and Product Orders
Toll-free: 1-800-854-3570
Visit Medtronic at
www.perfusionsystems.com