Medtronic

MYOthem XP™

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MYOthem XP™

Cardioplegia Delivery System with Cortiva™ BioActive Surface

1. Description

The MYOthem XP™ cardioplegia delivery system with Cortiva™ bioactive surface is designed to mix arterial blood from an oxygenator with asanguineous cardioplegia solution in specific ratios, depending upon the system chosen. The specifications provided by each standard system and their corresponding product designations are presented below in Table 1. The patient delivery line is packaged separately to allow for ease of sterile transfer to the operative field prior to connection to the delivery outlet of the cardioplegia system. A pressure monitoring line, three-way stopcock, gauge protector, and a pressure relief line with pressure relief valve are preconnected. Products coated with Cortiva bioactive surface include a “CB” prefix in the model number. Additional information regarding Cortiva bioactive surface is available in Section 7.

Table 1.

| MYOthem XP system with Cortiva bioactive surface model | Delivery ratio ^a | Tubing inner diameter dimensions | | Pump type |
|--|-----------------------------|----------------------------------|-------------------|----------------------|
| | | Blood line | Cardioplegia line | |
| CBXP41, CBXP41B | 4:1 | 0.64 cm (1/4 in) | 0.32 cm (1/8 in) | Standard roller pump |

^a Note: CBXP41B (Bridge) series provides the ability to administer 100% blood as prescribed by physician.

1.1. Specifications

| | |
|---|---|
| Recommended blood flow rate | ≤500 mL/min |
| Filter screen | 150 μm |
| Heat exchanger nominal volume | 47 mL |
| Heat exchanger blood contact surface area | 0.06 m ² (95 in ²) |
| Size | Adult |
| Inlet connector | 0.64 cm (1/4 in) |
| Outlet connector | 0.48 cm (3/16 in) |
| Pressure monitor/pressure relief ports | Standard female luer lock |
| Waterports | Quick coupling connectors [1.3 cm (1/2 in)] |
| Water phase pressure limit | 275.8 kPa (40 psi) |
| Blood phase pressure limit (postprime) | 750 mmHg |
| Maximum water temperature | 42°C (107°F) |
| Maximum transit temperature | 40°C (104°F) |

Store the product at room temperature.

1.2. Accessories available separately

- Temperature probe TP
- MYOthem XP™ cardioplegia holder

2. Indications for use

The MYOthem XP with Cortiva bioactive surface is a device intended for the mixing, warming/cooling and delivery of oxygenated blood/cardioplegia solution in a predetermined ratio during routine cardiopulmonary bypass (CPB) procedures up to 6 hours in duration. Blood/cardioplegia solution is delivered to the patient through the cardioplegia delivery system and appropriate cannula by the operation of a single occlusive roller pump.

3. Contraindications

Use the device only as indicated.

4. Warnings

Read all warnings, precautions, and instructions for use carefully prior to use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**

- Only persons thoroughly trained in CPB procedures should use this device. Operation of each device requires constant supervision by qualified personnel for patient safety. Perfusion must be carefully and constantly monitored.
- Each device is intended for single use only. Do not reuse or resterilize. Sterilized using ethylene oxide.
- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use the device if the package is opened or damaged or the protective caps are not in place.
- Attach tubing in such a manner as to prevent kinks or restrictions that may alter blood or water flow.

- Do not allow alcohol, alcohol-based fluids, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come into contact with the device as they may jeopardize the structural integrity.
- If air bubbles and/or leaks are observed during priming and/or operation, these conditions may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the device if these conditions are observed.
- Clear all gas emboli from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- This system is intended for use with the MYOthem XP cardioplegia holder to allow for secure mounting to a vertical pole during operation. Inspection and maintenance of this holder should be performed routinely.
- Never restrict the water outflow of the heat exchanger.
- Direction of water flow through heat exchanger should be countercurrent with blood flow for optimal thermal efficiency.
- It is recommended that appropriate clinical laboratory procedures prior to use of this system examine the presence of cold agglutinins in the patient's blood. The temperature of the blood/cardioplegia solution should not be less than the lowest temperature at which negative cold agglutinins were determined.
- Never exceed water inlet pressure of 275.8 kPa (40 psi).
- Do not exceed blood side pressures of 750 mmHg postprime.
- The heat exchanger must be in a 10° inclined position, as provided by the standard MYOthem XP holder position, at all times to help avoid the introduction of air into the delivery line.
- Follow institution CPB protocol when using disinfectants in the heater/cooler during bypass. The integrity of the water path has been verified with hydrogen peroxide (330 ppm). Contact Medtronic for information regarding the use of additional disinfectants.
- The heat exchanger is extremely efficient. It is recommended that water temperature never be allowed to exceed 42°C (107°F) as thermal blood damage may result.
- The luer cap on the monitoring port of this device is vented.
- The cardioplegia pump must be fully occlusive.
- A fully occlusive roller head must be used while operating this system. Failure to fully occlude both tubing segments within the roller head may result in backflow of cardioplegia solution into the oxygenator and inaccurate mixing ratios of blood to asanguineous cardioplegia solution.
- Do not start the cardioplegia pump unless the arterial pump is on.
- Do not stop the arterial pump unless the cardioplegia pump is off.
- Arterial flow must always exceed cardioplegia flow.
- Ensure that the system outflow line (patient delivery line) is not clamped while the pump is on.
- When using with membrane oxygenators the arterial pump must deliver flow through the oxygenator at a flow greater than that through the cardioplegia circuit. This will help prevent the introduction of air into the perfusion circuit.
- Use of a bottle for cardioplegia solution requires the presence of a functional standard air vent in the bottle, which extends at least 2.5 cm (1 in) above the top of the cardioplegia solution fluid level.
- The asanguineous cardioplegia solution should be maintained at an appropriate level within the bag or bottle at all times to avoid air entry into the system. Air entry into the system may result in air embolization to the patient.
- The concentrations of additives to the asanguineous cardioplegia solution should be adjusted to achieve the desired concentrations in the blood/cardioplegia solution being delivered to the patient.
- If a stopcock has been inserted between the patient end of the delivery line and the cardioplegia cannula, care should be taken to ensure that the stopcock is correctly oriented. If the stopcock becomes closed during delivery of cardioplegia solution, back pressure could build up resulting in bursting of the heat exchanger housing or tubing connections. Test the stopcock to ensure that acceptable resistance is generated at desired flow rate when used with a specific cannula. Monitor system pressure during operation to determine how the pressure drop generated by the cannula used in the procedure affects system performance.

Note: Additional warnings and cautions applicable to specific procedures can be found at appropriate places in these instructions for use.

5. Precautions

- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize any component of the system. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Use aseptic technique in all procedures.
- Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

- It is the responsibility of the user to dispose of the devices in accordance with local regulations and hospital procedures.
- This product contains phthalates.

6. Adverse effects

The following known adverse effects are associated with the use of the product: coagulopathy; excessive blood component activation or thrombogenicity; hemolysis; hemolytic anemia; hypotension; infection; ischemia; neurological dysfunction; and organ dysfunction.

7. Information for products with Cortiva bioactive surface

The primary blood-contacting surfaces of the product are coated with Cortiva bioactive surface. This coated surface enhances blood compatibility and provides a blood-contacting surface that is thromboresistant. Cortiva bioactive surface contains nonleaching heparin derived from porcine intestinal mucosa¹.

Caution: Product coated with Cortiva bioactive surface is intended for single use only. Reesterilization may adversely affect the Cortiva bioactive surface.

Caution: Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

8. Instructions for use

8.1. System setup (refer to figure 2 and figure 3)

1. Carefully remove the MYOthem XP system with Cortiva bioactive surface components from the packaging to ensure sterile fluid pathway.

Warning: Ensure aseptic technique is used during all stages of setup and use of this system.

Warning: Before removing the MYOthem XP system with Cortiva bioactive surface from the packaging, inspect the packaging and product for damage. If the package or product is damaged, do not use as sterility of the device may have been compromised and/or performance may be affected.

2. Securely mount the MYOthem XP holder to an appropriate mast near the pump head used. Verify that this attachment is secure and the mast is capable of safely supporting the system's weight. Orient the mounting bracket so the outlet chamber is in full view during operation.

Warning: If the holder is mounted to anything other than a vertical mast, care must be taken to ensure the MYOthem XP system with Cortiva bioactive surface is at the proper angle (10°) to maintain optimal air handling capability and avoid air embolism.

3. Slide the heat exchanger into the heat exchanger holder arms. The holder arms will securely seat the heat exchanger when inserted correctly.
4. Secure the quick coupling connectors from the water source to the inlet and outlet ports of the heat exchanger as labeled.
5. Perform a water leakage check by circulating water through the heat exchanger while examining the blood path chamber for evidence of water. This procedure must be performed prior to priming this system. Should there be any question of system integrity, do not use.

Warning:

- a. Do not restrict the water outflow of the heat exchanger.
 - b. Never exceed water inlet pressure of 275.8 kPa (40 psi).
 - c. As this heat exchanger is extremely efficient, it is recommended that water temperature never be allowed to exceed 42°C (107°F) as thermal blood damage may result.
 - d. Do not exceed blood side pressures of 750 mmHg postprime.
6. Install tubing in the pump head (using the appropriate head inserts) with the asanguineous cardioplegia solution tubing on top of the blood line within the roller pump. Adjust the pump head occlusion to achieve full occlusion of both tubing segments.

Warning: A fully occlusive roller head must be used while operating this system. Failure to fully occlude both tubing segments within the roller head may result in cardioplegia solution entering the oxygenator circuit and inaccurate mixing ratios of blood to cardioplegia solution.

7. Confirm secure connection of the male luer lock fitting of the pressure monitoring line to the pressure monitoring site at the top of the heat exchanger. Rotate the three-way stopcock to measure pressure from the heat exchanger.

Note: If the pressure monitoring site on the heat exchanger is not used, place a nonvented cap on this port.

8. Confirm secure connection of the pressure relief line to any available luer port on a vented cardiotomy reservoir. The line connected to the pressure monitoring port may be connected to a port on a cardiotomy reservoir for priming/recirculation purposes.

Note: The pressure relief and pressure monitoring lines are preattached to the heat exchanger.

¹ Nonleaching is defined as heparin at a level ≤ 0.1 IU/mL, as measured through clinically relevant extraction conditions.

9. Attach the oxygenated blood line to the extracorporeal circuit and clamp securely before priming the extracorporeal circuit.
10. Securely close both clamps, and insert the first IV spike of the cardioplegia line into the appropriate bag/ bottle containing the cardioplegia/priming solution. If only one spike is to be utilized, the second line clamp should be moved down as close as possible to the “Y” connector and securely clamped. If both spiked lines are to be used, the second line clamp should be moved up and securely clamped next to the spike to facilitate the prime.

Warning: Use of a bottle for cardioplegia solution requires the presence of a functional standard air vent in the bottle that extends at least 2.5 cm (1 in) above the top of the cardioplegia solution fluid level.
11. Attach the temperature probe to the temperature monitoring adaptor at the outlet port of the heat exchanger.

8.2. Priming

1. Prior to priming the cardioplegia system, the oxygenator and extracorporeal tubing circuit should be primed, recirculated, and debubbled.
2. Priming of the cardioplegia system may be completed either before or immediately after the commencement of cardiopulmonary bypass. If final priming of the cardioplegia circuit is to be deferred until after the commencement of bypass, the system must be primed to a point beyond the roller pump. This avoids inadvertent introduction of air into the oxygenator through the blood line from the cardioplegia system.

Note: When utilizing both spiked lines of the MYOthem XP system with Cortiva bioactive surface: using the primed line, retrograde prime the unprimed line until all air is removed, then clamp; insert the remaining spike into the second cardioplegia solution bag/bottle; gently tap the spike and the tubing to remove any bubbles; unclamp the second line only when ready to administer.

On bridge units only:

Ensure that all bridge clamps are opened momentarily to allow tubing bridge to be properly primed.

Note: Complete priming of this system with asanguineous pump prime will result in the delivery of approximately 190 mL of asanguineous cardioplegia solution before the direct delivery of blood/cardioplegia solution.

3. To prime the heat exchanger system:
 - a. Ensure that the stopcock on the pressure monitoring line is turned to the off position.

Note: The extracorporeal circuit must be actively recirculated during priming of the cardioplegia system.
 - b. Release the clamp from the blood line; rotate the pump slowly; release the clamp from the cardioplegia line. Fill the system to a maximum flow of 100 mL/min.

Note: Ensure that all bubbles are removed from both lines proximal to the heat exchanger before wetting the heat exchanger.
 - c. Ensure that the system pressure monitoring line and the pressure relief line are securely connected.
 - d. Stop the pump when the heat exchanger is appropriately filled. To facilitate complete debubbling of the outlet filter media, it is recommended that the heat exchanger be removed from the holder, inverted, and inspected for trapped air. Gently tap the heat exchanger to dislodge any remaining bubbles and pass them to the outlet chamber for removal.

Do not use clamps or other tools to tap on the device to dislodge bubbles to avoid damaging the seals and housing.
 - e. Turn the cardioplegia pump on at a low (<100 mL) flow, and momentarily clamp outlet line to activate the pressure relief valve until all visible air is purged through the pressure relief line, and then stop the pump. Open the stopcock to vent any remaining air in the heat exchanger.
 - f. Securely attach the pressure gauge to the transducer protector.
 - g. After priming the system and checking the system to ensure it contains no air bubbles, open the pressure monitoring line to the manometer by turning the stopcock to the appropriate position. Ensure that the monitoring line is filled with fluid to the stopcock. If this monitoring line is not used, turn the stopcock on the pressure monitoring line to the off position.
4. Remove the tubing clamp from the system outflow line.

Caution: If the outlet line is not unclamped, pressure buildup >750 mmHg can occur causing the pressure relief valve to activate.
5. Prime the extension/delivery line by slowly advancing the system pump head. Ensure that the delivery line is connected securely to the outflow line of the MYOthem XP system with Cortiva bioactive surface and that this connection is bubble-free.

Warning:

 - a. This system must be free of air bubbles before administration of cardioplegia to the patient.
 - b. The heat exchanger must be in a 10° inclined position, as provided by the standard MYOthem XP holder position, at all times to help avoid introduction of air into the delivery line.

- c. The various components and connections of the system should be checked for leaks during priming and operation. Leakage of any connections or components may result in air embolism and/or solution loss.
6. After priming the complete system, verify that the pump has been properly occluded and all clamps have been removed before operation.

8.3. System operating procedure

1. Cold water (2°C to 4°C) recirculation through the heat exchanger should be initiated 1 to 2 minutes before the delivery of cold blood/cardioplegia solution to the patient.
Warning: If cold agglutinins have been demonstrated in the patient's blood, the blood temperature must be maintained at a level greater than that at which agglutination was demonstrated.
Warning: When used with membrane oxygenators, the arterial pump must deliver flow through the oxygenator at a flow greater than that through the cardioplegia circuit. This will help prevent the introduction of air into the perfusion circuit.
2. To initiate blood/cardioplegia solution delivery to the patient, remove any tubing clamps from the delivery line, start the system pump, and adjust to the desired rate of infusion.
3. Stop the system pump to terminate the infusion of solution after achieving desired volume administration.
Note: A clamp should always be placed on the delivery line distal to the outlet chamber at the conclusion of each cardioplegia administration.
4. Cold water recirculation through the heat exchanger should be discontinued when the system pump has been stopped to avoid excessive cooling of the residual contents of the heat exchanger.
5. If additional blood/cardioplegia solution delivery is required, repeat steps 1 through 4.
Warning: The asanguineous cardioplegia solution should be maintained at an appropriate level within the bag or bottle at all times to avoid air entry into the system. Air entry into the system may result in air embolism to the patient.
6. Should replacement of the asanguineous cardioplegia solution bag (or bottle) become necessary, perform the following steps:
 - a. Ensure that the system administration pump is not running.
 - b. Clamp the asanguineous cardioplegia solution line between the supply bag/bottle and the pump head.
 - c. Replace the bag/bottle with a full bag/bottle.
 - d. Remove the tubing clamp from the asanguineous cardioplegia solution supply line.
 - e. Clear the solution supply line of any entrapped air by forcing it back into the bag/bottle before restarting system pump flow.

Warning: If a stopcock has been inserted between the patient end of the delivery line and the cardioplegia cannula, care should be taken to ensure that the stopcock is correctly oriented. If the stopcock becomes closed during delivery of cardioplegia solution, back pressure could build up resulting in bursting of the heat exchanger housing or tubing connections. Test the stopcock to ensure that acceptable resistance is generated at desired flow rate when used with a specific cannula. Monitor system pressure during operation to analyze how the pressure drop generated by the cannula used in the procedure affects system performance. To prevent device damage and subsequent leakage, do not exceed a pressure of 750 mmHg in the fluid path.

9. Additional Information

The following information and data are available upon request: test methods used to determine coverage, leaching, and bioactivity.

10. Important notice—limited warranty (US customers only)

- A. This **LIMITED WARRANTY** provides the following assurance to the patient who receives a MYOtherm XP system with Cortiva bioactive surface, hereafter referred to as the "Product."
 - (1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its Use-By Date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the Purchase Price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
 - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
 - (1) The Product must be used prior to its Use-by date.
 - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
 - (3) The Product must not have been altered or subjected to misuse, abuse, or accident.
 - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.

- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) Except as expressly provided by this **LIMITED WARRANTY**, Medtronic is not responsible for any direct, incidental, or consequential damages based on any defect, failure, or malfunction of the product, whether the claim is based on warranty, contract, tort, or otherwise.
 - (2) THIS **LIMITED WARRANTY** IS MADE ONLY TO THE PATIENT IN WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS **LIMITED WARRANTY** SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
 - (3) The exclusions and limitation set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid. This **LIMITED WARRANTY** gives the patient specific legal rights. The patient may also have other rights which vary from state to state.
 - (4) No person has any authority to bind Medtronic to any representation, condition, or warranty except this **LIMITED WARRANTY**.

Note: This **LIMITED WARRANTY** is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for exact terms of the **LIMITED WARRANTY**.

11. Important notice—limited warranty (for countries outside the US)

Important: This limited warranty does not apply in Australia.

- A. This **LIMITED WARRANTY** provides assurance for the purchaser who receives a MYOtherm XP system with Cortiva bioactive surface, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.

The Warnings contained in the product labeling are considered an integral part of this **LIMITED WARRANTY**. Contact your local Medtronic representative to obtain information on how to process a claim under this **LIMITED WARRANTY**.

- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its Use-by date.
 - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
 - (3) The Product may not have been used for any other patient.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation, or material alteration of the replaced Product.
 - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

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© 2018 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
Medtronic B.V.
Medtronic of Canada Ltd