
Quantu 
PureFlow







Standard heat exchanger

Spectrum  Medical

EN	ENGLISH	Instructions for use	1
IT	ITALIANO	Istruzioni per l'uso	8
FR	FRANÇAIS	Mode d'emploi	15
DE	DEUTSCH	Gebrauchsanleitung	22
ES	ESPAÑOL	Instrucciones de uso	29
NL	NEDERLANDS	Gebruiksaanwijzing	36
SV	SVENSKA	Bruksanvisning	43
DA	DANSK	Brugsanvisning	50
FI	SUOMI	Käyttöohjeet	57
NO	NORSK	Bruksanvisning	64
CS	ČEŠTINA	Návod k použití	71
PL	POLSKI	Instrukcje dotyczące użytkowania	78

EN - Heat Exchanger Performance	IT - Performace dello scambiatore di calore
FR - Performances de l'échangeur thermique	DE - Leistung des Wärmetauschers
ES - Rendimiento del intercambiador de calor	NL - Prestaties van de warmtewisselaar
SV - Prestanda för värmeväxlare	DA - Varmevekslerens ydeevne
FI - Lämmönvaihtimen suorituskyky	NO - Ytelse til varmeveksler
CS - Výkon tepelného výměníku	PL - Wydajność wymiennika ciepła

Quantum PureFlow Standard Heat Exchangers

HX-S

REF.	DESCRIPTION	CODE	DESCRIPTION
HX55V-S0	Standard Heat Exchanger High Flow 3/8	HX11V-S0	Standard Heat Exchanger Low Flow 3/16
HX22V-S0	Standard Heat Exchanger Medium Flow 1/4	HX11V-S1	Standard Heat Exchanger Low Flow 1/4

Standard Heat Exchanger

CONTENTS

- A. WARNINGS
- B. CONTRAINDICATIONS AND ADVERSE REACTIONS
- C. DESCRIPTION
- D. TECHNICAL CHARACTERISTICS
- E. INTENDED USE
- F. DEVICE PREPARATION AND SET-UP
- G. DEVICE PRIMING PROCEDURE
- H. DEVICE USE AND REMOVAL
- I. DEVICE REPLACEMENT
- J. MEDICAL DEVICES TO BE USED WITH THE DEVICE
- K. RETURN OF USED DEVICES AND NOTIFICATIONS
- L. LIMITED WARRANTY
- M. SYMBOLS USED

A. WARNINGS

CAUTION

- Carefully read the instructions before use.
- The device must be used in accordance with these instructions for use. QURA cannot be held responsible for damages resulting from incorrect use; the user is fully responsible for improper use or use different from that indicated.
- The device may only be used by adequately trained and qualified personnel.
- The device is sterilized by ethylene oxide and is sold STERILE. Sterility is guaranteed only when the package is intact. Do not use the device if the package is damaged, tampered with, open or wet.
- Visually inspect and carefully check the device before use. Inappropriate transport and/or storage conditions not compliant with the instructions may have damaged the product.
- Fragile, handle with care.
- Keep dry. Store in a dry place between 0°C/32°F and 60°C/140°F.
- Single-use product to be used only once and on one single patient. Reuse (on the same person or on other patients) may result in cross-contamination, infection and sepsis after the device comes into contact with blood/body fluids or any substances infused into the patient during the procedure. Reusing the product could also compromise its characteristics in terms of performance and safety.
- Do not resterilize, subject to further processing or clean for reuse.
- The device is not made with natural rubber latex and does not contain phthalates.
- The device is non-toxic and non-pyrogenic.
- Check the expiry date printed on the label. Do not use the device beyond this date.
- In order to avoid product contamination, the device must be used immediately after opening the sterile packaging and handled aseptically.
- The device must be used only in combination with the medical devices listed in section J "Medical devices to be used with the device" in these instructions for use. Carefully read the user manuals of the medical devices listed in section J before using the device.
- Do not allow alcohol, alcohol-based liquids, ethers, halogenated liquids (such as Halothane), or corrosive solvents (such as acetone) to come into contact with the device. This would cause damage such as to compromise the integrity and functionality of the device.
- The connection lines must be properly connected to prevent tube kinking or occlusion which could lead to a reduction in blood flow and/or an increase in circuit pressure.
- It is advisable to externally secure the tubes on the herringbone connectors of the device using cable-ties.
- Continuously and carefully monitor perfusion.
- A spare device must always be available during the procedure.
- Follow an anticoagulation protocol and monitor anticoagulation throughout the procedure. The benefit from using an extracorporeal support must however be evaluated against the risk of systemic anticoagulation and must be supported by a medical prescription.
- Frequently monitor both the patient and the device; do not leave the device unattended while it is running.
- After use, dispose of the device in accordance with the applicable regulations in force in the country of use.
- In the event of a malfunction, always keep the defective device and follow the instructions given in section J "Return of used devices and notifications" if it needs to be returned.
- For further information and/or in case of complaint, contact Qura S.r.l. or the distributor of the Qura product.
- For US market - Caution: Federal law restricts this device to sale by or on the order of a physician.

B. CONTRAINDICATIONS AND ADVERSE REACTIONS

The device may not be used for any purpose other than its intended use. The user is responsible for any use different from that intended. The following adverse reactions are possible: allergic reaction, haemorrhage, blood toxicity, embolism, excessive activation of blood components, bleeding, fever, hemolysis, hypercapnia and acidosis, hypotension, hypovolemia, hypoxia, compromised patient temperature control, infection, neurological dysfunction and stroke, organ dysfunction, renal failure, reproductive harm in male neonatal and paediatric patients, respiratory and pulmonary problems, sepsis, sterility, thrombosis, toxemia, carcinogenicity, mutagenicity, sensitization, cytotoxicity, vascular damage, death.

C. DESCRIPTION

The heat exchanger (Fig. 1) allows heating/cooling the blood during extracorporeal circulation. The device is available in the versions listed in Table 1 – section D. The device is sold sterile in a blister pack in a multiple box of 4 units.

All the device surfaces in contact with blood are treated with a phosphorylcholine-based biocompatible coating. Contact with blood for longer than 6 hours is not permitted. The device has been designed to be powered by heaters coolers that use a glycol-based heat transfer fluid (HTF).

The device must be used in combination with the medical devices listed in section J in these instructions for use.


Fig. 1

References:

- | | |
|-------------------------|---------------------------|
| 1. HTF inlet connector | 3. Blood inlet connector |
| 2. HTF outlet connector | 4. Blood outlet connector |

D. TECHNICAL CHARACTERISTICS

Characteristics	CODE				
	HX55V-S0	HX22V-S0	HX11V-S0	HX11V-S1	
	High Flow	Medium Flow	Low Flow	Low Flow	
Blood contact materials	Polycarbonate AISI 304 stainless steel sheet Epoxy resin Coating: phosphorylcholine				
Maximum operating pressure - blood side	750 mmHg / 100 kPa / 1 bar / 14.5 psi				
Maximum operating pressure - HTF side	1500 mmHg / 200 kPa / 2 bar / 29 psi				
Maximum HTF compartment supply temperature	42°C				
Static priming volume (ml)	55	22	11	11	
Exchange surface area (m ²)	0.12	0.055	0.028	0.028	
Maximum blood flow (l/min)	8.0	3.0	1.0	1.0	
Connections - blood side					
	IN	3/8" (9.525 mm)	1/4" (6.35 mm)	3/16" (4.76 mm)	1/4" (6.35 mm)
	OUT	3/8" (9.525 mm)	1/4" (6.35 mm)	3/16" (4.76 mm)	1/4" (6.35 mm)
Connections - HTF side (inlet/outlet)	CPC HFC profile	CPC HFC profile	CPC APC profile	CPC APC profile	

Table 1
E. INTENDED USE

The Standard Heat Exchanger has been designed to heat/cool blood during surgical procedures that require extracorporeal circulation and cardiopulmonary bypass for periods of maximum 6 hours.

The device is also indicated to be used in extracorporeal support systems (for periods of not more than 6 hours) that do not require cardiopulmonary bypass (e.g. valvuloplasty, surgery of the vena cava and the aorta, liver transplants).

F. DEVICE PREPARATION AND SET-UP

CAUTION

- A spare device must always be available during extracorporeal circulation.
- Check that the O-Ring is present on the HTF connectors (ref. 3c and 3d in Fig.2).
- The device has been designed to be used only with the dedicated holder provided by Spectrum Medical.
- During circuit set-up and before starting extracorporeal circulation, check that the spare medical devices and accessories indicated in the manual of the medical devices listed in section J are available.
- Be careful not to damage the connectors when connecting the tubes.
- Follow a suitable blood anticoagulation protocol before and during the procedure.
- The ACT (activated coagulation time) value must always be such as to guarantee anticoagulation suited to the procedure.

F.1. PREPARING THE HEATER-COOLER

CAUTION

- The pressure of the HTF solution inside the device must not exceed the value indicated in Table 1 - section D.
- The heat exchanger temperature must not exceed the value indicated in Table 1 - section D.

1. Prepare the quick-coupling tubes provided with the Spectrum Medical Quantum Heater-Cooler System.

F.2. INSTALLING THE HEAT EXCHANGER

1. Remove the device from the sterile packaging.
2. Fit the device on the dedicated holder provided by Spectrum Medical. For more details, refer to the user manual of the Spectrum Medical Quantum HLM heart-lung machine (Quantum Pump Console User Guide).
3. Make sure the device is firmly attached to the holder.
4. Remove the protective caps from the connectors on the HTF side of the device (ref. 3c and 3d in Fig.2).
5. Connect the tubes of the Spectrum Medical Quantum Heater-Cooler System (ref. 4 in Fig. 2) to the device (ref. 3c and 3d in Fig. 2) using the female quick-coupling connectors.

F.3. CHECKING THE INTEGRITY OF THE DEVICE

1. Start priming of the HTF circuit following the instructions given in the user manual of the Spectrum Medical Quantum Heater-Cooler System.
2. The integrity of the device is guaranteed if:
 - there are no HTF solution leaks from its compartment to the outside environment;
 - there are no HTF solution leaks in the blood compartment.

F.4. CONNECTING THE HEAT EXCHANGER TO THE EXTRACORPOREAL CIRCUIT

CAUTION

The maximum pressure of the blood compartment must not exceed the value indicated in Table 1 - section D.

Only use tubes with dimensions compatible with the blood inlet/outlet connectors of the device (see Table 1).

1. Connect the device to the extracorporeal circuit as shown in Figure 2.
2. Place the device between the arterial pump (ref. 2 in Fig. 2) and the oxygenator inlet (ref. 5 in Fig. 2) as shown in Fig. 2.
3. Remove the protective cap from the blood inlet connector (ref. 3a in Fig. 2).
4. Connect the appropriate tube to the inlet port of the device using an aseptic technique in order to maintain sterility.
5. Remove the protective cap from the blood outlet connector (ref. 3b in Fig. 2).
6. Connect the appropriate tube to the outlet port of the device using an aseptic technique in order to maintain sterility.


Fig. 2

References:

- | | | |
|----------------------------|---------------------------------|----------------------------|
| 1. Venous reservoir | 3c. HTF inlet connector | 5b. Blood outlet connector |
| 2. Arterial pump | 3d. HTF outlet connector | 5c. Water inlet connector |
| 3. Standard heat exchanger | 4. Quantum Heater-Cooler System | 5d. Water outlet connector |
| 3a. Blood inlet connector | 5. Oxygenator | 6. Arterial line |
| 3b. Blood outlet connector | 5a. Blood inlet connector | |

G. **DEVICE PRIMING PROCEDURE**

G.1. PRIMING THE CIRCUIT

CAUTION

1. Begin priming the device at a constant flow of 40-50 ml/min.
2. As soon as the fluid flows from the outlet connector, the priming flow can be increased up to a value compatible with the size of the device used (see Table 1 - section D).

G.2. REMOVING AIR BUBBLES

CAUTION

- The most effective technique to remove air bubbles is to recirculate fluid at a high flow in a circuit that allows this operation.
- During recirculation, the device must be kept in vertical position (i.e. Blood IN at the bottom and Blood OUT at the top).
- Do not use inappropriate tools to remove any air bubbles, as they may damage the device and alter its operation.

1. Facilitate air removal directing it towards the outlet and always hold the outlet facing up.
2. This operation must be done for a few minutes. The operation is successful when no more air bubbles come out of the blood outlet connector.

H. DEVICE USE AND REMOVAL

CAUTION

- To reduce the risk of air reaching the patient, make sure that the device and the circuit are free of air and check that they have been filled properly before starting the procedure. It is advisable to use an arterial filter downstream of the device.
- Visually check that all the air has been removed from the circuit.
- Carefully check for the absence of potential occlusions in the circuit before starting extracorporeal circulation.
- In case of leaks or other faults, remove the device in accordance with what is indicated in these instructions for use.
- Prolonged use of the device with stationary flow may cause the formation of clots.

H.1. USING THE DEVICE DURING THE PROCEDURE

CAUTION

- Monitor the blood temperature during the procedure.
- Monitor the HTF solution temperature during the procedure. For more details, refer to the user manual of the Spectrum Medical Quantum Heater-Cooler System.

H.2. REMOVING THE DEVICE

CAUTION

- Remove and dispose of the device after each procedure in accordance with what is indicated in section A. Follow the instructions given in section K if used devices need to be returned.
1. Empty the HTF circuit following the instructions reported in the user manual of the Spectrum Medical Quantum Heater-Cooler System.
 2. Disconnect the tubes of the Spectrum Medical Quantum Heater-Cooler System.
 3. Remove the device from the dedicated holder provided by Spectrum Medical.
 4. Remove the device from the circuit.

I. DEVICE REPLACEMENT

CAUTION

After six hours of use or if situations occur that, based on the judgement of the person responsible for perfusion, compromise the safety conditions for the patient (insufficient performance, loss of integrity of the HTF compartment or the blood compartment, etc.), proceed as follows to replace the device.

- If the device needs to be replaced during circulatory support, evaluate the potential risks for the patient before replacing the device.
 - Never disconnect the lines if the circuit is pressurized.
 - Visually inspect and carefully check the device before replacing it.
 - Clamp the venous and arterial lines before proceeding with device replacement.
 - Follow the instructions given in section K if the device needs to be returned.
1. Empty the HTF circuit following the instructions reported in the user manual of the Spectrum Medical Quantum Heater-Cooler System.
 2. Disconnect the tubes of the Spectrum Medical Quantum Heater-Cooler System.
 3. Remove the new device from the outer packaging and from the sterile packaging checking its integrity.
 4. Stop the pump following the instructions in the user manual of the Spectrum Medical Quantum HLM heart-lung machine (Quantum Pump Console User Guide).
 5. Place two clamps on both the blood inlet and outlet line of the device.
 6. Remove the device from the dedicated Spectrum Medical holder.
 7. Cut the device inlet and outlet lines at the point between the two clamps taking care to keep adequate lengths for reconnection.
 8. Position the new device on the dedicated Spectrum Medical holder.
 9. Fill the HTF compartment of the device as described in section F.2 and check its integrity as described in section F.3.
 10. Start the Spectrum Medical Quantum Heater-Cooler System.
 11. Connect the blood inlet line to the blood inlet connector and the blood outlet line to the blood outlet connector of the new device.
 12. Use cable-ties to externally lock the inlet and outlet tubes on the respective connectors of the device.
 13. Fill the device.
 14. Check that there is no air in the device.
 15. If there is air in the device, follow the instructions given in section G.2.
 16. Open the outlet line only after having checked that there is no air in the device.
 17. Open the venous and arterial lines.
 18. Resume perfusion.

J. MEDICAL DEVICES TO BE USED WITH THE DEVICE

CAUTION

- The user must respect the warnings and precautions and follow the instructions for use provided with the device.
- The blood circuit connections must be made using tubes of a diameter compatible with the dimensions of the connectors on the device in accordance with Table 1 - section D.

The device is qualified to be used exclusively with:

- Spectrum Medical Quantum Heater-Cooler System.
- Spectrum Medical Quantum HLM heart-lung machine and accessories.

K. RETURN OF USED DEVICES AND NOTIFICATIONS

CAUTION

- Always keep any defective medical devices. Follow the instructions below if the devices need to be sent for investigation.
- Healthcare institutions are responsible for preparing and identifying the products for return shipments. If the product to be returned has been in contact with the patient's blood, it must be treated, packed and handled in compliance with the provisions of the legislation in force in the country where the product has been used. Do not send products that have been contaminated by infectious pathogenic agents.
- Qura stresses the importance of reporting any incident that occurred regarding the device to Qura and to the competent authority of the Member State where the product was used. Should the user be dissatisfied with anything related to the quality of the product, the user may notify Qura or the distributor of the Qura product. All notifications considered critical by the user must be reported with particular care and urgency. Below is the minimum information to be provided:
 - thorough description of the event and, if pertinent, the conditions of the patient;
 - identification of the product involved;
 - batch number of the product involved;
 - availability of the product involved;
 - all the indications the user considers useful in order to understand the origin of the dissatisfaction.
- Qura reserves the right to authorize the return of the product involved in the notification for assessments.

L. LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the purchaser pursuant to applicable law. Qura guarantees that all reasonable precautions have been taken during the manufacture of this medical device, as required by the nature of the device and its intended use. Qura guarantees that its medical device is capable of functioning as indicated in these Instructions For Use when used in accordance with that specified in these Instructions by qualified personnel and within the expiry date indicated on the packaging. However, Qura cannot guarantee that the user will use the device correctly, nor that incorrect diagnosis or therapy and/or the particular physical and biological characteristics of an individual patient do not affect the performance and effectiveness of the device or cause damaging consequences for the patient, even though the specified Instructions for Use have been respected. Whilst emphasising the need to adhere strictly to the Instructions for Use and to adopt all the precautions necessary for correct use of the device, Qura therefore cannot assume any responsibility for any incidents, harm, loss, expense or consequences arising directly or indirectly from improper use of, tampering with or modifications made to this device. Qura will replace the medical device if it is defective at the time it is put on the market or while being shipped by Qura up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties, explicit or implicit, written or verbal, including warranties of marketability and/or functionality. No person, including any representative, agent, dealer, distributor or intermediary of Qura or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. Qura disclaims any variations to this Limited Warranty and to the information/Instructions for Use as explicitly contained herein. The purchaser pledges to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with Qura, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The relationship between the parties regarding the contract (even if not stipulated in writing) to whom this warranty is issued, as well as any dispute related to it or in any way connected with it, as well as any relationship or dispute concerning this warranty, its interpretation and execution, nothing excluded and/or reserved, are governed exclusively by Italian law and jurisdiction. The court appointed in the event of the above-said disputes shall exclusively be the Court of Modena (Italy).

Qura S.r.l.
Via di Mezzo, 23
41037 Mirandola (MO), Italy
Tel. +39 0535-1803050
www.quramed.com

M. SYMBOLS USED

	Manufacturer		Temperature limit
	Use-by Date		Non-pyrogenic
	Batch code		For US Audiences Only
	Catalogue number		Not made with natural rubber latex
	Sterilized using ethylene oxide		Serial number
	Do not re-sterilize		This way up
	Do not use if the packaging is damaged		Fragile, handle with care
	Do not re-use		Keep dry
	Consult the instructions for use		Keep away from sunlight
	Caution		Open here
	Compliance with Directive 93/42/EC		Packaging unit

Quantum PureFlow Standard Heat Exchangers

HX-S

REF.	DESCRIPTION	CODE	DESCRIPTION
HX55V-S0	Standard Heat Exchanger High Flow 3/8	HX11V-S0	Standard Heat Exchanger Low Flow 3/16
HX22V-S0	Standard Heat Exchanger Medium Flow 1/4	HX11V-S1	Standard Heat Exchanger Low Flow 1/4

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G.2. REMOVING AIR BUBBLES

CAUTION

- The most effective technique to remove air bubbles is to recirculate fluid at a high flow in a circuit that allows this operation.
 - During recirculation, the device must be kept in vertical position (i.e. Blood IN at the bottom and Blood OUT at the top).
 - Do not use inappropriate tools to remove any air bubbles, as they may damage the device and alter its operation.
1. Agevolare l'evacuazione dell'aria indirizzandola verso l'uscita e tenendo quest'ultima sempre orientata verso l'alto.
 2. Tale operazione deve protrarsi per alcuni minuti. L'indicazione della riuscita dell'operazione è data dalla cessata fuoriuscita di bolle d'aria dal connettore di uscita sangue.

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Medtronic Australasia Pty Ltd
Qura S.r.l. - Via di Mezzo, 23 - 41037 Mirandola (MO) Italy