# Medtronic

## Affinity Fusion®

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**USA** Caution: Federal law (USA) restricts this device to sale by  
or on the order of a physician.

Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup>, and Affinity Fusion<sup>®</sup> are registered trademarks of Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> et Affinity Fusion<sup>®</sup> sont des marques commerciales déposées de Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> und Affinity Fusion<sup>®</sup> sind eingetragene Marken von Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> y Affinity Fusion<sup>®</sup> son marcas comerciales registradas de Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> en Affinity Fusion<sup>®</sup> zijn geregistreerde handelsmerken van Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> e Affinity Fusion<sup>®</sup> sono marchi registrati della Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> og Affinity Fusion<sup>®</sup> er registrerte varemærker for Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> ja Affinity Fusion<sup>®</sup> ovat Medtronic, Inc:n rekisteröityjä tavaramerkkejä.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> och Affinity Fusion<sup>®</sup> är registrerade varumärken som tillhör Medtronic, Inc.  
Οι επωνυμίες Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> και Affinity Fusion<sup>®</sup> είναι σήματα κατατεθέντα της Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> og Affinity Fusion<sup>®</sup> er registrerede varemærker tilhørende Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> i Affinity Fusion<sup>®</sup> są zastrzeżonymi znakami towarowymi firmy Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> e Affinity Fusion<sup>®</sup> são marcas comerciais registadas da Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> ve Affinity Fusion<sup>®</sup> Medtronic, Inc. şirketinin tescilli ticari markalarıdır.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> a Affinity Fusion<sup>®</sup> jsou registrované ochranné známky společnosti Medtronic, Inc.  
Az Affinity<sup>®</sup>, a Balance<sup>®</sup>, az Affinity Orbit<sup>®</sup> és az Affinity Fusion<sup>®</sup> a Medtronic, Inc. bejegyzett védjegye.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> a Affinity Fusion<sup>®</sup> sú registrované ochranné známky spoločnosti Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> ja Affinity Fusion<sup>®</sup> on ettevõtte Medtronic, Inc. registreeritud kaubamärgid.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> и Affinity Fusion<sup>®</sup> са регистрирани търговски марки на Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> i Affinity Fusion<sup>®</sup> registrirani su zaštitni znakovi tvrtke Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> un Affinity Fusion<sup>®</sup> ir Medtronic, Inc reģistrētas preču zīmes.  
„Affinity<sup>®</sup>“, „Balance<sup>®</sup>“, „Affinity Orbit<sup>®</sup>“ ir „Affinity Fusion<sup>®</sup>“ yra registruotieji „Medtronic, Inc.“ prekių ženklai.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> și Affinity Fusion<sup>®</sup> sunt mărci comerciale înregistrate ale Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> и Affinity Fusion<sup>®</sup> являются зарегистрированными товарными знаками Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> i Affinity Fusion<sup>®</sup> su registrovani žigovi preduzeća Medtronic, Inc.  
Affinity<sup>®</sup>, Balance<sup>®</sup>, Affinity Orbit<sup>®</sup> in Affinity Fusion<sup>®</sup> so registrirane blagovne znamke družbe Medtronic, Inc.

**Explanation of symbols on package labeling / Explication des symboles des étiquettes sur l'emballage / Erläuterung der Symbole auf dem Verpackungsetikett / Explicación de los símbolos que aparecen en el etiquetado del envase / Verklaring van de symbolen op de verpakkinglabels / Spiegazione dei simboli sulle etichette della confezione / Forklaring av symboler på emballasjen / Pakkausetiketin symbolien selitykset / Förklaring av symboler på förpackningsetiketten / Επεξήγηση των συμβόλων στη σήμανση της συσκευασίας / Forklaring til symbolerne på emballagens mærkater / Objasnienie symboli znajdujących się na etykietach opakowania / Explicação dos símbolos na documentação da embalagem / Ambalaj etiketi üzerindeki sembollerin açıklaması / Vysvětlení symbolů uvedených na štítcích balení / A csomagoláson látható szimbólumok jelentése / Vysvetlivky k symbolom na označení balenia / Pakendi siltidel olevate tingmärkide tähendus / Объяснение на символите от етикетите на опаковката / Objašnjenje simbola s etikete na ambalaži / Uz iepakojuma etiķetes attēloto simbolu skaidrojums / Pakuotės etiketės nurodytų simbolių reikšmės / Legenda simbolurilor cuprinse în etichetele de pe ambalaj / Объяснение символов на этикетке упаковки / Objašnjenje simbola na oznakama na pakovanju / Razlaga simbolov na embalaži**

Refer to the outer package label to see which symbols apply to this product / Se référer à l'étiquette sur l'emballage extérieur pour savoir quels symboles s'appliquent à ce produit / Welche Symbole für dieses Produkt zutreffen, entnehmen Sie bitte dem Etikett auf der Verpackung. / Consulte la etiqueta del envase exterior para comprobar qué símbolos son aplicables a este producto / Controleer het label op de buitenverpakking om te zien welke symbolen op dit product van toepassing zijn / Fare riferimento all'etichetta sulla confezione esterna per sapere quali simboli si applicano a questo prodotto. / Se på etiketten utenpå pakningen for å fastslå hvilke symboler som gjelder for dette produktet / Katso ulkopakkauksen etiketistä, mitkä symbolit koskevat tätä tuotetta / Se etiketten på ytterförpackningen för de symboler som gäller denna produkt / Ανατρέξτε στην ετικέτα της εξωτερικής συσκευασίας για να δείτε ποια σύμβολα αφορούν αυτό το προϊόν / Se mærkaten på den ydre emballage for de symboler, der gælder for dette produkt / Należy zapoznać się z etykietą na zewnętrznym opakowaniu, aby określić, które symbole dotyczą niniejszego produktu / Consultar a etiqueta exterior da embalagem para ver que símbolos se aplicam a este produto / Bu üründe hangi sembollerin geçerli olduğunu görmek için dış ambalaj etiketine başvurun / Symboly, které se týkají tohoto produktu, naleznete na štítku vnějšího obalu / A termékre vonatkozó szimbólumok a külső csomagoláson találhatóak. / Prezrite si označenie na vonkajšom obale a zistite, ktoré symboly sa vzťahujú na tento produkt / Välispakendi sildilt näete, millised tingmärgid selle toote puhul kehtivad / Погледнете етикета на външната опаковка, за да видите какви символи се отнасят за този продукт / Na etiketama na ambalaži pogledajte koji se simboli primjenjuju na ovaj proizvod / Skatiet ārējā iepakojuma marķējumu, lai redzētu simbolus, kuri pielietojami šim izstrādājumam / Norādami sužinoti, koku simboliu pažymėtas šis gaminys, žiūrėkite į išorinę pakuotės etiketę / Consultați eticheta de pe ambalajul exterior pentru simbolurile aplicabile acestui produs / См. на этикетке внешней упаковки, какие символы применимы к данному продукту / Pogledajte spoljnu oznaku na pakovanju da biste videli koji se simboli primenjuju na ovaj proizvod / Opis simbolov, ki se nanašajo na izdelek, je na zunanji nalepki na ovojnini.

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
## Affinity Fusion®

### Oxygenator with Integrated Arterial Filter and Cardiotomy/Venous Reservoir with Balance® Biosurface

#### 1. Model

BB841	Oxygenator with Integrated Arterial Filter and Cardiotomy/Venous Reservoir with Balance® Biosurface <sup>1</sup>
-------	--

#### Accessories Sold Separately

ATP210	Affinity® Temperature Probe
AUH2093	Affinity Orbit® Holder System
AMH2014	Affinity® Manifold Holder
RCL841	Affinity Fusion® Recirculation Line

#### 2. Product Description

This product is single-use, nontoxic, nonpyrogenic, supplied STERILE in individual packaging. Sterilized by ethylene oxide.

##### 2.1. Oxygenator with Integrated Arterial Filter

The Affinity Fusion Oxygenator with Integrated Arterial Filter (OXY) and Balance Biosurface is a single-use, microporous, hollow-fiber, gas-exchange device with plasma-resistant fiber and integrated heat exchanger and filter. The OXY is bonded on its primary blood contacting surfaces with a nonleaching biocompatible surface to reduce platelet activation and adhesion and preserve platelet function.

##### 2.2. Cardiotomy/Venous Reservoir

The Affinity Fusion Cardiotomy/Venous Reservoir (CVR) with Balance Biosurface is a single-use device designed to collect and store blood during extracorporeal circulation. Cardiotomy blood is collected, filtered, and defoamed before mixing with the filtered venous blood. With the exception of the filter, the CVR is bonded on its primary venous blood contacting surfaces with a nonleaching biocompatible surface to reduce platelet activation and adhesion and preserve platelet function.

#### 3. Specifications

##### Oxygenator with Integrated Arterial Filter

Oxygenator Membrane Type	Microporous polypropylene hollow fiber
Maximum Gas:Blood ratio	2:1
Heat Exchanger	Polyethylene terephthalate (PET)
Oxygenator Membrane Surface Area	2.5 m <sup>2</sup>
Static Priming Volume	260 mL
Recommended Blood Flow Range	1.0–7.0 L/min
Maximum Water Pressure	206 kPa (1550 mm Hg)
Maximum Rated Blood Pressure	100 kPa (750 mm Hg)
Filtration	25 µm
Storage Temperature	–30°C–50°C (–22°F–122°F)
Storage Humidity	30–75% noncondensing

##### Cardiotomy/Venous Reservoir

Reservoir Volume Capacity	4500 mL
Recommended Blood Flow Range	1.0–7.0 L/min
Maximum Cardiotomy Flow Rate	6.0 L/min
Minimum Operating Level at 7.0 L/min	200 mL
Cardiotomy Filter	30 µm depth filter
Venous Inlet Screen	105 µm
Maximum Rated Pressure	+30 mm Hg / –150 mm Hg
Positive/Negative Pressure Relief Valve	<5 mm Hg positive / >100 mm Hg vacuum, average
Volume Marking Accuracy	±10%

<sup>1</sup> Technology licensed under agreement from BioInteractions, Limited, United Kingdom.

Dynamic Holdup Volume at 7.0 L/min	27 mL at 500 mL operating level
Storage Temperature	-30°C–50°C (-22°F–122°F)
Storage Humidity	30–75% noncondensing

## 4. Indications for Use

### 4.1. Oxygenator with Integrated Arterial Filter

The Affinity Fusion Oxygenator with Integrated Arterial Filter and Balance Biosurface is intended to be used in an extracorporeal perfusion circuit to oxygenate and remove carbon dioxide from the blood and to cool or warm the blood during routine cardiopulmonary bypass procedures up to 6 hours in duration.

The Affinity Fusion Oxygenator with Integrated Arterial Filter and Balance Biosurface is designed to filter from the circuit microemboli larger than the specified micron size for periods up to six hours during cardiopulmonary bypass surgery. The Affinity Fusion Oxygenator with Integrated Arterial Filter and Balance Biosurface is suitable for transfer of the volatile anesthetics isoflurane and sevoflurane by administration through the oxygenator gas inlet by means of a suitable gas vaporizer.

### 4.2. Cardiotomy/Venous Reservoir

The Affinity Fusion Cardiotomy/Venous Reservoir with Balance Biosurface is intended to be used in an extracorporeal perfusion circuit to collect venous and cardiotomy suctioned blood during routine cardiopulmonary procedures up to 6 hours in duration. The CVR is also intended for use during vacuum assisted venous drainage (VAVD) procedures.

The Affinity Fusion Cardiotomy/Venous Reservoir with Balance Biosurface is also intended for use after open heart surgery to collect autologous blood from the chest and to aseptically return the blood to the patient for blood volume replacement.

## 5. Contraindications

Do not use this device for any purpose other than indicated.

Do not use if air leaks are observed during priming and/or operation; this may result in air embolism to the patient and/or fluid loss.

The Affinity Fusion Cardiotomy/Venous Reservoir is contraindicated for use in postoperative chest drainage and autotransfusion procedures when:

- There is an air leak in the lung or gross perforations to the chest wall exist.
- Pericardial, mediastinal, pulmonary or systemic infection or malignancy is present.
- Gross contamination or a lymphatic failure is present or suspected.
- Suctioned blood is obtained from a site where a topical hemostatic agent has been used.
- The chest is open and vacuum is applied.
- Protamine has been administered prior to the reservoir being removed from the bypass circuit.
- The patient is returned to surgery for any reason.
- Vented chest tubes not incorporating vent flow regulation, such as a stopcock, are used.

**Caution:** An assessment should be made of the quality and suitability of the blood that has been collected before re-infusion begins.

## 6. Warnings

### 6.1. General Warnings

Read all warnings, precautions, and Instructions for Use carefully prior to use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**

- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device which could result in patient injury, illness, or death.
- Only persons thoroughly trained in cardiopulmonary bypass procedures should use this device. Operation of each device requires constant supervision by qualified personnel for patient safety.
- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use if the package is opened or damaged, the device is damaged, or the protective caps are not in place as the sterility of the device may have been compromised and/or performance may be affected. (Note: The OXY water inlet, water outlet, and gas outlet ports do not have protective caps.)
- All gas emboli must be cleared from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.

- In any extracorporeal circuit, the use of a prebypass and arterial line filter or OXY with integrated arterial filter is recommended.
- Monitoring of circuit pressure is recommended.
- The OXY and CVR should not be operated outside the recommended parameter ranges.
- Do not use beyond the specified duration (6 hours).
- After initiation of bypass, blood flow should circulate through the OXY and CVR at all times within the recommended blood flow range during cardiopulmonary bypass, except during emergency OXY replacement.
- Disinfectants must not be used in the heater/cooler system while the heat exchanger is in use; if disinfectants have been used in the heater/cooler system, the system must be thoroughly flushed prior to use.
- It is recommended to band all blood tubing connectors for added protection against disconnection.
- Ensure that unused ports remain capped to prevent contamination and that the caps are tightened to avoid leakage.
- Tubing should be attached in such a manner as to prevent kinks or restrictions that may alter the blood, water, or gas flow.
- Do not allow alcohol, alcohol-based fluid, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come in contact with the device as they may jeopardize its structural integrity.

## 6.2. Oxygenator-specific Warnings

- The blood phase pressure should be greater than the gas phase pressure at all times in the OXY.
- Do not exceed blood pressure of 100 kPa (750 mm Hg) in the blood pathway of the OXY.
- Once an OXY is primed with blood, adequate heparinization should be maintained per institution cardiopulmonary bypass (CPB) protocol and the blood pathway should be constantly recirculated within the recommended blood flow range.
- Heat exchanger temperature should not exceed 42°C.
- A replacement OXY should be readily available during CPB.
- Do not obstruct gas outlet vents to prevent gas-side pressures from exceeding blood-side pressures.
- Gas transfer characteristics of the OXY are not significantly affected by concentrations of isoflurane anesthetic gas up to 1.3% and sevoflurane anesthetic gas up to 2.6%. Above these concentrations, adjustments of FiO<sub>2</sub> and gas flow rates may be needed to achieve desired gas transfer performance.
- Gas transfer rates may change over time and adjustments of FiO<sub>2</sub> and gas flow rates may be needed to achieve desired gas transfer performance.
- Ensure that the venous reservoir outlet is always positioned above the highest point of the membrane compartment of the OXY.

## 6.3. CVR-specific Warnings

- Total combined flow into the cardiotomy filter should not exceed 6.0 L/min.
- Do not occlude or obstruct the vent port of the cardiotomy/venous reservoir during operation (except when using VAVD).
- The use of a blood level sensing mechanism during operation of this device is recommended.

## 7. Precautions

- Refer to package labeling for storage temperature requirements.
- Use aseptic technique in all procedures.
- Use of a gas scavenging system during transfer of anesthetic gas through the OXY should be considered.
- A strict anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician. Adequate heparinization should be maintained per institution CPB protocol.
- Dispose of the devices in accordance with local hospital, administrative, and/or other government policies.
- Do not use the purge line to monitor patient pressure.

## 8. Adverse Effects

### 8.1. OXY

- Possible side effects include, but are not limited to, hypoxia or hypercarbia, air embolism and thromboembolic phenomena, blood loss, hemolysis, platelet dysfunction, renal compromise, circulatory compromise, hypothermia or hyperthermia, activation (coagulation/complement), hypotension, and death.

### 8.2. CVR

- Possible side effects include, but are not limited to, air embolism, thromboembolic phenomena, hemolysis, platelet dysfunction, activation (coagulation/complement), circulatory compromise, hematoma, renal compromise, hypotension, blood loss, and death.

## 9. Additional Information for Products with Balance Biosurface

The blood-contacting surfaces of the device are coated with Balance Biosurface to reduce platelet activation and adhesion and preserve platelet function. See Figure 9.

**Caution:** Appropriate blood anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.

## 10. Instructions for Use

### 10.1. Instructions for Use for General CPB

**Note:** For instructions specific to Vacuum Assisted Venous Drainage, see Section 11. For instructions specific to Chest Drainage, see Section 12.

**Note:** Refer to the CPB pump speed controller Instructions for Use for details related to system set-up, priming, operating parameters, and emergency back-up equipment.

#### 10.1.1. System Setup

1. Remove the devices from the packaging carefully to ensure a sterile fluid pathway.  
**Warning:** Ensure aseptic technique is used during all stages of setup and use of this product.
2. Attach the holder mast clamp of the Affinity Orbit Holder System (sold separately) to the heart-lung machine mast at the desired height and tighten the mast clamp knob.
3. Turn the locking lever perpendicular to the holder arm to unlock it.
4. Slide the Orbit oxygenator arm onto the post if it is not already in place.
5. Slide the OXY holder connector onto the bottom of the Orbit oxygenator arm. A click will indicate the OXY is firmly seated.  
**Note:** To remove the OXY, press the holder release tab and pull the OXY away from the connector arm.
6. Slide the CVR holder sleeve onto the holder post.
7. Position the OXY and CVR as desired for the procedure. The OXY may be rotated clockwise or counterclockwise and the Orbit oxygenator arm may be moved laterally. The CVR may be rotated clockwise or counterclockwise.
8. Turn the locking lever parallel to the holder arm to lock components in place.
9. Turn the holder arm knob to lock the holder arm in place. See Figure 4 for final system assembly.
10. Remove the obturator flag from the CVR positive/negative pressure relief valve.
11. Connect all blood, gas, and suction lines to their appropriate connections per institution CPB protocol. The venous inlet comes with a 1.0 cm (3/8 in) adapter attached to the venous inlet. If using 1.27 cm (0.5 in) tubing, remove the 1.0 cm (3/8 in) adapter from the CVR blood inlet port.
12. Connect the premembrane purge line on the OXY to a luer port on the CVR. Using the luer ports of the venous inlet is recommended.
13. The sampling manifold should be attached so that the tubing segment with the one-way duckbill valve is positioned between the arterial sampling location and the manifold. Arterial blood should flow through the duckbill valve toward the manifold.
14. If using a recirculation line, connect the recirculation line to the recirculation port on the OXY. Connect the other end of the recirculation line to the recirculation port on the CVR.
15. If using the cardioplegia port, connect the blood access line prior to priming.
16. Connect the quick-connect 1.3 cm (1/2 in) water lines to the inlet and outlet ports on the OXY. Start the water flow and check for leaks from the water compartment to the blood compartment prior to priming.

**Caution:** Do not use the OXY if water is present in the blood compartment.

17. Ensure the CVR is vented to atmosphere by removing the vent port cap. If using the CVR for VAVD (see Section 11), connect the vacuum line to the vent port.

### 10.1.2. Priming & Recirculation

**Note:** A CO<sub>2</sub> flush of the circuit prior to priming may improve ease of priming.

1. Fill the CVR with the appropriate volume of priming solution per institution CPB protocol.
2. Remove all bubbles in the section of tubing (roller pump head/boot) or centrifugal pump (if used) from the CVR outlet to the OXY inlet.
3. Gravity prime the OXY or prime by forward pump flow per institution CPB protocol.
4. Gradually begin recirculation flow using a prebypass filter per institution CPB protocol.
5. Ensure the OXY is free of bubbles including unused ports.

**Warning:** Unused ports should be sealed tightly with the provided caps.

6. Add additional priming solution, if needed, to prime the rest of the extracorporeal circuit.
7. Upon completion of priming and debubbling, gradually reduce flow and stop pump, close all purge lines, and clamp arterial, venous, and recirculation lines.
8. Ensure the entire extracorporeal circuit is free of bubbles prior to initiating bypass.

**Note:** Priming solution may be prewarmed using a heat exchanger prior to initiating bypass.

### 10.1.3. Initiation of Bypass

**Warning:** The blood phase pressure should always be greater than the gas phase pressure.

1. Ensure the gas outlet port is not obstructed.
2. Check for adequate levels of anticoagulation prior to initiating bypass.
3. Remove arterial and venous clamps and gradually increase blood flow. If using the integrated arterial filter of the oxygenator, ensure the premembrane air purge line is open. Then initiate gas flow with appropriate gas settings per institution CPB protocol and clinical judgment specific to the patient and procedure.
4. Adjust water temperature to meet clinical demands.

### 10.1.4. Operation During Bypass

1. Arterial PO<sub>2</sub> is controlled by varying the percent concentration of oxygen present in the ventilating gas.
  - To DECREASE the PO<sub>2</sub>, DECREASE the percent of oxygen in the ventilating gas by adjusting the FiO<sub>2</sub> on your gas blender.
  - To INCREASE the PO<sub>2</sub>, INCREASE the percent of oxygen in the ventilating gas by adjusting the FiO<sub>2</sub> on your oxygen blender.
2. PCO<sub>2</sub> is controlled primarily by varying the total gas flow rate.
  - To DECREASE PCO<sub>2</sub>, INCREASE the total gas flow rate to increase the amount of CO<sub>2</sub> removed.
  - To INCREASE PCO<sub>2</sub>, DECREASE the total gas flow rate to decrease the amount of CO<sub>2</sub> removed.
3. Patient temperature is controlled by regulating the temperature of the water flow in the heat exchanger.

**Warning:** After initiation of bypass, blood flow should circulate through OXY and CVR at all times within blood flow range during cardiopulmonary bypass, except during emergency OXY or CVR replacement.

**Note:** The premembrane air purge line must remain open to ensure proper air removal when using the integrated filter function of the OXY.

### 10.1.5. Termination of Bypass

1. Terminate bypass per institution CPB protocol as individual case and patient condition indicates.
2. Turn off water flow to the heat exchanger, then clamp and remove the water lines.

### 10.1.6. Emergency Oxygenator Replacement

A standby OXY should always be available during cardiopulmonary bypass.

1. Turn off water flow to the heat exchanger, then clamp and remove the water lines.
2. Clamp the venous line at the CVR. Turn off arterial roller pump and clamp arterial line. (If using a centrifugal pump, clamp the arterial line first, prior to stopping the pump.)

3. Remove oxygen line from the gas inlet port.
4. Clamp all OXY blood lines. Leave adequate lengths for reconnections.
5. Remove any monitoring/sampling lines from the OXY.
6. Detach the OXY from the holder by pressing the release tab.
7. Attach the replacement OXY onto the holder.
8. Connect all blood lines to the replacement OXY. Ensure all connections are secure.
9. Reconnect the oxygen line to the gas inlet port.
10. Connect and unclamp the water lines, then turn on water source and inspect for water leaks.
 

**Warning:** Do not unclamp arterial outlet line or venous return line at this time.
11. With sufficient volume in the CVR, turn on the blood pump or gravity prime, and slowly fill the OXY.
12. Increase blood flow through the recirculation line.
13. Stop the pump and clamp the recirculation line.
 

**Note:** Ensure the OXY inlet is deaired.
14. Ensure the entire system is free of leaks and gas bubbles.
15. Remove the clamps from venous and arterial lines and reinitiate bypass.

### 10.1.7. Emergency CVR Replacement

A standby CVR should always be available during cardiopulmonary bypass.

1. If blood is to be salvaged, clamp the venous inlet on the CVR and pump the blood into the patient. If this is not possible, then pump the blood from the CVR into a storage bag or reservoir, or maintain in removed CVR for later reinfusion.
 

**Note:** The responsibility for the use of blood salvaged in all cases belongs solely to the physician ordering its use.
2. Clamp the venous inlet on the CVR and turn off the arterial roller pump. Clamp the CVR outlet (if using a centrifugal pump, clamp the pump outlet line prior to stopping the pump).
3. Clamp and disconnect the CVR inlet and outlet lines, OXY purge, and OXY recirculation lines, at a minimum. Detach the CVR from the holder.
4. Attach the new CVR onto the holder.
5. Reconnect the CVR inlet and outlet lines, OXY purge, and OXY recirculation lines, at a minimum.
6. Add volume if necessary to reinstitute bypass.
7. If air was introduced into the OXY, remove the CVR outlet clamp and reprime the OXY per steps 11 through 14 in Section 10.1.6.
8. Remove the clamps from venous and arterial lines and reinitiate bypass.
9. Connect other lines (eg, suction) as necessary, if they were not already connected.

## 11. Cardiotomy/Venous Reservoir, Vacuum Assisted Venous Drainage (VAVD)

### 11.1. Warnings

- Do not use a vacuum regulator with maximum negative vacuum pressure greater than  $-20$  kPa ( $-150$  mm Hg). Use of excessive negative pressure may lead to hemolysis.
- Do not occlude or obstruct the vent port of the cardiotomy/venous reservoir during operation (except when using VAVD).
- Draw only intended drug amounts into the syringe that will be used for drug administration into the CVR during VAVD. Negative pressure in the CVR could draw more drug from the syringe into the reservoir than intended and lead to excess drug delivery.
- All A/V shunt lines must be closed to atmosphere prior to stopping the pump or using low blood flow rates during VAVD. Failure to do so may allow air to be drawn into the blood side of the OXY from the fiber.
- Appropriate occlusion of the arterial roller pump head is necessary when using VAVD. Failure to appropriately occlude the pump may allow air to be drawn into the blood side of the OXY from the fiber.
- During use of VAVD techniques, arterial filtration should be used. The type of arterial filtration used, either the integrated filter in the Affinity Fusion oxygenator with purge line open or an additional stand-alone filter, is determined by the clinician.
- Gradually return to atmospheric pressure when ending use of VAVD. A sudden change in pressure may result in turbulent blood flow inside the CVR.
- The line between the OXY and centrifugal pump (if used) must be clamped before stopping the pump during use of VAVD. Failure to clamp off the arterial line may result in air drawn into the blood

side of the OXY from the fiber. Use of an arterial line one-way valve between the OXY and centrifugal pump is recommended.

- A controlled vacuum regulator is required during VAVD.
- A positive/negative pressure relief valve is required during VAVD (included, see Figure 3).

### 11.2. Precautions

- Use of a pressure measurement device on the CVR and an additional negative pressure relief valve that operates at  $-13$  kPa ( $-100$  mm Hg) should be considered.
- A vapor trap is required during use of VAVD due to condensation development.
- Do not allow the vapor trap to become completely filled during use. This may allow fluid to enter the vacuum controller, or prevent the CVR from being vented to atmosphere when no vacuum is applied.
- Do not apply vacuum to the venous reservoir when there is no forward blood flow through the OXY. This applies to arterial centrifugal and roller pumps (roller may not be occlusive in all positions). This will prevent air from being pulled across the membrane into the blood path by the reservoir vacuum.

### 11.3. VAVD Instructions for Use

1. Follow the Instructions for Use in Section 10 for system setup with the following modifications as listed in steps 2 through 3 below.
2. Attach a calibrated means of monitoring pressures to either the venous line or venous reservoir.
3. All manifold sampling ports must have nonvented caps (included) when using VAVD.
4. Follow the Instructions for Use in Section 10 for priming with the following additions as listed in steps 5 through 8 below.
5. Prior to initiating bypass, prepare the vacuum portion of the setup. Attach a vented vacuum line to the CVR vent port and to the vacuum regulator with vapor trap.

**Caution:** Clamp quick-prime line prior to application.

**Caution:** Do not allow the vapor trap to become completely filled during use. This may allow fluid to enter the vacuum controller, or prevent the CVR from being vented to atmosphere when no vacuum is applied.

6. To vent the vacuum line, ensure the Y-connector with side tubing is easily accessible and open to atmosphere. This may be clamped and unclamped regularly throughout the procedure to apply or stop vacuum.
7. Ensure all access ports on CVR and manifold are properly sealed and are secure before initiating bypass.

**Caution:** Vacuum will result in excess atmospheric air being entrained into the CVR and potentially into the blood pathway with opened blood pathway ports. Vacuum may result in increased flow through the sampling manifold during use.

8. Use vacuum assistance as needed to institute and maintain cardiopulmonary bypass.

**Caution:** Do not apply vacuum to the venous reservoir when there is no forward blood flow through the OXY. This applies to arterial centrifugal and roller pumps (roller may not be occlusive in all positions). This will prevent air from being pulled across the membrane into the blood path by the reservoir vacuum.

**Caution:** Do not allow the CVR to become pressurized as this could obstruct venous drainage, force air retrograde into the patient or into the blood path of the OXY.

## 12. Cardiotomy/Venous Reservoir, Postoperative Chest Drainage

### 12.1. Warnings

When using the Affinity Fusion Cardiotomy/Venous Reservoir for post-operative chest drainage:

- Salvaged blood should be reinfused on an hourly basis unless fewer than 50 mL of blood/fluid are collected hourly.
- A minimum reservoir level of 25 mL of blood should always be maintained to prevent passing an air embolus during reinfusion.
- Blood that has resided in the Affinity Fusion Cardiotomy/Venous Reservoir for 4 hours or more should not be transfused.
- It is not recommended that autotransfusion continue for more than 18 hours following surgery.
- If the filter occludes and chest drainage volume is excessive, unfiltered blood may cross the external shunt line to the outlet side of the reservoir. Any salvaged blood which has passed through the shunt line must be filtered prior to reinfusion.

- Follow all accepted medical and nursing care routines for patients undergoing chest drainage.
- If a vacuum source is being used, vacuum pressure should not exceed 20 cm H<sub>2</sub>O for chest drainage.
- The use of salvaged blood from this device may be contraindicated (for example, in the presence of sepsis or malignancy). The responsibility for the use of this device in all cases belongs solely to the physician ordering its use.
- Actual performance results may vary depending on many in-use variables.

## 12.2. Instructions for Use

**Caution:** The safe length of time that blood or blood products may remain in the appliance or plastic disposables is dependent on collection and storage methods.

1. Carefully remove the tubing and/or vented caps from the suction inlet ports, recirculation inlet port, and venous blood inlet port of the reservoir. Remove the stopcock sampling manifold and the attached coiled manifold sampling lines. Plug all 0.6 cm (1/4 in) suction inlet ports with sealing plugs. Plug the quick prime and recirculation ports.
2. Remove the 1.0 cm (3/8 in) adapter if used.
3. Plug the 1.3 cm (1/2 in) venous blood inlet port with a sealing plug.
4. Connect the chest drainage tube to the 1.0 cm (3/8 in) suction inlet port.
5. Cap all luer ports with nonvented luer caps.
6. Tilt the reservoir back such that any blood that is present is not at the outlet port but at the back of the reservoir. Remove the 1.0 cm (3/8 in) tubing from the venous blood outlet port and replace with 0.3 cm to 1.0 cm (1/8 in to 3/8 in) infusion adapter line. Clamp the line off and leave the end cap in place until autotransfusion is set up.

**Note:** Do not tilt the reservoir back so far as to cause blood to leak out of the pressure relief valve.

7. Place the cardiotomy/venous reservoir in a holder subthoracic to the patient. Maintain the reservoir below the thorax to facilitate chest drainage.
8. Using sterile tubing, connect the horizontal vent port on the top of the unit to a vacuum source with approximately 15 to 20 cm of H<sub>2</sub>O regulated vacuum.

**Note:** Use of a water seal device is recommended in all chest drainage cases.

9. Record the amount of drainage. Continually monitor drainage and elapsed time to determine the hourly drainage rate.
10. As soon as the reservoir has collected 50 mL of drainage, set up and prime the infusion pump tubing per standard hospital procedure.

**Note:** It is recommended that a transfusion filter be used in all reinfusion procedures.

**Note:** It is also recommended that an infusion pump be used with an air detection sensor for all reinfusion procedures.

11. Autotransfuse the fluid collected per hospital protocol after assessing the drainage rate, and hemodynamic status of the patient.

## 13. Additional Information

The following additional information and data are available upon request: sterilization method, materials of the blood pathway, plasma leakage across the semipermeable membrane, blood cell damage, particle release, air-handling capability, anti-foam characteristics, break-through volume, filtration efficiency, relevant tolerances for specifications in the Instructions for Use, and protocol summaries.

## 14. Important Notice — Limited Warranty (U.S. Customers Only)

- A. This **LIMITED WARRANTY** provides the following assurances to the patient who receives an Affinity Fusion<sup>®</sup> Oxygenator with Integrated Arterial Filter and Cardiotomy/Venous Reservoir with Balance<sup>®</sup> Biosurface, hereafter referred to as "The Product".
  - (1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its "Use By" date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the Purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
  - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
  - (1) The Product must be used prior to its "Use By" date.
  - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.

- (3) The Product must not have been altered or subjected to misuse, abuse or accident.
  - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) Except as expressly provided by this **LIMITED WARRANTY**, Medtronic is not responsible for any direct, incidental or consequential damages based on any defect, failure or malfunction of the product, whether the claim is based on warranty, contract, tort or otherwise.
  - (2) THIS **LIMITED WARRANTY** IS MADE ONLY TO THE PATIENT ON WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS **LIMITED WARRANTY** SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
  - (3) The exclusions and limitation set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid. This **LIMITED WARRANTY** gives the patient specific legal rights. The patient may also have other rights which vary from state to state.
  - (4) No person has any authority to bind Medtronic to any representation, condition or warranty except this **LIMITED WARRANTY**.

**Note:** This **LIMITED WARRANTY** is provided by Medtronic, Inc. 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States.

#### **15. Important Notice — Limited Warranty (For Countries Outside the U.S.)**

- A. This **LIMITED WARRANTY** provides assurance for the purchaser who receives an Affinity Fusion® Oxygenator with Integrated Arterial Filter and Cardiotomy/Venous Reservoir with Balance® Biosurface, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.
- The Warnings contained in the product labeling are considered an integral part of this **LIMITED WARRANTY**. Contact your local Medtronic representative to obtain information on how to process a claim under this **LIMITED WARRANTY**.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its "Use By" date.
  - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
  - (3) The Product may not have been used for any other patient.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
  - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

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© 2012–2013 Medtronic, Inc.