Medtronic

Affinity NT™

Oxygenator with Cortiva™ BioActive Surface
Oxygenator med Cortiva™ BioActive-overflade
Oxygenator mit bioaktiver Cortiva™ Oberfläche
Οξυγονωτής με βιοενεργή επίστρωση Cortiva™
Oxigenador con superficie biocompatible Cortiva™
Bioaktiivisella Cortiva™-pinnoitteella päällystetty hapetin
Oxygénateur avec surface bioactive Cortiva™
Ossigenatore con superficie bioattiva Cortiva™
Oxygenator met bioactieve Cortiva™-oppervlaktelaag
Oksygenator med Cortiva™ bioaktiv overflate
Oxigenador com superfície bioativa Cortiva™
Oxygenator med Cortiva™ bioaktiv yta
CB511

Instructions for Use • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης •
Instrucciones de uso • Käyttöohjeet • Mode d'emploi • Istruzioni per l'uso • Gebruiksaanwijzing •
Bruksanvisning • Instruções de utilização • Bruksanvisning

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
Affinity NT™

Oxygenator with Cortiva™ BioActive Surface

1. Model

CB511

2. Description

The Affinity NT™ oxygenator with Cortiva bioactive surface is a single-use gas exchange device with blood outside the fiber design and includes an integral heat exchanger.

Products coated with Cortiva bioactive surface include a “CB” prefix in the model number. Additional information regarding Cortiva bioactive surface is available in Section 8.

2.1. Specifications

Oxygenator

Membrane type	Microporous polypropylene hollow fibers
Membrane surface area	2.5 m ²
Static priming volume	270 mL
Residual blood volume ^a	144 mL
Minimum gas:blood ratio	0.5:1
Maximum gas:blood ratio	2:1
Recommended blood flow rate	1 to 7 L/min
Maximum water side pressure	206.8 kPa (30 psi)
Arterial outlet port	1.0 cm (3/8 in)
Venous inlet port	1.0 cm (3/8 in)
Access port	Female luer port
Arterial sample port	Female luer port
Recirculation port	0.6 cm (1/4 in)
Gas inlet port	0.6 cm (1/4 in)
Gas outlet port	1.0 cm (3/8 in) nonbarbed
Water ports	1.3 cm (1/2 in) quick disconnects
Size	Adult
Maximum transit temperature	40°C (104°F)

^a Tested with saline

Store the product at room temperature.

Upon written request, a list of all materials used in the construction of the oxygenator will be made available.

2.2. Accessories available separately

Model number	Description
TP	Temperature probe
61399401072	Affinity™ oxygenator holder

3. Indications for use

The Affinity NT oxygenator with Cortiva bioactive surface is intended to be used in an extracorporeal perfusion circuit to oxygenate and remove carbon dioxide from the blood and to cool or warm the blood during routine cardiopulmonary bypass (CPB) procedures up to 6 hours in duration.

4. Contraindications

Use the device only as indicated.

5. Warnings

Read all warnings, precautions, and instructions for use carefully prior to use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**

- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use the device if the package is opened or damaged, the device is damaged, or the protective caps are not in place.
- If air bubbles and/or leaks are observed during priming and/or operation, these conditions may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the device if these conditions are observed.
- Clear all gas emboli from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- In any extracorporeal circuit, the use of a prebypass and arterial line filter should be considered.

- Monitoring of preoxygenator circuit pressure should be considered.
- The blood phase pressure must be greater than the gas phase pressure at all times in the oxygenator.
- Do not exceed blood phase pressure of 750 mmHg.
- Do not exceed heat exchanger water phase pressure of 206.8 kPa (30 psi).
- Heat exchanger water temperature shall not exceed 42°C (107°F).
- Follow institution CPB protocol when using disinfectants in the heater/cooler during bypass. The integrity of the water path has been verified with hydrogen peroxide (330 ppm). Contact Medtronic for information regarding the use of additional disinfectants.
- Band all tubing connections for added protection against high fluid pressures.
- A replacement oxygenator should be readily available during perfusion.
- Do not obstruct gas outlet vents to prevent the gas side pressures from exceeding blood side pressures.
- Only persons thoroughly trained in CPB procedures should use this device. Operation of each device requires constant supervision by qualified personnel for patient safety. Perfusion must be carefully and constantly monitored.
- Each device is intended for single use only. Do not reuse or resterilize. Sterilized using ethylene oxide.
- Attach tubing in such a manner as to prevent kinks or restrictions that may alter blood, water, or gas flow.
- Do not allow alcohol, alcohol-based fluids, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come into contact with the device as they may jeopardize the structural integrity.

6. Precautions

- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize any component of the system. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Use aseptic technique in all procedures.
- Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures. Adequate anticoagulation must be maintained before and during bypass.
- Additional warnings and precautions applicable to specific procedures can be found at appropriate places in the instructions for use.
- It is the responsibility of the user to dispose of the devices in accordance with local regulations and hospital procedures.

7. Adverse effects

The following known adverse effects are associated with the use of the product: blood loss, coagulopathy, death, embolism, excessive blood component activation or thrombogenicity, exsanguination, hemolysis, hemolytic anemia, hypercarbia, hypotension, hypovolemia, hypoxia, infection, ischemia, neurological dysfunction, and organ dysfunction.

8. Information for products with Cortiva bioactive surface

The primary blood-contacting surfaces of the product are coated with Cortiva bioactive surface. This coated surface enhances blood compatibility and provides a blood-contacting surface that is thromboresistant. Cortiva bioactive surface contains nonleaching heparin derived from porcine intestinal mucosa¹.


Caution: Product coated with Cortiva bioactive surface is intended for single use only. Resterilization may adversely affect the Cortiva bioactive surface.

Caution: Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

¹ Nonleaching is defined as heparin at a level ≤ 0.1 IU/mL, as measured through clinically relevant extraction conditions.

9. Instructions for use

9.1. System setup

1. Filtered cardiomy reservoir
2. Venous reservoir
3. Arterial filter
4. Oxygenator
5. Prebypass filter

Figure 6.

1. Remove the device from the packaging carefully to ensure sterile fluid pathway.

Warning: Ensure that aseptic technique is used during all stages of setup and use of this system.

Warning: Before removing the device from the packaging, inspect the packaging and product for damage. If the package or product is damaged, do not use as sterility of the device may have been compromised and/or performance may be affected.

2. Place the oxygenator securely into its holder. Ensure that the venous reservoir blood outlet is higher than the oxygenator.
3. Set up the oxygenator, cardiomy reservoir, venous reservoir, blood lines, and oxygen lines. Do not occlude tubing in pump raceway.
4. Connect the 1.3 cm (1/2 in) water lines to the inlet and outlet ports on the heat exchanger. Start the heat exchanger water flow, and check for leaks from the water compartment to the blood compartment. Discard the oxygenator if water is present in blood compartment, and set up the standby oxygenator.
5. Clamp the tubing at cardiomy outlet, venous reservoir inlet, venous reservoir outlet, oxygenator arterial outlet, and oxygenator recirculation line outlet.

Note: A CO₂ flush of the oxygenator is not required, but is recommended before priming.

9.2. Priming


Figure 7.

1. Fill the clamped cardiomy reservoir with 1.5 to 2 L of priming solution.
2. With the venous reservoir at minimum volume setting, open air purge on the venous reservoir and unclamp the cardiomy reservoir outlet line. Allow the solution to slowly gravity prime the venous reservoir while carefully venting any air from purge port.

Note: Venous reservoir fluid level should always remain above the level of the oxygenator.

3. Open the venous reservoir outlet clamp and allow solution to gravity prime the oxygenator heat exchanger and fiber bundle. Remove any remaining air from venous reservoir purge port.

Note: Albumin has been reported to reduce the rare instance of transient high inlet pressure phenomenon in membrane oxygenators. Therefore, albumin in the prime solution should be considered.

9.3. Recirculation


Figure 8.

1. After placing tubing into a properly occluded pump raceway, remove clamp from recirculation outlet per Figure 8 to allow recirculation between the oxygenator and venous reservoir.
2. Begin recirculation flow at 500 mL/min and gradually increase flow to 5 L/min.
3. Ensure that oxygenator is bubble-free.
4. Carefully purge any remaining air from the venous reservoir.
5. Gradually reduce flow, stop the pump, and clamp the recirculation line.

9.4. Priming A/V circuit


Figure 9.

1. Unclamp the venous reservoir inlet line.
2. Unclamp the arterial outlet; slowly pump prime solution into the remaining circuit while carefully removing air from the venous reservoir purge port. Purge air from the arterial sample port and access port at this time.
3. The arterial filter should now be primed per manufacturer's instructions.
4. Inspect the entire circuit and all components for leaks.
5. Upon completion of priming and debubbling, gradually reduce blood flow and stop the pump, close all purge lines, and clamp the arterial, venous, and recirculation lines.
6. Ensure that the entire extracorporeal circuit is free of bubbles before initiating bypass.

9.5. Initiation of bypass

Warning: The blood phase pressure must be greater than the gas phase pressure at all times.

Note: Priming solution may be prewarmed through the heat exchanger before initiation of bypass.

1. Check for adequate levels of anticoagulation before initiation of bypass.
2. Remove the arterial and venous clamps and gradually increase blood flow. Then initiate gas to blood flow ratio of 1:1 with a minimum FiO_2 of 80%.
3. Initiate water flow through the heat exchanger.

Warning: Do not allow air to enter the circuit.

9.6. Operation during bypass

1. Arterial PO_2 is controlled by varying the percent concentration of oxygen present in the ventilating gas.
 - a. To **decrease** the PO_2 , **decrease** the amount of oxygen in the ventilating gas by decreasing the FiO_2 on your oxygen blender.
 - b. To **increase** the PO_2 , **increase** the amount of oxygen in the ventilating gas by increasing the FiO_2 on your oxygen blender.
2. PCO_2 is controlled primarily by varying total gas flow rate.
 - a. To **decrease** PCO_2 , **increase** the total gas flow rate to increase the amount of CO_2 removed.

- b. To **increase** PCO₂, **decrease** the total gas flow rate to decrease the amount of CO₂ removed.
 3. Patient temperature is controlled by regulating the temperature of water flow into the heat exchanger.
- Note:** Care should be taken to minimize temperature gradient during cooling and rewarming phases of bypass.

9.7. Termination of bypass

1. Terminate bypass as individual case and patient condition indicate.
2. Turn off gas flow before turning off blood flow.
3. Turn off water flow to heat exchanger.
4. After terminating bypass, continue to slowly recirculate through oxygenator using recirculation line and/or reconnected arterial venous loop.

9.8. Emergency oxygenator replacement

A standby Affinity NT oxygenator with Cortiva bioactive surface should always be available during perfusion.

1. With the cardiotomy reservoir outlet line clamped, add volume if required to prime replacement oxygenator.
2. Turn off water flow to the heat exchanger; clamp and remove water lines. Remove oxygen line from gas inlet port.
3. Clamp venous line at reservoir. Turn off blood pump and clamp arterial outlet.
4. Remove any monitoring/sample line from oxygenator.
5. Double clamp recirculation line, oxygenator inlet, and outlet tubing, and cut tubing between clamps, leaving adequate lengths for reconnection.
6. Remove old oxygenator from holder.
7. Place replacement oxygenator in holder. Connect and band all blood lines. Connect oxygen line.
8. Connect water lines, unclamp, turn on water source, and inspect for water leaks.
9. Remove clamps from cardiotomy line, recirculation line, and the oxygenator inlet.

Warning: Do not unclamp arterial line or venous return line at this time.

10. Turn on blood pump and slowly fill the oxygenator.
11. Increase blood flow to 5 L/min through recirculation line.
12. Vent air bubbles from venous reservoir.
13. Ensure that the entire system is free of leaks and gas bubbles.
14. Remove all clamps from venous and arterial lines, clamp recirculation line, and reinstate bypass.

10. Additional information

The following additional information and data are available upon request: sterilization method, materials of the blood pathway, plasma leakage across the semipermeable membrane, blood and gas pathway pressure drop at the maximum blood and gas flow rates, blood cell damage, particle release, relevant tolerances for data in the Instructions for Use, and test methods used to determine coverage, leaching, and bioactivity.

11. Important notice — limited warranty (US customers only)

- A. This **LIMITED WARRANTY** provides the following assurance to the patient who receives an Affinity NT oxygenator with Cortiva bioactive surface, hereafter referred to as the "Product."
 - (1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its Use-by date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the Purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
 - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
 - (1) The Product must be used prior to its Use-by date.
 - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
 - (3) The Product must not have been altered or subjected to misuse, abuse, or accident.
 - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
 - (1) Except as expressly provided by this **LIMITED WARRANTY**, Medtronic is not responsible for any direct, incidental, or consequential damages based on any defect, failure or malfunction of the product, whether the claim is based on warranty, contract, tort, or otherwise.
 - (2) THIS **LIMITED WARRANTY** IS MADE ONLY TO THE PATIENT IN WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED,

INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS **LIMITED WARRANTY** SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.

- (3) The exclusions and limitation set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid. This **LIMITED WARRANTY** gives the patient specific legal rights. The patient may also have other rights which vary from state to state.
- (4) No person has any authority to bind Medtronic to any representation, condition or warranty except this **LIMITED WARRANTY**.

Note: This **LIMITED WARRANTY** is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for exact terms of the **LIMITED WARRANTY**.

12. Important notice — limited warranty (for countries outside the US)

Important: This limited warranty does not apply in Australia.

- A. This **LIMITED WARRANTY** provides assurance for the purchaser who receives an Affinity NT oxygenator with Cortiva bioactive surface, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.

The Warnings contained in the product labeling are considered an integral part of this **LIMITED WARRANTY**. Contact your local Medtronic representative to obtain information on how to process a claim under this **LIMITED WARRANTY**.

- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its Use-by date.
 - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
 - (3) The Product may not have been used for any other patient.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
 - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

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© 2018 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
Medtronic B.V.
Medtronic of Canada Ltd