Hollow Fiber Oxygenator with Plasma Resistant Fiber
Oxygénateur à fibre creuse plasma-résistante
Oxygenator mit plasmaresistenter Hohlfaser
Oxigenador de fibra hueca resistente al plasma
Hollow-fiber-oxygenator met plasmabestendige vezel (PRF)
Ossigenatore a fibre cave plasmaresistenti
Hullfiberoksygenator med plasmaresistent fiber
Plasmaresistentistä kuidusta valmistettu ontelokuituhapetin
Hålfiberoksygenator med plasmaresistent fiber
Οξυγονωτής τριχοειδών κοίλων ινών με ίνες ανθεκτικές στο πλάσμα
Hollow-fiber oxygenator med plasmaresistent fiber
Oksygenator membranowy HFM z włóknami odpornymi na działanie osocza
Oxigenador de fibra oca resistente ao plasma
Plazmaya Dirençli Fiberli Hollow Fiber Oksijenatör
Oxygenátor s dutými vlákny odolnými vůči plazmě
Membránoxigenátor plazmarezisztens kapillárisal
Oxygenátor s dutými vlákny rezistentnými voči plazme
Plasmaresistentsete õðneskiududega fiiberoksügenaator
Оксигенатор тип кyxо влакно с плазма резистентно платно
Oksygenator sa šupljim vlaknima otpornim na plazmu
Dobškieðru oksygenators ar plazmu neuzsūcošām šķiedrām
Tuščiaviduris skaidulų oksygenatorius su plazmai atspariomis skaidulomis
Oxigenator cu fibre concave și rezistente la plasmă
Половолоконный оксигенатор с плазморезистентным волокном
Oksygenator sa šupljim vlaknima otpornim na plazmu
Oksygenator z votlimi vlakni, odpornimi na plazmo

Instructions For Use ▪ Mode d'emploi ▪ Gebrauchsanweisung ▪ Instrucciones de uso
Gebruiksaanwijzing ▪ Istruzioni per l'uso ▪ Bruksanvisning ▪ Käyttöohjeet ▪ Bruksanvisning
Οδηγίες χρήσης ▪ Brugsanvisning ▪ Instrukcja użytkowania ▪ Instruções de utilização
Kullanım Talimatları ▪ Návod k použití ▪ Használati utasítás ▪ Pokyny pre používanie
Kasutusjuhend ▪ Инструкции за употреба ▪ Upute za upotrebu ▪ Lietošanas pamācība
Naudojimo instrukcijos ▪ Instrucțiuni de utilizare ▪ Инструкция по эксплуатации
Uputstva za upotrebu ▪ Navodila za uporabo

Explanation of symbols on package labeling / Explicación des symboles des étiquettes sur l'emballage / Erklärung der Symbole auf dem Verpackungsetikett / Explicación de los símbolos que aparecen en el etiquetado del envase / Verklaring van de symbolen op de verpakingslabels / Spiegazione dei simboli sulle etichette della confezione / Forklaring av symboler på emballasjen / Pakkauksen merkintöjen selitykset / Förklaring av symboler på förpackningsetiketten / Επεξήγηση των συμβόλων στην επικετα της συσκευασίας / Forklaring af symboler på emballagens mærkater / Objašnjenja simboli najdujućacych się na etykietach opakowania / Explicação dos símbolos nas etiquetas da embalagem / Ambalaj etiketindeki sembollerin açıklamaları / Vysvětlení symbolů uvedených na štítcích obalu / A csomagoláson látható szimbólumok jelentése / Vysvětlivky k symbolom na označení balenia / Pakendi siltidel olevate tingmärkide tähendus / Обяснение на символите от етикетите на опаковката / Objašnjenje simbola na etiketi na ambalaži / Iepakojuma etiķetē attēlotu simbolu skaidrojums / Simboli, esančiji ant pakuočių, paaiškinimas / Explicarea simbolurilor cuprinse în etichetele de pe ambalaj / Объяснение обозначений на наклейках продукта / Objašnjenje simbola na oznakama na pakovanju / Razlaga znakov na ovojnini

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




AFFINITY® NT

Hollow Fiber Oxygenator with Plasma Resistant Fiber

Description

AFFINITY® NT Hollow Fiber Oxygenator with Plasma Resistant Fiber is a single-use gas exchange device with blood outside the fiber design and includes an integral stainless steel heat exchanger.

Specifications

Oxygenator

Membrane type:	Microporous Polypropylene Hollow Fibers
Membrane Surface Area:	2.5 m ²
Static Priming Volume:	270 mL
Recommended Blood Flow Rate:	1 - 7 liters/minute
Maximum Water Side Pressure:	30 psi
Arterial Outlet Port:	1.0 cm (3/8 in)
Venous Inlet Port:	1.0 cm (3/8 in)
Access Port:	Female Luer Port
Arterial Sample Port:	Female Luer Port
Recirculation Port:	0.6 cm (1/4 in)
Gas Inlet Port:	0.6 cm (1/4 in)
Gas Outlet Port:	1.0 cm (3/8 in) non-barbed
Water Ports:	1.3 cm (1/2 in) quick disconnects
Size:	Adult

Upon written request, a list of all materials used in the construction of the oxygenator will be made available.

Indications for Use

The AFFINITY® NT Hollow Fiber Oxygenator with Plasma Resistant Fiber is intended to be used in an extracorporeal perfusion circuit to oxygenate and remove carbon dioxide from the blood and to cool or warm the blood during routine cardiopulmonary bypass procedures up to 6 hours in duration.

Contraindications

This device used for any other purposes than for the indicated intended use is the responsibility of the user.

Warnings

Read all Warnings, Precautions, and Instructions for Use carefully prior to use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**

- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use if the package is opened or damaged or if the device is damaged or the protective caps are not in place.
- If air bubbles and/or leaks are observed during priming and/or operation, these conditions may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the device if these conditions are observed.
- All gas emboli must be cleared from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- In any extracorporeal circuit, the use of a pre-bypass and arterial line filter should be considered.
- Monitoring of pre-oxygenator circuit pressure should be considered.
- The blood phase pressure must be greater than the gas phase pressure at all times in the oxygenator.
- Do not exceed blood phase pressure of 750 mm Hg.
- Do not exceed heat exchanger water phase pressure of 30 psi.
- Heat exchanger water temperature shall not exceed 42°C.
- Disinfectants must not be used in the heater/cooler system during the time the heat exchanger is in use; if disinfectants have been used in the heater/cooler system, the system must be thoroughly flushed prior to use.
- All blood tubing connections should be banded for added protection against high fluid pressures.
- A replacement oxygenator should be readily available during perfusion.
- Do not obstruct gas outlet vents to prevent the gas side pressures from exceeding blood side pressures.
- Only persons thoroughly trained in cardiopulmonary bypass procedures should use this device. Operation of each device requires constant supervision by qualified personnel for patient safety.
- Each device has been carefully manufactured, tested and packaged; however, the state of the art has not been developed to the point that Medtronic is able to ensure that the device will not leak, crack, or fail during use. Perfusion must be carefully and constantly monitored.
- Each device is intended for single use only. Do not reuse or resterilize. EO sterilized.
- Tubing should be attached in such a manner as to prevent kinks or restrictions that may alter blood, water or gas flow.

- Do not allow alcohol, alcohol-based fluids, anesthetic fluids (such as isoflurane), or corrosive solvents (such as acetone) to come into contact with the device as they may jeopardize the structural integrity.
- Possible side effects include, but are not limited to infections, mechanical failure, hemolysis, air embolism, blood loss, circulatory compromise and thromboembolic phenomena. These are potential side effects with all extracorporeal blood systems.
- Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.**

Precautions

- Refer to package labeling for storage temperature requirements.
- Use aseptic technique in all procedures.
- A strict anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician. Adequate heparinization must be maintained before and during bypass.
- Additional warnings and precautions applicable to specific procedures can be found at appropriate places in the Instructions for Use.

Instructions for Use

System Setup (Refer to Figure 1)


Figure 1.

Warning: The setup and use of the device(s) is the responsibility of the attending clinician.

1. Remove the device(s) from the packaging carefully to ensure sterile fluid pathway.

Warning: Ensure aseptic technique is used during all stages of setup and use of this system.

Warning: Before removing the device(s) from the packaging, inspect the packaging and product for damage. If the package or product is damaged, do not use as sterility of the device may have been compromised and/or performance may be affected.

2. Place oxygenator securely into its holder. Ensure that venous reservoir blood outlet is higher than the oxygenator.
3. Set up oxygenator, cardiotomy reservoir, venous reservoir, blood lines and oxygen lines. Do not occlude tubing in pump raceway.
4. Connect 1.3 cm (1/2 in) water lines to inlet and outlet ports on the heat exchanger. Start the heat exchanger water flow and check for leaks from the water compartment to the blood compartment. Discard oxygenator if water is present in blood compartment and set up standby oxygenator.
5. Clamp tubing at cardiotomy outlet, venous reservoir inlet, venous reservoir outlet, oxygenator arterial outlet and oxygenator recirculation line outlet.

Note: A CO₂ flush of the oxygenator is not required, but recommended prior to priming.

Priming (Refer to Figure 2)


Figure 2.

1. Fill clamped cardiomy reservoir with 1.5 to 2 liters of priming solution.
2. With the venous reservoir at minimum volume setting, open air purge on venous reservoir and unclamp cardiomy reservoir outlet line. Allow solution to slowly gravity prime venous reservoir while carefully venting any air from purge port.

Note: Venous reservoir fluid level should always remain above the level of the oxygenator.

3. Open the venous reservoir outlet clamp and allow solution to gravity prime the oxygenator heat exchanger and fiber bundle. Remove any remaining air from venous reservoir purge port.

Note: Albumin has been reported to reduce the rare instance of transient high inlet pressure phenomenon in membrane oxygenators. Therefore albumin in the prime solution should be considered.

Recirculation (Refer to Figure 3)


Figure 3.

1. After placing tubing into a properly occluded pump raceway, remove clamp from recirculation outlet per Figure 3 to allow recirculation between the oxygenator and venous reservoir.
2. Begin recirculation flow at 500 mL/minute and gradually increase flow to 5 liters/minute.
3. Ensure that oxygenator is bubble free.
4. Carefully purge any remaining air from venous reservoir.
5. Gradually reduce flow, stop pump and clamp recirculation line.

Priming A/V Circuit (Refer to Figure 4)


Figure 4.

1. Unclamp venous reservoir inlet line.
2. Unclamp arterial outlet, slowly pump prime solution into remaining circuit while carefully removing air from venous reservoir purge port. Purge air from arterial sample port and access port at this time.
3. The arterial filter should now be primed per manufacturer's instructions.
4. Inspect entire circuit and all components for leaks.

5. Upon completion of priming and debubbling, gradually reduce blood flow and stop pump, close all purge lines and clamp arterial, venous and recirculation lines.
6. Ensure the entire extracorporeal circuit is free of bubbles prior to initiating bypass.

Initiation of Bypass

Warning: The blood phase pressure must be greater than gas phase pressure at all times.

Note: Priming solution may be pre-warmed through heat exchanger prior to initiation of bypass.

1. Check for adequate levels of anticoagulation prior to initiation of bypass.
2. Remove arterial and venous clamps and gradually increase blood flow. Then initiate gas to blood flow ratio of 1:1 with a minimum FiO_2 of 80%.
3. Initiate water flow through heat exchanger.

Warning: Do not allow air to enter the circuit.

Operation During Bypass

1. Arterial PO_2 is controlled by varying the percent concentration of oxygen present in the ventilating gas.
 - To DECREASE the PO_2 , DECREASE the amount of oxygen in the ventilating gas by decreasing the FiO_2 on your oxygen blender.
 - To INCREASE the PO_2 , INCREASE the amount of oxygen in the ventilating gas by increasing the FiO_2 on your oxygen blender.
2. PCO_2 is controlled primarily by varying total gas flow rate.
 - To DECREASE PCO_2 , INCREASE the total gas flow rate to increase the amount of CO_2 removed.
 - To INCREASE PCO_2 , DECREASE the total gas flow rate to decrease the amount of CO_2 removed.
3. Patient temperature is controlled by regulating the temperature of water flow into the heat exchanger.

Note: Care should be taken to minimize temperature gradient during cooling and rewarming phases of bypass.

Termination of Bypass

1. Terminate bypass as individual case and patient condition indicates.
2. Turn gas flow off prior to turning off blood flow.
3. Turn off water flow to heat exchanger.
4. After terminating bypass, continue to slowly recirculate through oxygenator using recirculation line and/or reconnected arterial venous loop.

Emergency Oxygenator Replacement

A standby AFFINITY® NT Hollow Fiber Oxygenator should always be available during perfusion.

1. With the cardiotomy reservoir outlet line clamped, add volume if required to prime replacement oxygenator.
2. Turn off water flow to the heat exchanger, clamp and remove water lines. Remove oxygen line from gas inlet port.
3. Clamp venous return line at venous reservoir. Turn off blood pump and clamp arterial outlet.
4. Remove any monitoring/sample line from oxygenator.
5. Double clamp recirculation line, oxygenator inlet and outlet tubing and cut between clamps, leaving adequate lengths for reconnection.
6. Remove old oxygenator from holder.
7. Place replacement oxygenator in holder. Connect and band all blood lines. Connect oxygen line.
8. Connect water lines, unclamp, turn on water source and inspect for water leaks.
9. Remove clamps from cardiotomy line, recirculation line, and the oxygenator inlet.

Warning: Do not unclamp arterial outlet line or venous return line at this time.

10. Turn on blood pump and slowly fill the oxygenator.
11. Increase blood flow to 5 liters/minute through recirculation line.
12. Vent air bubbles from venous reservoir.
13. Ensure entire system is free of leaks and gas bubbles.
14. Remove all clamps from venous and arterial lines, clamp recirculation line and reinitiate bypass.

Important Notice – Limited Warranty (U.S. Customers Only)

- A. This **LIMITED WARRANTY** provides the following assurance to the patient who receives a Medtronic® AFFINITY® NT Hollow Fiber Oxygenator with Plasma Resistant Fiber, hereafter referred to as the "Product."
- (1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its "Use By" date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the Purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
 - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its "Use By" date.
 - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
 - (3) The Product must not have been altered or subjected to misuse, abuse or accident.
 - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) Except as expressly provided by this **LIMITED WARRANTY**, Medtronic is not responsible for any direct, incidental or consequential damages based on any defect, failure or malfunction of the product, whether the claim is based on warranty, contract, tort or otherwise.
 - (2) **THIS LIMITED WARRANTY IS MADE ONLY TO THE PATIENT IN WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.**
 - (3) The exclusions and limitation set out above are not intended to, and should not be constructed so as to contravene mandatory provisions of applicable law. If any part or term of this Limited Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid. This **LIMITED WARRANTY** gives the patient specific legal rights. The patient may also have other rights which vary from state to state.
 - (4) No person has any authority to bind Medtronic to any representation, condition or warranty except this **LIMITED WARRANTY**.
Note: *This **LIMITED WARRANTY** is provided by Medtronic, Inc. 710 Medtronic Parkway, Minneapolis, MN 55432-5604. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for exact terms of the **LIMITED WARRANTY**.*

Important Notice – Limited Warranty (For Countries Outside the U.S.)

- A. This **LIMITED WARRANTY** provides assurance for the purchaser who receives a Medtronic® AFFINITY® NT Hollow Fiber Oxygenator with Plasma Resistant Fiber, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.
 The Warnings contained in the product labeling are considered an integral part of this **LIMITED WARRANTY**. Contact your local Medtronic representative to obtain information on how to process a claim under this **LIMITED WARRANTY**.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its "Use By" date.
 - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
 - (3) The Product may not have been used for any other patient.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
 - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

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Alleviating Pain · Restoring Health · Extending Life

Europe

Europe/Africa/Middle East Headquarters

Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale 84
CH - 1131 Tolochenaz
Switzerland
Internet: www.medtronic.co.uk
Tel. 41-21-802-7000
Fax 41-21-802-7900

Medtronic E.C. Authorized Representative/Distributed by

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. 31-45-566-8000
Fax 31-45-566-8668

Asia-Pacific

Japan

Medtronic Japan
Comodio Shiodome 5F
2-14-1 Higashi-Shimbashi, Minato-ku
Tokyo 105-0021
Japan
Tel. 3-6430-2011
Fax 3-6430-7140

Australia

Medtronic Australasia Pty. Ltd.
Business Address:
97 Waterloo Road
North Ryde NSW 2113
Mailing Address:
PO Box 945
North Ryde NSW 1670
Australia
Tel. 61-2-9857-9000
Fax 61-2-9878-5100

Asia

Medtronic International Ltd.
Suite 1602 16/F, Manulife Plaza
The Lee Gardens, 33 Hysan Avenue
Causeway Bay
Hong Kong
Tel. 852-2891-4068
Fax 852-2591-0313

Americas

Latin America

Medtronic Latin America
3750 NW 87th Avenue
Suite 700
Miami, FL 33178
Tel. 305-500-9328
Fax 786-709-4244

Canada

Medtronic of Canada Ltd.
6733 Kitimat Road, Mississauga
Ontario L5N 1W3
Canada
Tel. 905-826-6020
Fax 905-826-6620
Toll-free in Canada:
1-800-268-5346

United States

Manufacturer:

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Internet: www.medtronic.com or
www.perfusionsystems.com
Tel. 763-514-4000
Fax 763-391-9100
Toll-free in the USA:
1-800-328-2518
(24-hour consultation service)

Medtronic Perfusion Systems

7611 Northland Drive
Minneapolis, MN 55428-9947
Tel. 763-391-9000
Fax 763-391-9100
Customer Service and Product Orders
Toll-free: 1-800-854-3570
Visit Medtronic at
www.perfusionsystems.com

* M 9 3 0 6 2 3 A 0 0 1 *

© Medtronic, Inc. 2005, 2007
All Rights Reserved

M930623A001 Rev 1.0