# Medtronic

## Affinity Pixie™

Oxygenation System with Cortiva™ BioActive Surface  
Система за оксигенация с биоактивна повърхност Cortiva™ BioActive  
Oksyličovací systém s bioaktivním povrchem Cortiva™  
Oxygenationssystem med Cortiva™ BioActive-overflade  
Oxygenator mit bioaktiver Cortiva™ Oberfläche  
Σύστημα οξυγόνωσης με βιοενεργή επίστρωση Cortiva™  
Sistema de oxigenación con superficie biocompatible Cortiva™  
Cortiva™ bioaktiivse pinnakattega oksügenisatsioonisüsteem  
Bioaktiivisella Cortiva™-pinnoitteella päällystetty hapetinjärjestelmä  
Système d'oxygénation avec surface bioactive Cortiva™  
Oksigenacijski sustav s bioaktivnim premazom Cortiva™  
Oxygenátor-rendszer Cortiva™ bioaktív felülettel  
Sistema di ossigenazione con superficie bioattiva Cortiva™  
Oksigenacijos sistema su „Cortiva™“ bioaktyviaja danga  
Oksigenācijas sistēma ar Cortiva™ bioloģiski aktīvo pārklājumu  
Oxygenatiesysteem met bioactieve Cortiva™-oppervlaktelaag  
Oxygeneringssystem med Cortiva™ bioaktiv overflate  
System oksygenacji z powłoką bioaktywną Cortiva™  
Sistema de oxigenação com superfície bioativa Cortiva™  
Sistem de oxigenare cu suprafață bioactivă Cortiva™  
Система оксигенации с биоактивным покрытием Cortiva™  
Oxygenačný systém s bioaktívnou povrchovou úpravou Cortiva™  
Oksigenacijski sistem z biološko aktivno površino Cortiva™  
Sistem za oksigenizaciju sa Cortiva™ BioActive površinom  
Oxygeneringssystem med Cortiva™ bioaktiv yta  
Cortiva™ Biyo Aktif Yüzeyli Oksijenasyon Sistemi  
CBP211, CBP241

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning •  
Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet •  
Mode d'emploi • Upute za upotrebu • Használati útmutató • Istruzioni per l'uso •  
Naudojimo instrukcijos • Lietošanas pamācība • Gebruiksaanwijzing • Bruksanvisning •  
Instrukcja użytkowania • Instruções de utilização • Instrucțiuni de utilizare •  
инструкция по применению • Pokyny na používanie • Navodila za uporabo •  
Uputstva za upotrebu • Bruksanvisning • Kullanim Talimatları

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# Affinity Pixie™

## Oxygenation System with Cortiva™ BioActive Surface

### 1. Model

|        |  |
|--------|--|
| CBP211 | Hollow fiber oxygenator with Cortiva bioactive surface                                 |
| CBP241 | Hollow fiber oxygenator and cardiotomy/venous reservoir with Cortiva bioactive surface |

### Accessories available separately

|         |  |
|---------|--|
| ATP210  | Affinity™ temperature probe                |
| AUH2093 | Affinity Orbit™ holder system              |
| AUH2070 | Affinity Orbit™ oxygenator arm replacement |
| AMH2014 | Affinity™ manifold holder                  |

### 2. Description

This product is a single-use device, which is nontoxic, nonpyrogenic, and supplied STERILE in individual packaging. It is sterilized by ethylene oxide.

#### 2.1. Hollow fiber oxygenator

The Affinity Pixie hollow fiber oxygenator with Cortiva bioactive surface is a single-use, microporous, hollow-fiber, gas-exchange device with plasma-resistant fiber and integrated heat exchanger. Products coated with Cortiva bioactive surface include a "CB" prefix in the model number. Additional information regarding Cortiva bioactive surface is available in Section 9.

The following accessories are packaged for use with the oxygenator (Figure 3):

- Two 0.5 cm (3/16 in) adapters for blood inlet and outlet ports
- One recirculation line — 0.5 cm (3/16 in)

#### 2.2. Cardiotomy/venous reservoir

The Affinity Pixie cardiotomy/venous reservoir (CVR) with Cortiva bioactive surface is a single-use device designed to collect and store blood during extracorporeal circulation. Cardiotomy blood is collected, filtered, and defoamed before mixing with the venous blood, which is filtered. Products coated with Cortiva bioactive surface include a "CB" prefix in the model number. Additional information regarding Cortiva bioactive surface is available in Section 9. Coverage is verified on all primary blood-contacting surfaces and bioactivity is verified on the venous flow path surfaces.

The following accessories are packaged for use with the CVR (Figure 3):

- Two 0.5 cm (3/16 in) adapters for venous return inlet and outlet ports
- Two 1.0 cm (3/8 in) adapters for venous return inlet and outlet ports
- Two flexible luer lock adapters
- One sampling manifold assembly

### 3. Specifications

#### Oxygenator

|                                  |  |
|----------------------------------|--|
| Oxygenator membrane type         | Microporous polypropylene hollow fiber |
| Maximum gas:blood ratio          | 2:1                                    |
| Heat exchanger                   | Polyethylene terephthalate (PET)       |
| Oxygenator membrane surface area | 0.67 m <sup>2</sup>                    |
| Static priming volume            | 48 mL                                  |
| Recommended blood flow range     | 0.1 to 2.0 L/min                       |
| Maximum water pressure           | 206 kPa (1550 mmHg)                    |
| Maximum rated blood pressure     | 100 kPa (750 mmHg)                     |
| Maximum transit temperature      | 40°C (104°F)                           |

Store the product at room temperature.

#### Cardiotomy/venous reservoir

|                              |                                      |
|------------------------------|--------------------------------------|
| Reservoir volume capacity    | 1200 mL                              |
| Recommended blood flow range | 0.1 to 2.0 L/min                     |
| Maximum cardiotomy flow rate | 2.0 L/min                            |
| Minimum operating level      | 20 mL                                |
| Cardiotomy filter            | 30 µm nominal polyester depth filter |
| Venous inlet screen          | 64 µm                                |
| Maximum rated pressure       | +20 mmHg / -100 mmHg                 |

## Cardiotomy/venous reservoir

Positive/negative pressure relief valve

<5 mmHg positive / >60 mmHg vacuum

Maximum transit temperature

40°C (104°F)

Store the product at room temperature.

## 4. Indications for use

The Affinity Pixie hollow fiber oxygenator and cardiotomy/venous reservoir with Cortiva bioactive surface are indicated for use for neonate, infant, and small pediatric patients undergoing cardiopulmonary bypass (CPB) procedures requiring a blood flow rate up to 2.0 L/min.

The Affinity Pixie hollow fiber oxygenator is intended to be used in an extracorporeal perfusion circuit to oxygenate and remove carbon dioxide from the blood and to cool or warm the blood during routine CPB procedures up to 6 hours in duration. The Affinity Pixie hollow fiber oxygenator is suitable for transfer of the volatile anesthetics isoflurane, sevoflurane, desflurane, and enflurane by administration through the oxygenator gas inlet by means of a suitable gas vaporizer.

The Affinity Pixie cardiotomy/venous reservoir is intended to be used in an extracorporeal perfusion circuit to collect venous and cardiotomy suctioned blood during routine cardiopulmonary procedures up to 6 hours in duration. The CVR is also intended for use during vacuum-assisted venous drainage (VAVD) procedures.

## 5. Contraindications

Use the device only as indicated.

## 6. Warnings

- Read all warnings, precautions, and instructions for use carefully before use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**
- Only persons thoroughly trained in CPB procedures should use this device. Operation of each device requires constant supervision by qualified personnel for patient safety.
- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device before use. Do not use the device if the package is opened or damaged, the device is damaged, or the protective caps are not in place. **Note:** The oxygenator water inlet, water outlet, and gas outlet ports do not have protective caps.
- Clear all gas emboli from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- If air leaks are observed during priming or operation, this may result in air embolism to the patient or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the device if these conditions are observed.
- In any extracorporeal circuit, the use of a prebypass and arterial line filter is recommended.
- Monitoring of circuit pressure is recommended.
- The blood phase pressure should be greater than the gas phase pressure at all times in the oxygenator.
- The oxygenator and CVR should not be operated outside the recommended parameter ranges.
- Total combined flow into the cardiotomy filter should not exceed 2.0 L/min.
- Do not exceed blood pressure of 100 kPa (750 mmHg) in the blood pathway.
- Once an oxygenator is primed with blood, adequate anticoagulation should be maintained per institution CPB protocol and the blood pathway should be constantly recirculated at a low flow rate.
- After initiation of bypass, blood flow should circulate through oxygenator and CVR at all times within the recommended blood flow range during CPB, except during emergency oxygenator replacement.
- Heat exchanger temperature should not exceed 42°C (107°F).
- Follow institution CPB protocol when using disinfectants in the heater/cooler during bypass. The integrity of the water path has been verified with hydrogen peroxide (330 ppm). Contact Medtronic for information regarding the use of additional disinfectants.
- Band all tubing connections for added protection against disconnection.
- Ensure that unused ports remain capped to prevent contamination and that the caps are tightened to avoid leakage.
- A replacement oxygenator should be readily available during CPB.
- Do not obstruct the gas outlet vents to prevent gas side pressures from exceeding blood side pressures.
- Do not occlude or obstruct the vent port of the CVR during operation (except when using VAVD).
- Gas transfer characteristics of the oxygenator are not significantly affected by concentrations of anesthetic gas up to 2%. Above 2% anesthetic gas concentrations, adjustments of FiO<sub>2</sub> and gas flow rates may be needed to achieve desired gas transfer performance.
- Gas transfer rates may change over time and adjustments of FiO<sub>2</sub> and gas flow rates may be needed to achieve desired gas transfer performance.

- The use of a blood level sensing mechanism during operation of this device is recommended.
- Ensure that the venous reservoir outlet is always positioned above the highest point in the membrane compartment of the oxygenator.
- CPB must be carefully and constantly monitored.
- Attach tubing in such a manner as to prevent kinks or restrictions that may alter the blood, water, or gas flow.
- Do not allow alcohol, alcohol-based fluid, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come in contact with the device as they may jeopardize the structural integrity.

## 7. Precautions

- Each device is sterilized using ethylene oxide.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device which could result in patient injury, illness, or death.
- Use aseptic technique in all procedures.
- Use of a gas scavenging system during transfer of anesthetic gas through the oxygenator should be considered.
- Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures. Adequate anticoagulation should be maintained per institution CPB protocol.
- It is the responsibility of the user to dispose of the devices in accordance with local regulations and hospital procedures.

## 8. Adverse effects

The following known adverse effects are associated with the use of the product: blood loss, coagulopathy, death, embolism, excessive blood component activation or thrombogenicity, exsanguination, hemolysis, hemolytic anemia, hypercarbia, hypotension, hypovolemia, hypoxia, infection, ischemia, neurological dysfunction, and organ dysfunction.

## 9. Information for products with Cortiva bioactive surface

The primary blood-contacting surfaces of the product are coated with Cortiva bioactive surface. This coated surface enhances blood compatibility and provides a blood-contacting surface that is thromboresistant. Cortiva bioactive surface contains nonleaching heparin derived from porcine intestinal mucosa<sup>1</sup>.

**Caution:** Product coated with Cortiva bioactive surface is intended for single use only. Resterilization may adversely affect the Cortiva bioactive surface.

**Caution:** Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

## 10. Instructions for use

### 10.1. System setup

1. Remove the devices from the packaging carefully to ensure sterile fluid pathway.
  - Warning:** Use aseptic technique during all stages of setup and use of this product.
  - Warning:** Before removing the devices from the packaging, inspect the packaging and product for damage. If the packaging or product is damaged, do not use; sterility may have been compromised or performance may be affected.
2. If using the Affinity Orbit holder system (sold separately), attach the holder mast clamp to the heart-lung machine mast at the desired height and tighten the mast clamp knob.
3. Turn the locking lever perpendicular to the holder arm to unlock it.
4. Slide the oxygenator arm onto the post if it is not already in place.
5. Slide the oxygenator holder connector onto the bottom of the oxygenator arm. A click will indicate that the oxygenator is firmly seated.
  - Note:** To remove the oxygenator, press the holder release tab and pull the oxygenator away from the connector arm.
6. If using the Affinity Pixie CVR, slide the CVR holder sleeve onto the holder post.
7. Position the oxygenator and CVR as desired for procedure. The oxygenator may be rotated clockwise or counterclockwise and the oxygenator arm may be moved laterally. The CVR may be rotated clockwise or counterclockwise.
8. Turn the locking lever parallel to the holder arm to lock components in place.
9. Turn the holder arm knob to lock the holder arm in place. See Figure 4 for final system assembly.

<sup>1</sup> Nonleaching is defined as heparin at a level  $\leq 0.1$  IU/mL, as measured through clinically relevant extraction conditions.

10. Remove the obturator cap from the CVR positive/negative pressure relief valve.
11. Connect all blood, gas, and suction lines to their appropriate connections per institution CPB protocol. If necessary, attach an appropriately sized adapter to the blood inlet or outlet port to accommodate 0.5 cm (3/16 in) or 1.0 cm (3/8 in) tubing.
12. The sampling manifold should be attached so that the tubing segment with the one-way duckbill valve is positioned between the oxygenator arterial sample port and manifold and that arterial blood flows through the duckbill valve toward the manifold.
13. Connect the recirculation line to the recirculation port on the oxygenator. Connect the other end of the recirculation line to a luer port on the CVR.  
**Warning:** If not used, the recirculation port should be sealed tightly with a luer cap.
14. Connect quick-connect 1.2 cm (1/2 in) water lines to inlet and outlet ports on the oxygenator. Start the water flow and check for leaks from the water compartment to the blood compartment. Discard oxygenator if water is present in the blood compartment and set up the standby oxygenator.
15. Ensure that the CVR is vented to atmosphere by removing the vent port cap.  
**Note:** A CO<sub>2</sub> flush of the oxygenator is recommended before priming.

### 10.2. Priming and recirculation

1. Fill the CVR with the appropriate volume of priming solution per institution CPB protocol.
2. Remove all bubbles in the section of tubing (roller pump head/boot) or centrifugal pump (if used) from the CVR outlet to the oxygenator inlet.
3. Gravity prime the oxygenator or prime by forward pump flow per institution CPB protocol.
4. Gradually begin recirculation flow using a prebypass filter per institution CPB protocol.
5. Ensure that the oxygenator is free of bubbles.
6. Add additional priming solution, if needed, to prime the rest of the extracorporeal circuit.
7. Upon completion of priming and debubbling, gradually reduce blood flow and stop pump, close all purge lines, and clamp arterial, venous, and recirculation lines.
8. Ensure that the entire extracorporeal circuit is free of bubbles before initiating bypass.  
**Note:** Priming solution may be prewarmed through heat exchanger before initiating bypass.

### 10.3. Initiation of bypass

**Warning:** The blood phase pressure should be greater than the gas phase pressure at all times.

1. Ensure that the gas outlet port is not obstructed.
2. Check for adequate levels of anticoagulation before initiating bypass.
3. Remove arterial and venous clamps and gradually increase blood flow. Then initiate gas flow with appropriate gas settings per institution CPB protocol and clinical judgment specific to the patient and procedure.
4. Adjust water temperature to meet clinical demands.

### 10.4. Operation during bypass

1. Arterial PO<sub>2</sub> is controlled by varying the percent concentration of oxygen present in the ventilating gas.
  - To **decrease** the PO<sub>2</sub>, **decrease** the percent of oxygen in the ventilating gas by adjusting the FiO<sub>2</sub> on your gas blender.
  - To **increase** the PO<sub>2</sub>, **increase** the percent of oxygen in the ventilating gas by adjusting the FiO<sub>2</sub> on your gas blender.
2. PCO<sub>2</sub> is controlled primarily by varying the total gas flow rate.
  - To **decrease** PCO<sub>2</sub>, **increase** the total gas flow rate to increase the amount of CO<sub>2</sub> removed.
  - To **increase** PCO<sub>2</sub>, **decrease** the total gas flow rate to decrease the amount of CO<sub>2</sub> removed.
3. Patient temperature is controlled by regulating the temperature of the water flow in the heat exchanger.  
**Warning:** After initiation of bypass, blood flow should circulate through oxygenator and CVR at all times within blood flow range during CPB, except during emergency oxygenator replacement.  
**Note:** Care should be taken to minimize temperature gradient during cooling and warming phases of bypass.

### 10.5. Termination of bypass

Terminate bypass per institution CPB protocol as individual case and patient condition indicates.

### 10.6. Emergency oxygenator replacement

A standby oxygenator and CVR should always be available during CPB.

1. Turn off water flow to the heat exchanger, then clamp and remove the water lines.

2. Clamp the venous line at the CVR. Turn off the arterial roller pump and clamp the arterial line. (If using a centrifugal pump, clamp the arterial line first, before stopping the pump.)
3. Remove the oxygen line from the gas inlet port.
4. Remove any monitoring/sampling lines from the oxygenator.
5. Double clamp the recirculation line, oxygenator inlet, and outlet tubing. Cut between the clamps leaving adequate lengths for reconnections.
6. Detach the oxygenator from the holder by pressing the release tab.
7. Attach the replacement oxygenator onto the holder.
8. Reconnect the recirculation, oxygenator, and outlet tubing. Ensure that all connections are secure.
9. Reconnect the oxygen line to the gas inlet port.
10. Connect and unclamp the water lines, then turn on water source and inspect for water leaks.

**Warning:** Do not unclamp the arterial line or venous return line at this time.

11. With sufficient volume in the CVR, turn on the blood pump and slowly fill the oxygenator.
12. Increase blood flow through the recirculation line.
13. Stop the pump and clamp the recirculation line.
14. Ensure that the entire system is free of leaks and gas bubbles.
15. Remove the clamps from the venous and arterial lines and reinstate bypass.

## 11. Cardiomy/venous reservoir, vacuum-assisted venous drainage (VAVD)

### 11.1. Warnings

- Do not use a vacuum regulator with maximum negative vacuum pressure greater than  $-20$  kPa ( $-150$  mmHg). Use of excessive negative pressure may lead to hemolysis.
- Do not occlude or obstruct the vent port of the CVR during operation (except when using VAVD).
- Draw only intended drug amounts into the syringe that will be used for drug administration into the CVR during VAVD. Negative pressure in the CVR could draw more drug from the syringe into the reservoir than intended and lead to excess drug delivery.
- All A/V shunt lines must be closed to atmosphere before stopping the pump or using low blood flow rates during VAVD. Failure to do so may allow air to be drawn into the blood side of the oxygenator from the fiber.
- Appropriate occlusion of the arterial roller pump head is necessary when using VAVD. Failure to appropriately occlude the pump may allow air to be drawn into the blood side of the oxygenator from the fiber.
- An arterial filter should be used during VAVD techniques to minimize potential delivery of air emboli.
- Gradually return to atmospheric pressure when ending use of VAVD. A sudden change in pressure may result in turbulent blood flow inside the CVR.
- The line between the oxygenator and centrifugal pump (if used) must be clamped before stopping the pump during use of VAVD. Failure to clamp off the arterial line may result in air drawn into the blood side of the oxygenator from the fiber. Use of an arterial line one-way valve between the oxygenator and centrifugal pump is recommended.

### 11.2. Precautions

- A controlled vacuum regulator is required during VAVD.
- A positive/negative pressure relief valve is required during VAVD (included, see Figure 2).
- Use of a pressure measurement device on the CVR and an additional negative pressure relief valve that operates at  $-20$  kPa ( $-150$  mm Hg) should be considered.
- A vapor trap is required during use of VAVD due to condensation development.

### 11.3. VAVD instructions for use

1. Follow previously listed instructions for use for system setup with the following modifications as listed in steps 2 through 3.
2. Attach a calibrated means of monitoring pressures to either the venous line or venous reservoir.
3. All manifold sampling ports must have nonvented caps when using VAVD.
4. Follow previously listed instructions for use for priming with the following additions as listed in steps 5 through 8 below.
5. Before initiating bypass, prepare the vacuum portion of the setup. Attach a vented vacuum line to the CVR vent port and to the vacuum regulator with vapor trap.

**Warning:** Clamp the quick-prime line before application.

**Warning:** Do not allow the vapor trap to become completely filled during use. This may allow fluid to enter the vacuum controller, or may prevent the CVR from being vented to atmosphere when no vacuum is applied.

6. To vent the vacuum line, ensure that the Y-connector with side tubing is easily accessible and open to atmosphere. This may be clamped and unclamped regularly throughout the procedure to apply or stop vacuum.  
**Caution:** Do not exceed vacuum levels greater than  $-6.67$  kPa ( $-50$  mmHg) as measured in the venous line.
7. Ensure that all access ports on the CVR and manifold are properly sealed and are secure before initiating bypass.  
**Caution:** Vacuum will result in excess atmospheric air being entrained into the CVR and potentially into the blood pathway with opened blood pathway ports. Vacuum may result in increased flow through the sampling manifold during use.
8. Use vacuum assistance as needed to institute and maintain CPB.  
**Warning:** Do not apply vacuum to the venous reservoir when there is no forward blood flow through the oxygenator. This applies to arterial centrifugal and roller pumps (roller may not be occlusive in all positions). This will prevent air from being pulled across the membrane into the blood path by the reservoir vacuum.  
**Warning:** Do not allow the CVR to become pressurized as this could obstruct venous drainage, force air retrograde into the patient or into the blood path of the oxygenator.

## 12. Additional information

The following additional information and data are available upon request: sterilization method, materials of the blood pathway, plasma leakage across the semipermeable membrane, blood and gas pathway pressure drop at the maximum blood and gas flow rates, blood cell damage, particle release, air-handling capability, antifoam characteristics, breakthrough volume, filtration efficiency, relevant tolerances for data in the instructions for use, protocol summaries, and test methods used to determine coverage, leaching, and bioactivity.

## 13. Important notice — limited warranty (for countries outside the US)

**Important: This limited warranty does not apply in Australia.**

- A. This **LIMITED WARRANTY** provides assurance for the purchaser who receives an Affinity Pixie hollow fiber oxygenator and/or cardiotomy/venous reservoir with Cortiva bioactive surface, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.  
The Warnings contained in the product labeling are considered an integral part of this **LIMITED WARRANTY**. Contact your local Medtronic representative to obtain information on how to process a claim under this **LIMITED WARRANTY**.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
  - (1) The Product must be used prior to its Use-by date.
  - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
  - (3) The Product may not have been used for any other patient.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
  - (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
  - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

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© 2018 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. ™\* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
Medtronic B.V.
Medtronic of Canada Ltd