Medtronic

Affinity™ AF100

Arterial Filter with Cortiva™ BioActive Surface
Артериален филтър с биоактивна повърхност Cortiva™ BioActive
Arteriální filtr s bioaktivním povrchem Cortiva™
Arteriefilter med Cortiva™ BioActive-overflade
Arterieller Filter mit bioaktiver Cortiva™ Oberfläche
Αρτηριακό φίλτρο με βιοενεργή επίστρωση Cortiva™
Filtro arterial con superficie biocompatible Cortiva™
Cortiva™ bioaktiivse pinnakattega arteriaalne filter
Bioaktiivisella Cortiva™-pinnoitteella päällystetty valtimoverisuodatin
Filtre artériel avec surface bioactive Cortiva™
Arterijski filter s bioaktivnim premazom Cortiva™
Artériás szűrő Cortiva™ bioaktív felülettel
Filtro arterioso con superficie bioactiva Cortiva™
Arterinis filtrs su „Cortiva™“ bioaktyviaja danga
Arteriālais filtrs ar Cortiva™ bioloģiski aktīvo pārklājumu
Arteriële filter met bioactieve Cortiva™-oppervlaktelaag
Arteriefilter med Cortiva™ bioaktiv overflate
Filtr tętniczy z powłoką bioaktywną Cortiva™
Filtro arterial com superficie bioativa Cortiva™
Filtru arterial cu suprafață bioactivă Cortiva™
Артериальный фильтр с биоактивным покрытием Cortiva™
Arteriálny filter s bioaktívnuou povrchovou úpravou Cortiva™
Filter za arterijsko kri z biološko aktivno površino Cortiva™
Arterijski filter sa Cortiva™ BioActive površinom
Artärfilter med Cortiva™ bioaktiv yta
Cortiva™ Biyo Aktif Yüzeyli Arteriyel Filtre

CB851

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Istruzioni per l'uso • Naudojimo instrukcijos • Lietošanas pamācība • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по эксплуатации • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları

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Affinity™ AF100

Arterial Filter with Cortiva™ BioActive Surface

1. Model

CB851 30 µm arterial filter with Cortiva bioactive surface

2. Product description

The AF100 is designed to filter from the circuit microemboli larger than the specified micron size for periods up to six hours during cardiopulmonary bypass surgery.

Products coated with Cortiva bioactive surface include a "CB" prefix in the model number. Additional information regarding Cortiva bioactive surface is available in Section 7.

The device is single-use, nontoxic, nonpyrogenic, and is supplied sterile in individual packaging. Sterilized by ethylene oxide.

Note: The Affinity arterial filter holders (HR4014 or AUH4014) are sold separately.

3. Specifications

Recommended blood flow range	1–7 L/min
Screen size	30 µm nominal
Prime volume	100 mL
Size	Adult
Inlet connector	1.0 cm (3/8 in)
Outlet connector	1.0 cm (3/8 in)
Vent port	Standard female luer lock
Storage temperature	–30°C–40°C (–22°F–105°F)
Storage humidity	30–75% noncondensing
Maximum operating pressure	750 mmHg (14.5 psi)

4. Indications for use

The AF100 is indicated for use in cardiopulmonary bypass procedures up to 6 hours in duration for the removal of particulate and gaseous microemboli.

5. Contraindications

Use the device only as indicated.

6. Warnings

Read all warnings, precautions, and instructions for use carefully prior to use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**

- Only persons thoroughly trained in cardiopulmonary bypass procedures should use this device. Operation of each device requires supervision by qualified personnel for patient safety. Perfusion must be carefully and constantly monitored.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Device leakage during operation may result in air embolism and/or fluid loss.
- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use if the package is opened or damaged or the protective caps are not in place, as the sterility of the device may have been compromised and/or performance may be affected.
- Tubing should be attached in such a manner as to prevent kinks or restrictions that may alter blood flow.
- Use aseptic technique in all procedures.
- Do not use at flow rates outside of the specification (1–7 L/min).
- Do not use beyond the specified duration (6 hours).
- Do not allow alcohol, alcohol-based fluids, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come in contact with the product because they may jeopardize the structural integrity.
- Clear all gas emboli from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.

- Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

7. Precautions

- Refer to the package labeling for storage temperature requirements.
- If a device leak is observed during priming, replace the device.
- If thrombus is observed in the device during use, replace the device.
- Avoid agitating, inverting, or tapping the device once cardiopulmonary bypass is initiated.
- The female luer lock port must be closed when the pump is stopped to prevent backflow of blood and possible exsanguination of the patient.
- All blood tubing connections should be banded for added protection against high fluid pressures.
- Before blood is introduced into any extracorporeal circuit, the use of a prebypass filter should be considered.
- The use of a bubble detector downstream from the arterial line filter is strongly recommended.
- Dispose of the devices in accordance with local hospital, administrative, and/or other government policies.

8. Adverse effects

Possible side effects include, but are not limited to, air embolism, thromboembolic phenomena, platelet dysfunction, blood loss, circulatory compromise, renal compromise, hypotension, activation (coagulation/complement), hemolysis, and death.

9. Information for products with Cortiva bioactive surface

The primary blood-contacting surfaces of the product are coated with Cortiva bioactive surface. This coated surface enhances blood compatibility and provides a blood-contacting surface that is thromboresistant. Cortiva bioactive surface contains nonleaching heparin derived from porcine intestinal mucosa¹.

Caution: Product coated with Cortiva bioactive surface is intended for single use only. Resterilization may adversely affect the Cortiva bioactive surface.

Caution: Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

Warning: Do not store product coated with Cortiva bioactive surface above 40 °C (104 °F).

10. Instructions for use

10.1. Setup

Note: Refer to the CPB pump speed controller Instructions for Use for details related to system set-up, priming, operating parameters, and emergency back-up equipment.

1. Remove the device(s) from the packaging carefully to ensure a sterile fluid pathway.
Caution: Ensure aseptic technique is used during all stages of setup and use of this device.
2. Insert the arterial filter into the arterial line. The use of tie bands is recommended. Ensure proper direction of fluid flow.
3. Position the holder disk of the AF100 in line with the clamp of the filter holder (Figure 3) and slide in the AF100. Ensure the AF100 is securely in its holder with the blood outlet port pointing down.
4. A filter bypass line may be used to facilitate priming and to provide an alternate blood flow path around the filter (Figure 2).
5. A purge line with stopcock and one-way valve may be used to vent the arterial filter.
Caution: Always purge the arterial filter through the female luer lock port to a low pressure port to prevent back flow of air through the filter purge line.

10.2. Priming

Note: To aid priming, 100% medical grade CO₂ may be used to displace air and aid in debubbling.

1. Prime the filter per institutional protocol.
2. The filter may be inverted and tapped gently to aid in debubbling.
3. Visually inspect the filter and ensure that the filter is completely debubbled prior to initiating bypass.

¹ Nonleaching is defined as heparin at a level ≤ 0.1 IU/mL, as measured through clinically relevant extraction conditions.

4. Double check all connections and visually inspect for leaks.

Caution: Before blood is introduced into any extracorporeal circuit, the use of a prebypass filter should be considered.

10.3. Operating

Warning: Ensure the filter is in an upright position. If a one-way valve is employed, ensure proper orientation to allow purging of air.

1. Check for adequate levels of anticoagulation before initiating bypass.
2. During extracorporeal circulation, the stopcock should be operated in the open position.

Caution: The stopcock of the arterial purge line must be closed when the pump is stopped to prevent retrograde blood flow from the patient through the purge line.

11. Additional information

The following information and data are available upon request: a list of all materials used in the construction of the arterial line filter, blood pathway materials, blood cell damage, air handling capacity, blood pathway pressure drop, and relevant tolerances for specifications in the instructions for use.

12. Important notice - limited warranty (US customers only)

- A. This Limited Warranty provides the following assurance to the patient who receives an Affinity AF100 arterial filter, hereafter referred to as the "Product".
 - (1) Should the product fail to function within normal tolerances due to a defect in materials or workmanship prior to its Use-by date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
 - (2) As used herein, purchase price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the Limited Warranty, these conditions must be met:
 - (1) The Product must be used prior to its Use-by date.
 - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
 - (3) The Product must not have been altered or subjected to misuse, abuse or accident.
 - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.
- C. This Limited Warranty is limited to its express terms. In particular:
 - (1) Except as expressly provided by this LIMITED WARRANTY, Medtronic is not responsible for any direct, incidental or consequential damages based on any defect, failure, or malfunction of the product, whether the claim is based on warranty, contract, tort, or otherwise.
 - (2) THIS LIMITED WARRANTY IS MADE ONLY TO THE PATIENT FOR WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER, OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
 - (3) The exclusions and limitations set out above are not intended to, and should not be constructed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid. This LIMITED WARRANTY gives the patient specific legal rights. The patient may also have other rights which vary from state to state.
 - (4) No person has any authority to bind Medtronic to any representation, condition, or warranty except this LIMITED WARRANTY

Note: This LIMITED WARRANTY is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for exact terms of the LIMITED WARRANTY.

13. Important notice - limited warranty (for countries outside the US)

- A. This LIMITED WARRANTY provides assurance for the purchaser who receives an Affinity AF100 arterial filter, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.
- The Warnings contained in the product labeling are considered an integral part of this LIMITED WARRANTY. Contact your local Medtronic representative to obtain information on how to process a claim under this LIMITED WARRANTY.
- B. To qualify for the LIMITED WARRANTY, these conditions must be met:
- (1) The Product must be used prior to its Use-by date.
 - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
 - (3) The Product may not have been used for any other patient.
- C. This LIMITED WARRANTY is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence or improper handling, improper implantation, or material alteration of the replaced Product.
 - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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© 2012, 2016 Medtronic
Trademarks may be registered and are the property of their respective owners.
Medtronic Australasia Pty Ltd
Medtronic of Canada Ltd