Medtronic

Affinity® AF100

Arterial Filter with Balance® Biosurface
Filtre artériel avec biosurface Balance®
Arterieller Filter mit Balance® Biosurface
Filtro arterial con superficie biocompatible Balance®
Arteriële filter met Balance®-bio-oppervlaktelaag
Filtro arterioso con superficie biopassiva Balance®
Arteriefilter med Balance® biooverflate
Balance®-biopinnoitteella päällystetty valtimoverisuodatin
Artärfilter med Balance® Biosurface
Αρτηριακό φίλτρο με βιοεπιφάνεια Balance®
Arteriefilter med Balance® Bio-overflade
Filtr tętniczy z biopowłoką Balance®
Filtro arterial com bio-superfície Balance®
Balance® Biyo Yüzeyli Arteriyel Filtre
Arteriální filtr s biopovrchem Balance®
Artériás szűrő Balance® biofelülettel
Arteriálny filter s biologickou povrchovou úpravou Balance® Biosurface
Balance® biokattega arteriaalne filter
Αρτηριαλєν φίλτρє с биоповърхност Balance®
Arterijski filter s biopovršinom Balance®
Arteriālais filtrs ar Balance® Biosurface virsmas pārklājumu
Arterinis filtras su „Balance®“ biologine danga
Filtru arterial cu biosuprafață Balance®
Αρτηριαλєν φίλτρє с биопокрытием Balance®
Arterijski filter sa Balance® Biosurface površinom
Filter za arterijsko kri z biološko površino Balance®
BB851

Instructions for Use • Mode d'emploi • Gebrauchsanweisung
Instrucciones de uso • Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning
Käyttöohjeet • Bruksanvisning • Οδηγίες χρήσης • Brugsanvisning
Instrukcja użytkowania • Instruções de utilização • Kullanım Talimatları
Návod k použiti • Használati útmutató • Pokyny na používanie • Kasutusjuhend
Инструкции за употреба • Upute za upotrebu • Lietošanas pamācība
Naudojimo instrukcijos • Instrucțiuni de utilizare • Инструкция по эксплуатации
Uputstva za upotrebu • Navodila za uporabo

USA **Caution:** Federal law (USA) restricts this device to sale by
or on the order of a physician.

Affinity® and Balance® are registered trademarks of Medtronic, Inc.
Affinity® et Balance® sont des marques déposées de Medtronic, Inc.
Affinity® und Balance® sind eingetragene Marken von Medtronic, Inc.
Affinity® y Balance® son marcas comerciales registradas de Medtronic, Inc.
Affinity® en Balance® zijn geregistreerde handelsmerken van Medtronic, Inc.
Affinity® e Balance® sono marchi registrati della Medtronic, Inc.
Affinity® og Balance® er registrerte varemerker for Medtronic, Inc.
Affinity® ja Balance® ovat Medtronic, Inc:n rekisteröityjä tavaramerkkejä.
Affinity® och Balance® är registrerade varumärken som tillhör Medtronic, Inc.
Τα Affinity® και Balance® είναι σήματα κατατεθέντα της Medtronic, Inc.
Affinity® og Balance® er registrerede varemærker tilhørende Medtronic, Inc.
Affinity® i Balance® są zastrzeżonymi znakami towarowymi firmy Medtronic, Inc.
Affinity® e Balance® são marcas comerciais registadas da Medtronic, Inc.
Affinity® ve Balance®, Medtronic, Inc. şirketinin tescilli ticari markalarıdır.
Affinity® a Balance® jsou registrované ochranné známky společnosti Medtronic, Inc.
Az Affinity® és a Balance® a Medtronic, Inc. bejegyzett védjegye.
Affinity® a Balance® sú registrovanými ochrannými známkami spoločnosti Medtronic, Inc.
Affinity® ja Balance® on ettevõtte Medtronic, Inc. registreeritud kaubamärgid.
Affinity® и Balance® са регистрирани търговски марки на Medtronic, Inc.
Affinity® i Balance® registrirani su zaštitni znakovi tvrtke Medtronic, Inc.
Affinity® un Balance® ir Medtronic, Inc. reģistrētas preču zīmes.
„Affinity”™ ir „Balance”™ yra registruotieji „Medtronic, Inc.” prekių ženklai.
Affinity® și Balance® sunt mărci comerciale înregistrate ale Medtronic, Inc.
Affinity® и Balance® являются зарегистрированными товарными знаками корпорации Medtronic, Inc.
Affinity® i Balance® su registrovani žigovi preduzeća Medtronic, Inc.
Affinity® in Balance® sta registrirani blagovni znamki družbe Medtronic, Inc.

Explanation of symbols on product or package / Explication des symboles indiqués sur le produit et l'emballage / Erläuterung der auf den Produkten und der Verpackung befindlichen Symbole / Explicación de los símbolos que aparecen en el producto o en el envase / Verklaring van de symbolen op het product of de verpakking / Spiegazione dei simboli del prodotto o della confezione / Forklaring av symboler på produktet eller pakningen / Tuotteen tai pakkauksen symbolien selitykset / Förklaring av symboler på produkt och förpackning / Επεξήγηση των συμβόλων στο προϊόν ή στη συσκευασία / Forklaring på symboler på produktet og emballagen / Objaśnienia symboli na produkcie lub opakowaniu / Explicação dos símbolos no produto ou na embalagem / Ürün veya ambalaj üzerindeki sembollerin açıklamaları / Vysvětlení symbolů na produktu nebo obalu / A terméken vagy annak csomagolásán lévő szimbólumok magyarázata / Vysvetlivky k symbolom na produkte alebo balení / Tootel või pakendil olevate sümbolite tähendus / Объяснение на символите върху продукта или опаковката / Objašnjenje simbola na proizvodu ili ambalaži / Uz produkta vai iepakojuma atrodamo simbolu skaidrojums / Simbolių ant produkto arba pakuotės paaiškinimai / Legenda simbolurilor de pe produs sau de pe ambalajul produsului / Объяснение символов на продукте или на упаковке / Objašnjenje simbola na proizvodu ili pakovanju / Razlaga simbolov na izdelku ali embalaži

Refer to the device labeling to see which symbols apply to this product / Se référer à la documentation sur le produit pour savoir quels symboles s'appliquent à ce produit / Welche Symbole für dieses Produkt zutreffen, entnehmen Sie bitte dem Etikett auf dem Gerät / Consulte la documentación del dispositivo para comprobar qué símbolos se utilizan con este producto / Controleer de verpakking om te zien welke symbolen op dit product van toepassing zijn / Fare riferimento alle etichette sul dispositivo per sapere quali simboli si applicano a questo prodotto / Se på pakningen for å fastslå hvilke symboler som gjelder for dette produktet / Katso laitteen merkinnöistä, mitkä symbolit koskevat tätä tuotetta / Se dokumentationen för enheten för de symboler som gäller denna produkt / Ανατρέξτε στη σήμανση του προϊόντος για να δείτε ποια σύμβολα ισχύουν για το παρόν προϊόν / Se mærkaten på enheden for de symboler, der gælder for dette produkt / Aby określić, które symbole dotyczą niniejszego produktu, należy zapoznać się z etykietą na opakowaniu produktu. / Consultar a documentação do dispositivo para ver que símbolos se aplicam a este produto / Bu üründe hangi sembollerin geçerli olduğunu görmek için cihazın etiketlerine bakın / Symboly, které se vztahují k tomuto produktu, najdete na štítcích obalu / A termékre vonatkozó szimbólumok a csomagoláson találhatók / Prezrite si označenje zariadenia, aby ste zistili, ktoré symboly sa vzťahujú na tento produkt / Toote puhul kehtivad sümbolid leiate seadme siltidel / Погледнете етикета на устройството, за да видите какви символи се отнасят за този продукт / Na etiketama na uređaju pogledajte koji su simboli primjenjivi na ovaj proizvod / Skatiet marķējumu uz ierīces, lai redzētu kādus simbolus piemēro šim izstrādājumam / Norēdami sužinoti, koki simbolii pažymėtas gaminys, žiūrėkite į pakuotės etiketę / Consultați eticheta dispozitivului pentru simbolurile aplicabile acestui produs / См. на етикетках устройства, какви символи применимы к данному продукту / Pogledajte spoljnu oznaku na pakovanju da biste videli koji se simboli primenjuju na ovaj proizvod / Simboli, ki veljajo za ta izdelek, so na nalepkah na izdelku

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Affinity® AF100

Arterial Filter with Balance® Biosurface

1. Model

BB851 30 µm Arterial Filter with Balance® Biosurface¹

2. Product Description

The AF100 is designed to filter from the circuit microemboli larger than the specified micron size for periods up to six hours during cardiopulmonary bypass surgery.

The AF100 with Balance Biosurface (BB851) is coated with a nonleaching biocompatible surface to reduce platelet activation and adhesion and preserve platelet function.

The device is single-use, nontoxic, nonpyrogenic, supplied STERILE in individual packaging. Sterilized by ethylene oxide.

Note: The Arterial Filter Holders (HR4014 or AUH4014) are sold separately.

3. Specifications

Recommended Blood Flow Range	1–7 L/min
Screen Size	30 µm nominal
Prime Volume	100 mL
Size	Adult
Inlet Connector	1.0 cm (3/8 in)
Outlet Connector	1.0 cm (3/8 in)
Vent Port	Standard Female Luer Lock
Storage Temperature	–30°C–50°C (–22°F–122°F)
Storage Humidity	30–75% noncondensing
Maximum Operating Pressure	750 mm Hg (14.5 psi)

4. Indications for Use

The AF100 is indicated for use in cardiopulmonary bypass procedures up to 6 hours in duration for the removal of particulate and gaseous microemboli.

5. Contraindications

Do not use this device for any purpose other than indicated.

6. Warnings

Read all warnings, precautions, and Instructions for Use carefully prior to use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**

- Only persons thoroughly trained in cardiopulmonary bypass procedures should use this device. Operation of each device requires supervision by qualified personnel for patient safety. Perfusion must be carefully and constantly monitored.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device which could result in patient injury, illness, or death.
- Device leakage during operation may result in air embolism and/or fluid loss.
- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use if the package is opened or damaged or the protective caps are not in place as the sterility of the device may have been compromised and/or performance may be affected.
- Tubing should be attached in such a manner as to prevent kinks or restrictions that may alter blood flow.
- Use aseptic technique in all procedures.
- Do not use at flow rates outside of the specification (1–7 L/min).
- Do not use beyond the specified duration (6 hours).
- Do not allow alcohol, alcohol-based liquids, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come into contact with the device as they may jeopardize its structural integrity.

¹ Technology licensed under agreement from BiolInteractions, Limited, United Kingdom.

- All gas emboli must be cleared from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- Appropriate anticoagulation protocol should be followed per institutional CPB protocol and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.

7. Precautions

- Refer to package labeling for storage temperature requirements.
- If a device leak is observed during priming, replace the device.
- If thrombus is observed in the device during use, replace the device.
- Avoid agitating, inverting, or tapping the device once cardiopulmonary bypass is initiated.
- The female luer lock port must be closed when the pump is stopped to prevent backflow of blood and possible exsanguination of the patient.
- All blood tubing connections should be banded for added protection against high fluid pressures.
- Before blood is introduced into any extracorporeal circuit, the use of a prebypass filter should be considered.
- The use of a bubble detector downstream from the arterial line filter is strongly recommended.
- Dispose of the devices in accordance with local hospital, administrative, and/or other government policies.

8. Adverse Effects

Possible side effects include, but are not limited to, air embolism, thromboembolic phenomena, platelet dysfunction, blood loss, circulatory compromise, renal compromise, hypotension, activation (coagulation/complement), hemolysis, and death.

9. Additional Information for Products with Balance Biosurface

The blood-contacting surfaces of the device are coated with Balance Biosurface to reduce platelet activation and adhesion and preserve platelet function. See Figure 4.

Caution: Appropriate blood anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.

10. Instructions for Use

10.1. Setup

Note: Refer to the CPB pump speed controller Instructions for Use for details related to system set-up, priming, operating parameters, and emergency back-up equipment.

1. Remove the device(s) from the packaging carefully to ensure a sterile fluid pathway.
Caution: Ensure aseptic technique is used during all stages of setup and use of this device.
2. Insert the arterial filter into the arterial line. The use of tie bands is recommended. Ensure proper direction of fluid flow.
3. Position the holder disk of the AF100 in line with the clamp of the filter holder (Figure 3) and slide in the AF100. Ensure the AF100 is securely in its holder with the blood outlet port pointing down.
4. A filter bypass line may be used to facilitate priming and to provide an alternate blood flow path around the filter (Figure 2).
5. A purge line with stopcock and one-way valve may be used to vent the arterial filter.
Caution: Always purge the arterial filter through the female luer lock port to a low pressure port to prevent back flow of air through the filter purge line.

10.2. Priming

Note: To aid priming, 100% medical grade CO₂ may be used to displace air and aid in debubbling.

1. Prime the filter per institutional protocol.
2. The filter may be inverted and tapped gently to aid in debubbling.
3. Visually inspect the filter and ensure that the filter is completely debubbled prior to the initiation of bypass.
4. Double check all connections and visually inspect for leaks.
Caution: Before blood is introduced into any extracorporeal circuit, the use of a prebypass filter should be considered.

10.3. Operating

Caution: Ensure the filter is in an upright position. If a one-way valve is employed, ensure proper orientation to allow purging of air.

1. Check for adequate levels of anticoagulation before initiating bypass.
2. During extracorporeal circulation, the stopcock should be operated in the open position.

Caution: The stopcock of the arterial purge line must be closed when the pump is stopped to prevent retrograde blood flow from the patient through the purge line.

11. Additional Information

The following information and data are available upon request: a list of all materials used in the construction of the arterial line filter, blood pathway materials, blood cell damage, air handling capacity, blood pathway pressure drop, and relevant tolerances for specifications in the Instructions for Use.

12. Important Notice - Limited Warranty (U.S. Customers Only)

- A. This Limited Warranty provides the following assurance to the patient who receives an Affinity[®] AF100 Arterial Filter, hereafter referred to as the "Product".
- (1) Should the product fail to function within normal tolerances due to a defect in materials or workmanship prior to its "Use By" date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the Purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
 - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the Limited Warranty, these conditions must be met:
- (1) The Product must be used prior to its "Use By" date.
 - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
 - (3) The Product must not have been altered or subjected to misuse, abuse or accident.
 - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.
- C. This Limited Warranty is limited to its express terms. In particular:
- (1) Except as expressly provided by this LIMITED WARRANTY, Medtronic is not responsible for any direct, incidental or consequential damages based on any defect, failure or malfunction of the product, whether the claim is based on warranty, contract, tort or otherwise.
 - (2) THIS LIMITED WARRANTY IS MADE ONLY TO THE PATIENT FOR WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
 - (3) The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid. This LIMITED WARRANTY gives the patient specific legal rights. The patient may also have other rights which vary from state to state.
 - (4) No person has any authority to bind Medtronic to any representation, condition or warranty except this LIMITED WARRANTY

Note: This LIMITED WARRANTY is provided by Medtronic, Inc. 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for exact terms of the LIMITED WARRANTY.

13. Important Notice - Limited Warranty (For Countries Outside the U.S.)

- A. This LIMITED WARRANTY provides assurance for the purchaser who receives an Affinity[®] AF100 Arterial Filter, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.

The Warnings contained in the product labeling are considered an integral part of this LIMITED WARRANTY. Contact your local Medtronic representative to obtain information on how to process a claim under this LIMITED WARRANTY.

- B. To qualify for the LIMITED WARRANTY, these conditions must be met:
 - (1) The Product must be used prior to its "Use By" date.
 - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
 - (3) The Product may not have been used for any other patient.
- C. This LIMITED WARRANTY is limited to its express terms. In particular:
 - (1) In no event shall any replacement credit be granted where there is evidence or improper handling, improper implantation or material alteration of the replaced Product.
 - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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© 2012 Medtronic