Medtronic

Affinity™ CP AP40AST

Adapter • Адаптер • Adaptér • Adapter • Adapter • Προσαρμογέας • Adaptador • Adapter • Sovitin • Adaptateur • Prilagodnik • Adapter • Adaptor • Adattatore • Adapteris • Adapteris • Адаптер • Adapter • Adapter • Adapter • Adaptador • Adaptador • Adaptor • Адаптер • Adaptér • Adapter • Adapter • Adapter • Adaptör • Адаптер • Bộ đổi nguồn

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Petunjuk Penggunaan • Istruzioni per l'uso • Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по эксплуатации • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з використання • Hướng dẫn sử dụng

! USA

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
1 Description

The Medtronic Affinity CP adapter Model AP40AST (Figure 1) is intended to be used to enable the Medtronic Affinity CP centrifugal blood pumps (Model AP40, Model BBAP40, and Model CBAP40) to be operated with the Stöckert™ and Sorin™ centrifugal pump systems or the Sarns™ and Terumo™ centrifugal systems.

Compatible systems:

- Stöckert and Sorin centrifugal pump systems (with drive unit Model 60-01-04 and Model 60-01-00 and emergency drive unit Model 60-01-35 and Model 60-01-50)
- Sarns and Terumo centrifugal systems (with drive unit Model 164267 and manual drive unit Model 164268)

The adapter uses magnetic coupling to transfer drive motor torque to the disposable centrifugal blood pump.

The adapter is supplied nonsterile.

The use of the adapter does not require modification of the Stöckert and Sorin centrifugal pump systems or the Sarns and Terumo centrifugal systems.

The adapter is designed for multiple uses with different patients. The adapter may be reused until the Use-by date on the device label.

For additional information on operating the Affinity CP centrifugal blood pumps, refer to the instructions for use provided with the pumps.

2 Intended purpose

The Affinity CP adapter has no direct effect on the body. It enables the use of a Medtronic Affinity CP blood pump with compatible systems (see Chapter 1).

2.1 Indications for use

Refer to the Affinity CP centrifugal blood pump IFU for indications for use.

2.2 Contraindications

Refer to the Affinity CP centrifugal blood pump IFU for contraindications.

2.3 Intended users

This device should be used only by clinicians thoroughly trained in extracorporeal circulation procedures.

2.4 Patient target group

Refer to the Affinity CP centrifugal blood pump IFU for patient target groups.

2.5 Expected clinical benefits

Refer to the Affinity CP centrifugal blood pump IFU for expected clinical benefits.

3 Warnings

Read all warnings, precautions, and instructions for use carefully before use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.

- Always have a spare Affinity CP adapter and backup equipment available with appropriate protocol for replacement. If the adapter needs to be replaced during circulatory support, the potential risk to the patient should be evaluated before changing the adapter.
- Ensure that the adapter is properly attached to the drive motor. Verify that the adapter is locked onto the drive motor. Improper attachment may adversely affect pump performance.
- Ensure that the pump is properly attached to the adapter. The pump should be secure under the retaining ledge and the retaining latch of the adapter. Verify that the pump is locked securely and in close contact with the adapter receptacle. Improper attachment may adversely affect pump performance.
- The Affinity CP centrifugal blood pump and the Affinity CP adapter are not compatible with operation in a pulsatile mode.
- Do not use the adapter at speeds greater than 3600 rpm. The adapter is not qualified for speeds greater than 3600 rpm.
- Do not sterilize the adapter. Sterilization may damage the adapter.
- The Medtronic Affinity CP adapter is to be used only with the Medtronic Affinity CP centrifugal blood pumps. Use the Medtronic Affinity CP centrifugal blood pump as described in the instructions for use accompanying each device.
- Use the Stöckert and Sorin centrifugal pump systems and the Sarns and Terumo centrifugal systems as described in the instructions for use accompanying each device.
- When using the Affinity CP centrifugal blood pump and the Affinity CP adapter with a compatible system that displays inlet pressure, use a measured value for inlet pressure rather than a calculated value. A

calculated value is based on characteristic curves that are not representative of the Affinity CP centrifugal blood pump and may result in inaccurate readings.

- The device is not serviceable. Do not disassemble the device because the performance of the device may be adversely affected. If the device has been tampered with, do not use it.

4 Precautions

- This device should be used only by clinicians thoroughly trained in extracorporeal circulation procedures.
- Inspect the adapter for physical damage before use. Do not use the adapter if it appears to be damaged. Return the damaged product to a Medtronic representative.
- The adapter contains magnets. Keep the adapter away from metals, magnetic tapes, and magnets.
- Do not strike or drop the device. A shock may damage the device.
- Do not insert or remove the adapter or the pump until the drive motor has completely stopped.
- To reduce the risk of the adapter and pump detaching from the drive motor or any tubing connections, do not bump or strike the adapter or the circuit components during use. If detachment occurs, refer to Chapter 8.
- Do not store at extremes of temperature and humidity. Avoid storing the adapter in direct sunlight.
- For complete installation information, refer to the instructions for use provided with the Stöckert and Sorin centrifugal pump systems or the Sarns and Terumo centrifugal systems.
- Do not use alcohol-based solutions, other organic solvents, or concentrated acids or bases on the adapter. The adapter may be damaged by these solutions.
- Do not use the adapter after the Use-by date on the device label.
- Refer to the package labeling or instructions for use for storage temperature requirements.
- Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

5 Adverse effects

The following known adverse effects are associated with the use of the device: abrasion, blood loss, bruise, burn, death, exsanguination, hemolysis, hypovolemia, irritation, ischemia, laceration, and organ dysfunction.

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

6 Instructions for use

This section describes how to use the adapter with the Stöckert and Sorin centrifugal pump systems or the Sarns and Terumo centrifugal systems. Instructions are also provided for inserting the pump on the adapter.

6.1 Inserting the Affinity CP adapter on the Stöckert™ and Sorin™ centrifugal pump systems

Caution: For complete installation information, refer to the instructions for use provided with the Stöckert and Sorin centrifugal pump systems. Failure to refer to the instructions for use may result in patient injury or death.

Note: Before inserting the adapter on the drive motor, inspect the adapter for visible biocontamination, such as blood or foreign material. If necessary, clean the adapter or drive motor receptacle according to the instructions in Chapter 9.

Note: These instructions also apply when the adapter is used with a compatible handcrank from the Stöckert and Sorin centrifugal pump systems.

Insert the Affinity CP adapter on the drive motor of the Stöckert and Sorin centrifugal pump systems by performing the following steps:

1. Seat the bottom of the adapter into the receptacle of the drive motor. Push down the retaining key on the drive motor during this process.
2. Completely insert the adapter into the drive motor receptacle (Figure 2).
3. To rotate the adapter to the desired location, press the retaining key, turn the adapter, and ensure that the retaining key locks into the notch on the underside of the adapter.

Warning: Ensure that the adapter is properly attached to the drive motor. Verify that the adapter is locked onto the drive motor. Improper attachment may result in patient injury or death.

Caution: When inserting the adapter on the drive motor, ensure that the underside of the adapter is in close contact with the receptacle of the drive

motor and that the retaining key is in the locked position. Improper attachment may adversely affect adapter performance.

6.2 Inserting the Affinity CP adapter on the Sarns™ and Terumo™ centrifugal systems

Caution: For complete installation information, refer to the instructions for use provided with the Sarns and Terumo centrifugal systems. Failure to refer to the instructions for use may result in patient injury or death.

Caution: During use, do not bump or strike the adapter or the circuit components to reduce the risk of the adapter and pump detaching from the drive motor or any tubing connections. If detachment occurs, refer to Chapter 8.

Note: Before inserting the adapter on the drive motor, inspect the adapter for visible biocontamination, such as blood or foreign material. If necessary, clean the adapter or drive motor receptacle according to the instructions in Chapter 9.

Note: These instructions also apply when the adapter is used with a compatible handcrank from the Sarns and Terumo centrifugal systems.

Insert the Affinity CP adapter on the drive motor of the Sarns and Terumo centrifugal systems by performing the following steps:

1. Open the receptacle latch of the drive motor and remove the protective cover from the drive motor receptacle.
2. Position the bottom of the adapter over the drive motor receptacle. Slip the adapter under the receptacle clip (Figure 3).
3. Open the receptacle latch and press the adapter into place. Release the latch to secure the adapter in place.

Warning: Ensure that the adapter is properly attached to the drive motor. Verify that the adapter is locked onto the drive motor. Improper attachment may adversely affect adapter performance.

Caution: When inserting the adapter on the drive motor, ensure that the underside of the adapter is in close contact with the receptacle of the drive motor and that the receptacle latch is in the locked position.

6.3 Inserting the Affinity CP centrifugal blood pump on the Affinity CP adapter

Insert the Affinity CP centrifugal blood pump on the Affinity CP adapter by performing the following steps:

1. Position the bottom of the pump over the adapter receptacle. Slip the pump under the retaining ledge.
2. Open the retaining latch and press the pump into place (Figure 4). Release the retaining latch and adjust the pump until the retaining latch locks.

Warning: Ensure that the pump is properly attached to the adapter. The pump should be secure under the retaining ledge and the retaining latch of the adapter. Verify that the pump is locked securely and in close contact with the adapter receptacle. Improper attachment may adversely affect pump performance.

3. To rotate the pump to the desired location, press the retaining latch, turn the pump, and release the retaining latch. Adjust the pump until the retaining latch locks into place.

Caution: Prime the pump before starting circulation according to the directions in the pump's instructions for use. Operate the pump, coupled to the drive motor with the adapter, at maximum RPM to check for vibration, unusual noises, lack of flow, and other anomalies. If an anomaly is noted, refer to Chapter 8.

4. Proceed with extracorporeal circulation. Refer to the instructions for use supplied with each device.

Warning: Always have a spare Affinity CP adapter and backup equipment available with the appropriate protocol for replacement. If the adapter needs to be replaced during circulatory support, the potential risk to the patient should be evaluated before changing the adapter.

Warning: The Affinity CP centrifugal blood pump and the Affinity CP adapter are not compatible with operation in a pulsatile mode. Operating the adapter in pulsatile mode may result in patient injury or death.

Warning: Do not use the adapter at speeds greater than 3600 rpm. The adapter is not qualified for speeds greater than 3600 rpm and may result in patient injury or death.

7 Emergency backup equipment

Caution: Verify the Use-by date on the device label of the adapter.

Ensure that the following equipment is available:

- Backup power
- An appropriate handcrank and drive unit from a compatible system
- An appropriate controller from a compatible system
- A spare Affinity CP adapter

8 Troubleshooting guide

This section covers 2 situations that may occur when using the adapter with a centrifugal pump. Follow the clinical protocol specified by the hospital and consider the following troubleshooting suggestions.

8.1 Lack of flow

If the pump does not rotate or there is no flow, consider the following troubleshooting suggestions:

- Verify that the pump is properly connected to the adapter. Ensure that the pump is secure under the retaining ledge and the retaining latch of the adapter.
- Verify that the adapter is properly connected to the drive motor.
- Reduce the pump speed to zero and restart the drive motor.
- Follow hospital protocol to reestablish flow.

8.2 Unusual noise

If the pump or the adapter vibrate or make a grinding noise, consider the following troubleshooting suggestions:

- Verify that the pump is properly connected to the adapter. Ensure that the pump is secure under the retaining ledge and the retaining latch of the adapter.
- Verify that the adapter is properly connected to the drive motor.
- Reduce the pump speed to zero and restart the drive motor.
- Follow hospital protocol to reestablish flow.
- Consider replacing the pump or the adapter if the vibration or grinding noise continues.

9 Storage and cleaning of the Affinity CP adapter

Warning: Do not sterilize the adapter. Sterilization may damage the adapter.

Caution: Do not store at extremes of temperature and humidity. Avoid storing the adapter in direct sunlight.

Caution: Do not use alcohol-based solutions, other organic solvents, or concentrated acids or bases on the adapter. The adapter may be damaged by these solutions.

Follow your hospital's protocol to clean or disinfect the adapter and consider the following actions:

- Use only hospital-approved disinfectants to clean the adapter. The following solutions are acceptable for cleaning or disinfecting the adapter: water, Betadine™, hydrogen peroxide, and mild detergent. After cleaning, wipe the adapter with a sponge or soft cloth and water.
- Immediately remove any contamination from the adapter by using a sponge or soft cloth and water.
- If contamination cannot be removed with a sponge or soft cloth and water, the adapter may be immersed in water for up to 10 minutes. Do not immerse the adapter in any other liquids.

10 Service

Warning: The adapter is not serviceable. Do not disassemble the adapter because the performance of the adapter may be adversely affected. If the adapter has been tampered with, do not use it.

11 Specifications for the Affinity CP adapter

Part number:	AP40AST
Compatible pump:	<ul style="list-style-type: none"> • Affinity CP centrifugal blood pump Model AP40 • Affinity CP centrifugal blood pump Model BBAP40 • Affinity CP centrifugal blood pump Model CBAP40
Compatible systems:	<ul style="list-style-type: none"> • Stöckert and Sorin centrifugal pump systems (with drive unit Model 60-01-04 and Model 60-01-00 and emergency drive unit Model 60-01-35 and Model 60-01-50) • Sarns and Terumo centrifugal systems (with drive unit Model 164267 and manual drive unit Model 164268)
Maximum adapter speed:	3600 rpm
IP rating:	IP47

11.1 Gross dimensions

Weight:	315 g (0.69 lb)
Height:	28 mm (1.1 in)
Diameter:	94 mm (3.7 in)

11.2 Storage and transport environment

Temperature:	–30°C to 70°C (–22°F to 158°F)
Humidity:	15% to 93% noncondensing
Atmospheric pressure:	70 kPa to 106.3 kPa

11.3 Operating environment

Temperature:	10°C to 40°C (50°F to 104°F)
Humidity:	20% to 75% noncondensing
Atmospheric pressure:	70 kPa to 106.3 kPa

12 Limited warranty

The following LIMITED WARRANTY applies to all customers with the exception of Canada:

- A. This LIMITED WARRANTY is available solely to the original purchaser of the Affinity™ CP Adapter Model AP40AST (hereinafter, referred to as “Equipment”).
- (1) In the event that the Equipment fails to function within normal tolerances due to a defect in materials or workmanship within a period of one (1) year, commencing with the delivery of the Equipment to the purchaser, Medtronic will at its sole option: (a) replace the Equipment; (b) issue a credit to the purchaser equal to the Purchase Price (as defined in Subsection A(2) below) against the purchase price of the replacement Equipment; or (c) provide a functionally comparable replacement Equipment at no charge.
 - (2) As used herein, Purchase Price shall mean the lesser of (i) the net invoiced price of the original Equipment, (ii) the price of currently functional comparable equipment, or (iii) the price of replacement equipment.
- B. To qualify for the replacement or credit set forth in Section A above, the following conditions must be met:
- (1) The Equipment must be returned to Medtronic within sixty (60) days after discovery of the defect;
 - (2) The Equipment must not have been: (a) used other than in strict accordance with the Instructions for Use, (b) altered in any way which, in the judgment of Medtronic, affects the stability or reliability

of the Equipment, or (c) subjected to misuse, abuse or accident;
and

(3) The “use-by date” for the Equipment must not have passed.

- C. THE LIMITED WARRANTY OFFERED HEREBY IS EXTENDED SOLELY TO THE ORIGINAL PURCHASER OF THE EQUIPMENT, ONLY FOR THE LIMITED PERIOD DESCRIBED, AND OTHERWISE LIMITED TO THE EXPRESS TERMS SET FORTH HEREIN. EXCEPT AS EXPRESSLY SET FORTH IN THIS Chapter 12, MEDTRONIC EXPRESSLY DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, WHETHER AS TO MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR ANY OTHER MATTER, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOM, OR OTHERWISE. THE REMEDIES SET FORTH IN THIS LIMITED WARRANTY ARE THE EXCLUSIVE REMEDIES AVAILABLE FOR BREACH OF THIS LIMITED WARRANTY AND MEDTRONIC SHALL HAVE NO OTHER LIABILITY FOR DIRECT, INDIRECT, INCIDENTAL, CONSEQUENTIAL, OR PUNITIVE DAMAGES OF ANY DESCRIPTION BASED ON ANY DEFECT, FAILURE OR MALFUNCTION OF THE EQUIPMENT, WHETHER ARISING OUT OF WARRANTY, OTHER CONTRACT, TORT, OR OTHERWISE.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.
- E. No person has any authority to bind Medtronic to any representation, condition or warranty except this LIMITED WARRANTY.

13 Limited warranty¹

The following LIMITED WARRANTY applies to customers in Canada:

- A. This LIMITED WARRANTY provides the following assurances to the health care institution (“the purchaser”) that purchases any model of the Medtronic Affinity™ CP Adapter AP40AST, hereafter referred to as the “Equipment,” that should the Equipment fail to function within normal

¹ Warranty only applies to models licensed and purchased in Canada.

tolerances due to a defect in materials or workmanship within a period of one (1) year, commencing with the delivery of the Equipment to the purchaser, Medtronic will at its option:

- (1) Replace the Equipment;
- (2) Issue a credit equal to the original Equipment purchase price (but not to exceed the value of the replacement Equipment), against the purchase of replacement Equipment; or
- (3) Provide functionally comparable replacement Equipment at no charge.

B. To qualify for this replacement or credit, the following conditions must be met:

- (1) The Equipment must be returned to Medtronic within sixty (60) days after discovery of the defect.
- (2) The Equipment must not have been altered by anyone in any way.
- (3) The Equipment must not have been subjected to misuse, abuse or accident.

This LIMITED WARRANTY is provided by Medtronic of Canada Ltd. This LIMITED WARRANTY only applies to models licensed and purchased in Canada.

C. This LIMITED WARRANTY is limited to its express terms. In particular:

- (1) Except as expressly provided by this LIMITED WARRANTY, MEDTRONIC IS NOT RESPONSIBLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE, OR MALFUNCTION OF THE EQUIPMENT, IN A MANNER CONSISTENT WITH ITS INTENDED OPERATION AND PERFORMANCE DUE TO THE QUALITY OF MATERIALS AND WORKMANSHIP, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.
- (2) This LIMITED WARRANTY is made only to the original purchaser of the Equipment. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY OR CONDITION, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED TERM OR CONDITION OF SALE OR WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE WHETHER ARISING FROM ANY LAWS AND REGULATIONS, COMMON LAW, CUSTOM, OR OTHERWISE. NO SUCH EXPRESS OR IMPLIED WARRANTY OR CONDITION TO THE PURCHASER SHALL EXTEND BEYOND THE MODEL SPECIFIED ABOVE.

THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.

- (3) To the extent permitted by applicable law, the terms, conditions and warranties and the remedies provided in this LIMITED WARRANTY set out the sole and exclusive terms, conditions, warranties and remedies made. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid. This LIMITED WARRANTY gives the purchaser specific legal rights. The purchaser may also have other rights, which may vary by province.
- (4) No person has any authority to bind Medtronic to any representation, condition, or warranty except this LIMITED WARRANTY.
- (5) This LIMITED WARRANTY is not applicable to any device used with this Equipment.

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© 2019 Medtronic
The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.
Medtronic B.V.