# Medtronic

## Affinity™ CP CBAP40

Centrifugal Blood Pump with Cortiva™ BioActive Surface  
Центрофужна кръвна помпа с биоактивна повърхност Cortiva™  
Odstředivá krevní pumpa s bioaktivním povrchem Cortiva™  
Centrifugalblodpumpe med Cortiva™ BioActive-overflade  
Zentrifugalblutpumpe mit Cortiva™ BioActive Surface  
Φυγόκεντρη αντλία αίματος με βιοενεργή επίστρωση Cortiva™  
Bomba centrífuga de sangre con superficie biocompatible Cortiva™  
Cortiva™ bioaktiivse pinnakattega tsentrifugaalverepump  
Bioaktiivisella Cortiva™-pinnoitteella päällystetty keskipakoveripumppu  
Pompe sanguine centrifuge avec surface bioactive Cortiva™  
Centrifugalna pumpa za krv s bioaktivnim premazom Cortiva™  
Centrifugális pumpa Cortiva™ bioaktív felülettel  
Biopompa centrifuga con superficie bioattiva Cortiva™  
Centrifuginė kraujo pompa su „Cortiva™“ bioaktyvuoju paviršiumi  
Centrbēdzes asins sūkņis ar Cortiva™ bioloģiski aktīvo pārklājumu  
Centrifugaalpompe met bioactieve Cortiva™-oppervlaktelaag  
Sentrifugalpumpe med Cortiva™ bioaktiv overflate  
Wirnikowa pompa krwi z powłoką bioaktywną Cortiva™  
Bomba centrífuga de sangue com superfície bioativa Cortiva™  
Pompă de sânge centrifugală cu suprafață bioactivă Cortiva™  
Центробежный насос для крови с биоактивным покрытием Cortiva™  
Odstředivá krvná pumpa s bioaktivným povrchem Cortiva™  
Centrifugalna črpalka za kri z biološko aktivno površino Cortiva™  
Centrifugalna pumpa za krv sa Cortiva™ BioActive površinom  
Centrifugalblodpump med Cortiva™ bioaktiv yta  
Cortiva™ Biyo Aktif Yüzeyli Santrifüjli Kan Pompası

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης •  
Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Istruzioni per  
l'uso • Naudojimo instrukcijos • Lietošanas pamācība • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções  
de utilização • Instrucțiuni de utilizare • Инструкция по эксплуатации • Pokyny na používanie • Navodila za uporabo • Uputstva  
za upotrebu • Bruksanvisning • Kullanım Talimatları

! USA

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## 1 Description

The Affinity CP centrifugal blood pump with Cortiva bioactive surface (CBAP40) is intended for use in extracorporeal circulation procedures. It is designed to move blood by centrifugal force generated by a combination of a smooth rotating cone and low-profile impeller fins. Energy is transferred from the pump in the form of pressure and velocity as the blood is driven toward the outlet port of the pump. To limit friction and heat generation, the Affinity CP centrifugal blood pump with Cortiva bioactive surface utilizes a pivot-bearing design on a dual ceramic pivot.

The Affinity CP centrifugal blood pump with Cortiva bioactive surface has been sterilized using ethylene oxide.

Products coated with Cortiva bioactive surface include a “CB” prefix in the model number. Additional information regarding Cortiva bioactive surface is available in Chapter 7.


The Affinity CP centrifugal blood pump with Cortiva bioactive surface couples to a remote magnetic drive unit called the external drive motor model 560A, which interfaces with a Medtronic pump speed controller. It also couples with the emergency handcrank model HC150A in the event of controller or power failure. Refer to the appropriate Medtronic pump speed controller operator’s manual for additional information about the external drive motor and the emergency handcrank.

The Affinity CP centrifugal blood pump with Cortiva bioactive surface also couples with the Affinity CP adapter to enable the pump to be operated with the Stöckert™ and Sorin™ centrifugal pump systems or the Sarns™ and Terumo™ centrifugal systems.

The adapter may not be available in your country. Contact your Medtronic representative to confirm if the adapter is available.

For additional information on using the Affinity CP adapter or operating the Medtronic pump speed controller or accessories, refer to the instructions for use or operator’s manual.

**Figure 1.** Affinity CP centrifugal blood pump with Cortiva bioactive surface

- 1 Inlet port
- 2 Outlet port

## 2 Indications for use

The Affinity CP centrifugal blood pump with Cortiva bioactive surface is used to pump blood through the extracorporeal bypass circuit for extracorporeal circulatory support for periods appropriate to cardiopulmonary bypass (up to 6 hours).

It is also indicated for use in extracorporeal support systems (for periods up to 6 hours) not requiring complete cardiopulmonary bypass (for example, valvuloplasty, circulatory support during mitral valve reoperation, surgery of the vena cava or aorta, liver transplants).

The Affinity CP centrifugal blood pump with Cortiva bioactive surface is driven by the external drive motor or the emergency handcrank.

The Affinity CP centrifugal blood pump with Cortiva bioactive surface is intended for use with Medtronic controllers or may be used with the Stöckert and Sorin centrifugal pump systems or the Sarns and Terumo centrifugal systems by attaching the Affinity CP adapter.

## 3 Contraindications

The Affinity CP centrifugal blood pump with Cortiva bioactive surface is contraindicated for use as a cardiotomy suction device. Use the device only as indicated.

## 4 Warnings

Read all warnings, precautions, and instructions for use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.

## 4.1 Affinity CP centrifugal blood pump

- The centrifugal blood pump has not been qualified through in vitro, in vivo, or clinical studies for long-term use (longer than 6 hours) as a bridge-to-transplant or for pending recovery of the natural heart.
- The Affinity CP centrifugal blood pump and the Affinity CP adapter are not compatible with operation in a pulsatile mode.
- Do not use the centrifugal blood pump with an external drive motor or adapter that it is not compatible with because this may result in reduced performance, decoupling, or damage.
- Ensure that the centrifugal blood pump and circuit have been debubbled and primed properly before initiating bypass to minimize the risk of air reaching the patient. The use of an arterial filter is recommended. You must visually confirm that all air has been removed from the bearing assembly.
- Never clamp or restrict flow proximally to the pump; doing so may result in pump failure. If flow is occluded at the inlet, the pump may run at high rpms and possibly generate a high vacuum greater than 500 mmHg.
- Do not operate the centrifugal blood pump unprimed; damage to the internal components will occur.
- Massive air entry into the centrifugal blood pump deprimes the pump and stops blood flow. Stop the pump and remove air before resuming circulation.
- Never operate the pump without fluid.
- Do not operate the centrifugal blood pump at rotational speeds greater than the specified maximum (4000 rpm). Doing so may result in excessive wear to the pivot bearings.
- Do not expose the centrifugal blood pump to liquid chemical agents, because they may affect the integrity of this device. Anesthetic liquids, such as isoflurane, are known to degrade polycarbonate plastics. Avoid contact of these solutions with the centrifugal blood pump.
- To prevent retrograde blood flow when the centrifugal blood pump outlet tubing is open, establish and maintain a minimum pump speed that overcomes line and patient resistance. If positive flow cannot be maintained, the pump outlet tubing must be clamped to isolate flow from the patient. Failure to do so could allow retrograde flow and exsanguinate the patient.
- To prevent retrograde flow, arterial lines must always be clamped distally when the centrifugal blood pump is at rest.
- Always have a spare Affinity CP centrifugal blood pump and backup equipment available with appropriate protocol for replacement. If the centrifugal blood pump needs to be replaced during circulatory support, the potential risk to the patient should be evaluated before changing the pump.
- Frequent patient and device monitoring is recommended; do not leave the centrifugal blood pump unattended while in use. Monitor flow rate carefully for signs of occlusion throughout the circuit.
- Do not operate the centrifugal blood pump for more than 30 seconds in the absence of flow. The temperature within the pump could rise, and increased cellular damage may result.
- The centrifugal blood pump must be handled using aseptic technique.

## 4.2 Affinity CP adapter

- Always have a spare Affinity CP adapter and backup equipment available with appropriate protocol for replacement. If the adapter needs to be replaced during circulatory support, the potential risk to the patient should be evaluated before changing the adapter.
- Ensure that the adapter is properly attached to the drive motor. Verify that the adapter is locked onto the drive motor. Improper attachment may adversely affect pump performance.
- Ensure that the pump is properly attached to the adapter. The pump should be secure under the retaining ledge and the retaining latch of the adapter. Verify that the pump is locked securely and in close contact with the adapter receptacle. Improper attachment may adversely affect pump performance.
- The Affinity CP centrifugal blood pump and the Affinity CP adapter are not compatible with operation in a pulsatile mode.
- Do not use the adapter at speeds greater than 3600 rpm. The adapter is not qualified for speeds greater than 3600 rpm.
- Do not sterilize the adapter. Sterilization may damage the adapter.
- The adapter is not serviceable. Do not disassemble the adapter because the performance of the adapter may be adversely affected. If the adapter has been tampered with, do not use it.

## 5 Precautions

**Caution:** This device should only be used by clinicians thoroughly trained in extracorporeal circulation procedures.

### 5.1 Affinity CP centrifugal blood pump

- Do not forcefully strike or drop the centrifugal blood pump. Shock may cause damage to the device, which may cause the device to malfunction.

- The fluid pathway of the centrifugal blood pump is sterile and nonpyrogenic. Inspect the device and package carefully before use. Do not use if the unit package or the product has been damaged or soiled, or if caps are not in place.
- The centrifugal blood pump was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Follow hospital protocol for maintaining adequate heparinization for a period appropriate to cardiopulmonary bypass (up to 6 hours). Formation of thrombus in the circuit may increase the risk of damage to the perfusion system equipment.
- Attach tubing in a manner that prevents kinks or restrictions that may alter flow.
- Do not use excessive force to install tubing on the pump because damage to the pump might occur.
- Take care to prevent damage to connectors when repositioning or attaching tubing.
- Do not place the centrifugal blood pump near items adversely affected by magnetic fields.
- Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures. Adequate heparinization should be maintained per institution CPB protocol.
- Refer to the package labeling or instructions for use for storage temperature requirements.
- It is the user's responsibility to dispose of the devices in accordance with local regulations and hospital procedures.

## 5.2 Affinity CP adapter

- Use the Affinity CP adapter as described in the instructions for use provided with the adapter.
- Use the Medtronic Affinity CP adapter only with the Medtronic Affinity CP centrifugal blood pumps.
- Use the Stöckert and Sorin centrifugal pump systems and the Sarns and Terumo centrifugal systems as described in the instructions for use accompanying each device.
- Inspect the adapter for physical damage before use. Do not use the adapter if it appears to be damaged. Return the damaged product to a Medtronic representative.
- The adapter contains magnets. Keep the adapter away from metals, magnetic tapes, and magnets.
- Do not strike or drop the adapter. A shock may damage the adapter.
- Do not insert or remove the adapter or the pump until the drive motor has completely stopped.
- When inserting the adapter on the drive motor, ensure that the underside of the adapter is in close contact with the receptacle of the drive motor and that the receptacle latch is in the locked position.
- During use, do not bump or strike the adapter or the circuit components to reduce the risk of the adapter and pump detaching from the drive motor or any tubing connections. If detachment occurs, refer to Chapter 9.
- Do not store at extremes of temperature and humidity. Avoid storing the adapter in direct sunlight.
- For complete installation information, refer to the instructions for use provided with the Stöckert and Sorin centrifugal pump systems or the Sarns and Terumo centrifugal systems.
- Prime the pump before starting circulation. Operate the pump, coupled to the drive motor with the adapter, at maximum RPM to check for vibration, unusual noises, lack of flow, and other anomalies. If an anomaly is noted, refer to Chapter 9.
- When using the Affinity CP centrifugal blood pump and the Affinity CP adapter with a compatible system that displays inlet pressure, use a measured value for inlet pressure rather than a calculated value. A calculated value is based on characteristic curves that are not representative of the Affinity CP centrifugal blood pump and may result in inaccurate readings.
- Do not use alcohol-based solutions, other organic solvents, or concentrated acids or bases on the adapter. The adapter may be damaged by these solutions.
- Do not use the adapter after the Use-by date on the device label.
- Refer to the package labeling or instructions for use for storage temperature requirements.
- It is the user's responsibility to dispose of the devices in accordance with local regulations and hospital procedures.

## 6 Adverse effects

The following known adverse effects are associated with the use of the Affinity CP centrifugal blood pump with Cortiva bioactive surface (CBAP40):

- blood loss
- coagulopathy
- death
- embolism
- excessive blood component activation or thrombogenicity

- exsanguination
- hemolysis
- hemolytic anemia
- infection
- ischemia
- neurological dysfunction
- organ dysfunction

## 7 Information for products with Cortiva bioactive surface

The primary blood-contacting surfaces of the product are coated with Cortiva bioactive surface. This coated surface enhances blood compatibility and provides a blood-contacting surface that is thromboresistant. Cortiva bioactive surface contains nonleaching heparin derived from porcine intestinal mucosa<sup>1</sup>.

**Caution:** Product coated with Cortiva bioactive surface is intended for single use only. Resterilization may adversely affect the Cortiva bioactive surface.

**Caution:** Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

**Warning:** Do not store product coated with Cortiva bioactive surface above 40°C (104°F).

## 8 Instructions for use

If using the Affinity CP centrifugal blood pump without the Affinity CP adapter, refer to Section 8.1, Section 8.2, Section 8.7, and Section 8.8 of these instructions for use.


If using the Affinity CP centrifugal blood pump with the Affinity CP adapter, refer to Section 8.3, Section 8.4, Section 8.5, Section 8.6, Section 8.7, and Section 8.8 of these instructions for use.

### 8.1 Installing the Affinity CP centrifugal blood pump on the Medtronic external drive motor model 560A

Connect the external drive motor model 560A to the appropriate Medtronic controller. For complete installation information, refer to the appropriate Medtronic controller operator's manual.

1. Slide the Affinity CP centrifugal blood pump into the external drive motor model 560A (Figure 2) until the thumb latch pin engages.

**Figure 2.** Installing and removing the Affinity CP centrifugal blood pump

1 Thumb latch with locking pin on external drive motor model 560A

**Note:** Whenever possible, install the pump at the lowest point in the circuit, directly below the reservoir.

2. Prime the centrifugal blood pump according to the procedures described in Section 8.7 of these instructions for use.

**Caution:** Operating the centrifugal blood pump without priming the pump may cause damage to the centrifugal blood pump.

3. Turn the rpm knob to the clicked-off position. As the knob is turned to zero, the user will feel a slight pressure and hear a click. The black button on the knob serves as the pointer to the rpm settings.

**Warning:** Use of the Affinity CP centrifugal blood pump beyond the labeled recommendations may result in failure of the centrifugal blood pump, reduced pumping capacity, leaks, excessive blood trauma, or degradation or corrosion of blood contact materials that may pass through the blood to the patient.

**Caution:** Installing the centrifugal blood pump with the drive motor revolving may harm the pump.

4. Proceed with extracorporeal circulation. Refer to the instructions for use supplied with each device.

<sup>1</sup> Nonleaching is defined as heparin at a level  $\leq 0.1$  IU/mL, as measured through clinically relevant extraction conditions.

## 8.2 Removing the Affinity CP centrifugal blood pump from the Medtronic external drive motor model 560A


Remove and discard the pump after each procedure. Discard the pump according to local and hospital protocols.

1. Push the drive motor thumb latch to release the centrifugal blood pump (Figure 2).
2. Slide the centrifugal blood pump out of the drive motor receptacle (Figure 2).

## 8.3 Inserting the Affinity CP adapter

The Medtronic Affinity CP adapter (Figure 3) is intended to be used to enable the Medtronic Affinity CP centrifugal blood pump (model CBAP40) to be operated with the Stöckert and Sorin centrifugal pump systems or the Sarns and Terumo centrifugal systems.

**Figure 3.** Affinity CP adapter

1 Retaining ledge

3 Retaining latch

2 Adapter receptacle

The adapter uses magnetic coupling to transfer drive motor torque to the disposable centrifugal blood pump.

The adapter is supplied nonsterile.

The use of the adapter does not require modification of the Stöckert and Sorin centrifugal pump systems or the Sarns and Terumo centrifugal systems.

The adapter is designed for multiple uses with different patients. The adapter may be reused until the Use-by Date on the device label.

The adapter may not be available in your country. Contact your Medtronic representative to confirm if the adapter is available.

## 8.4 Inserting the Affinity CP adapter on the Stöckert™ and Sorin™ centrifugal pump systems

**Caution:** For complete installation information, refer to the instructions for use provided with the Stöckert and Sorin centrifugal pump systems.


**Note:** Before inserting the adapter on the drive motor, inspect the adapter for visible biocontamination, such as blood or foreign material. If necessary, clean the adapter or drive motor receptacle according to the instructions in Chapter 10.

**Note:** These instructions also apply when the adapter is used with a compatible handcrank from the Stöckert and Sorin centrifugal pump systems.

Insert the Affinity CP adapter on the drive motor of the Stöckert and Sorin centrifugal pump systems by performing the following steps:

1. Seat the bottom of the adapter into the receptacle of the drive motor. Push down the retaining key on the drive motor during this process.
2. Completely insert the adapter into the drive motor receptacle (Figure 4).

**Figure 4.** Inserting the Affinity CP adapter into the Stöckert and Sorin centrifugal pump system's drive motor

- 1 Affinity CP adapter
- 2 Retaining key

- 3 Drive motor
- 4 Drive motor receptacle

3. To rotate the adapter to the desired location, press the retaining key, turn the adapter, and ensure that the retaining key locks into the notch on the underside of the adapter.

**Warning:** Ensure that the adapter is properly attached to the drive motor. Verify that the adapter is locked onto the drive motor. Improper attachment may adversely affect adapter performance.

**Caution:** When inserting the adapter on the drive motor, ensure that the underside of the adapter is in close contact with the receptacle of the drive motor and that the retaining key is in the locked position.

### 8.5 Inserting the Affinity CP adapter on the Sarns™ and Terumo™ centrifugal systems

**Caution:** For complete installation information, refer to the instructions for use provided with the Sarns and Terumo centrifugal systems.

**Caution:** During use, do not bump or strike the adapter or the circuit components to reduce the risk of the adapter and pump detaching from the drive motor or any tubing connections. If detachment occurs, refer to Chapter 9.


**Note:** Before inserting the adapter on the drive motor, inspect the adapter for visible biocontamination, such as blood or foreign material. If necessary, clean the adapter or drive motor receptacle according to the instructions in Chapter 10.

**Note:** These instructions also apply when the adapter is used with a compatible handcrank from the Sarns and Terumo centrifugal systems.

Insert the Affinity CP adapter on the drive motor of the Sarns and Terumo centrifugal systems by performing the following steps:

1. Open the receptacle latch of the drive motor and remove the protective cover from the drive motor receptacle.
2. Position the bottom of the adapter over the drive motor receptacle. Slip the adapter under the receptacle clip (Figure 5).

**Figure 5.** Inserting the Affinity CP adapter into the Sarns and Terumo centrifugal system's drive motor

- |                          |                    |
|--------------------------|--------------------|
| 1 Affinity CP adapter    | 4 Receptacle latch |
| 2 Drive motor receptacle | 5 Drive motor      |
| 3 Receptacle clip        |                    |

- Open the receptacle latch and press the adapter into place. Release the latch to secure the adapter in place.

**Warning:** Ensure that the adapter is properly attached to the drive motor. Verify that the adapter is locked onto the drive motor. Improper attachment may adversely affect adapter performance.


**Caution:** When inserting the adapter on the drive motor, ensure that the underside of the adapter is in close contact with the receptacle of the drive motor and that the receptacle latch is in the locked position.

## 8.6 Inserting the Affinity CP centrifugal blood pump on the Affinity CP adapter

Insert the Affinity CP centrifugal blood pump on the Affinity CP adapter by performing the following steps:

- Position the bottom of the pump over the adapter receptacle. Slip the pump under the retaining ledge.
- Open the retaining latch and press the pump into place (Figure 6). Release the retaining latch and adjust the pump until the retaining latch locks.

**Figure 6.** Inserting the Affinity CP centrifugal blood pump onto the Affinity CP adapter

- 1 Affinity CP centrifugal blood pump  
2 Affinity CP adapter

- 3 Drive motor

**Warning:** Ensure that the pump is properly attached to the adapter. The pump should be secure under the retaining ledge and the retaining latch of the adapter. Verify that the pump is locked securely and in close contact with the adapter receptacle. Improper attachment may adversely affect pump performance.

- To rotate the pump to the desired location, press the retaining latch, turn the pump, and release the retaining latch. Adjust the pump until the retaining latch locks into place.

**Caution:** Prime the pump before starting circulation according to the directions in the pump's instructions for use. Operate the pump, coupled to the drive motor with the adapter, at maximum RPM to check for vibration, unusual noises, lack of flow, and other anomalies. If an anomaly is noted, refer to Chapter 9.

- Proceed with extracorporeal circulation. Refer to the instructions for use supplied with each device.

**Warning:** Always have a spare Affinity CP adapter and backup equipment available with the appropriate protocol for replacement. If the adapter needs to be replaced during circulatory support, the potential risk to the patient should be evaluated before changing the adapter.

**Warning:** The Affinity CP centrifugal blood pump and the Affinity CP adapter are not compatible with operation in a pulsatile mode.

**Warning:** Do not use the adapter at speeds greater than 3600 rpm. The adapter is not qualified for speeds greater than 3600 rpm.

## 8.7 Priming the Affinity CP centrifugal blood pump circuit

- Using accepted aseptic technique to maintain a sterile fluid pathway, attach the appropriate tubing to the inlet and outlet ports of the centrifugal blood pump.
- If possible, flush the circuit and the pump with carbon dioxide (CO<sub>2</sub>).
- Using gravity, fill the centrifugal blood pump with priming solution to a point beyond the pump outlet port and clamp. Remove air from the outlet tubing.

**Warning:** Ensure that the centrifugal blood pump and circuit have been debubbled and primed properly before initiating bypass to minimize the risk of air reaching the patient. The use of an arterial filter is recommended.

**Warning:** Massive air entry into the centrifugal blood pump causes the pump to deprime and blood flow to stop. Stop the pump and remove air before resuming circulation.

**Caution:** Do not strike the centrifugal blood pump with instruments. Shock may cause damage to the device, which may cause device malfunction.

4. With the outlet tubing clamped, turn on the Medtronic controller power.
5. Connect the flow probe to the circuit tubing according to the appropriate compatible system operator's manual. If using a Medtronic controller, zero the flow probe according to the appropriate Medtronic controller operator's manual.
6. With the outlet tubing still clamped, turn the rpm knob of the Medtronic controller to the desired level to adequately prime and debubble the circuit. Observe the pump for leakage or other anomalies.

**Warning:** Do not operate the centrifugal blood pump for more than 30 seconds in the absence of flow. The temperature within the pump could rise, and increased cellular damage may result.

7. Set the rpm knob to zero and recheck the integrity of the centrifugal blood pump.

**Warning:** If leaks are found, remove the centrifugal blood pump and replace with a new, sterile pump. If other anomalies are found, refer to Chapter 9. Repeat Step 1 to Step 7 to prime.

8. If no anomalies are noted, continue priming the circuit.
9. Verify all connections and the integrity and flow of the circuit before use.

**Warning:** Do not operate the centrifugal blood pump with its inlet tubing clamped, as negative pressure would be generated in the pump and air bubbles may be formed in the blood.

## 8.8 Affinity CP centrifugal blood pump checklist

Follow the protocol specified by your clinical team. Information from this blood pump bypass checklist can be used to revise, upgrade, or expand existing checklists appropriately.

### 1. Equipment assembly

- Mount external drive motor model 560A or drive unit from a compatible system.
- Check that all electrical connections are secure.
- Test control module power and display.
- Check date and integrity of sterile blood pump (and disposable probe if applicable) package(s).
- If using the adapter, check the Use-by Date on the device label.
- Check that flow transducer/sensor/disposable probe are sized properly.
- Assemble perfusion circuit in a sterile manner.
- Allow sufficient tubing length for standby drive unit or handcrank.
- Connect flow transducer/sensor/disposable probe to circuit in correct location and flow direction.

### 2. Prime pump and circuit

- Check the circuit to ensure that no kinks or occlusions are present.
- CO<sub>2</sub> flush pump and circuit, if indicated; turn off CO<sub>2</sub>.
- Gravity prime and debubble pump and perfusion circuit.
- Check pump for leaks, irregular motion, and noise.
- Check circuit for visible air.
- Check that all tubing connections are secure.
- Clamp pump outlet tubing completely.
- Clamp venous return tubing completely.

### 3. Operating parameters

Perform the following, per the appropriate compatible system operator's manual:

- Calibrate transducers/sensors according to the manufacturer's instructions.
- Set and verify low/high flow alarms and other available safety systems.

### 4. Emergency backup equipment

Ensure that the following equipment is available:

- Backup power
- Appropriate handcrank and drive unit from a compatible system
- Spare Affinity CP centrifugal blood pump and disposable probe
- If using the adapter, have a spare Affinity CP adapter available
- Appropriate controller from a compatible system

## 5. Perfusion

- Achieve minimum pump flow before unclamping lines.
- Monitor control module for messages and alarms.
- Monitor perfusion circuit for visible air and kinks in the lines.
- Maintain minimum pump flow before clamping lines.

## 6. Cleanup

- Turn off power.
- Properly discard disposable components.
- Clean motor and flow transducer/sensor.
- If using the adapter, clean it according to the instructions in the Storage and Cleaning section.

## 7. Check equipment

- Inspect and verify that equipment is operational.
- Maintain indicated preventive maintenance schedule.
- Ensure that batteries are charged to full capacity.

# 9 Troubleshooting guide

This section covers 3 situations that may occur when using a centrifugal blood pump. Follow the clinical protocol specified by the hospital and consider the following troubleshooting suggestions.

## 9.1 Pump failure

An unusual or high-pitched sound indicates that pump failure is occurring. Replace the Affinity CP centrifugal blood pump immediately.

## 9.2 Lack of flow

If the pump does not rotate or there is no flow, consider the following:

- Verify that the pump and the drive motor are properly connected.

If using the adapter, consider the following:

- Verify that the pump is properly connected to the adapter. Ensure that the pump is secure under the retaining ledge and the retaining latch of the adapter.
- Verify that the adapter is properly connected to the drive motor.
- Reduce the pump speed to zero and restart the drive motor.
- Follow hospital protocol to reestablish flow.

## 9.3 Unusual noise

If the pump vibrates or makes a grinding noise, consider the following:

- Consider replacing the pump if the vibration or grinding noise continues.

If using the adapter, consider the following:

- Verify that the pump is properly connected to the adapter. Ensure that the pump is secure under the retaining ledge and the retaining latch of the adapter.
- Verify that the adapter is properly connected to the drive motor.
- Reduce the pump speed to zero and restart the drive motor.
- Follow hospital protocol to reestablish flow.
- Consider replacing the pump or the adapter if the vibration or grinding noise continues.

# 10 Storage and cleaning of the Affinity CP adapter

**Warning:** Do not sterilize the adapter. Sterilization may damage the adapter.

**Caution:** Do not store at extremes of temperature and humidity. Avoid storing the adapter in direct sunlight.

**Caution:** Do not use alcohol-based solutions, other organic solvents, or concentrated acids or bases on the adapter. The adapter may be damaged by these solutions.

Follow your hospital's protocol to clean or disinfect the adapter and consider the following:

- Use only hospital-approved disinfectants to clean the adapter. The following solutions are acceptable for cleaning or disinfecting the adapter: water, Betadine™, hydrogen peroxide, and mild detergent. After cleaning, wipe the adapter with a sponge or soft cloth and water.
- Immediately remove any contamination from the adapter by using a sponge or soft cloth and water.
- If contamination cannot be removed with a sponge or soft cloth and water, the adapter may be immersed in water for up to 10 minutes. Do not immerse the adapter in any other liquids.

## 11 Service

**Warning:** The adapter is not serviceable. Do not disassemble the adapter because the performance of the adapter may be adversely affected. If the adapter has been tampered with, do not use it.

## 12 Specifications for the Affinity CP centrifugal blood pump with Cortiva bioactive surface

Part number:	CBAP40
Priming volume:	40 mL (approximate)
Inlet/Outlet ID:	9.5 mm (3/8 in)
Maximum operating pressure:	760 mmHg (101 kPa)
Maximum flow rate:	10 L/min
Maximum outlet pressure:	700 mmHg (93.3 kPa)
Maximum pump speed:	<ul style="list-style-type: none"> <li>• 4000 rpm (with Medtronic external motor drive 560A)</li> <li>• 3600 rpm (with Affinity CP adapter on Stöckert and Sorin centrifugal pump systems or Sarns and Terumo centrifugal systems)</li> </ul>
Compatible systems:	Medtronic speed controller systems (with external motor drive model 560A, emergency hand-crank model HC150A)

Store the product at room temperature.

### 12.1 Transport range

Temperature:	-30°C to 40°C (-22°F to 105°F)
Humidity:	20% to 85% noncondensing

## 13 Specifications for the Affinity CP adapter

Part number:	AP40AST
Compatible pump:	Affinity CP centrifugal blood pump with Cortiva bioactive surface (CBAP40)
Compatible systems:	<ul style="list-style-type: none"> <li>• Stöckert and Sorin centrifugal pump systems (with drive unit models 60-01-04 and 60-01-00 and emergency drive unit models 60-01-35 and 60-01-50)</li> <li>• Sarns and Terumo centrifugal systems (drive unit model 164267 and manual drive unit model 164268)</li> </ul>
Maximum adapter speed:	3600 rpm
IP rating:	IP47

### 13.1 Gross dimensions

Weight:	315 g (0.69 lb)
Height:	28 mm (1.1 in)
Diameter:	94 mm (3.7 in)

### 13.2 Storage and transport environment

Temperature:	-30°C to 70°C (-22°F to 158°F)
Humidity:	15% to 93% noncondensing
Atmospheric pressure:	70 kPa to 106.3 kPa

### 13.3 Operating environment


Temperature:	10°C to 40°C (50°F to 104°F)
Humidity:	20% to 75% noncondensing
Atmospheric pressure:	70 kPa to 106.3 kPa

### 14 Flow rate of Affinity CP centrifugal blood pump

**Note:** The actual obtainable flow is dependent on afterload of the pump which results from the extracorporeal circuit components and the patient's arterial resistance.

**Note:** Pressure drop is the difference between the outlet of the pump and the inlet of the pump.

**Figure 7.** Hydraulic performance of the Affinity CP centrifugal blood pump

1 Pressure drop (mmHg)

2 Flow rate (L/min)

3 <sup>a</sup> = RPM

### 15 Limited warranty

The following LIMITED WARRANTY applies to United States customers only:

- A. This LIMITED WARRANTY is available solely to the original purchaser of the Affinity CP centrifugal blood pump with Cortiva bioactive surface, Model CBAP40 (hereinafter, referred to as "Product").
- (1) In the event that the Product fails to function within normal tolerances due to a defect in materials or workmanship within a period of one (1) year, commencing with the delivery of the Product to the purchaser, Medtronic will at its sole option: (a) repair or replace any defective part or parts of the Product; (b) issue a credit to the purchaser equal to the Purchase Price (as defined in Subsection A(2) below) against the purchase price of the replacement Product; or (c) provide a functionally comparable replacement Product at no charge.
  - (2) As used herein, Purchase Price shall mean the lesser of (i) the net invoiced price of the original Product, (ii) the price of currently functional comparable equipment, or (iii) the price of replacement equipment.
- B. To qualify for the repair, replacement or credit set forth in Section A above, the following conditions must be met:
- (1) The Product must be returned to Medtronic within sixty (60) days after discovery of the defect;
  - (2) The Product must not have been used other than in strict accordance with the Instructions for Use, must not have been used on more than a single patient, and must not have been altered or repaired in any way which, in the judgment of Medtronic, affects the stability and reliability of the Product. The Product must not have been subjected to misuse, abuse or accident; and
  - (3) The Use-by date for the Product must not have passed.

- C. This LIMITED WARRANTY is limited to its express terms. In particular:
- (1) Except as expressly provided by this LIMITED WARRANTY, MEDTRONIC IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.
  - (2) This LIMITED WARRANTY is made only to the purchaser of the Equipment. AS TO ALL OTHERS, MEDTRONIC MAKES NO OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOM OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PURCHASER SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.
- E. No person has any authority to bind Medtronic to any representation, condition or warranty except this LIMITED WARRANTY.

## 16 Limited warranty<sup>2</sup>

The following Limited Warranty applies to customers outside the United States:

- A. This LIMITED WARRANTY is available only to the original purchaser of the Medtronic Affinity CP centrifugal blood pump with Cortiva bioactive surface, Model CBAP40 (hereinafter referred to as "Product"). In the event that the Product fails to function to the specifications set forth in the Instructions for Use, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient. THE WARNINGS CONTAINED IN THE PRODUCT LABELING ARE CONSIDERED AN INTEGRAL PART OF THIS LIMITED WARRANTY. CONTACT YOUR LOCAL MEDTRONIC REPRESENTATIVE TO OBTAIN INFORMATION ON HOW TO PROCESS A CLAIM UNDER THIS LIMITED WARRANTY.
- B. To qualify for the LIMITED WARRANTY, these conditions must be met:
- (1) The Product must be used prior to its Use-by date.
  - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
  - (3) The Product must not have been used for any other patient, must not have been misused, altered, abused or subject to accident, and the Product must not have been altered or repaired in any way which, in the judgment of Medtronic, affects its stability and reliability.
- C. This LIMITED WARRANTY is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
  - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

---

<sup>2</sup> This LIMITED WARRANTY is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only outside the United States.

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© 2017 Medtronic
The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.
Medtronic B.V.
Medtronic Australasia Pty Ltd
Medtronic of Canada Ltd