NAUTILUS™

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




1. GENERAL STATEMENT

Federal (U.S.A.) law restricts this device to sale, distribution and use by or on the order of a physician. Any medical techniques described in this document are meant for informational purposes only and do not replace the expertise of the physician and/or medical protocols.

2. INDICATION FOR USE

The Nautilus™ ECMO Oxygenator with integrated heat exchanger is intended to provide assisted extracorporeal circulation and physiologic gas exchange of the patient's blood for up to 14 days in patients with acute respiratory failure or acute cardiopulmonary failure, where other available treatment options have failed, and continued clinical deterioration is expected or the risk of death is imminent. The integrated heat exchanger is intended to heat or cool the blood as needed during use.

3. CONTRAINDICATION

This device is not designed, sold, or intended for use except as indicated.

4. DEVICE DESCRIPTION

The Nautilus ECMO diffusion polymethylpentene (PMP) membrane oxygenator with Balance™ Biosurface is used in extracorporeal life support procedures to oxygenate blood, remove carbon dioxide and regulate blood temperature. Blood enters the device and passes through both the heat exchange membrane, where temperature is adjusted, and the gas transfer membrane, where oxygen is added and carbon dioxide is removed.

The device is bonded on its primary blood contacting surfaces with a nonleaching biocompatible surface to reduce platelet activation and adhesion and preserve platelet function.

The device is single-use, nontoxic, nonpyrogenic, and not made from natural latex rubber materials.

5. ACCESSORIES SOLD SEPARATELY

To ensure a proper connection between the device and holder, the Nautilus ECMO Oxygenator is to be used with the Nautilus ECMO Oxygenator Holder only. This holder is sold separately.

6. GENERAL WARNINGS AND PRECAUTIONS

Warning: Only use this device in accordance with these instructions for use.

Warning: Read all warnings, precautions, and instructions for use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.

Warning: Only clinicians thoroughly trained in extracorporeal life support procedures should use this device.

Warning: No modification of this equipment is allowed. This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device which could result in patient injury, illness, or death.

Warning: The polymethylpentene (PMP) oxygenator fiber is not suitable for use with volatile anesthetics (anesthetic gases).

Warning: Do not use disinfectants in the heater/cooler while the heat exchanger is in use.

Note: Device is compatible for use with 330 ppm hydrogen peroxide in the heat exchanger during extracorporeal circulation.

Warning: This device should not be operated outside the recommended parameter ranges.

Warning: The following patient parameters must be monitored by an external system: core body temperature, arterial pressure, venous pressure, and blood oxygen saturation.

Warning: Gas transfer rates may change over time and adjustments of FiO₂ and gas flow rates may be needed to achieve desired gas transfer performance.

Warning: Ensure supply of oxygen is sufficient.

Warning: Avoid pulling or putting strain on the tubing.

Warning: Only use a medical gas supply with dry air and oxygen.

Warning: The use of PVC pump segments with a roller pump causes an electrostatic charge to build up on the blood side. Upon discharge, the integrity of the heat exchanger may be compromised.

Warning: Tubing should be attached in such a manner as to prevent kinks or restrictions that may alter the blood, water, or gas flow.

Warning: If air leaks are observed during priming and/or operation, this may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the device if these conditions are observed.

Warning: Positive pressure must be maintained on the blood side at all times during perfusion.

Warning: The blood path pressure should be greater than the gas path pressure at all times.

Warning: Do not obstruct gas outlet port and vents.

Warning: Use proper technique when using blood access ports.

Warning: Ensure that adequate anticoagulation procedures are employed during use.

Warning: Once flow to the patient is initiated, blood flow should circulate through the oxygenator at all times within the recommended blood flow range, except during emergency oxygenator change out.

Warning: Circuit access is the responsibility of the attending physician.

Caution: Do not allow anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come in contact with the device or accessories as they may jeopardize structural integrity.

Caution: Gas transfer rates decrease as a result of lower atmospheric pressures.

Caution: The device is not recommended for use with a pulsatile pump.

Caution: Exposing the oxygenator to extreme cold may compromise its mechanical integrity.

Caution: A replacement oxygenator should be readily available during ECLS.

Caution: Avoid mechanical impacts to the device while carrying or during operation.

7. ADDITIONAL INFORMATION FOR PRODUCTS WITH BALANCE BIOSURFACE

The blood-contacting surfaces of the device are coated with Balance Biosurface to reduce platelet activation and adhesion and preserve platelet function.

8. PACKAGING AND STORAGE

The fluid pathway is delivered sterile and nonpyrogenic. Sterilized by ethylene oxide. Verify that the packaging is not opened or damaged and the use-by date has not expired prior to use.

Warning: Do not use device if expired.

Warning: Do not use device if package is opened or damaged or if the device is damaged as the sterility of the device may have been compromised and/or performance may be affected.

9. INSTRUCTIONS FOR USE

9.1. Set Up of Holder

Remove holder from packaging, and inspect for damage prior to use. Attach the holder clamp to the ECLS cart at the desired height, and tighten the clamp knob. Adjust the remaining holder arms as desired, and tighten the holder arm knobs to secure in place.

9.2. Set Up of Device

Warning: Aseptic technique is to be used during all stages of setup and use of the device.

Warning: Setup and use of the device is the responsibility of the attending clinician.

Carefully remove the device from the sterile packaging. Slide the oxygenator holder connector onto the top of the oxygenator arm. A click indicates the oxygenator is correctly in place.

Note: To remove the oxygenator, press up on the holder release tab and pull the oxygenator away from the mounting arm.

Warning: Sterility is not ensured after device packaging is open.

Warning: Ensure the device is positioned lower than the patient or at the same level, but no higher.

Caution: The device must be securely mounted in an upright position on its intended holder.

9.3. Priming the Device

Caution: Tubes and luer lock connections must be fastened by hand only. Using tools or other aids may damage the ports.

Warning: Ensure that unused ports remain capped to prevent contamination and that the caps are tightened to avoid leakage.

Connect the water lines to the water ports on the oxygenator. Run water through the heat exchanger, and check for leaks from the water compartment to the blood compartment.

Warning: Do not use the oxygenator if water is present in the blood compartment.

Connect all blood and gas lines to their appropriate connections under aseptic conditions per institution protocol.

Note: Banding of all blood tubing connectors is recommended.

Fill the priming bag with the appropriate volume of priming solution per institution protocol. Remove the yellow luer lock cap on the de-airing membrane located on the blood inlet side of the device. Gravity prime the oxygenator or prime by forward pump flow per institution protocol. Gradually begin recirculation flow using a pre-bypass filter per institution protocol.

Warning: All gas emboli must be cleared from the extracorporeal circuit before initiating flow to the patient. Gas emboli are hazardous to the patient.

Add additional priming solution as required to prime the rest of the circuit. Upon completion of priming, re-cap the de-airing membrane. Ensure the entire system, including ports, is fully de-aired prior to beginning extracorporeal circulation.

Note: Device integrity has been verified by MC3 after 30 days of primed storage with phosphate buffered saline solution.

9.4. Extracorporeal Circulation

Caution: Check all lines are properly connected prior to use.

Caution: The de-airing membrane must be closed during operation.

Ensure there is no air in the circuit. Check for adequate levels of anticoagulation prior to and during extracorporeal circulation. Remove arterial and venous clamps and gradually increase blood flow. Start the gas flow with the appropriate gas settings per institution protocol. Adjust water temperature per institution protocol.

Note: Monitoring circuit pressure is recommended.

Caution: Do not exceed 42°C degrees in the heat exchanger.

9.5. Oxygenator Change out

The oxygenator should be evaluated for change out if the clinician determines during the extracorporeal procedure that adequate patient support is not being achieved. If a change out is deemed necessary, obtain a primed replacement device (see section 9.3) before discontinuing circulation through the original device. Have materials needed for the change available before starting the change, including clamps and sterile scissors. Always use aseptic technique during change out.

1. Double clamp the tubing on the inlet side and outlet side of the oxygenator that is to be replaced, and stop the pump. Clamp far enough away from the oxygenator blood inlet/ outlet to allow adequate room for reconnection to the replacement device.
2. Detach the oxygenator from the holder by pressing the release tab.
3. Attach the primed replacement oxygenator to the holder.
4. Cut the venous and arterial tubing between the clamps connected to the original oxygenator.
5. Using air-free technique, connect the arterial line to the blood outlet and the venous line to the blood inlet of the primed replacement oxygenator.
6. Open the clamp on the blood inlet side and turn the pump on at a low speed.
7. Connect the gas line to the gas inlet of the replacement device.
8. Confirm that the system is free of bubbles before removing the outlet clamp, and gradually increase the pump flow to the necessary value.
9. Connect the water lines to the replacement oxygenator.

10. SPECIFICATIONS

Oxygenation Fiber Type	Polymethylpentene (PMP)		
Fiber Surface Area – Oxygen Membrane	1.8 m ²		
Heat Exchanger	Polyethylene terephthalate (PET)		
Fiber Surface Area – Heat Exchange	0.3 m ²		
Prime Volume	226 mL		
Residual Blood Volume	189 mL		
Operating Blood Volume	226 mL		
Blood Flow Range	0.5 – 7 L/min		
Maximum Rated Blood Pressure	750 mmHg		
Gas Flow Rate (Gas:Blood)	0.5:1 – 3:1		
Maximum Gas Path Pressure	100 mmHg (1.9 psi)		
Maximum Water Path Pressure	1125 mmHg (21.8 psi)		
Environmental	Operating	Storage	Transportation
Temperature	15 – 30 °C	10 – 30 °C	-30 – 50 °C
Relative Humidity (non-condensing)	30 – 75 %	5 – 95 %	5 – 95 %

The following information and data are available upon request: sterilization method, materials of the blood pathway, blood pathway pressure drop, gas pathway pressure drop, blood cell damage, particle release, and relevant tolerances for specification in the Instructions for Use.

11. DISPOSAL

Warning: Dispose of the device in accordance with local hospital, administrative, and/or other government policies.

The device is a single-use only disposable. Dispose of the oxygenator according to institution protocol and local regulations.

12. CLEANING

If the outside of the oxygenator becomes soiled, wipe down surfaces using a cloth moistened with suitable cleaning solution. Do NOT spray cleaner directly on the device. The oxygenator is compatible with the following cleaner types: 3% hydrogen peroxide, bleach, isopropyl alcohol, aldehydes, and quaternary ammonium compounds.

Cleaning of the holder is recommended after each use. All surfaces of the holder can be cleaned and disinfected for blood, saline, or contaminate spills using normal medical equipment cleaners and disinfectants.





13. POTENTIAL COMPLICATIONS

Potential complications, including those normally associated with ECLS procedures and anticoagulation, may necessitate physician intervention. Some potential complications related to the use of ECLS oxygenators include but are not limited to hypoxia or hypercarbia due to inadequate gas exchange, air embolism, thrombus formation leading to deep-vein thrombosis or pulmonary embolism, hemolysis, thrombocytopenia, mechanical failure of ports and connections, infection, hypothermia or hyperthermia due to inadequate heat transfer, hemorrhage related to bleeding, and inadequate blood flow related to clotting or hypovolemia resulting in inadequate gas exchange. There are risks and adverse events related to all ECLS procedures and anticoagulation including heart, vessel, or lung damage, hypoxia, anemia, infection, hemorrhage, liver or kidney failure, stroke, or death.

14. GENERAL INFORMATION

ECLS General Information. This product is one component of an ECLS circuit. A typical ECLS circuit includes tubing, connectors, a gas exchange device and a blood pump. Configuration of the circuit and selection of individual component types may vary according to each specific use and treatment center. Each component of the circuit is commercially available. The extracorporeal circuit is typically managed with a systemic anticoagulant (usually heparin). The anticoagulation is managed to titrate either the active whole blood clotting time or activated partial thromboplastin time to a level above the normal ranges (typically 1-2 times normal). The principle pathway of the ECLS circuit consists of blood being drained from the drainage lumen of the catheter, passing through a blood pump which disperses the blood through a gas exchange device and back into the patient through the cannula's reinfusion lumen. The systemic temperature is controlled by a heat exchanger integrated into the gas exchange device. The components of the circuit are selected to meet the blood flow and gas exchange needs of each patient based on their metabolic requirements. The circuit has been widely described in the literature and is summarized in each iteration of the textbook published by the Extracorporeal Life Support Organization.¹ Maintenance of the ECLS circuit is managed by personnel who have been trained on the intricacies of extracorporeal support including physicians, perfusionists, nurses, respiratory therapists and other medical specialists. These individuals are qualified to manage the circuit based on training and continuing education per specific institutional guidelines.

¹Toomasian JM, Lawson DS, and Harris WH. ECLS Equipment. Chapter 8 in ECMO: Extracorporeal Cardiopulmonary Support in Critical Care, Fourth edition, Annich G, Bartlett RH, Lynch WR, MacLaren G, Wilson JM (eds). Extracorporeal Life Support Organization, Ann Arbor, 2012. pp 107-32.

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¹ Toomasian JM, Lawson DS, and Harris WH. ECLS Equipment. Chapter 8 in ECMO: Extracorporeal Cardiopulmonary Support in Critical Care, Fourth edition, Annich G, Bartlett RH, Lynch WR, MacLaren G, Wilson JM (eds). Extracorporeal Life Support Organization, Ann Arbor, 2012. pp 107-32.

MC3, Inc. Disclaimer of Warranty

MC3, INC. MANUFACTURES ITS PRODUCTS (REFERRED TO AS THE "PRODUCTS") UNDER CONDITIONS THAT MEET APPLICABLE REGULATORY REQUIREMENTS, BUT HAS NO CONTROL OVER HOW THE PRODUCTS ARE ACTUALLY USED OR THE CONDITIONS UNDER WHICH THEY ARE USED. AS A RESULT, MC3, INC. DISCLAIMS ALL PRODUCT WARRANTIES INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MC3, INC. SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE.

Limitation on Liability

IN NO EVENT WILL MC3, INC.'S TOTAL LIABILITY TO PURCHASER OR ANY THIRD PARTY ARISING FROM THE SALE OF PRODUCTS EXCEED THE PURCHASE PRICE OF THE PRODUCTS SUBJECT TO THE CLAIM REGARDLESS OF THE THEORY OF LIABILITY OF THE CLAIM.

Limitation on Damages

MC3, INC. IS NOT RESPONSIBLE TO PURCHASER OR ANY THIRD PARTY FOR ANY DIRECT, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY CLAIMED BREACH OF WARRANTY, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT. Applicable law in your jurisdiction may limit the ability of a manufacturer to disclaim some warranties or exclude or limit damages. While the above limitations are an integral part of the business relationship between MC3, Inc. and the purchasers of its Products, nothing in the above is intended to contravene these provisions of applicable law. If any part or term of this Disclaimer of Warranty or the Limitations on Liability or Damages is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of all remaining provisions shall not be affected and shall be construed and enforced as if the Disclaimer of Warranty or the Limitations on Liability or Damages did not contain the particular part or term held to be invalid.

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MC3, Inc. -takuun vastuuvapauslauseke
MC3, Inc. Garantifriskrivning
MC3, Inc. Garanti Feragati