Medtronic

Autotransfusion Blood Collection Reservoir EL2 Series
R servoir de r cup ration du sang pour autotransfusion s rie EL2
Blutsammelreservoir f r Autotransfusion der Serie EL2
Reservorio de recogida de sangre de autotransfusi n Serie EL2
Autotransfusiereservoir voor het verzamelen van bloed, EL2-serie
Serbatoio di raccolta del sangue per autotrasfusione Serie EL2
Reservoir for oppsamling av blod til autotransfusjon, EL2-serien
Autotransfuusioon veris illi pussi, EL2-sarja
Blodinsamlingsreservoar f r autotransfusion i EL2-serien
 εξαμεν  συλλογ ς α ματος αυτομετ γγισης Σειρά EL2
EL2-serien af blodopsamlingsbeholdere til autotransfusion
Zbiornik na krew pobieran  do autotransfuzji, seria EL2
Reservat rio de colheita de sangue para autotransfus o da s rie EL2
Ototransf zyon Kan Toplama Rezervuarı EL2 Serisi
Z sobn k na krev pro autotransfuzi řady EL2
EL2-es sorozat  autotranszf zi s v rgy jtt tart ly
Autotransf zny zbern  rezervo r krvi s ria EL2
Seeria EL2 autotransfusiooni vere kogumise mahuti
Резервоар за събиране на кръв за автотрансфузия серия EL2
Spremnik za skupljanje krvi kod autotranfuzija serije EL2
Autotransf zijas asins sav kšanas rezervu rs, EL2 s rija
Automatinio perpylimo kraujo surinkimo rezervuaras, serija EL2
Rezervor pentru colectarea s ngelui pentru autotransfuzie, seria EL2
Резервуар для сбора крови во время аутогемотрансфузии серии EL2
Rezervoar za prikupljanje krvi za autotransfuziju serije EL2
Zbiralnik krvi za avtotransfuziju serije EL2
EL240, EL2120

Instructions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso
Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • K ytt ohjeet • Bruksanvisning
Οδηγίες χρ σης • Brugsanvisning • Instrukcja użytowania • Instru  es de utiliza  o
Kullanım Talimatları • N vod k pou it  • Haszn lati  tmutat  • Pokyny na pou ivanie
Kasutusjuhend • Инструкци за употреба • Upute za upotrebu • Lietošanas pam c ba
Naudojimo instrukcijos • Instruc iuni de utilizare • Инструкция по эксплуатации
Uputstva za upotrebu • Navodila za uporabo

USA **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

Medtronic® is a registered trademark of Medtronic, Inc.
Medtronic® est une marque déposée de Medtronic, Inc.
Medtronic® ist eine eingetragene Marke von Medtronic, Inc.
Medtronic® es una marca registrada de Medtronic, Inc.
Medtronic® is een geregistreerd handelsmerk van Medtronic, Inc.
Medtronic® è un marchio registrato di Medtronic, Inc.
Medtronic® er et registrert varemerke for Medtronic, Inc.
Medtronic® on Medtronic, Inc.-n rekisteröity tavamerkki.
Medtronic® är ett registrerat varumärke som tillhör Medtronic, Inc.
Το Medtronic® αποτελεί σήμα κατατεθέν της Medtronic, Inc.
Medtronic® er et registreret varemærke tilhørende Medtronic, Inc.
Medtronic® jest zastrzeżonym znakiem towarowym firmy Medtronic, Inc.
Medtronic® é uma marca comercial registada da Medtronic, Inc.
Medtronic®, Medtronic, Inc. şirketinin tescilli ticari markasıdır
Medtronic® je registrovaná ochranná známka spoločnosti Medtronic, Inc.
A Medtronic® a Medtronic, Inc. bejegyzett védjegye.
Medtronic® je registrovaná ochranná známka spoločnosti Medtronic, Inc.
Medtronic® on ettevõtte Medtronic, Inc. registreeritud kaubamärk
Medtronic® e регистрирана търговска марка на Medtronic, Inc.
Medtronic® je registrirani zaštitni znak tvrtke Medtronic, Inc.
Medtronic® ir Medtronic, Inc. reģistrēta preču zīme.
„Medtronic™” yra registruotas „Medtronic, Inc.” prekės ženklas.
Medtronic® este o marcă comercială înregistrată a Medtronic, Inc.
Medtronic® является зарегистрированным товарным знаком Medtronic, Inc.
Medtronic® je registrovani žig preduzeća Medtronic, Inc.
Medtronic® je registrirana blagovna znamka družbe Medtronic, Inc.

Explanation of symbols on package labeling / Explication des symboles des étiquettes sur l'emballage / Erläuterung der Symbole auf dem Verpackungsetikett / Explicación de los símbolos que aparecen en el etiquetado del envase / Verklaring van de symbolen op de verpakkingslabels / Spiegazione dei simboli sulle etichette della confezione / Forklaring av symboler på emballasjen / Pakkauksen symbolien selitykset / Förklaring av symboler på förpackningsetiketten / Επεξήγηση των συμβόλων στην ετικέτα της συσκευασίας / Forklaring på symboler på emballagens mærkater / Objaśnienie symboli znajdujących się na etykietach opakowania / Explicação dos símbolos nas etiquetas da embalagem / Ambalaj etiketi üzerindeki sembollerin açıklaması / Vysvětlení symbolů uvedených na štítcích balení / A csomagoláson látható szimbólumok jelentése / Vysvetlivky k symbolom na označení balenia / Pakendi siltidel olevate tingmärkide tähendus / Объяснение на символите от етикетите на опаковката / Objašnjenje simbola na ambalažnoj etiketi / Iesaiņojuma marķējuma simbolu skaidrojums / Pakuotės etiketės nurodytų simbolių reikšmės / Explicarea simbolurilor cuprinse în etichetele de pe ambalaj / Объяснение символов, приведенных на этикетке упаковки / Objašnjenje simbola na oznakama na pakovanju / Razlaga simbolov na embalaži

Refer to the device labeling to see which symbols apply to this product / Se référer aux étiquettes sur le produit pour savoir quels symboles s'appliquent à ce produit / Welche Symbole für dieses Produkt zutreffen, entnehmen Sie bitte dem Etikett auf dem Gerät / Consulte la etiqueta del dispositivo para comprobar qué símbolos se utilizan con este producto / Controleer het verpakkingslabel om te zien welke symbolen op dit product van toepassing zijn / Vedere le etichette sul dispositivo per sapere quali simboli si applicano a questo prodotto / Se på pakningen for å fastslå hvilke symboler som gjelder for dette produktet / Katso laitteen merkinnöistä, mitkä symbolit koskevat tätä tuotetta / Se etiketten på enheten för de symboler som gäller denna produkt / Ανατρέξτε στην ετικέτα της συσκευής για να δείτε ποια σύμβολα ισχύουν για το προϊόν αυτό / Se mærkaterne på enheden for de symboler, der gælder for dette produkt / Aby określić, które symbole dotyczą niniejszego produktu, należy zapoznać się z etykietą na opakowaniu produktu / Consultar as etiquetas do dispositivo para ver que símbolos se aplicam a este produto / Bu üründe hangi sembollerin geçerli olduğunu görmek için cihazın etiketlerine bakın / Symboly, které se vztahují k tomuto produktu, najdete na štítcích obalu / A termékre vonatkozó szimbólumok a csomagoláson találhatóak / Prežite si označenie prístroja, aby ste zistili, ktoré symboly sa vzťahujú na tento výrobok / Toote puhul kehtivad tingmärgid leiata seadme siltidelt / Погледнете етикета на устройството, за да видите какви символи се отнасят за този продукт / Na etiketama na uređaju pogledajte koji su simboli primjenjivi na ovaj proizvod / Skatiet ierīces marķējumu, lai redzētu simbolus, kuri attiecas uz šo izstrādājumu / Norėdami sužinoti, kokių simbolių pažymėtas gaminys, žiūrėkite pakuotės etiketę / Consultați eticheta dispozitivului pentru simbolurile aplicabile acestui produs / См. на этикетке устройства, какие символы применимы к данному продукту / Pogledajte spoljnu oznaku na pakovanju da biste videli koji se simboli primenjuju na ovaj proizvod / Opis simbolov, ki se nanašajo na izdelek, je na oznakah na izdelku

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Autotransfusion Blood Collection Reservoir EL2 Series

Models

EL240, EL2120

Contents

2-port blood collection reservoir

Indications for Use

The EL2 Series Autotransfusion Blood Collection Reservoir is indicated for use with patients undergoing autotransfusion procedures. The unit has a capacity of 4000 mL and contains either a 40 or 120 micron integral microfilter. The reservoir has been designed to collect, filter, and store blood aspirated during the perioperative period of surgery or trauma prior to subsequent cell washing.

Contraindications

- The use of reinfused blood from the product may be contraindicated (eg, in the presence of sepsis or malignancy). This product used for any other purposes than for the indicated intended use is the responsibility of the user.
- The safe length of time that blood may remain in this reservoir is a medical decision and may vary with each case where intraoperative salvage is used.

Warnings

- Read all Warnings, Precautions, and Instructions for Use carefully prior to use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**
- Only persons thoroughly trained in cardiopulmonary bypass or autotransfusion procedures should use this product. Operation of each product requires constant supervision by qualified personnel for patient safety.
- Each product has been carefully manufactured, tested, and packaged; however, the state of the art has not been developed to the point that Medtronic is able to ensure that the product will not leak, crack, or fail during use. Perfusion must be carefully and constantly monitored.
- The product must be used immediately after the removal of the protective packaging.
- The fluid pathway is sterile and nonpyrogenic. Inspect each package and product prior to use. Do not use if the package is opened or damaged or if the protective caps are not in place.
- Use aseptic technique in all procedures.
- Tubing should be attached in such a manner as to prevent kinks or restrictions that may alter blood or water flow.
- Do not allow alcohol, alcohol-based fluids, anesthetic fluids (such as isoflurane), or corrosive solvents (such as acetone) to come into contact with the product as they may jeopardize the structural integrity.
- Polycarbonate may react adversely upon direct contact with paregoric or halothane. Avoid contact between these agents and the reservoir.
- If air bubbles and/or leaks are observed during priming and/or operation, these conditions may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the product if these conditions are observed.
- All gas emboli must be cleared from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- The fluid level in the venous reservoir bag should be maintained above the oxygenator at all times.
- Monitor the blood level in the venous reservoir bag at all times during perfusion.
- Never agitate or invert the product once cardiopulmonary bypass is initiated.
- Should any evidence of damage to components be found during setup visual inspection, do not use the product.
- Never exceed 150 mm Hg negative pressure when aspirating blood, as excessive vacuum may create accelerated blood trauma.
- Failure to verify the security of the reservoir holder and reservoir position may result in inadvertent overturn of the system and/or improper operation.
- This product is not intended for use as a cardiotomy reservoir.
- **Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

Precautions

- Store in a dry place away from extremes of temperature. Refer to package labeling for storage temperature requirements.

- A strict anticoagulation protocol should be followed, and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.
- “Perioperative products intended for transfusion shall be transfused through a filter designed to retain particles that are potentially harmful to the patient.”¹
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Adverse Effects

None known.

Instructions for Use

Warning: The setup and use of this product is the responsibility of the attending clinician.

1. Remove the product(s) from the packaging carefully to ensure sterile fluid pathway.

Warning: Ensure aseptic technique is used during all stages of setup and use of this product.

Warning: Before removing the product from the packaging, inspect the packaging and product for damage. If the packaging or product is damaged, do not use as sterility of the product may have been compromised and/or performance may be affected.

2. Mount holder by attaching it firmly to the IV pole.
3. Remove the reservoir from its sterile wrapper and place it in the holder. Tighten luer cap(s).
4. Remove the blue cap from the inlet connector and connect the suction/anticoagulant line.
5. Remove the yellow cap from the vacuum port on the reservoir. Attach the vacuum line.
6. Connect an outlet drain line; either a 1.0 to 0.6 cm (3/8 in to 1/4 in) reducer or a 1.0 to 0.6 × 0.6 cm (3/8 in to 1/4 in × 1/4 in) “Y” adapter, to the 1.0 cm (3/8 in) outlet port on the bottom of the reservoir. Close the clamp(s) on the leg(s) of the drain line.
7. Turn on the vacuum from 80 to 120 mm Hg. Do not exceed 150 mm Hg.
8. Prime the blood collection reservoir with a minimum of 100 to 200 mL of anticoagulant solution to completely wet the filter material.
9. Adjust the flow of anticoagulant to a ratio of approximately 15 mL of anticoagulant solution/100 mL blood.
10. If the blood collection is not going to start immediately (within 5 to 10 minutes), the anticoagulant flow may be stopped. Be sure to restart the anticoagulant solution prior to the collection of blood.
11. Blood collection may now begin.
12. If the blood is to be used for autotransfusion, follow the setup and operation instructions in the applicable manufacturer’s operator’s manual.

¹AABB. Standards for Perioperative Autologous Blood Collection and Administration. 4th ed. Bethesda, MD: AABB; 2009:19.

The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE AUTOTRANSFUSION BLOOD COLLECTION RESERVOIR, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE AUTOTRANSFUSION BLOOD COLLECTION RESERVOIR, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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© 2004, 2008, 2010 Medtronic
Medtronic of Canada Ltd.