# Medtronic

## Intersept™ 1351/CB1351

Filtered Cardiotomy Reservoir

Kardiotomireservoir med filter

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## 1 Description

Model Number	Description
1351	Filtered cardiotomy reservoir
CB1351	Filtered cardiotomy reservoir with Cortiva™ BioActive Surface

The Intersept filtered cardiotomy reservoir is a single-use device with a sterile, nonpyrogenic fluid path. The maximum capacity of each reservoir is 2,600 mL. The maximum recommended flow rate is 2 L/min.

Each reservoir has eight 0.6 cm (0.25 in) ID access (luer) ports:

- 4 suction access ports
- 2 predefoamer access ports
- 2 postdefoamer access ports. One of the postdefoamer access ports is intended for use as a vent.

**Note:** The predefoamer access port labeled VENT/VACUUM is fitted with a yellow vented cap. All other access ports are fitted with nonventing universal caps. To convert a universal cap to a female luer fitting, twist the top portion of the cap completely off (clockwise). To provide an access port for 0.6 cm (0.25 in) ID tubing, turn the entire universal cap counterclockwise and remove.

### 1.1 Model 1351/CB1351

Each Intersept filtered cardiotomy reservoir, with or without Cortiva bioactive surface, contains an open cell polyurethane defoamer with a 20-micron microaggregate filter covered with a polyester sleeve.

Products coated with Cortiva bioactive surface include a “CB” prefix in the model number. Additional information regarding Cortiva bioactive surface is available in Chapter 7.

## 2 Indications for use

### 2.1 Model 1351

The Intersept filtered cardiotomy reservoir is indicated for use in the cardiopulmonary bypass circuit during surgery for:

- an air/fluid separation chamber
- a temporary storage reservoir for priming solution and blood
- the filtration of particulates from bank blood and the storage and filtration of blood recovered from the field by suction
- the addition of medications or other fluids

The Intersept filtered cardiotomy reservoir is also indicated for use after open heart surgery for the postoperative collection of autologous blood from the chest and the aseptic return of blood to the patient for volume replacement.

### 2.2 Model CB1351

The Intersept filtered cardiotomy reservoir with Cortiva bioactive surface is indicated for use in the cardiopulmonary bypass circuit during surgery for:

- an air/fluid separation chamber
- a temporary storage reservoir for priming solution and blood
- the filtration of particulates from bank blood and the storage and filtration of blood recovered from the field by suction
- the addition of medications or other fluids

## 3 Contraindications

Use the device only as indicated.

### 3.1 Model 1351

The Intersept filtered cardiotomy reservoir (Model 1351 only) is contraindicated for use in postoperative chest drainage/autotransfusion procedures when:

- there is an air leak in the lung or gross perforations to the chest wall exist
- pericardial, mediastinal, pulmonary or systemic infection or malignancy is present
- gross contamination (eg, bone cement, lymphatic failure, perforated intestine) is present or suspected
- aspirating blood from a site containing topical hemostatic agents, bactericidal wound irritants or antibiotics not intended for parenteral administration
- the chest is open and vacuum is applied
- protamine has been administered prior to the reservoir being removed from the bypass circuit
- the patient is returned to surgery for any reason
- vented chest tubes not incorporating vent flow regulation such as stopcock are used

Complications that have been associated with autotransfusion include blood trauma, blood coagulation, coagulopathies, and particulate or air embolism.

## 4 Warnings

Read all warnings, precautions, and instructions for use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.

- Only persons thoroughly trained in cardiopulmonary bypass procedures should use this device. Perfusion must be carefully and constantly monitored.
- Each device is intended for single use only. Do not reuse or resterilize. Sterilized using ethylene oxide.
- To vent the reservoir, remove the yellow vented cap on the 0.6 cm (1/4 in) access port labeled VENT/VACUUM prior to use.
- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use if the package is opened or damaged or if the device is damaged or the protective caps are not in place.
- Use aseptic technique in all procedures.
- Attach the tubing in such a manner as to prevent kinks or restrictions that may alter blood flow.
- Do not allow alcohol, alcohol-based fluids, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come in contact with the product because they may jeopardize the structural integrity.
- If air bubbles and/or leaks are observed during priming and/or operation, these conditions may result in air embolism to the patient and/or fluid loss. Continually monitor the extracorporeal circuit. Do not use the product if these conditions are observed.
- Clear all gas emboli from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- A standby reservoir should always be available during perfusion.
- Pressurizing the cardiotomy reservoir could result in leakage or damage. The reservoir must be vented at all times when using a roller pump in the suction lines.
- If a vacuum source is being used, negative pressure should not exceed –200 mmHg.
- Blood that has been in the cardiotomy reservoir for 4 hours or more must not be reinfused.

#### 4.1 Model 1351 Only

- The safe reinfusion of collected fluids and the assurance that fluids collected are safe for reinfusion is the responsibility of the attending physician.
- Chest drainage should be autotransfused on an hourly basis unless fewer than 50 mL of blood/fluid are collected hourly.
- If the filter occludes, nonfiltered blood may pass through the external shunt and enter the reservoir on the filtered side (Figure 3, number 1). If chest drainage volume is excessive, nonfiltered blood may pass through the external shunt to the vent connector and back into the reservoir on the filtered side of the reservoir (Figure 3, number 2). Any drainage that has passed to the filtered side of the cardiotomy reservoir through the external shunt must be filtered prior to autotransfusion.
- A minimum reservoir level of 25 mL of blood should always be maintained to prevent passing an air embolus during reinfusion.
- It is not recommended that autotransfusion continue for more than 18 hours following surgery.
- Follow all accepted medical and nursing care routines for patients undergoing chest drainage.

## 5 Precautions

- Refer to the package labeling for storage temperature requirements.
- The maximum recommended flow rate is 2 L/min.
- Band all tubing connections.
- Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Dispose of the devices in accordance with local hospital, administrative, and other government policies.
- Additional warnings and precautions applicable to specific procedures can be found at appropriate places in the Instructions for Use.

## 6 Adverse effects

The following adverse effects are associated with the use of the product: allergic reaction, sensitization; blood loss, exsanguination, hypovolemia; capillary leak syndrome; carcinogenicity; circulatory compromise; cytotoxicity; death; fever; hemolysis; hypercarbia/acidosis; hypotension; hypoxia; impaired patient temperature management; infection; mutagenicity; neurological dysfunction/stroke, such as air embolism and post perfusion syndrome (pump head); organ dysfunction; platelet dysfunction, thrombocytopenia; renal failure; reproductive harm in male neonatal/pediatric patients; respiratory/pulmonary issues; sepsis, toxemia, blood toxicity; thrombotic and inflammatory responses, including excessive blood component activation; thrombus dislodgement; and vasculature damage.

## 7 Information for products with Cortiva bioactive surface

The primary blood-contacting surfaces of the product are coated with Cortiva bioactive surface. This coated surface enhances blood compatibility and provides a blood-contacting surface that is thromboresistant. Cortiva bioactive surface contains nonleaching heparin derived from porcine intestinal mucosa<sup>1</sup>.

**Caution:** Product coated with Cortiva bioactive surface is intended for single use only. Resterilization may adversely affect the Cortiva bioactive surface.

**Caution:** Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

**Warning:** Do not store product coated with Cortiva bioactive surface above 40°C (104°F).

## 8 Instructions for use

### 8.1 Operating room - cardiopulmonary bypass

When using **Model 1351/CB1351** for cardiopulmonary bypass in the operating room, follow this procedure:

1. Remove the device from the packaging carefully to ensure a sterile fluid pathway.

<sup>1</sup> Nonleaching is defined as heparin at a level ≤0.1 IU/mL, as measured through clinically relevant extraction conditions.

**Warning:** Ensure aseptic technique is used during all stages of setup and use of this device.

**Warning:** Before removing the device from the packaging, inspect the packaging and product for damage. If the packaging or product is damaged, do not use, as the sterility of the device may have been compromised and performance may be affected.

2. Place the cardiotomy reservoir in an Intersept cardiotomy reservoir holder.

**Note:** The cardiotomy reservoir should be mounted high enough (30.5 cm to 91.5 cm [12 in to 36 in]) above the oxygenator or venous reservoir to ensure adequate gravity flow. Adjust the height as necessary during priming and perfusion.

3. To vent the reservoir, remove the yellow vented cap on the 0.6 cm (0.25 in) access port labeled VENT/VACUUM prior to use. A spare yellow nonvented luer connecting cap has been placed on the side of the unit for your convenience. If the cap is used as a luer connection site, ensure that the opening is complete and free of debris after twisting the top off of the cap.

**Warning:** Be sure the vent port is open and vented. Failure to properly vent the port could result in air being introduced into the bypass circuit, causing serious injury or death to the patient.

**Note:** Do not open unused access ports.

4. Remove the protective cap from the outlet port at the bottom of the cardiotomy reservoir. Attach an appropriate length of 1.0 cm (0.39 in) ID tubing from the outlet port on the cardiotomy reservoir to the cardiotomy return port on the oxygenator inlet or venous reservoir.
5. Remove the universal cap from one of the suction access ports. Connect a 0.6 cm (0.25 in) ID cardiotomy suction line to the suction access port. Repeat this procedure for additional suction lines.

## 8.2 Operating room - postbypass chest drainage conversion

This procedure applies to **Model 1351 only**.

The Intersept filtered cardiotomy reservoir may be converted after bypass to an autotransfusion system to recover blood from the chest. Autotransfusion may be initiated after the patient reaches the intensive care unit, final setup is completed, and the hourly drainage rate is established.

To convert the Intersept filtered cardiotomy reservoir to an autotransfusion system:

1. Remove the device from the packaging carefully to ensure a sterile fluid pathway.

**Warning:** Ensure aseptic technique is used during all stages of setup and use of this device.

**Warning:** Before removing the device from the packaging, inspect the packaging and product for damage. If the packaging or product is damaged, do not use as sterility of the device may have been compromised and/or performance may be affected.

2. Place the reservoir in a holder on an IV pole, subthoracic to the patient, until ready for transfer to the ICU/CCU.
3. Using aseptic technique, carefully remove tubing from suction access ports (Figure 2, number 3 or Figure 3, number 3). Seal 3 of these ports with sterile (0.6 cm [0.25 in] ID) occluding plugs.
4. Remove tubing (if applicable) or the universal cap from the nonfiltered prime port (Figure 2 or Figure 3, number 9).
5. Attach sterile 0.6 cm (0.25 in) ID tubing to form a shunt between the 1 remaining unplugged suction access port and the nonfiltered prime port (Figure 2 or Figure 3, number 9).
6. Convert the 1.0 cm (0.39 in) ID outlet port to 0.3 cm (0.12 in) ID infusion tubing (Figure 2 or Figure 3, number 6) to permit connection to standard intravenous tubing for reinfusion.
7. Attach sterile 0.6 cm (0.25 in) ID tubing and a sterile 0.6 x 1.0 cm (0.25 x 0.39 in) connector to a filtered prime port (Figure 2 or Figure 3, number 1 or 2). Plug the unused port.
8. Attach the chest drainage blood tubing from the chest to the 0.6 x 1.0 cm (0.25 x 0.39 in) connector on the filtered prime port (Figure 2 or Figure 3, number 1 or 2).
9. Using sterile tubing, connect the vent/vacuum port (Figure 2 or Figure 3, number 7) on the top of the unit to a vacuum source with approximately 11 mmHg to 15 mmHg of regulated vacuum.

**Important:** Use of a water seal device is recommended in all chest drainage cases.

## 8.3 Intensive care unit

This section applies to **Model 1351 only**.

1. Remove the device from the packaging carefully to ensure a sterile fluid pathway.

**Warning:** Ensure aseptic technique is used during all stages of setup and use of this device.

**Warning:** Before removing the device from the packaging, inspect the packaging and product for damage. If the packaging or product is damaged, do not use as sterility of the device may have been compromised and/or performance may be affected.

2. Attach the cardiotomy reservoir holder to an IV pole and place the cardiotomy reservoir in the holder. Place the infusion pump (Figure 2, number 5 or Figure 3, number 5) on an IV pole.
3. Adjust the level of the cardiotomy reservoir to a subthoracic level. Maintain the cardiotomy reservoir at a subthoracic level to facilitate chest tube drainage.
4. Connect the chest tube drainage line (Figure 2 or Figure 3, number 8) to the chest drainage port.
5. Record the amount of drainage in the cardiotomy reservoir. Continually monitor drainage to determine the hourly drainage rate.
6. Using sterile tubing, connect the vent/vacuum port (Figure 2 or Figure 3, number 7) on the top of the unit to a vacuum source with approximately 11 mmHg to 15 mmHg of regulated vacuum.
7. Assure that all unused inlet ports on the cardiotomy reservoir are covered with airtight caps.
8. As soon as the cardiotomy reservoir contains 50 mL of chest tube drainage, set up and prime the infusion pump tubing per the standard procedure.
9. Monitor the hourly rate of drainage.
10. Refer to the infusion pump manufacturer's instructions for infusion of cellular fluids.
11. Autotransfuse drainage collected in the cardiotomy reservoir on an hourly basis after assessment of the hourly drainage rate and the hemodynamic status of the patient.

## 9 Important notice – limited warranty (US customers only)

- A. This **LIMITED WARRANTY** provides the following assurances to the patient who receives an Intersept filtered cardiotomy reservoir, hereafter referred to as the "Product":
- (1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its Use-by date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
  - (2) As used herein, purchase price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its Use-by date.
  - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
  - (3) The Product must not have been altered or subjected to misuse, abuse, or accident.
  - (4) The Product must be used in accordance with the labeling and Instructions for Use provided with the Product.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) Except as expressly provided by this **LIMITED WARRANTY**, Medtronic is not responsible for any direct, incidental, or consequential damages based on any defect, failure, or malfunction of the product, whether the claim is based on warranty, contract, tort, or otherwise.
  - (2) THIS **LIMITED WARRANTY** IS MADE ONLY TO THE PATIENT IN WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER, OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS **LIMITED WARRANTY** SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
  - (3) The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid. This **LIMITED WARRANTY** gives the patient specific legal rights. The patient may also have other rights which vary from state to state.
  - (4) No person has any authority to bind Medtronic to any representation, condition, or warranty except this **LIMITED WARRANTY**.

**Note:** This **LIMITED WARRANTY** is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for exact terms of the **LIMITED WARRANTY**.

## 10 Important notice – limited warranty (for countries outside the US)

- A. This **LIMITED WARRANTY** provides assurances for the purchaser who receives an Intersept filtered cardiotomy reservoir, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient. The Warnings contained in the product labeling are considered an integral part of this **LIMITED WARRANTY**. Contact your local Medtronic representative to obtain information on how to process a claim under this **LIMITED WARRANTY**.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its Use-by date.
  - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
  - (3) The Product may not have been used for any other patient.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation, or material alteration of the replaced Product.
  - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

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The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.
Medtronic Australasia Pty Ltd
Medtronic of Canada Ltd