
Soft-Flow[®] Aortic Cannula

-  EN Instructions for Use
Aortic Cannula
-  AR إرشادات الاستخدام
كانيولا الشريان الأورطي
-  DA Brugsanvisning
Aortakanyle
-  DE Gebrauchsanweisung
Aorta-Kanüle
-  ES Instrucciones de uso
Cánula Aórtica
-  FR Mode d'emploi
Canule Aortique
-  IT Istruzioni per l'uso
Cannula Aortica
-  JA 使用説明書
大動脈カニューレ
-  NL Gebruiksaanwijzing
Aortacanule
-  NO Bruksanvisning
Aortakanyle
-  PT Instruções de Uso
Cânula Aórtica
-  SV Bruksanvisning
Aortakanyl
-  ZH 使用说明
主动脉插管



0050

Symbols Glossary

	Manufacturer (ISO 15223-1:2016 5.1.1)		Do not use if package is damaged. (ISO 15223-1:2016 5.2.8)
	Authorized representative in the European Community (ISO 15223-1:2016 5.1.2)		Contains phthalates (EN15986:2011)
	Use-by date (ISO 15223-1:2016 5.1.4)		Single Use (ISO 15223-1:2016 5.4.2)
	Batch code (ISO 15223-1:2016 5.1.5)		Do not re-sterilize (ISO 15223-1:2016 5.2.6)
	Catalogue number (ISO 15223-1:2016 5.1.6)		Quantity
	Sterilized using ethylene oxide (ISO 15223-1:2016 5.2.3)		Keep dry (ISO 15223-1:2016 5.3.4)
	Non-pyrogenic (ISO 15223-1:2016 5.6.3)		Consult Instructions for Use (ISO 15223-1:2016 5.4.3)
	Caution (ISO 15223-1:2016 5.4.4)		

REF	1	2	3	4a	4b	5	6	7
7071	6.0 mm (18 Fr)		✓	✓			✓	✓
7072	6.0 mm (18 Fr)		✓	✓		✓	✓	✓
7080	6.0 mm (18 Fr)	✓		✓		✓	✓	✓
5761	7.0 mm (21 Fr)		✓	✓		✓	✓	✓
5767	7.0 mm (21 Fr)		✓	✓			✓	✓
5798	7.0 mm (21 Fr)		✓		✓	✓	✓	✓
5842	7.0 mm (21 Fr)		✓		✓		✓	✓
6392	7.0 mm (21 Fr)	✓		✓			✓	✓
6394	7.0 mm (21 Fr)	✓		✓		✓	✓	✓
5762	8.0 mm (24 Fr)		✓	✓		✓	✓	✓
5768	8.0 mm (24 Fr)		✓	✓			✓	✓
5841	8.0 mm (24 Fr)		✓		✓	✓	✓	✓
5843	8.0 mm (24 Fr)		✓		✓		✓	✓
6384	8.0 mm (24 Fr)	✓		✓			✓	✓
6386	8.0 mm (24 Fr)	✓		✓		✓	✓	✓

- 1. Tip outer diameter
- 2. Suture flange
- 3. Suture bulb
- 4a. Angled tip
- 4b. Straight tip
- 5. Luer port
- 6. Clamp site
- 7. Vented cap

Indications

Soft-Flow® Aortic Cannula are indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery for up to 6 hours.

Contraindications

This device is not designed, sold, or intended for use except as indicated.

Description

These cannulae feature a diffuser tip. They are designed with a closed tip and inner cone which direct flow through the vents, decreasing velocity while maintaining flow rate. All cannulae have a connector for 3/8 in. (0.95cm) ID tubing and wire-reinforced body.

Sterile and Non-pyrogenic.
Sterilized by ethylene oxide.

Blood Contact Material: polyvinyl chloride, polyethylene, and ABS (acrylonitrile-butadiene-styrene).
Not made with natural rubber latex.

Use

1. Inspect package and cannula for damage. Do not use if damaged.
2. Check that the vented cap and luer cap are secured before use.
3. Insert the tip portion of the cannula. Do not insert the cannula beyond the suture bulb/flange.
4. Holes in the vented cap allow air to escape and blood to fill the cannula.
5. When blood escapes from the vented cap, clamp tube where indicated. Pull tab on vented cap to remove from connector.
6. Complete priming and removal of air using conventional techniques.

Removal

1. Clamp cannula.
2. Remove cannula and close entry wound using standard surgical technique.
3. When the surgical procedure is completed, dispose of the cannula according to hospital policy. Do not reuse.

Cautions and Warnings

Proper surgical procedures and techniques are the responsibility of the medical professional.

These devices are designed for use in cardiovascular surgery by trained surgeons.

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Insertion beyond the suture bulb/flange could cause aortic damage.

Sterile and non-pyrogenic unless package is opened or damaged.

Only clamp on the clear tube as shown in the diagram. Do not clamp on the reinforced tubing areas or connector.

Do not use if package is damaged.

Adequate anticoagulation is required.

Phthalates

This device contains phthalates. Based on animal studies, significant exposure to phthalates may interfere with the normal development of the male reproductive tract. For children, pregnant and nursing women, alternative devices may be appropriate.

Single Use

For single use only. Do not reuse.

Do not re-sterilize. Do not reprocess.
Reprocessing may compromise the sterility, biocompatibility, and functional integrity of the device.

Adverse Effects

Cardiopulmonary bypass surgery includes clinical risks of adverse effects such as hemolysis, gas embolism, thrombosis, and device problems (leaks, component breaks, kinked tubing that reduces blood flow, foreign matter/particulates, and vessel injury due to surface irregularities during device insertion/removal) which may cause injury or death. Care and caution should be taken to avoid these clinical risks during cannulation or other cardiac surgery procedures.

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MC3, INC. Disclaimer of Warranty

MC3, INC. MANUFACTURES ITS PRODUCTS (REFERRED TO AS THE "PRODUCTS") UNDER CONDITIONS THAT MEET APPLICABLE REGULATORY REQUIREMENTS, BUT HAS NO CONTROL OVER HOW THE PRODUCTS ARE ACTUALLY USED OR THE CONDITIONS UNDER WHICH THEY ARE USED. AS A RESULT, MC3, INC. DISCLAIMS ALL PRODUCT WARRANTIES INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MC3, INC. SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE.

Limitation on Liability

IN NO EVENT WILL MC3, INC.'S TOTAL LIABILITY TO PURCHASER OR ANY THIRD PARTY ARISING FROM THE SALE OF PRODUCTS EXCEED THE PURCHASE PRICE OF THE PRODUCTS SUBJECT TO THE CLAIM REGARDLESS OF THE THEORY OF LIABILITY OF THE CLAIM.

Limitation on Damages

MC3, INC. IS NOT RESPONSIBLE TO PURCHASER OR ANY THIRD PARTY FOR ANY DIRECT, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY CLAIMED BREACH OF WARRANTY, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT.

Applicable law in your jurisdiction may limit the ability of a manufacturer to disclaim some warranties or exclude or limit damages. While the above limitations are an integral part of the business relationship between MC3, INC. and the purchasers of its Products, nothing in the above is intended to contravene these provisions of applicable law. If any part or term of this Disclaimer of Warranty or the Limitations on Liability or Damages is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of all remaining provisions shall not be affected and shall be construed and enforced as if the Disclaimer of Warranty or the Limitations on Liability or Damages did not contain the particular part or term held to be invalid.

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