NAUTILUS™ SMART

ECMO Module with Balance™ Biosurface

Instructions for Use

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




1. GENERAL STATEMENT

Federal (U.S.A.) law restricts this device to sale, distribution and use by or on the order of a physician. Any medical techniques described in this document are meant for informational purposes only and do not replace the expertise of the physician and/or medical protocols.

2. INDICATION FOR USE

The Nautilus™ Smart ECMO Module with integrated heat exchanger is intended to provide assisted extracorporeal circulation and physiologic gas exchange of the patient's blood for up to 14 days in patients with acute respiratory failure or acute cardiopulmonary failure, where other available treatment options have failed, and continued clinical deterioration is expected or the risk of death is imminent. The integrated heat exchanger is intended to heat or cool the blood as needed during use. Integrated fluid path pressure, temperature, and oxygen saturation monitoring is achieved by built-in sensor modules and display.

3. CONTRAINDICATION

This device is not designed, sold, or intended for use except as indicated.

4. DEVICE DESCRIPTION

The Nautilus Smart ECMO Module diffusion polymethylpentene (PMP) membrane oxygenator with Balance™ Biosurface (REF 48135E) is used in extracorporeal life support procedures to oxygenate blood, remove carbon dioxide and regulate blood temperature. Blood enters the device and passes through both the heat exchange membrane, where temperature is adjusted, and the gas transfer membrane, where oxygen is added and carbon dioxide is removed.

The device contains integrated sensors with an electronic touch screen display. The following measured parameters are visible on the electronic display: inlet pressure, inlet oxygen saturation, outlet pressure, outlet oxygen saturation, and outlet blood temperature. The difference between the inlet and outlet pressure, delta pressure, is calculated and also displayed on the screen.

The touchscreen display allows users to set alarm limits for all measured parameters. The device will alarm visually and audibly when the limits are exceeded.

The device is bonded on its primary blood contacting surfaces with a nonleaching biocompatible surface to reduce platelet activation and adhesion and preserve platelet function.

The device is single-use, nontoxic, nonpyrogenic, and not made from natural latex rubber materials.

5. ACCESSORIES SOLD SEPARATELY

To ensure a proper connection between the device and holder, the Nautilus Smart ECMO Module is to be used with the following accessories only: Nautilus ECMO Oxygenator Holder, Nautilus Smart ECMO Module Power Supply and country specific Nautilus Smart ECMO Module Power Cord.

6. GENERAL WARNINGS AND PRECAUTIONS

Warning: Only use this device in accordance with these instructions for use.

Warning: Read all warnings, precautions, and instructions for use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.

Warning: Only clinicians thoroughly trained in extracorporeal life support procedures should use this device.

Warning: No modification of this equipment is allowed. This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device which could result in patient injury, illness, or death.

Warning: The polymethylpentene (PMP) oxygenator fiber is not suitable for use with volatile anesthetics (anesthetic gases).

Warning: Do not use disinfectants in the heater/cooler while the heat exchanger is in use.

Note: Device is compatible for use with 330 ppm hydrogen peroxide in the heat exchanger during extracorporeal circulation.

Warning: This device should not be operated outside the recommended parameter ranges.

Warning: The following patient parameters must be monitored by an external system: core body temperature, arterial pressure, venous pressure, and blood oxygen saturation.

Warning: Gas transfer rates may change over time and adjustments of FiO₂ and gas flow rates may be needed to achieve desired gas transfer performance.

Warning: Ensure supply of power and oxygen are sufficient.

Warning: Avoid pulling or putting strain on the tubing or power cable.

Warning: Only use a medical gas supply with dry air and oxygen.

Warning: The use of PVC pump segments with a roller pump causes an electrostatic charge to build up on the blood side. Upon discharge, the integrity of the heat exchanger may be compromised.

Warning: Tubing should be attached in such a manner as to prevent kinks or restrictions that may alter the blood, water, or gas flow.

Warning: If air leaks are observed during priming and/or operation, this may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the device if these conditions are observed.

Warning: Positive pressure must be maintained on the blood side at all times during perfusion.

Warning: The blood path pressure should be greater than the gas path pressure at all times.

Warning: Do not obstruct gas outlet port and vents.

Warning: Use proper technique when using blood access ports.

Warning: Ensure that adequate anticoagulation procedures are employed during use.

Warning: Once flow to the patient is initiated, blood flow should circulate through the oxygenator at all times within the recommended blood flow range, except during emergency oxygenator change out.

Warning: In noisy environments, there is a risk that audible alarms may not be heard.

Warning: Circuit access is the responsibility of the attending physician.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Nautilus Smart ECMO Module, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Caution: Do not allow anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come in contact with the device or accessories as they may jeopardize structural integrity.

Caution: The displayed oxygen saturation measurements will be affected by Methylene blue, indocyanine green, Evans blue, fetal hemoglobin, and unusual levels of bilirubin, lipids, beta-carotene, carboxyhemoglobin, methemoglobin, and sulhemoglobin.

Caution: Gas transfer rates decrease as a result of lower atmospheric pressures.

Caution: The device is not recommended for use with a pulsatile pump.

Caution: Exposing the oxygenator to extreme cold may compromise its mechanical integrity.

Caution: A replacement oxygenator should be readily available during ECLS.

Caution: Avoid mechanical impacts to the device while carrying or during operation.

7. ADDITIONAL INFORMATION FOR PRODUCTS WITH BALANCE BIOSURFACE

The blood-containing surfaces of the device are coated with Balance Biosurface to reduce platelet activation and adhesion and preserve platelet function.

8. PACKAGING AND STORAGE

The fluid pathway is delivered sterile and nonpyrogenic. Sterilized by ethylene oxide. Verify the packaging is not opened or damaged and the use-by date has not expired prior to use.

Warning: Do not use device if expired.

Warning: Do not use device if package is opened or damaged, or if the device is damaged as the sterility of the device may have been compromised and/or performance may be affected.

9. INSTRUCTIONS FOR USE

9.1. Set Up of Holder, Power Cord, and Power Supply

Remove all accessories from packaging, and inspect for damage prior to use. Attach the holder clamp to the ECLS cart at the desired height, and tighten the clamp knob. Adjust the remaining holder arms as desired, and tighten the holder arm knobs to secure in place. Plug the power cord into the power supply and appropriate wall outlet.

9.2. Set Up of Device

Warning: Aseptic technique is to be used during all stages of setup and use of the device.

Warning: Setup and use of the device is the responsibility of the attending clinician.

Carefully remove the device from the sterile packaging. Slide the oxygenator holder connector onto the top of the oxygenator arm. A click indicates the device is correctly in place.

Warning: Sterility is not ensured after device packaging is open.

Warning: Ensure the device is positioned lower than the patient or at the same level, but no higher.

Caution: The device must be securely mounted in an upright position on its intended holder.

Note: To remove the device, press up on the holder release tab and pull the device away from the mounting arm.


Note: The device must be positioned such that the operator can see the display and indicator light and access the touch screen interface at all times.

9.3. Powering on the Device

Caution: Only use a Nautilus Smart ECMO Module power cord.

Remove the battery pull tab by gently pulling in the direction indicated by the arrow, and discard the pull tab. Insert the power cord into the device. Press the power switch to turn on the touchscreen. Verify the screen is on, the indicator light illuminates, and the device sounds during start-up. The display will show the Ready Screen. On the Ready Screen, measured parameter alarms are not yet active. Remain on the Ready Screen until ready to begin extracorporeal circulation.

Ready Screen

Note: The ECMO module comes with pre-zeroed pressure sensors. No further action is needed to ensure pressure sensors are zeroed prior to use.

Note: A press and hold of the [▶] button on the Ready Screen will activate the measured parameter alarms. Once measured parameter alarms are activated, there is no option to return to the Ready Screen.

Note: The device gas transfer and heat exchange functions DO NOT require power. The device functions without activation of the Touch Screen Display in the event that user does not wish to utilize the sensing capabilities of this device.

Note: The device can be powered via the power cord or 4 AA batteries (included).

Note: If the device is unplugged after it is powered on, the device will automatically switch to battery power.

9.4. Priming the Device

Caution: Tubes and luer lock connections must be fastened by hand only. Using tools or other aids may damage the ports.

Warning: Ensure that unused ports remain capped to prevent contamination and that the caps are tightened to avoid leakage.

Connect the water lines to the water ports on the oxygenator. Run water through the heat exchanger, and check for leaks from the water compartment to the blood compartment.

Warning: Do not use the oxygenator if water is present in the blood compartment.

Connect all blood and gas lines to their appropriate connections under aseptic conditions per institution protocol.

Note: Banding of all blood tubing connectors is recommended.

Fill the priming bag with the appropriate volume of priming solution per institution protocol. Remove the yellow luer lock cap on the de-airing membrane located on the blood inlet side of the device. Gravity prime the oxygenator or prime by forward pump flow per institution protocol. Gradually begin recirculation flow using a pre-bypass filter per institution protocol.

Warning: All gas emboli must be cleared from the extracorporeal circuit before initiating flow to the patient. Gas emboli are hazardous to the patient.

Add additional priming solution as required to prime the rest of the circuit. Upon completion of priming, re-cap the de-airing membrane. Ensure the entire system, including ports, is fully de-aired prior to beginning extracorporeal circulation.

Note: If the device is being stored after priming, it is recommended that the device be powered off prior to storage. When the device is powered back on, the touch screen will display the Ready Screen.

Note: Device integrity has been verified by MC3 after 30 days of primed storage with phosphate buffered saline solution.

Caution: Do not leave batteries in the device after batteries have been depleted.

9.5. Extracorporeal Circulation

Caution: Check all lines are properly connected prior to use.

Caution: The de-airing membrane must be closed during operation.


Ensure that there is no air in the circuit. Check for adequate levels of anticoagulation prior to and during extracorporeal circulation. Remove arterial and venous clamps and gradually increase blood flow. Start the gas flow with the appropriate gas settings per institution protocol. Adjust water temperature per institution protocol.


Caution: Do not exceed 42°C degrees in the heat exchanger.

9.6. Activation of Measured Parameter Alarms

Activate the measured parameter alarms by pressing and holding the [▶] button on the Ready Screen. The Main Screen is now displayed, and measured parameter alarms are active. After measured parameter alarms are activated, the module will give both audible and visual alarms when alarm limits are exceeded. The audible and visual signals will cease if the measured parameters return to within set alarm limits. See Section 12 Main Screen Tiles for a detailed explanation of items displayed on the Main Screen.

Ready Screen

Main Screen

9.7. Setting Alarm Limits


Caution: Setting alarm limits to extreme values can render the alarm system useless.

Caution: Where other equipment is used to perform the same or similar measurements as the Nautilus Smart ECMO Module in the same environment (e.g. intensive care unit or cardiac operating room), use of differing alarm limit settings may cause a hazard.

Press the tile of interest on the Main Screen to bring up its respective alarm limit setting screen. Adjust the lower alarm limit using [^] and [v] on the left side of the screen. Adjust the upper alarm limit using [^] and [v] on the right side of the screen. Press [v] to confirm the change or [X] to cancel. A checkmark screen appears when alarm limits are successfully saved.

Note: The Outlet SO₂ does not have an upper alarm limit.

Note: The alarm limit settings for each measured parameter are retained when the device is powered off.

Alarm Limit Screen


Alarm Limit Change Successful

9.8. SO₂ Calibration

Press the center SO₂ Save button on the Main Screen to initiate SO₂ calibration. Store current SO₂ readings by pressing the Save button on the SO₂ Save Screen. A checkmark screen appears when SO₂ values are successfully saved. At the same time as saving, take reference blood samples for blood gas and hematocrit analysis.


Main Screen

SO₂ Save Screen

SO₂ Save Successful

SO₂ History

Note: After saving SO₂ readings, the Main Screen will show the SO₂ Lab Values button in place of the SO₂ Save button.

Note: Tap the Clipboard before saving current values to review the previous SO₂ calibration. Press [←] to return to the SO₂ Save Screen.

Note: The History button will only appear if previous calibrations have been completed.


When blood sample results have been obtained, press the center SO₂ Lab Values button on the Main Screen. Select the patient hematocrit range obtained via the reference blood sample. If a hematocrit measurement was not gathered, leave the hematocrit at the default level of 27-32%. Press [→] to continue. Press [X] to cancel. Use [^] and [v] to enter the inlet (left) and outlet (right) laboratory SO₂ values. Press [✓] to confirm or [←] to return to the hematocrit selection screen. A checkmark screen appears when the lab values are successfully saved.


Main Screen

SO₂ Hct Screen

SO₂ Lab Values Screen

SO₂ Calibration Successful

Delete Saved Values

Note: Including a HCT value will increase the accuracy of the SO₂ calibration.

Note: SO₂ calibration can be done on the inlet, outlet or both. To adjust only the inlet SO₂ reading, leave the outlet lab value unchanged. To adjust only the outlet SO₂ reading, leave the inlet lab value unchanged.

Note: Lab values may only be entered once for each set of saved SO₂ readings.

Note: The SO₂ calibration process can be completed as many times as desired.

Note: The SO₂ offset calibrations are retained when the device is powered off.

Note: The device will only allow entry of SO₂ lab value within ±15% (absolute) of the saved value.

Note: To delete stored values and start over, press the SO₂ Delete button located top center on the SO₂ Hct or SO₂ Lab Values. Screen at any time. Press [✓] to confirm or [X] to cancel on the Delete Saved Values screen.


9.9. Power Notifications

The device will automatically switch between Mains and battery power when plugged or unplugged. When switching from Mains to battery power, the Mains to Battery Notification screen will be shown. Press [✓] to dismiss and return to the previous screen.


If the battery power drops to 25%, the 25% Percent Notification Screen will appear. Press [✓] to dismiss and return to the previous screen.

If on Mains power and battery power is critically low, the Battery Critical Notification Screen will be shown. Press [✓] to dismiss and return to the previous screen.


If on battery power and battery power is critically low, the Battery Critical Alarm screen will be shown. Plug in or replace the batteries to dismiss the screen.

Main to Battery Notification

25% Power Notification (Yellow fill)

Battery Critical Notification (Red fill)

Battery Critical Alarm (Red fill)

9.10. Replacing Batteries

Caution: Replace all batteries at one time. Do not mix old and new batteries.

Caution: Only use alkaline AA batteries.

Slide the battery door up to expose the battery compartment. Remove all 4 batteries and replace, using alkaline batteries only. Slide the battery door back into place, and snap it closed to secure.

9.11. Oxygenator Change out

The oxygenator should be evaluated for change out if the clinician determines during the extracorporeal procedure that adequate patient support is not being achieved. If a change out is deemed necessary, obtain a primed replacement device (see sections 9.3 – 9.6) before discontinuing circulation through the original device. Have materials needed for the change available before starting the change, including clamps and sterile scissors. Always use aseptic technique during change out.

1. Double clamp the tubing on the inlet side and outlet side of the oxygenator that is to be replaced, and stop the pump. Clamp far enough away from the oxygenator blood inlet/outlet to allow adequate room for reconnection to the replacement device.
2. Detach the oxygenator from the holder by pressing the release tab.
3. Attach the primed replacement oxygenator to the holder.
4. Cut the venous and arterial tubing between the clamps connected to the original oxygenator.
5. Using air-free technique, connect the arterial line to the blood outlet and the venous line to the blood inlet of the primed replacement oxygenator.
6. Open the clamp on the blood inlet side and turn the pump on at a low speed.
7. Connect the gas line to the gas inlet of the replacement device.
8. Confirm that the system is free of bubbles before removing the outlet clamp, and gradually increase the pump flow to the necessary value.
9. Connect the water lines to the replacement oxygenator.

10. SPECIFICATIONS

Oxygenation Fiber Type	Polymethylpentene (PMP)		
Fiber Surface Area – Oxygen Membrane	1.8 m ²		
Heat Exchanger	Polyethylene terephthalate (PET)		
Fiber Surface Area – Heat Exchange	0.3 m ²		
Prime Volume	226 mL		
Residual Blood Volume	189 mL		
Operating Blood Volume	226 mL		
Blood Flow Range	0.5 – 7 L/min		
Maximum Rated Blood Pressure	750 mmHg		
Gas Flow Rate (Gas:Blood)	0.5:1 – 3:1		
Maximum Gas Path Pressure	100 mmHg (1.9 psi)		
Maximum Water Path Pressure	1125 mmHg (21.8 psi)		
Power			
Nautilus Medical Grade Power Cord:			
- Mains voltage input	100 – 240 VAC		
- Frequency input	50 – 60 Hz		
- Output Voltage	12 VDC		
- Rated power	7.8 W		
- Output Current at rated power	0.65A		
Battery:			
- Type	4 AA Alkaline		
- Life	> 120 min		
Sensor Accuracy			
Parameter	Measuring Range	Resolution	Measuring Accuracy
Pressure: Pin, Pout	0 – 750 mmHg	1	0 – 249 mmHg: ± 10 mmHg 250 – 750 mmHg: ± 7 % of the measured value Maximum ± 15 mmHg drift over 30 days
Delta Pressure	0 – 750 mmHg	1	Calculated value
Oxygen Saturation: SO ₂ in, SO ₂ out	40 – 100 %	0.1	± 5% (absolute)
Temperature: Tout	0 – 42 °C	0.1	± 0.5 °C
Environmental			
	Operating	Storage	Transportation
Temperature	15 – 30 °C	10 – 30 °C	-30 – 50 °C
Relative Humidity (non-condensing)	30 - 75 %	5 - 95 %	5 - 95 %
Alarm Settings			
	Default Alarm Limits (Lower / Upper)	Minimum Setting	Maximum Setting
Inlet Pressure	50 / 300 mmHg	0 mmHg	750 mmHg
Delta Pressure	10 / 75 mmHg	0 mmHg	750 mmHg


Outlet Pressure	50 / 300 mmHg	0 mmHg	750 mmHg
Inlet SO ₂	50.0 / 90.0 %	40.0 %	100 %
Outlet SO ₂	95.0 / n/a %	40.0 %	100 %
Outlet Temperature	35.0 / 37.5 °C	20.0 °C	40.0 °C





Note: The upper alarm limit must always be greater than the lower alarm limit for any given parameter.

11. AUDIBLE AND VISUAL SIGNALS

Indicator Light	Volume	Signals	Status
Not illuminated	-	-	Device off or measured parameter alarms inactive
Green/Yellow Flash 1x	-	1x 2 beeps, 1x 3 beeps	Power on test sequence
Green	-	-	Device on, measured parameter alarms active, no alarm condition
Yellow, solid	60 - 68 db(A)	2 beeps repeating every 30 seconds	Low priority alarm
Yellow, flashing	66-71 db(A)	3 beeps repeating every 7.5 seconds	Medium priority alarm

12. MAIN SCREEN TILES

Display	Details
 Blue inlet pressure	Pressure at the oxygenator inlet barb
 Red outlet pressure	Pressure at the oxygenator outlet barb
 Gray delta pressure	Difference between the outlet and inlet membrane pressures
 Blue inlet oxygen saturation	Blood oxyhemoglobin saturation at the oxygenator inlet barb
 Red outlet oxygen saturation	Blood oxyhemoglobin saturation at the oxygenator outlet barb
 Red outlet temperature	Blood temperature at the oxygenator outlet barb
 Blue, white, or red color numerical value	Reading within alarm limits
 Black numerical value with yellow solid or yellow flashing background. Arrow up.	Reading exceeding upper alarm limit

	Black numerical value with yellow solid or yellow flashing background. Arrow down.	Reading exceeding lower alarm limit
	Black "Hi" with yellow solid or yellow flashing background. Arrow up.	Reading above valid measurement range
	Black "Lo" with yellow solid or yellow flashing background. Arrow down.	Reading below valid measurement range
	Yellow bell	To pause the audible alarm, press the button.
	Gray bell with dashed "X"	Alarms are temporarily paused for 1 minute. If a new alarm occurs, the device will un-pause the audible alarm.
	Gray bell	No alarm is currently active
	Light gray screen rotate	Press to rotate the screen 180°.
	Dark gray screen rotate	The screen rotate button is disabled in the alarm state.
	SO ₂ Save	Press to begin SO ₂ calibration process and store measured SO ₂ readings.
	SO ₂ Lab Values	Press to enter results from blood gas analysis. Displayed after saving SO ₂ readings.
	Light gray plug, dark gray full battery	Device running on Mains power. Battery full.
	Light gray plug, dark gray 75% battery	Device running on Mains power. Battery at 75%.
	Light gray plug, dark gray 50% battery	Device running on Mains power. Battery at 50%.
	Light gray plug, yellow fill 25% battery	Device running on Mains power. Battery at 25%.
	Light gray plug, red fill critical battery	Device running on Mains power. Battery level critically low.
	Dark gray plug with slash, green full battery	Device running on battery power. Battery full.
	Dark gray plug with slash, green 75% battery	Device running on battery power. Battery at 75%.
	Dark gray plug with slash, green 50% battery	Device running on battery power. Battery at 50%.
	Dark gray plug with slash, yellow fill 25% battery	Device running on battery power. Battery at 25%.

13. DISPOSAL

Warning: Dispose of the device in accordance with local hospital, administrative, and/or other government policies.

The device is a single-use only disposable. Turn off the device, unplug the power cord, and remove the batteries. Leave the power cord with the holder, and dispose of the device and batteries according to institution protocol and local regulations.

14. CLEANING

If the outside of the oxygenator becomes soiled, wipe down surfaces using a cloth moistened with suitable cleaning solution. Do NOT spray cleaner directly on the device. The device, excluding the touchscreen, is compatible with the following cleaner types: 3% hydrogen peroxide, bleach, isopropyl alcohol, aldehydes, and quaternary ammonium compounds. The touchscreen should not be cleaned with detergent based or corrosive cleaners.

Cleaning of the holder, power cord and power supply is recommended after each use or if contaminated. All surfaces of the holder can be cleaned and disinfected for blood, saline, or contaminate spills using normal medical equipment cleaners and disinfectants. Suitable cleaning solutions for sensitive medical devices should be used when cleaning the power cord and power supply.

Caution: Unplug the power cord and power supply from AC Mains prior to cleaning.

15. POTENTIAL COMPLICATIONS

Potential complications, including those normally associated with ECLS procedures and anticoagulation, may necessitate physician intervention. Some potential complications related to the use of ECLS oxygenators include but are not limited to hypoxia or hypercarbia due to inadequate gas exchange, air embolism, thrombus formation leading to deep-vein thrombosis or pulmonary embolism, hemolysis, thrombocytopenia, mechanical failure of ports and connections, infection, hypothermia or hyperthermia due to inadequate heat transfer, hemorrhage related to bleeding, and inadequate blood flow related to clotting or hypovolemia resulting in inadequate gas exchange. There are risks and adverse events related to all ECLS procedures and anticoagulation including heart, vessel, or lung damage, hypoxia, anemia, infection, hemorrhage, liver or kidney failure, stroke, or death.

16. ALARMS

16.1. Measured Parameter Alarms – Low Priority

Measured Parameter Alarms	Possible Causes
Outlet Temperature Above Upper Limit	Heater/cooler settings, incorrect measurement, incorrect limit set
Outlet Temperature Below Lower Limit	Heater/cooler settings, incorrect measurement, incorrect limit set
Outlet Temperature Outside of Measurement Range	Value outside of valid range; heater/cooler settings, incorrect measurement
Inlet SO ₂ Above Upper Limit	Gas input settings, incorrect measurement (including SO ₂ calibration), incorrect limit set
Inlet SO ₂ Below Lower Limit	Gas input settings, incorrect measurement (including SO ₂ calibration), incorrect limit set


Inlet SO ₂ Outside of Measurement Range	Value outside of valid range; gas input settings, incorrect measurement (including SO ₂ calibration)
Outlet Pressure Above Upper Limit	Occlusion post-oxygenator, flow setting, incorrect measurement, incorrect limit set
Outlet Pressure Below Lower Limit	Flow settings, incorrect measurement, incorrect limit set
Outlet Pressure Outside Measurement Range	Value outside of valid range; occlusion post-oxygenator, flow setting, incorrect measurement
Delta Pressure Above Upper Limit	Occlusion in oxygenator, flow settings, incorrect measurement, incorrect limit set
Delta Pressure Below Lower Limit	Mechanical failure, flow settings, incorrect measurement, incorrect limit set
Delta Pressure Outside Measurement Range	Value outside of valid range; occlusion in oxygenator, mechanical failure, incorrect measurement
Inlet Pressure Above Upper Limit	Occlusion in or post-oxygenator, flow settings, incorrect measurement, incorrect limit set
Inlet Pressure Below Lower Limit	Mechanical failure, incorrect measurement, incorrect limit set
Inlet Pressure Outside Measurement Range	Value outside of valid range; occlusion in or post-oxygenator, flow settings, mechanical failure, incorrect measurement

16.2. Measured Parameter Alarms – Medium Priority


Measured Parameter Alarms	Possible Causes
Outlet SO ₂ Below Lower Limit	Gas input settings, occlusion in oxygenator, incorrect measurement, incorrect limit set
Outlet SO ₂ Outside of Measurement Range	Value outside of valid range; gas input settings, occlusion in oxygenator, incorrect measurement

16.3. Technical Error – On Screen Display

In the event of a technical error condition, the Unrecoverable Error Screen or Recoverable Error Screen is displayed. Refer to the error code displayed on the screen for more information. To dismiss the Unrecoverable Error Screen, power the device off and back on. To dismiss the Recoverable Error Screen, consult the IFU, then press [✓] to continue.

Recoverable Error

Unrecoverable Error

16.4. Power On Self-Test Errors

Error Code (or code range)	Description	Possible Causes	Possible Actions
20005 – 20010 20014 – 20015 20027 20029 20031 20036	Internal memory fault	Technical malfunction	Power cycle device. Use alternate circuit measurement methods, if necessary.
51000 – 57999	Internal hardware fault	Technical malfunction	Power cycle device. Use alternate circuit measurement methods, if necessary.
30003 – 30011 40000 – 40159	Internal software fault	Technical malfunction	Power cycle device. Use alternate circuit measurement methods, if necessary.

16.5. Technical Alarms – Low Priority

Error Code (or code range)	Description	Possible Causes	Possible Actions
20028	Alarm limits invalid	Technical malfunction	Press [✓] to dismiss. NOTE: All alarm limits will be reset to factory defaults.
20030	SO ₂ calibration invalid	Technical malfunction	Press [✓] to dismiss. NOTE: SO ₂ calibration will be reset to factory default. Perform an SO ₂ calibration per the instructions for use.

16.6. Technical Alarms – Medium Priority

Error Code (or code range)	Description	Possible Causes	Possible Actions
10030 – 10031 51000 – 57999	Internal hardware fault	Technical malfunction	Power cycle device. Use alternate circuit measurement methods, if necessary.
10000 – 10018 30000 – 30008 40000 – 40159	Internal software fault	Technical malfunction	Power cycle device. Use alternate circuit measurement methods, if necessary.
10027 – 10029	Communications failure	Technical malfunction	Power cycle device. Use alternate circuit measurement methods, if necessary.
20003 – 20004 20011 – 20013 20016 – 20026 20032 – 20033	Internal memory fault	Technical malfunction	Power cycle device. Use alternate circuit measurement methods, if necessary.

17. TECHNICAL INFORMATION

There are no user serviceable parts or fuses.

Electric shock protection (IEC 60601-1)	Class II
Applied part classification (IEC 60601-1)	Integrated pressure sensors – Defibrillator-protected type CF
Degree of Ingress Protection (IEC 60529)	IPX2
Mode of Operation	Continuous

The following information and data are available upon request: sterilization method, materials of the blood pathway, blood pathway pressure drop, gas pathway pressure drop, blood cell damage, particle release, and relevant tolerances for specifications in the instructions for use.

The essential performance of the Nautilus Smart ECMO Module is to oxygenate blood, remove carbon dioxide, and heat or cool blood. The essential performance of the device is independent of all electronics. It will not be affected by electromagnetic interference.

If continued operation of the device is required even in the event of power supply interruptions, it is recommended the user run the device on battery power.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The Nautilus Smart ECMO Module is intended to be used in hospitals (except near active HF surgical equipment and the RF shielded room of magnetic resonance imaging, where the intensity of electromagnetic disturbances is high). The user of the device should ensure it is used in an appropriate environment. Mains power quality and the strength of power-frequency magnetic fields should be that of a typical professional healthcare facility.

Note: The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Emission Test	Compliance	Electromagnetic Environment Guidelines
RF emissions (CISPR11)	Group 1	The Nautilus Smart ECMO Module uses H/RF energy solely for its internal operation. Therefore, its RF emissions are very low and not likely to cause interference in nearby electrical equipment.
RF emissions (CISPR11)	Class A	The device is suitable for use in professional healthcare and industrial establishments not directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions (IEC 61000-3-2)	Class A	
Voltage fluctuation/flicker emissions (IEC 61000-3-3)	Complies	

The Nautilus Smart ECMO Module is intended to be used in the electromagnetic environment specified below. The user should ensure that it is used in such an environment.





Immunity Test	IEC 60601-1-2 Test Level
Electrostatic Discharge (IEC 61000-4-2)	±8 kV contact ±2, ±4, ±8, ±15 kV air
Electrical fast transients/bursts (IEC 61000-4-4)	Mains power: ±2 kV 100 kHz repetition frequency
Surges (IEC 61000-4-5)	±0.5, ±1kV line to line ±0.5, ±1kV, ±2 kV line to ground

Rated Power Frequency Magnetic Fields (IEC 61000-4-8)	30 A/m 50 Hz or 60 Hz
Conducted disturbances induced by RF fields (IEC 61000-4-6)	3V 0.15MHz – 80MHz 6V in ISM (and amateur radio bands for HHC) between 0.15MHz and 80MHz 80% AM at 1kHz
Radiated RF Electromagnetic Fields (IEC 61000-4-3)	3 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz
Voltage Dips (IEC 61000-4-11)	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
	0 % UT; 1 cycle 70 % UT; 25/30 cycles Single phase: at 0°
Voltage Interruptions (IEC 61000-4-11)	0 % UT; 250/300 cycle
Note: UT is the AC supply voltage to be used for application of the test rule.	

18. GENERAL INFORMATION

ECLS General Information. This product is one component of an ECLS circuit. A typical ECLS circuit includes tubing, connectors, a gas exchange device and a blood pump. Configuration of the circuit and selection of individual component types may vary according to each specific use and treatment center. Each component of the circuit is commercially available. The extracorporeal circuit is typically managed with a systemic anticoagulant (usually heparin). The anticoagulation is managed to titrate either the active whole blood clotting time or activated partial thromboplastin time to a level above the normal ranges (typically 1-2 times normal). The principle pathway of the ECLS circuit consists of blood being drained from the drainage lumen of the catheter, passing through a blood pump which disperses the blood through a gas exchange device and back into the patient through the cannula's reinfusion lumen. The systemic temperature is controlled by a heat exchanger integrated into the gas exchange device. The components of the circuit are selected to meet the blood flow and gas exchange needs of each patient based on their metabolic requirements. The circuit has been widely described in the literature and is summarized in each iteration of the textbook published by the Extracorporeal Life Support Organization.¹ Maintenance of the ECLS circuit is managed by personnel who have been trained on the intricacies of extracorporeal support including physicians, perfusionists, nurses, respiratory therapists and other medical specialists. These individuals are qualified to manage the circuit based on training and continuing education per specific institutional guidelines.

¹ Toomasian JM, Lawson DS, and Harris WH. ECLS Equipment. Chapter 8 in ECMO: Extracorporeal Cardiopulmonary Support in Critical Care, Fourth edition, Annich G, Bartlett RH, Lynch WR, MacLaren G, Wilson JM (eds). Extracorporeal Life Support Organization, Ann Arbor, 2012. pp 107-32.

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MC3, Inc. Disclaimer of Warranty

MC3, INC. MANUFACTURES ITS PRODUCTS (REFERRED TO AS THE "PRODUCTS") UNDER CONDITIONS THAT MEET APPLICABLE REGULATORY REQUIREMENTS, BUT HAS NO CONTROL OVER HOW THE PRODUCTS ARE ACTUALLY USED OR THE CONDITIONS UNDER WHICH THEY ARE USED. AS A RESULT, MC3, INC. DISCLAIMS ALL PRODUCT WARRANTIES INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MC3, INC. SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE.

Limitation on Liability

IN NO EVENT WILL MC3, INC.'S TOTAL LIABILITY TO PURCHASER OR ANY THIRD PARTY ARISING FROM THE SALE OF PRODUCTS EXCEED THE PURCHASE PRICE OF THE PRODUCTS SUBJECT TO THE CLAIM REGARDLESS OF THE THEORY OF LIABILITY OF THE CLAIM.

Limitation on Damages

MC3, INC. IS NOT RESPONSIBLE TO PURCHASER OR ANY THIRD PARTY FOR ANY DIRECT, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY CLAIMED BREACH OF WARRANTY, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT. Applicable law in your jurisdiction may limit the ability of a manufacturer to disclaim some warranties or exclude or limit damages. While the above limitations are an integral part of the business relationship between MC3, Inc. and the purchasers of its Products, nothing in the above is intended to contravene these provisions of applicable law. If any part or term of this Disclaimer of Warranty or the Limitations on Liability or Damages is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of all remaining provisions shall not be affected and shall be construed and enforced as if the Disclaimer of Warranty or the Limitations on Liability or Damages did not contain the particular part or term held to be invalid.

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MC3, Inc. -takuun vastuuvapauslauseke
MC3, Inc. Garantifriskrivning
MC3, Inc. Garanti Feragati