Medtronic

Affinity Pixie™

Oxygenation System with Balance™ Biosurface
Système d'oxygénation avec biosurface Balance™
Oxygenator mit Balance™ Biosurface
Sistema de oxigenación con superficie biocompatible Balance™
Oxygenatiesysteem met Balance™-bio-oppervlaktelaag
Sistema di ossigenazione con superficie biopassiva Balance™
Oksygeneringssystem med Balance™ biooverflate
Balance™-biopinnoitteella päällystetty hapetinjärjestelmä
Oxygeneringssystem med Balance™ Biosurface
Σύστημα οξυγόνωσης με βιοεπιφάνεια Balance™
Oxygenationssystem med Balance™ Bio-overflade
System oksygenacji z biopowłoką Balance™
Sistema de oxigenação com bio-superfície Balance™
Balance™ Biyo Yüzevli Oksijenasyon Sistemi
Oxygenační systém s biopovrchem Balance™
Oxygenátor-rendszer Balance™ biofelülettel
Système oxygénacie s biologickou povrchovou úpravou Balance™ Biosurface
Balance™ biokattega oksügenatsioonisüsteem
Система за оксигениране с биоповърхност Balance™
Sustav za oksigenaciju s biopovršinom Balance™
Oksigenācijas sistēma ar Balance™ Biosurface virsmas pārklājumu
Oksigenācijas sistema su „Balance™“ biologine danga
Oxygenator cu biosuprafață Balance™
Система оксигенации с биопокрытием Balance™
Sistem za oksigenizaciju sa Balance™ Biosurface površinom
Oksigenacijski sistem z biološko površino Balance™
BBP211, BBP241

Instructions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso
Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • Käyttöohjeet • Bruksanvisning
Οδηγίες χρήσης • Brugsanvisning • Instrukcja użytkowania • Instruções de utilização
Kullanım Talimatları • Návod k použití • Használati útmutató • Pokyny na používanie
Kasutusjuhend • Инструкции за употреба • Upute za upotrebu • Lietošanas pamācība
Naudojimo instrukcijos • Instrucțiuni de utilizare • Инструкция по эксплуатации
Uputstva za upotrebu • Navodila za uporabo

USA Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Affinity[®] is a registered trademark of Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™], and Balance[™] are trademarks of Medtronic, Inc. Affinity[®] est une marque déposée de Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] et Balance[™] sont des marques commerciales de Medtronic, Inc. Affinity[®] ist eine eingetragene Marke von Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] und Balance[™] sind Marken von Medtronic, Inc. Affinity[®] es una marca comercial registrada de Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] y Balance[™] son marcas comerciales de Medtronic, Inc. Affinity[®] is een geregistreerd handelsmerk van Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] en Balance[™] zijn handelsmerken van Medtronic, Inc. Affinity[®] è un marchio registrato della Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] e Balance[™] sono marchi di fabbrica della Medtronic, Inc. Affinity[®] er et registrert varemerke for Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] og Balance[™] er varemerker for Medtronic, Inc. Affinity[®] on Medtronic, Inc:n rekisteröity tavaramerkki. Affinity Pixie[™], Affinity Orbit[™] ja Balance[™] ovat Medtronic, Inc:n tavaramerkkejä. Affinity[®] är ett registrerat varumärke som tillhör Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] och Balance[™] är varumärken som tillhör Medtronic, Inc. Το Affinity[®] είναι σήμα κατατεθέν της Medtronic, Inc. Το Affinity Pixie[™], Affinity Orbit[™] και Balance[™] είναι εμπορικά σήματα της Medtronic, Inc. Affinity[®] er et registreret varemærke tilhørende Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] og Balance[™] er varemærker tilhørende Medtronic, Inc. Affinity[®] jest zastrzeżonym znakiem towarowym firmy Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] i Balance[™] są znakami towarowymi firmy Medtronic, Inc. Affinity[®] é uma marca comercial registada da Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] e Balance[™] são marcas comerciais da Medtronic, Inc. Affinity[®] Medtronic, Inc. şirketinin tescilli ticari markasıdır. Affinity Pixie[™], Affinity Orbit[™] ve Balance[™] Medtronic, Inc. şirketinin ticari markalarıdır. Affinity[®] je registrovaná ochranná známka spoločnosti Medtronic, Inc. Affinity Pixie[™] a Affinity Orbit[™] a Balance[™] jsou ochranné známky společnosti Medtronic, Inc. Az Affinity[®] a Medtronic, Inc. bejegyzett védjegye; az Affinity Pixie[™], az Affinity Orbit[™] és a Balance[™] a Medtronic, Inc. védjegye. Affinity[®] je registrovaná ochranná známka spoločnosti Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] a Balance[™] sú ochranné známky spoločnosti Medtronic, Inc. Affinity[®] on ettevõtte Medtronic, Inc. registreeritud kaubamärk. Affinity Pixie[™], Affinity Orbit[™] ja Balance[™] on ettevõtte Medtronic, Inc. kaubamärgid. Affinity[®] e регистрирана търговска марка на Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] и Balance[™] са търговски марки на Medtronic, Inc. Affinity[®] je zaštitni znak tvrtke Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] i Balance[™] zaštitni su znakovi tvrtke Medtronic, Inc. Affinity[®] ir registrēta Medtronic, Inc. preču zīme. Affinity Pixie[™], Affinity Orbit[™] un Balance[™] ir Medtronic, Inc. preču zīmes. „Affinity[™]” yra registruotas „Medtronic, Inc.” prekės ženklas. „Affinity Pixie[™]”, „Affinity Orbit[™]” ir „Balance[™]” yra „Medtronic, Inc.” prekių ženklai. Affinity[®] este marcă comercială înregistrată a Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] și Balance[™] sunt mărci comerciale ale Medtronic, Inc. Affinity[®] является зарегистрированным товарным знаком корпорации Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] и Balance[™] являются товарными знаками корпорации Medtronic, Inc. Affinity[®] je registrovani žig proizvedea Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] i Balance[™] su žigovi proizvedea Medtronic, Inc. Affinity[®] je registrirana blagovna znamka družbe Medtronic, Inc. Affinity Pixie[™], Affinity Orbit[™] in Balance[™] so blagovne znamke družbe Medtronic, Inc.

Explanation of symbols on package labeling / Explication des symboles des étiquettes sur l'emballage / Erläuterung der Symbole auf der Verpackung / Explicación de los símbolos que aparecen en el etiquetado del envase / Verklaring van de symbolen op de verpakkinglabels / Spiegazione dei simboli sulle etichette della confezione / Forklaring av symboler på produktet og pakningen / Pakkauksen symbolien selitykset / Förklaring av symboler på förpackningsetiketten / Επεξήγηση των συμβόλων στη σήμανση της συσκευασίας / Forklaring af symboler på emballagens mærkater / Objasnienie symboli zamieszczonych na etykietach opakowania / Explicação dos símbolos nas etiquetas da embalagem / Ambalaj etiketindeki sembollerin açıklamaları / Vysvětlení symbolů uvedených na štítcích obalu / A csomagoláson látható szimbólumok jelentése / Vysvetlivky k symbolom na označení balenia / Pakendi siltidel olevate tingmärkide tähendus / Объяснение на символите от этикетите на опаковката / Objašnjenje simbola na ambalažnoj etiketi / Iesaiņojuma marķējuma simbolu skaidrojums / Pakuotės etikečių simbolių reikšmės / Explicarea simbolurilor cuprinse în etichetele de pe ambalaj / Объяснение символов, приведенных на этикетке упаковки / Objašnjenje simbola na oznakama na pakovanju / Razlaga simbolov na embalaži

Refer to the outer package label to see which symbols apply to this product / Se référer à l'étiquette sur l'emballage extérieur pour savoir quels symboles s'appliquent à ce produit / Welche Symbole für dieses Produkt zutreffen, entnehmen Sie bitte dem Etikett auf der Verpackung / Consulte la etiqueta del envase exterior para comprobar qué símbolos son aplicables a este producto / Controleer het label op de buitenverpakking om te zien welke symbolen op dit product van toepassing zijn / Fare riferimento all'etichetta sulla confezione esterna per conoscere i simboli applicati a questo prodotto / Se på etiketten utenpå pakningen for å fastslå hvilke symboler som gjelder for dette produktet / Katso ulkopakkauksesta, mitkä symbolit koskevat tätä tuotetta / Se etiketten på ytterförpackningen för de symboler som gäller denna produkt / Ανατρέξτε στην ετικέτα της εξωτερικής συσκευασίας για να δείτε ποια σύμβολα αφορούν αυτό το προϊόν / Se mærkaten på den ydre emballage for de symboler, der gælder for dette produkt / Należy zapoznać się z etykietą na zewnętrznym opakowaniu, aby określić, które symbole dotyczą niniejszego produktu / Consultar as etiquetas exteriores da embalagem para ver que símbolos se aplicam a este produto / Bu üründe hangi sembollerin geçerli olduğunu görmek için dış ambalaj etiketine başvurun / Symbole, které se týkají tohoto produktu, najdete na štítcích obalu / A termékre vonatkozó szimbólumok a külső csomagoláson találhatóak / Prezrite si označenie na vonkajšom obale a zistíte, ktoré symboly sa vzťahujú na tento výrobok / Välispakendi sildilt näete, millised tingmärgid selle toote puhul kehtivad / Погледнете етикета на външната опаковка, за да видите какви символи се отнасят за този продукт / Na etiketama na ambalaži pogledajte koji se simboli primjenjuju na ovaj proizvod / Skatiet ārējā iepakojuma marķējumu, lai redzētu simbolus, kuri pielietojami šim izstrādājumam / Norēdami sužinoti, koki simbolii pažymėtas šis gaminy, žiūrėkite į išorinę pakuotės etiketę / Simbolurile aplicabile pentru acest produs sunt pe eticheta ambalajului / См. на внешней этикетке упаковки, какие символы применимы к данному продукту / Pogledajte spoljnu oznaku na pakovanju da biste videli koji se simboli primenjuju na ovaj proizvod / Opis simbolov, ki se nanašajo na izdelek, je na nalepkah na ovojnini

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
Affinity Pixie™

Oxygenation System with Balance™ Biosurface

Model

BBP211	Hollow Fiber Oxygenator with Balance™ Biosurface ¹
BBP241	Hollow Fiber Oxygenator and Cardiotomy/Venous Reservoir with Balance™ Biosurface

Accessories Sold Separately

ATP210	Affinity® Temperature Probe
AUH2093	Affinity Orbit™ Holder System
AMH2014	Affinity® Manifold Holder

Description

Hollow Fiber Oxygenator

The Affinity Pixie Hollow Fiber Oxygenator with Balance Biosurface is a single-use, microporous, hollow-fiber, gas-exchange device with plasma-resistant fiber and integrated heat exchanger. The primary blood-contacting surfaces of the oxygenator are coated with Balance Biosurface to enhance blood compatibility and provide thromboresistant blood-contacting surfaces.

The following accessories are packaged for use with the oxygenator (Figure 3):

- Two 0.5 cm (3/16 in) adapters for blood inlet and outlet ports
- One recirculation line — 0.5 cm (3/16 in)

Cardiotomy/Venous Reservoir

The Affinity Pixie Cardiotomy/Venous Reservoir (CVR) with Balance Biosurface is a single-use device designed to collect and store blood during extracorporeal circulation. Cardiotomy blood is collected, filtered, and defoamed before mixing with the venous blood, which is filtered. The primary blood-contacting surfaces of the CVR are coated with Balance Biosurface to enhance blood compatibility and provide thromboresistant blood-contacting surfaces.

The following accessories are packaged for use with the CVR (Figure 3):

- Two 0.5 cm (3/16 in) adapters for venous return inlet and outlet ports
- Two 1.0 cm (3/8 in) adapters for venous return inlet and outlet ports
- Two flexible luer lock adapters
- One sampling manifold assembly

Specifications

Oxygenator

Oxygenator Membrane Type	Microporous Polypropylene Hollow Fiber
Maximum Gas:Blood ratio	2:1
Heat Exchanger	Polyethylene Terephthalate (PET)
Oxygenator Membrane Surface Area	0.67 m ²
Static Priming Volume	48 mL
Recommended Blood Flow Range	0.1–2.0 L/min
Maximum Water Pressure	206 kPa (1550 mm Hg)
Maximum Rated Blood Pressure	100 kPa (750 mm Hg)

Cardiotomy/Venous Reservoir

Reservoir Volume Capacity	1200 mL
Recommended Blood Flow Range	0.1–2.0 L/min
Maximum Cardiotomy Flow Rate	2.0 L/min
Minimum Operating Level	20 mL
Cardiotomy Filter	30 µm nominal Polyester Depth Filter
Venous Inlet Screen	64 µm
Maximum Rated Pressure	+20 mm Hg / –100 mm Hg
Positive/Negative Pressure Relief Valve	<5 mm Hg positive / >60 mm Hg vacuum

¹ Technology licensed under agreement from BioInteractions, Limited, United Kingdom.

Indications for Use

The Affinity Pixie Hollow Fiber Oxygenator with Balance Biosurface is intended to be used in an extracorporeal perfusion circuit to oxygenate and remove carbon dioxide from the blood and to cool or warm the blood during routine cardiopulmonary bypass procedures up to 6 hours in duration. The Affinity Pixie Hollow Fiber Oxygenator is suitable for transfer of the volatile anesthetics isoflurane, sevoflurane, desflurane, and enflurane by administration through the oxygenator gas inlet by means of a suitable gas vaporizer.

The Affinity Pixie Cardiotomy/Venous Reservoir with Balance Biosurface is intended to be used in an extracorporeal perfusion circuit to collect venous and cardiotomy suctioned blood during routine cardiopulmonary procedures up to 6 hours in duration. The CVR is also intended for use during vacuum assisted venous drainage (VAVD) procedures.

Contraindications

This device used for any purpose other than the indicated use is the responsibility of the user.

Warnings

- Read all Warnings, Precautions, and Instructions for Use carefully prior to use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**
- Only persons thoroughly trained in cardiopulmonary bypass procedures should use this device. Operation of each device requires constant supervision by qualified personnel for patient safety.
- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use if the package is opened or damaged, the device is damaged, or the protective caps are not in place.
Note: The oxygenator water inlet, water outlet, and gas outlet ports do not have protective caps.
- All gas emboli must be cleared from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- If air leaks are observed during priming and/or operation, this may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the device if these conditions are observed.
- In any extracorporeal circuit, the use of a pre-bypass and arterial line filter is recommended.
- Monitoring of circuit pressure is recommended.
- The blood phase pressure should be greater than the gas phase pressure at all times in the oxygenator.
- The oxygenator and CVR should not be operated outside the recommended parameter ranges.
- Total combined flow into the cardiotomy filter should not exceed 2.0 L/min.
- Do not exceed blood pressure of 100 kPa (750 mm Hg) in the blood pathway.
- Once an oxygenator is primed with blood, adequate heparinization should be maintained per institution CPB protocol and the blood pathway should be constantly recirculated at a low flow rate.
- After initiation of bypass, blood flow should circulate through oxygenator and CVR at all times within the recommended blood flow range during cardiopulmonary bypass, except during emergency oxygenator replacement.
- Heat exchanger temperature should not exceed 42°C.
- Disinfectants must not be used in the heater/cooler system while the heat exchanger is in use; if disinfectants have been used in the heater/cooler system, the system must be thoroughly flushed prior to use.
- All blood tubing connectors should be banded for added protection against disconnection.
- Ensure that unused ports remain capped to prevent contamination and that the caps are tightened to avoid leakage.
- A replacement oxygenator should be readily available during cardiopulmonary bypass.
- Do not obstruct gas outlet vents to prevent gas side pressures from exceeding blood side pressures.
- Do not occlude or obstruct the vent port of the cardiotomy/venous reservoir during operation (except when using VAVD).
- Gas transfer characteristics of the oxygenator are not significantly affected by concentrations of anesthetic gas up to 2%. Above 2% anesthetic gas concentrations, adjustments of FiO₂ and gas flow rates may be needed to achieve desired gas transfer performance.
- Gas transfer rates may change over time and adjustments of FiO₂ and gas flow rates may be needed to achieve desired gas transfer performance.
- The use of a blood level sensing mechanism during operation of this device is recommended.
- Ensure that the venous reservoir outlet is always positioned above the highest point in the membrane compartment of the oxygenator.

- Each device has been carefully manufactured, tested, and packaged; however the state of the art has not been developed to the point that Medtronic is able to ensure that the device will not leak, crack, or fail during use. Cardiopulmonary bypass must be carefully and constantly monitored.
- Tubing should be attached in such a manner as to prevent kinks or restrictions that may alter the blood, water, or gas flow.
- Do not allow alcohol, alcohol-based fluid, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come in contact with the device as they may jeopardize the structural integrity.

Precautions

- Refer to package labeling for storage temperature requirements.
- Each device is sterilized using EtO.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device which could result in patient injury, illness, or death.
- Use aseptic technique in all procedures.
- Use of a gas scavenging system during transfer of anesthetic gas through the oxygenator should be considered.
- A strict anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician. Adequate anticoagulation should be maintained per institution CPB protocol.

Adverse Effects

- Possible side effects include, but are not limited to, hypoxia, hypercarbia, hypovolemia, infection, toxemia, allergenicity, mechanical failure, hemolysis, platelet dysfunction/excessive post-op bleeding, air embolism, blood loss, circulatory compromise, user exposed to patient blood, and thromboembolic phenomena. These are potential side effects with all extracorporeal blood circuits.

Instructions for Products with Balance Biosurface

The blood-contacting surfaces of the device are coated with Balance Biosurface to enhance blood compatibility and provide thromboresistant blood-contacting surfaces.

Warning: A product coated with Balance Biosurface is intended for single use only.

Caution: Appropriate blood anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.

Caution: Do not store a Balance Biosurface-coated product above 50°C.

Instructions for Use

System Setup

Warning: The setup and use of the device(s) is the responsibility of the attending clinician.

1. Remove the device(s) from the packaging carefully to ensure sterile fluid pathway.

Warning: Ensure aseptic technique is used during all stages of setup and use of this product.

Warning: Before removing the device(s) from the packaging, inspect the packaging and product for damage. If the packaging or product is damaged, do not use; sterility may have been compromised and/or performance may be affected.

2. If using the Affinity Orbit Holder System (sold separately), attach the holder mast clamp to the heart-lung machine mast at the desired height and tighten the mast clamp knob.
3. Turn the locking lever perpendicular to the holder arm to unlock it.
4. Slide the oxygenator arm onto the post if it is not already in place.
5. Slide the oxygenator holder connector onto the bottom of the oxygenator arm. A click will indicate the oxygenator is firmly seated.

Note: To remove the oxygenator, press the holder release tab and pull the oxygenator away from the connector arm.

6. If using the Affinity Pixie CVR, slide the CVR holder sleeve onto the holder post.
7. Position the oxygenator and CVR as desired for procedure. The oxygenator may be rotated clockwise or counterclockwise and the oxygenator arm may be moved laterally. The CVR may be rotated clockwise or counterclockwise.
8. Turn the locking lever parallel to the holder arm to lock components in place.
9. Turn the holder arm knob to lock the holder arm in place. See Figure 4 for final system assembly.

10. Remove the obturator cap from the CVR positive/negative pressure relief valve.
11. Connect all blood, gas, and suction lines to their appropriate connections per institution CPB protocol. If necessary, attach an appropriately sized adapter to the blood inlet or outlet port to accommodate 0.5 cm (3/16 in) or 1.0 cm (3/8 in) tubing.
12. The sampling manifold should be attached so that the tubing segment with the one-way duckbill valve is positioned between the oxygenator arterial sample port and manifold and that arterial blood flows through the duckbill valve toward the manifold.
13. Connect the recirculation line to the recirculation port on the oxygenator. Connect the other end of the recirculation line to a luer port on the CVR.
Warning: If not used, the recirculation port should be sealed tightly with a luer cap.
14. Connect quick-connect 1.3 cm (1/2 in) water lines to inlet and outlet ports on the oxygenator. Start the water flow and check for leaks from the water compartment to the blood compartment. Discard oxygenator if water is present in the blood compartment and set up the standby oxygenator.
15. Ensure the CVR is vented to atmosphere by removing the vent port cap.
Note: A CO₂ flush of the oxygenator is recommended prior to priming.

Priming & Recirculation

1. Fill the CVR with the appropriate volume of priming solution per institution CPB protocol.
2. Remove all bubbles in the section of tubing (roller pump head/boot) or centrifugal pump (if used) from the CVR outlet to the oxygenator inlet.
3. Gravity prime the oxygenator or prime by forward pump flow per institution CPB protocol.
4. Gradually begin recirculation flow using a pre-bypass filter per institution CPB protocol.
5. Ensure the oxygenator is bubble free.
6. Add additional priming solution, if needed, to prime the rest of the extracorporeal circuit.
7. Upon completion of priming and debubbling, gradually reduce blood flow and stop pump, close all purge lines, and clamp arterial, venous, and recirculation lines.
8. Ensure the entire extracorporeal circuit is free of bubbles prior to initiating bypass.
Note: Priming solution may be prewarmed through heat exchanger prior to initiating bypass.

Initiation of Bypass

Warning: The blood phase pressure should be greater than the gas phase pressure at all times.

1. Ensure the gas outlet port is not obstructed.
2. Check for adequate levels of anticoagulation prior to initiating bypass.
3. Remove arterial and venous clamps and gradually increase blood flow. Then initiate gas flow with appropriate gas settings per institution CPB protocol and clinical judgment specific to the patient and procedure.
4. Adjust water temperature to meet clinical demands.

Operation During Bypass

1. Arterial PO₂ is controlled by varying the percent concentration of oxygen present in the ventilating gas.
 - To DECREASE the PO₂, DECREASE the percent of oxygen in the ventilating gas by adjusting the FiO₂ on your gas blender.
 - To INCREASE the PO₂, INCREASE the percent of oxygen in the ventilating gas by adjusting the FiO₂ on your oxygen blender.
2. PCO₂ is controlled primarily by varying the total gas flow rate.
 - To DECREASE PCO₂, INCREASE the total gas flow rate to increase the amount of CO₂ removed.
 - To INCREASE PCO₂, DECREASE the total gas flow rate to decrease the amount of CO₂ removed.
3. Patient temperature is controlled by regulating the temperature of the water flow in the heat exchanger.

Warning: After initiation of bypass, blood flow should circulate through oxygenator and CVR at all times within blood flow range during cardiopulmonary bypass, except during emergency oxygenator replacement.

Note: Care should be taken to minimize temperature gradient during cooling and warming phases of bypass.

Termination of Bypass

Terminate bypass per institution CPB protocol as individual case and patient condition indicates.

Emergency Oxygenator Replacement

A standby oxygenator and CVR should always be available during cardiopulmonary bypass.

1. Turn off water flow to the heat exchanger, then clamp and remove the water lines.

2. Clamp the venous line at the CVR. Turn off arterial roller pump and clamp arterial line. (If using a centrifugal pump, clamp the arterial line first, prior to stopping the pump.)
3. Remove oxygen line from the gas inlet port.
4. Remove any monitoring/sampling lines from the oxygenator.
5. Double clamp the recirculation line, oxygenator inlet, and outlet tubing. Cut between clamps leaving adequate lengths for reconnections.
6. Detach the oxygenator from the holder by pressing the release tab.
7. Attach the replacement oxygenator onto the holder.
8. Reconnect the recirculation, oxygenator, and outlet tubing. Ensure all connections are secure.
9. Reconnect the oxygen line to the gas inlet port.
10. Connect and unclamp the water lines, then turn on water source and inspect for water leaks.

Warning: Do not unclamp arterial outlet line or venous return line at this time.

11. With sufficient volume in the CVR, turn on the blood pump and slowly fill the oxygenator.
12. Increase blood flow through the recirculation line.
13. Stop the pump and clamp the recirculation line.
14. Ensure the entire system is free of leaks and gas bubbles.
15. Remove the clamps from venous and arterial lines and reinstate bypass.

Cardiotomy/Venous Reservoir, Vacuum Assisted Venous Drainage (VAVD)

Warnings

- Do not use a vacuum regulator with maximum negative vacuum pressure greater than -20 kPa (-150 mm Hg). Use of excessive negative pressure may lead to hemolysis.
- Do not occlude or obstruct the vent port of the cardiotomy/venous reservoir during operation (except when using VAVD).
- Draw only intended drug amounts into the syringe that will be used for drug administration into the CVR during VAVD. Negative pressure in the CVR could draw more drug from the syringe into the reservoir than intended and lead to excess drug delivery.
- All A/V shunt lines must be closed to atmosphere prior to stopping the pump or using low blood flow rates during VAVD. Failure to do so may allow air to be drawn into the blood side of the oxygenator from the fiber.
- Appropriate occlusion of the arterial roller pump head is necessary when using VAVD. Failure to appropriately occlude the pump may allow air to be drawn into the blood side of the oxygenator from the fiber.
- An arterial filter should be used during VAVD techniques to minimize potential delivery of air emboli.
- Gradually return to atmospheric pressure when ending use of VAVD. A sudden change in pressure may result in turbulent blood flow inside the CVR.
- The line between the oxygenator and centrifugal pump (if used) must be clamped before stopping the pump during use of VAVD. Failure to clamp off the arterial line may result in air drawn into the blood side of the oxygenator from the fiber. Use of an arterial line one-way valve between the oxygenator and centrifugal pump is recommended.

Precautions

- A controlled vacuum regulator is required during VAVD.
- A positive/negative pressure relief valve is required during VAVD (included, see Figure 2).
- Use of a pressure measurement device on the CVR and an additional negative pressure relief valve that operates at -20 kPa (-150 mm Hg) should be considered.
- A vapor trap is required during use of VAVD due to condensation development.

VAVD Instructions for Use

1. Follow previously listed Instructions for Use for system setup with the following modifications as listed in steps 2 through 3.
2. Attach a calibrated means of monitoring pressures to either the venous line or venous reservoir.
3. All manifold sampling ports must have nonvented caps when using VAVD.
4. Follow previously listed Instructions for Use for priming with the following additions as listed in steps 5 through 10 below.
5. Prior to initiating bypass, prepare the vacuum portion of the setup. Attach a vented vacuum line to the CVR vent port and to the vacuum regulator with vapor trap.

Warning: Clamp quick-prime line prior to application.

Warning: Do not allow the vapor trap to become completely filled during use. This may allow fluid to enter the vacuum controller, or may prevent the CVR from being vented to atmosphere when no vacuum is applied.

6. To vent the vacuum line, ensure the Y-connector with side tubing is easily accessible and open to atmosphere. This may be clamped and unclamped regularly throughout the procedure to apply or stop vacuum.

Caution: Do not exceed vacuum levels greater than -6.67 kPa (-50 mm Hg) as measured in the venous line.

7. Ensure all access ports on CVR and manifold are properly sealed and are secure before initiating bypass.

Caution: Vacuum will result in excess atmospheric air being entrained into the CVR and potentially into the blood pathway with opened blood pathway ports. Vacuum may result in increased flow through the sampling manifold during use.

8. Use vacuum assistance as needed to institute and maintain cardiopulmonary bypass.

Warning: Do not apply vacuum to the venous reservoir when there is no forward blood flow through the oxygenator. This applies to arterial centrifugal and roller pumps (roller may not be occlusive in all positions). This will prevent air from being pulled across the membrane into the blood path by the reservoir vacuum.

Warning: Do not allow the CVR to become pressurized as this could obstruct venous drainage, force air retrograde into the patient or into the blood path of the oxygenator.

Important Notice — Limited Warranty (U.S. Customers Only)

- A. This **LIMITED WARRANTY** provides the following assurances to the patient who receives an Affinity Pixie™ Hollow Fiber Oxygenator and/or Cardiotomy/Venous Reservoir with Balance™ Biosurface, hereafter referred to as “The Product”.
- (1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its “Use By” date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the Purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
 - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its “Use By” date.
 - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
 - (3) The Product must not have been altered or subjected to misuse, abuse or accident.
 - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) Except as expressly provided by this **LIMITED WARRANTY**, Medtronic is not responsible for any direct, incidental or consequential damages based on any defect, failure or malfunction of the product, whether the claim is based on warranty, contract, tort or otherwise.
 - (2) THIS **LIMITED WARRANTY** IS MADE ONLY TO THE PATIENT ON WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS **LIMITED WARRANTY** SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
 - (3) The exclusions and limitation set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid. This **LIMITED WARRANTY** gives the patient specific legal rights. The patient may also have other rights which vary from state to state.
 - (4) No person has any authority to bind Medtronic to any representation, condition or warranty except this **LIMITED WARRANTY**.

Note: This **LIMITED WARRANTY** is provided by Medtronic, Inc. 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States.

Important Notice — Limited Warranty (For Countries Outside the U.S.)

- A. This **LIMITED WARRANTY** provides assurance for the purchaser who receives an Affinity Pixie™ Hollow Fiber Oxygenator and/or Cardiotomy/Venous Reservoir with Balance™ Biosurface, hereafter referred to as the “Product”, that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.
- The Warnings contained in the product labeling are considered an integral part of this **LIMITED WARRANTY**. Contact your local Medtronic representative to obtain information on how to process a claim under this **LIMITED WARRANTY**.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its “Use By” date.
 - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
 - (3) The Product may not have been used for any other patient.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
 - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall

be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

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© 2010 Medtronic