VitalFlow™ Caddy 58120  
Контейнер VitalFlow™ 58120  
Úložný modul VitalFlow™ 58120  
VitalFlow™-ramme 58120  
VitalFlow™ Traggestell 58120  
Πλάσιο VitalFlow™ 58120  
Soporte VitalFlow™ 58120  
Kandehoidik VitalFlow™ 58120  
VitalFlow™-kehikko 58120  
Caddie VitalFlow™ 58120  
Nosač VitalFlow™ 58120  
VitalFlow™ 58120-as típusú tároló  
Kedi VitalFlow™ 58120  
Contentitore VitalFlow™ 58120  
VitalFlow™ 58120 ұстағышы  
VitalFlow™ কেടി 58120  
„VitalFlow™“ korpusas 58120  
VitalFlow™ kārba (58120. modelis)  
Носач VitalFlow™ 58120  
VitalFlow™-houder Model 58120  
VitalFlow™-modul 58120  
Stojak VitalFlow™ 58120  
Estrutura de apoio VitalFlow™ 58120  
Base de suporte VitalFlow™ 58120  
Cutie VitalFlow™ 58120  
Контейнер VitalFlow™ 58120  
Rám VitalFlow™ 58120  
Nosilec VitalFlow™ 58120  
VitalFlow™ nosiljka 58120  
VitalFlow™ hållare 58120  
VitalFlow™ Taşıyıcı 58120  
Контейнер VitalFlow™ 58120  
Khung gán VitalFlow™ 58120

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning •  
Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet •  
Mode d'emploi • Upute za upotrebu • Használati útmutató • Petunjuk Penggunaan • Istruzioni  
per l'uso • Пайдалану нұсқаулары • 사용 지침 • Naudojimo instrukcija • Lietošanas pamācība •  
Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania •  
Instruções de uso • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по  
применению • Pokyny na používanie • Navodila za uporabo • Uputstvo za upotrebu •  
Bruksanvisning • Kullanım Talimatları • Інструкція з використання • Hướng dẫn sử dụng

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## English

### 1 Description

The VitalFlow caddy is nonsterile and reusable. It is intended to hold the VitalFlow motor drive (Model 58102), VitalFlow emergency handcrank (Model 58103), and Nautilus™ VF ECMO oxygenator with Balance™ Biosurface (Model 58110). Refer to the instructions for use for these devices for indications, contraindications, intended users, and patient target groups.

Verify that the packaging is not opened or damaged before use.

Contents:

- Caddy
- Hook and loop strap for cord management
- Motor drive 5/8 in stud wrench
- Motor drive hex wrench
- Motor drive mounting screw
- Emergency handcrank hex wrench
- Emergency handcrank bracket mounting screws (2)

#### 1.1 Specifications

Temperature limitation	-20 °C to 50 °C
Humidity limitation	15% to 95%

### 2 Intended purpose

The VitalFlow Caddy has no direct effect on the body. It is intended to hold the VitalFlow motor drive, VitalFlow emergency handcrank, and Nautilus VF ECMO oxygenator with Balance Biosurface.

### 3 Adverse effects


The following known adverse effects are associated with the use of this product: bruise, fracture, hypotension, ischemia, neurological dysfunction, and organ dysfunction.

If a serious incident related to the device occurs, immediately report the incident to the manufacturer and the applicable competent authority or regulatory body.

### 4 Attaching the motor drive


1. Use the motor drive 5/8 in stud wrench to remove the motor drive from the motor drive stem.

**Figure 1.**

2. Slide the motor drive into the back of the caddy with the connector cord facing out. Align the hole in the caddy with the hole on the motor drive.
3. Insert the motor drive mounting screw into the hole and tighten the screw with the motor drive hex wrench.


Figure 2.

## 5 Attaching the emergency handcrank


1. Remove the emergency handcrank caddy bracket from the caddy by pulling the pin out and lifting the bracket up.
2. Use the emergency handcrank hex wrench to unscrew and remove the pole mount adapter from the emergency handcrank.

Figure 3.


3. Use the emergency handcrank hex wrench and emergency handcrank bracket mounting screws to attach the caddy bracket onto the emergency handcrank. The angled corner on the caddy bracket should face the pump end of the handcrank.

Figure 4.

4. Place the emergency handcrank onto the caddy by sliding the bracket into the holder. Push the pin attached to the caddy into the hole to secure.

Figure 5.

## 6 Attaching the oxygenator

Place the oxygenator holder connector onto the caddy with the blood inlet facing in toward the motor drive.


## 7 Transporting the caddy

1. Inspect the caddy for damage and wear. Ensure that all devices are securely attached to the caddy before transport and use. Do not use the caddy if it is damaged or if the devices cannot be secured.
2. Carefully place the caddy on a cart or another flat and stable surface.

**Warning:** Use care when moving the caddy to avoid damage to the system or detaching of components.

**Warning:** Do not kink the tubing when moving the caddy. Kinked tubing could result in loss of blood flow to the patient.

Figure 6.

## 8 Cleaning the caddy

Clean the surfaces of the caddy with a suitable disinfectant wipe or damp cloth per hospital protocol. The following disinfectants may be used on the caddy:

- Commonly used hospital wipes
- Alcohol (ethanol 96% and isopropyl 90%)
- Hydrogen peroxide (7.5%)
- Bleach (1:10 bleach to water)

Immediately remove any contamination from the caddy by using a sponge or soft cloth and water.

## 9 Reuse limitations

- The material and construction of these devices is suited to handle multiple cleaning cycles without impacting safe use or functionality.
- End of life is determined by excessive wear or damage from normal use.
- Complete an examination to determine if the device is at the end of its useful life.
- Devices must always be examined before and after cleaning and before use.
- Examination should be thorough and must include a visual inspection of working surfaces. Never use devices that are damaged or nonfunctional.
- Ensure that the following are true:
  - No cracks, deformation, corrosion, or distortion are present on any part of the device.
  - All device markings are legible, including the device name and size.
  - All parts are present and free of damage and deterioration.
- If the device passes the visual inspection as described in these Instructions for Use, the device can be expected to perform as intended.

## 10 Disposal

**Warning:** Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

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© 2025
**Medtronic, Inc.**
Medtronic B.V.