VitalFlow™

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
1 Intended purpose

The intended use for the VitalFlow Centrifugal Pump is to provide cardiac and/or respiratory support in an ECLS circuit by moving blood via centrifugal force.

1.1 Indications for use

The VitalFlow Centrifugal Pump with Balance Biosurface¹ is intended to pump blood through the extracorporeal circuit for circulatory support up to 14 days, in adult patients with acute respiratory failure or acute cardiopulmonary failure, where other available treatment options have failed, and continued clinical deterioration is expected or the risk of death is imminent.

The VitalFlow Centrifugal Pump is driven by the VitalFlow™ Console and VitalFlow™ Motor Drive or the Emergency Hand Crank.

1.2 Contraindications

The VitalFlow Centrifugal Pump should not be used in destination therapy. This device, used for any other purposes than the intended use, is the responsibility of the user.

1.3 Intended users

The device is to be operated by trained medical clinicians that are familiar with ECMO procedures. Potential roles of the medical staff include physicians, perfusionists, nurses, respiratory therapists, and other trained ECMO professionals.

1.4 Expected clinical benefits

The device functions as part of an extracorporeal circuit during ECLS procedures. The device has no direct clinical benefit beyond the support of ECLS procedures.

2 Device description

The VitalFlow Centrifugal Pump with Balance Biosurface is intended for use in extracorporeal circulation procedures. It is designed to move blood by centrifugal force generated by a combination of a smooth rotating cone and low-profile impeller fins. Energy is transferred from the pump in the form of pressure and velocity as the blood is driven toward the outlet port of the pump. To limit friction and heat generation, the VitalFlow Centrifugal Pump utilizes a pivot bearing design on a dual ceramic pivot.

The pump's blood contacting surface is coated with Balance Biosurface which is a non-heparin biocompatible coating that reduces platelet adhesion. The VitalFlow Centrifugal Pump couples to a remote magnetic drive unit called the VitalFlow Motor Drive, which interfaces with a pump speed controller. It also couples with the VitalFlow Emergency Hand Crank in the event of controller or power failure. The VitalFlow Centrifugal Pump has been sterilized using ethylene oxide gas.

For additional information on operating the pump speed controller, refer to the VitalFlow Console operator's manual.

3 General warnings and precautions

Warning: Read all warnings, cautions, and instructions for use carefully prior to use. Failure to read and follow instructions or failure to observe the stated warnings could cause serious injury or death to the patient.

Warning: The VitalFlow Centrifugal Pump is qualified to be used only with the VitalFlow Motor Drive, and the VitalFlow Emergency Hand Crank in case of an emergency. Do not use the pump with a different motor drive; this may result in reduced performance, decoupling, or damage.

Warning: Ensure that the centrifugal blood pump and circuit have been debubbled and primed properly prior to initiating blood flow to minimize the risk of air reaching the patient. You must visually confirm that all air has been removed from the pump head.

Warning: Never clamp or restrict flow proximal to the pump; doing so may result in failure. If flow is occluded at the inlet, the pump may run at high rpms and possibly generate a high negative pressure.

Warning: Do not operate the centrifugal blood pump unprimed or without fluid; damage to the internal components will occur.

Warning: Massive air entry into the centrifugal blood pump deprimes the pump and stops blood flow. Stop the pump and remove air prior to resuming circulation.

Warning: Do not operate the centrifugal blood pump at rotational speeds greater than the specified maximum (3800 rpm). This may result in excessive wear to the pivot bearings.

Warning: Do not expose the centrifugal blood pump to liquid chemical agents, as they may affect the integrity of this device. Adhesive removers, alcohol, and anesthetic liquids, (such as isoflurane) and corrosive solvents (such as acetone) are known to degrade polycarbonate plastics. Avoid contact of these solutions with the centrifugal blood pump.

Warning: To prevent retrograde flow, the tubing at the outlet of the pump must always be clamped distal to the pump when the centrifugal blood pump is at rest.

Warning: To prevent retrograde blood flow when the centrifugal blood pump outlet tubing is unclamped, establish and maintain a minimum pump speed that overcomes line and patient resistance. If positive flow cannot be maintained, the pump outlet tubing must be clamped to isolate flow from the patient.

Warning: Always have a spare VitalFlow Centrifugal Pump and backup equipment available with appropriate protocol for change out. If the centrifugal blood pump needs replacement during circulatory support, the potential risk to the patient should be evaluated prior to changing the pump.

Warning: Frequent patient and device monitoring is recommended; do not leave the device unattended while in use. Monitor flow rate carefully for signs of malfunction throughout the circuit.

Warning: Do not operate the centrifugal blood pump in the absence of blood flow. The blood temperature within the pump could rise, and increased cellular damage may result.

Warning: The centrifugal blood pump must be handled using aseptic technique.

¹ Technology licensed under agreement from BioInteractions, Ltd., UK

Warning: The VitalFlow Centrifugal Pump is MR Unsafe. Keep it outside of the MRI scanner room.

Caution: This device should be used only by persons thoroughly trained in extracorporeal circulation procedures.

Caution: Do not forcefully hit or drop the centrifugal blood pump. Shock may cause damage to the device, which may cause device malfunction.

Caution: The centrifugal blood pump is sterile and nonpyrogenic. Inspect device and package carefully prior to use. Do not use if the unit package or the product has been damaged or soiled, or if caps are not in place.

Caution: This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Caution: Attach tubing in a manner that prevents kinks or restrictions that may alter flow.

Caution: Do not use excessive force to install tubing on pump as damage to pump might occur.

Caution: Take care to prevent damage to connectors when repositioning or attaching tubing.

Caution: Do not place the centrifugal blood pump near items adversely affected by magnetic fields.

Caution: A strict anticoagulation protocol should be followed, and anticoagulation should be routinely monitored during all procedures. The prescribing physician must weigh the benefits of extracorporeal support against the risk of systemic anticoagulation. Adequate heparinization should be maintained per institution ECLS protocol.

Caution: Refer to package labeling for storage temperature requirements.

Caution: It is the responsibility of the user to dispose of the devices in accordance with local regulations and hospital procedures.

4 Instructions for products with Balance Biosurface

The blood-contacting surfaces of the device are coated with Balance Biosurface to reduce platelet adhesion.

Warning: A product coated with Balance Biosurface is intended for single use only.

Caution: Appropriate blood anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.

Caution: Do not store product coated with Balance Biosurface above 50 °C (122 °F).

5 Storage environment

Storage temperature	-30 °C to 50 °C (-22 °F to 122 °F)
Storage humidity	20% to 85% noncondensing

6 Specifications

Priming volume	(approximate) 40 mL
Inlet/Outlet ID	9.5 mm (3/8 in)
Maximum operating pressure	760 mmHg (101 kPa)
Maximum flow rate	7 L/min
Maximum outlet pressure	700 mmHg (93.3 kPa)
Maximum pump speed	3800 rpm

7 Instructions for use

Verify the packaging is not opened or damaged and the use-by date has not expired prior to use.

7.1 Connection to controller

Connect the VitalFlow Motor Drive cable to the appropriate controller. For complete setup information, refer to the appropriate controller operator's manual.

7.2 Installing and removing the VitalFlow centrifugal pump

7.2.1 Installing the VitalFlow centrifugal pump

1. Slide the VitalFlow centrifugal pump into the VitalFlow motor drive (*Figure 2*) until the thumb latch pin engages.

Caution: To ensure optimal air handling and pump longevity, mount the pump with the outlet port angled as shown in *Figure 2*.

Note: Whenever possible, install the pump at or below the level of the patient's heart.

2. Prime the centrifugal blood pump according to the procedures described below.

Caution: Operating the centrifugal blood pump without being primed will cause damage to the pump.

Caution: Installing the centrifugal blood pump with the motor revolving may cause damage to the pump.

Warning: Use of the VitalFlow centrifugal pump beyond the labeled recommendations may result in failure of the centrifugal blood pump, reduced pumping capacity, leaks, excessive blood trauma, or degradation or corrosion of blood contact materials that may pass through the blood to the patient.

7.2.1.1 Removing the VitalFlow centrifugal pump

Set the pump speed to zero RPM prior to removing the pump.

Remove and discard the pump after each procedure. Discard the pump according to local and hospital protocols.

Push the pump motor thumb latch to release the centrifugal blood pump (*Figure 2*).

Slide the centrifugal blood pump out of the pump motor receptacle (*Figure 2*).

7.3 Recommended procedures for priming the VitalFlow centrifugal pump

1. Using accepted aseptic technique to maintain a sterile fluid pathway, attach the appropriate tubing to the port(s) of the centrifugal blood pump.
2. Using gravity, fill the centrifugal blood pump with priming solution to a point beyond the pump outlet port and clamp the tubing.
Caution: Do not hit or strike the centrifugal blood pump with instruments. Shock may cause damage to the device, which may cause device malfunction.
Warning: During extracorporeal life support procedures ensure that the centrifugal blood pump and circuit have been de-aired and primed properly prior to initiating support to minimize the risk of air reaching the patient.
3. Turn on the pump console and increase RPMs using the RPM knob for further priming. Circulate priming fluid to remove residual air.
4. Observe the pump for leakage or other anomalies.
5. Perform functional checks to the flow/bubble sensor and zero the flow probe in accordance to the applicable operator's manual.
6. Verify a leak-free, air-free circuit with correct flow path; patient drainage line connects to pump inlet, return line connects to pump outlet. Check that all tubing connections are secure.
Note: Check circuit to ensure no kinks or occlusions are present.
Warning: If leaks or other anomalies are found, remove the centrifugal blood pump and replace with a new, sterile pump; repeat steps 1 to 6 to prime.
7. To initiate circulation, turn pump RPM to minimum forward flow, release all clamps from circuit and increase pump RPM to achieve clinical target flow rate.
Warning: Massive air entry into the centrifugal blood pump may cause the pump to deprime and blood flow to stop. Stop the pump and remove air prior to resuming circulation.
Warning: Do not operate the centrifugal blood pump with its inlet tubing clamped, as negative pressure would be generated in the pump and air bubbles may be formed in the blood.

7.4 VitalFlow centrifugal pump checklist

Follow the protocol specified by your clinical program. Information from this blood pump checklist can be used to revise, upgrade, or expand existing checklists appropriately.

1. **Equipment assembly**
 - Mount VitalFlow motor drive so that the pump is mounted as shown in *Figure 2*.
 - Check that all electrical connections are secure.
 - Test control module power and display.
 - Check date and integrity of sterile blood pump package.
 - Check that flow sensor is sized properly.
 - Assemble extracorporeal circuit in a sterile manner.
 - Allow sufficient tubing length for standby pumping unit.
 - Connect flow sensor to circuit in correct location and flow direction.
2. **Prime pump and ECLS circuit**
 - Check circuit to ensure no kinks or occlusions are present.
 - Orient the pump properly with outlet port at highest point and parallel to the ground.
 - Prime and de-air pump and extracorporeal circuit.
 - Check pump for leaks, irregular motion, and noise.
 - Check circuit for visible air.
 - Check that all tubing connections are secure.
 - Check for adequate levels of anticoagulation prior to extracorporeal circulation.
3. **Operating parameters**

Perform the following, per the appropriate VitalFlow console operator's manual:

 - Calibrate sensors according to the console's instructions.
 - Set and verify alarm parameters.
4. **Emergency backup equipment**

The following backup equipment is available:

 - Backup power
 - VitalFlow emergency hand crank
 - Spare VitalFlow centrifugal pump
 - VitalFlow console and accessories
5. **During use**
 - Achieve minimum pump flow prior to unclamping lines.
 - Monitor the console for messages and alarms.
 - Monitor ECLS circuit for visible air and kinks in the lines.
 - Check for adequate levels of anticoagulation during extracorporeal circulation.

8 Emergency procedures

Pump failure

An unusual grinding or chattering sound can indicate pump failure is occurring. Replace the VitalFlow centrifugal pump according to the change out protocol.

Mechanical failure

An unusual or high-pitched sound indicates that drive failure is occurring. Evaluate the VitalFlow motor drive and proceed to change out protocol if necessary.

Power failure

If the AC/DC and battery power fails, or the VitalFlow motor drive or controller stops functioning, connect the VitalFlow centrifugal pump to the VitalFlow emergency hand crank. The VitalFlow console operator's manual describes hardware emergency procedures.

9 Pump change out

The blood pump is indicated for 14 days of use and should be changed out at or before 14 days, if there are clinical signs or device performance indications that would require a pump change out. Examples of device performance related indications that a change out is needed may include excessive pump noise or severe hemolysis.

The blood pump should be evaluated for change out if the clinician determines during the extracorporeal procedure that adequate patient support is not being achieved. The clinician should weigh the risks and benefits of potential pump replacement. If a change out is deemed necessary, obtain a replacement device before discontinuing circulation through the original device.

Note: Have materials needed for the change available before starting the change, including clamps and sterile scissors. Always use aseptic technique during change out.

Note: Allow sufficient tubing length for hand crank accessibility and changeout.

1. Reduce pump speed to minimum forward flow and clamp the pump outlet and pump inlet lines.
2. Double clamp both the pump outlet tubing followed by the pump inlet tubing. Clamp far enough away from the pump outlet/inlet to allow adequate room for reconnection to the replacement device.
3. Decrease RPM down to zero to stop the pump.
4. Remove the pump from the motor drive.
5. Using aseptic technique, cut the pump inlet and outlet tubing between the clamps connected to the original blood pump.
6. Connect new device using air-free technique. Connect the inlet tubing to the pump inlet and the pump outlet tubing to the pump outlet of the replacement blood pump. Confirm the tubing and pump are free of air.
7. Secure pump inlet and outlet tubing connections with tie bands per institution protocol.
8. Insert new pump back into motor drive and secure.
9. To initiate circulation, turn pump RPM to minimum forward flow, release all clamps from circuit and increase pump RPM to achieve clinical target flow rate.

10 Potential complications

Potential complications, including those normally associated with ECLS procedures and anticoagulation, may necessitate physician intervention. Some potential complications related to the use of ECLS pumps include but are not limited to air embolism or entrainment, thrombus formation leading to deep-vein thrombosis or pulmonary embolism, hemolysis, thrombocytopenia, mechanical failure of pump or connections, infection, hypothermia or hyperthermia due to inadequate heat transfer, hemorrhage related to bleeding, and inadequate blood flow related to clotting or hypovolemia resulting in inadequate gas exchange. There are risks and adverse events related to all ECLS procedures and anticoagulation including heart, vessel or lung damage, hypoxia or hypercarbia due to inadequate gas exchange, anemia, infection, hemorrhage, liver or kidney failure, stroke, or death.

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

11 ECLS general information

This product is one component of an ECLS circuit. A typical ECLS circuit includes tubing, connectors, a gas exchange device and a blood pump. Configuration of the circuit and selection of individual component types may vary according to each specific use and treatment center. Each component of the circuit is commercially available. The VitalFlow centrifugal pump is compatible with circuits using 1 cm (3/8 in) tubing in the range of flow and pressure (see *Chapter 6*) and used with anticoagulation typically used during ECLS procedures.

The extracorporeal circuit is typically managed with a systemic anticoagulant (usually heparin). The anticoagulation is managed to titrate either the active whole blood clotting time or activated partial thromboplastin time to a level above the normal ranges (typically 1-2x normal).

The principal pathway of the ECLS circuit consists of blood being drained from the drainage lumen of the catheter, passing through a blood pump which disperses the blood through a gas exchange device and back into the patient through the catheter's reinfusion lumen.

The systemic temperature may be controlled by a heat exchanger integrated into the gas exchange device. The components of the circuit are selected to meet the blood flow and gas exchange needs of each patient based on their metabolic requirements. The circuit has been widely described in the literature and is summarized in each iteration of the textbook published by the Extracorporeal Life Support Organization. Maintenance of the ECLS circuit is managed by personnel who have been trained on the intricacies of extracorporeal support including physicians, perfusionists, nurses, respiratory therapists and other medical specialists. These individuals are qualified to manage the circuit based on training and continuing education per specific institutional guidelines.

12 Disclaimer of warranty

MC3, INC. MANUFACTURES ITS PRODUCTS (REFERRED TO AS THE "PRODUCTS") UNDER CONDITIONS THAT MEET APPLICABLE REGULATORY REQUIREMENTS, BUT HAS NO CONTROL OVER HOW THE PRODUCTS ARE ACTUALLY USED OR THE CONDITIONS UNDER WHICH THEY ARE USED. AS A RESULT, MC3, INC. DISCLAIMS ALL PRODUCT WARRANTIES INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MC3, INC. SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE.

Limitation on Liability

IN NO EVENT WILL MC3, INC.'S TOTAL LIABILITY TO PURCHASER OR ANY THIRD PARTY ARISING FROM THE SALE OF PRODUCTS EXCEED THE PURCHASE PRICE OF THE PRODUCTS SUBJECT TO THE CLAIM REGARDLESS OF THE THEORY OF LIABILITY OF THE CLAIM.

Limitation on Damages

MC3, INC. IS NOT RESPONSIBLE TO PURCHASER OR ANY THIRD PARTY FOR ANY DIRECT, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY CLAIMED BREACH OF WARRANTY, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT.

Applicable law in your jurisdiction may limit the ability of a manufacturer to disclaim some warranties or exclude or limit damages. While the above limitations are an integral part of the business relationship between MC3, Inc. and the purchasers of its Products, nothing in the above is intended to contravene these provisions of applicable law. If any part or term of this Disclaimer of Warranty or the Limitations on Liability or Damages is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of all remaining provisions shall not be affected and shall be construed and enforced as if the Disclaimer of Warranty or the Limitations on Liability or Damages did not contain the particular part or term held to be invalid.

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© 2025 Medtronic
Medtronic and Medtronic logo are trademarks of Medtronic. TM* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. MC3 Cardiopulmonary is a Medtronic company.
Medtronic, Inc.
Medtronic B.V.