Nautilus™ VF

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English

1 Indication for use

The Nautilus VF ECMO Oxygenator with integrated heat exchanger is intended to provide assisted extracorporeal circulation and physiologic gas exchange (oxygenation and carbon dioxide removal) of the patient's blood for up to 14 days in adult and pediatric adolescent patients with acute respiratory failure or acute cardiopulmonary failure, where other available treatment options have failed, and continued clinical deterioration is expected or the risk of death is imminent. The integrated heat exchanger is intended to heat or cool the blood as needed during use. Integrated fluid path pressure, temperature, and oxygen saturation monitoring is achieved by built-in sensor modules.

2 Contraindications

This device is not designed, sold, or intended for use except as indicated.

3 Device description

The Nautilus VF ECMO Oxygenator, a diffusion polymethylpentene (PMP) fiber membrane oxygenator with Balance™ Biosurface¹ (REF 58110), is used in extracorporeal life support procedures to oxygenate blood, remove carbon dioxide and regulate blood temperature. Blood enters the device and passes through both the heat exchange membrane, where temperature is adjusted, and the gas transfer membrane, where oxygen is added and carbon dioxide is removed.

The device contains integrated sensors. The following measured parameters can be transmitted to the VitalFlow ECMO Console through the Nautilus VF Oxygenator Cable: inlet pressure, inlet oxygen saturation, outlet pressure, outlet oxygen saturation, and outlet temperature. The difference between the inlet and outlet pressure, delta pressure, is calculated and also transmitted.

The device is bonded on its primary blood contacting surfaces with a nonleaching biocompatible surface to reduce platelet adhesion.

The device is single-use, nontoxic, nonpyrogenic, and not made from natural latex rubber materials.

3.1 Intended use

The intended use for the Nautilus VF ECMO Oxygenator is to provide cardiac and/or respiratory oxygen support in an ECMO circuit.

3.2 Intended users

The device is to be operated by trained medical clinicians that are familiar with ECMO procedures. Potential roles of the medical staff include physicians, perfusionists, nurses, respiratory therapists, and other trained ECMO professionals.

4 Accessories sold separately

To ensure a proper connection between the device and holder, the Nautilus VF ECMO Oxygenator is to be used with a Nautilus oxygenator holder. These holders are sold separately. This device should only be used with the Nautilus VF Oxygenator Cable and VitalFlow ECMO Console.

5 General warnings and precautions

Warning: Only use this device in accordance with these instructions for use.

Warning: Read all warnings, precautions, and instructions for use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.

Warning: Only clinicians thoroughly trained in extracorporeal life support procedures should use this device.

Warning: Use of this device for other patient populations is the responsibility of the clinician. MC3 has not yet demonstrated clinical performance of this device in the young pediatric population.

Warning: No modification of this equipment is allowed. This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device which could result in patient injury, illness, or death.

Warning: The polymethylpentene (PMP) oxygenator fiber is not suitable for use with volatile anesthetics (anesthetic gases).

Warning: The use of all possible disinfectants in heater/coolers while the heat exchanger is in use has not been fully evaluated. H₂O₂ at low concentrations (≤ 300 ppm) is acceptable for use with PET heat exchangers, but higher concentrations should be avoided. Other disinfectant use has not been evaluated.

Warning: This device should not be operated outside the recommended parameter ranges.

Warning: The following patient parameters must be monitored by an external system: core body temperature, arterial pressure, venous pressure, and blood oxygen saturation.

Warning: Gas transfer rates may change over time and adjustments of FiO₂ and gas flow rates may be needed to achieve desired gas transfer performance.

Warning: Ensure supply of power and oxygen are sufficient.

Warning: Avoid pulling or putting strain on the tubing or power/data cable.

Warning: Only use a medical gas supply with dry air and oxygen.

Warning: The use of PVC pump segments with a roller pump causes an electrostatic charge to build up on the blood side. Upon discharge, the integrity of the heat exchanger may be compromised.

Warning: Tubing should be attached in such a manner as to prevent kinks or restrictions that may alter the blood, water, or gas flow.

Warning: If air leaks are observed during priming and/or operation, this may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the device if these conditions are observed.

Warning: Positive pressure must be maintained on the blood side at all times during perfusion.

Warning: The blood path pressure should be greater than the gas path pressure at all times.

Warning: Do not obstruct gas outlet port and vents.

Warning: Use proper technique when using blood access ports.

Warning: Ensure that adequate anticoagulation procedures are employed during use.

Warning: Once flow to the patient is initiated, blood flow should circulate through the oxygenator at all times within the recommended blood flow range, except during emergency oxygenator change out.

Warning: Circuit access is the responsibility of the attending physician.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

¹ Technology licensed under agreement from BioInteractions, Ltd., UK

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the Nautilus VF ECMO Oxygenator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Caution: Do not allow anesthetic liquids (such as isoflurane) or corrosive solvents (such as acetone) to come in contact with the device or accessories as they may jeopardize structural integrity.

Caution: The oxygen saturation measurements will be affected by Methylene blue, indocyanine green, Evans blue, fetal hemoglobin, and unusual levels of bilirubin, lipids, beta-carotene, carboxyhemoglobin, methemoglobin, and sulfhemoglobin.

Caution: Gas transfer rates decrease as a result of lower atmospheric pressures.

Caution: The device is not recommended for use with a pulsatile pump.

Caution: Exposing the oxygenator to extreme cold may compromise its mechanical integrity.

Caution: A replacement oxygenator should be readily available during ECLS.

Caution: Avoid mechanical impacts to the device while carrying or during operation.

6 Balance Biosurface

The blood-contacting surfaces of the device are coated with Balance Biosurface to reduce platelet adhesion.

7 Packaging and storage

The fluid pathway is delivered sterile and nonpyrogenic. Sterilized by ethylene oxide. Verify the packaging is not opened or damaged and the use-by date has not expired prior to use.

Warning: Do not use device if expired.

Warning: Do not use device if package is opened or damaged, or if the device is damaged as the sterility of the device may have been compromised and/or performance may be affected.

8 Instructions for use

8.1 Set up of holder

Remove all accessories from packaging, and inspect for damage prior to use. Attach the holder clamp to the ECLS cart at the desired height, and tighten the clamp knob. Adjust the remaining holder arms as desired, and tighten the holder arm knobs to secure in place.

8.2 Set up of device

Warning: Aseptic technique is to be used during all stages of setup and use of the device.

Warning: Setup and use of the device is the responsibility of the attending clinician.

Carefully remove the device from the sterile packaging. Slide the oxygenator holder connector onto the top of the oxygenator arm. A click indicates the device is correctly in place.

Warning: Sterility is not ensured after device packaging is open.

Warning: Ensure the device is positioned lower than the patient or at the same level, but no higher.

Caution: The device must be securely mounted in an upright position on its intended holder.

Note: To remove the device, press up on the holder release tab and pull the device away from the mounting arm.

8.3 Connecting the device

Caution: Only use a Nautilus VF Oxygenator Cable.

Caution: Only connect oxygenator to the VitalFlow ECMO Console.

Plug one end of the power/data cable into the Nautilus VF ECMO Oxygenator and the other end into the VitalFlow ECMO Console.

Note: The Nautilus VF ECMO Oxygenator comes with pre-zeroed pressure sensors. No further action is needed to ensure pressure sensors are zeroed prior to use.

Note: The device gas transfer and heat exchange functions DO NOT require power. The device functions without connection of the Nautilus VF Oxygenator Cable in the event that user does not wish to utilize the sensing capabilities of this device.

Note: The device is powered by the VitalFlow ECMO Console via the Nautilus VF Oxygenator Cable only.

Note: If the device is unplugged the sensing capabilities will terminate.

8.4 Priming the device

Caution: Tubes and luer lock connections must be fastened by hand only. Using tools or other aids may damage the ports.

Warning: Ensure that unused ports remain capped to prevent contamination and that the caps are tightened to avoid leakage.

See section 8.7 Extracorporeal Circuit Component Compatibility

Connect the water lines to the water ports on the oxygenator. Run water through the heat exchanger, and check for leaks from the water compartment to the blood compartment.

Warning: Do not use the oxygenator if water is present in the blood compartment.

Connect all blood and gas lines to their appropriate connections under aseptic conditions per institution protocol.

Note: Banding of all blood tubing connectors is recommended.

Fill the priming bag with the appropriate volume of priming solution per institution protocol. Remove the yellow luer lock cap on the de-airing membrane located on the blood inlet side of the device. Gravity prime the oxygenator or prime by forward pump flow per institution protocol, avoiding fluid exposure to the sensor area.

Warning: All gas emboli must be cleared from the extracorporeal circuit before initiating flow to the patient. Gas emboli are hazardous to the patient.

Add additional priming solution as required to prime the rest of the circuit. Upon completion of priming, re-cap the de-airing membrane. Ensure the entire system, including ports, is fully de-aired prior to beginning extracorporeal circulation.

Note: Device integrity and performance has been verified by MC3 after 30 days of primed storage with phosphate buffered saline solution.

8.5 Extracorporeal circulation

Caution: Check all lines are properly connected prior to use.

Caution: The de-airing membrane must be closed during operation.

Ensure that there is no air in the circuit. Check for adequate levels of anticoagulation prior to and during extracorporeal circulation. Remove arterial and venous clamps and gradually increase blood flow. Start the gas flow with the appropriate gas settings per institution protocol. Adjust water temperature per institution protocol.

Caution: Do not exceed 42 °C degrees in the heat exchanger.

8.6 Activation of measured parameters

Once the Nautilus VF Oxygenator Cable has been connected and the VitalFlow ECMO Console is powered on, the sensing capabilities of the Nautilus VF ECMO Oxygenator are active.

8.7 Extracorporeal circuit component compatibility

The Nautilus VF ECMO Oxygenator is to be used in conjunction with other components to make up the extracorporeal circuit. The components used with the Nautilus device are not specified, rather it is left to the user to select devices that are compatible. The sensing capabilities of the Nautilus VF ECMO Oxygenator are only intended to be used in conjunction with the Nautilus VF Oxygenator Cable and the VitalFlow ECMO Console.

Tubing circuit: The device blood ports are designed for use with 1 cm (3/8 in) ID tubing.

Pump: The device is designed for use with a constant-flow blood pump (e.g., centrifugal pump) with flows greater than or equal to 500 mL/min but less than or equal to 7 L/min.

For reference: Pump must provide sufficient pressure to overcome the resistance to flow through all circuit components. This pressure will depend on patient conditions and specific circuit components selected. The user should reference the pressure vs. flow plot provided in this IFU, as well as information provided with the selected catheter and pump to determine whether the pump output is adequate.

Catheters: Access to major veins and/or arteries is needed to facilitate extracorporeal circulation. The device is compatible with all catheters that allow at least 500 mL/min flow.

Heater/cooler: The device water ports are designed for use with 1.3 cm (1/2 in) quick-connects.

Gas supply: The device inlet gas port is designed for use with 0.6 cm (1/4 in) ID tubing. Other equipment should be used, as needed, to regulate the gas pressure and gas flow.

Other components: The device is compatible with other common extracorporeal circuit components, including hemofilters, sampling/monitoring ports, and tubing configurations. The device also provides standard luer connectors, which are compatible with standard sampling, monitoring, and infusion devices.

General ELSO guidelines: For further information on extracorporeal circuit component selection and use, refer to ELSO guidelines.

8.8 Circuit change out general information

The ECMO circuit, including the oxygenator, pump, catheter, and tubing set, should be monitored throughout extracorporeal circulation. The ECMO circuit components should be evaluated for change out if the clinician determines during the extracorporeal procedure that adequate patient support is not being achieved.

If a change out is deemed necessary, obtain a replacement device before discontinuing circulation through the original device and/or circuit.

Reasons for change outs may include the following:

- General: observed leaks, cracks, air entrainment, patient infection, severe or increasing hemolysis trends, or visible clotting within the circuit
- Catheter: inadequate flow
- Pump: excessive pump noise or reduction of flow performance
- Oxygenator: reduced gas exchange or heat exchange performance that no longer provides adequate patient support or an increase of pressure gradient across the device.

8.9 Change out

The oxygenator should be evaluated for change out if the clinician determines during the extracorporeal procedure that adequate patient support is not being achieved. If a change out is deemed necessary, obtain a primed replacement device (see section 8.4) before discontinuing circulation through the original device. Have materials needed for the change available before starting the change, including clamps and sterile scissors. Always use aseptic technique during change out.

1. Double clamp the tubing on the inlet side and outlet side of the oxygenator that is to be replaced and stop the pump. Clamp the tubing far enough away from the oxygenator blood inlet/ outlet to allow adequate tubing length for reconnection to the replacement device.
2. Detach the oxygenator from the holder by pressing the release tab.
3. Attach the primed replacement oxygenator to the holder.
4. Cut the venous and arterial tubing between the clamps connected to the original oxygenator.
5. Using air-free technique, connect the arterial line to the blood outlet and the venous line to the blood inlet of the primed replacement oxygenator.
6. Open the clamp on the blood inlet side and turn the pump on at a low speed.
7. Connect the gas line to the gas inlet of the replacement device.

8. Confirm that the system is free of bubbles before removing the outlet tubing clamp, and gradually increase the pump flow to the necessary value.
9. Connect the water lines to the replacement oxygenator.

9 Specifications

Oxygenation Fiber Type	Polymethylpentene (PMP)		
Fiber Surface Area – Oxygen Membrane	1.8 m ²		
Heat Exchanger	Polyethylene terephthalate (PET)		
Fiber Surface Area – Heat Exchange	0.3 m ²		
Prime Volume	226 mL		
Residual Blood Volume	189 mL		
Operating Blood Volume	226 mL		
Blood Flow Range	0.5 – 7 L/min		
Maximum Rated Blood Pressure	750 mmHg		
Gas Flow Rate (Gas:Blood)	0.5:1 – 3:1		
Maximum Gas Path Pressure	100 mmHg (1.9 psi)		
Maximum Water Path Pressure	1125 mmHg (21.8 psi)		
Power			
Rated input power Note: ECMO Oxygenator to be powered by Nautilus VF Oxygenator Cable	6 V 500 mA		
Sensor Accuracy			
Parameter	Measuring Range	Resolution	Measuring Accuracy
Pressure: P _{in} , P _{out}	0 – 750 mmHg	1	0 – 249 mmHg: ± 10 mmHg 250 – 750 mmHg: ± 7% of the measured value Maximum ± 15 mmHg drift over 30 days
Delta Pressure	0 – 750 mmHg	1	Calculated value
Oxygen Saturation: SO ₂ _{in} , SO ₂ _{out}	40 – 100%	0.1	± 5% (absolute)
Temperature: T _{out}	20 – 42 °C	0.1	± 1 °C
Environmental			
	Operating	Storage	Transportation
Temperature	15 – 30 °C	10 – 30 °C	-30 – 50 °C
Relative Humidity (non-condensing)	30 – 75%	5 – 95%	5 – 95%

10 Disposal

Warning: Dispose of the device in accordance with local hospital, administrative, and/or other government policies.

The device is a single-use only disposable. Unplug the power/data cable prior to disposal. Leave the power/data cable with the console, and dispose of the device according to institution protocol and local regulations.

11 Cleaning

If the outside of the oxygenator becomes soiled, wipe down surfaces using a cloth moistened with suitable cleaning solution. Do NOT spray cleaner directly on the device. The device is compatible with the following cleaner types: 3% hydrogen peroxide, bleach, isopropyl alcohol, aldehydes, and quaternary ammonium compounds.

Cleaning of the holder and power/data cable is recommended after each use or if contaminated. All surfaces of the holder can be cleaned and disinfected for blood, saline, or contaminate spills using normal medical equipment cleaners and disinfectants. Suitable cleaning solutions for sensitive medical devices should be used when cleaning the power/data cable.

Caution: Unplug the power/data cable prior to cleaning.

12 Potential complications

Potential complications, including those normally associated with ECLS procedures and anticoagulation, may necessitate physician intervention. Some potential complications related to the use of ECLS oxygenators include but are not limited to hypoxia or hypercarbia due to inadequate gas exchange, air embolism, thrombus formation leading to deep-vein thrombosis or pulmonary embolism, hemolysis, thrombocytopenia, mechanical failure of ports and connections, infection, hypothermia or hyperthermia due to inadequate heat transfer, hemorrhage related to bleeding, and inadequate blood flow related to clotting or hypovolemia resulting in inadequate gas exchange. There are risks and adverse events related to all ECLS procedures and anticoagulation including heart, vessel, or lung damage, hypoxia, anemia, infection, hemorrhage, liver or kidney failure, stroke, or death.

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

13 Technical information

There are no user serviceable parts or fuses.

Electric shock protection (IEC 60601-1)	Nautilus VF powered by VitalFlow Console. Console is Class II or Internally Powered.
---	--

Applied part classification (IEC 60601-1)	Integrated pressure sensors – Defibrillator-protected type CF
Degree of Ingress Protection (IEC 60529)	IPX2
Mode of Operation	Continuous

The following information and data are available upon request: sterilization method, materials of the blood pathway, blood pathway pressure drop, gas pathway pressure drop, blood cell damage, particle release, and relevant tolerances for specifications in the instructions for use.

The essential performance of the Nautilus VF ECMO Oxygenator is to oxygenate blood, remove carbon dioxide, and heat or cool blood. The essential performance of the device is independent of all electronics. It will not be affected by electromagnetic interference.

Guidance and Manufacturer's Declaration – Electromagnetic Compatibility

Warning: Common RF emitters including diathermy, lithotripsy, electrocautery, RFID, electromagnetic anti-theft systems, and metal detectors may interfere with the electronic functions of the Nautilus VF ECMO Oxygenator. Avoid exposure to these RF sources as the cause of interference may not be obvious.

The Nautilus VF ECMO Oxygenator is intended to be used in hospitals (except near active HF surgical equipment and the RF shielded room of magnetic resonance imaging, where the intensity of electromagnetic disturbances is high). The user of the device should ensure it is used in an appropriate environment.

⊗ The Nautilus VF ECMO Oxygenator is MR Unsafe.

Emission Test	Compliance	Electromagnetic environment guidelines
RF emissions (CISPR11)	Group 1	The Nautilus VF ECMO Oxygenator uses H/RF energy solely for its internal operation. Therefore, its RF emissions are very low and not likely to cause interference in nearby electrical equipment.
RF emissions (CISPR11)	Class B	The device complies to CISPR11 Class B emissions for increased assurance. However, the device is suitable for use in professional healthcare and industrial establishments not directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions (IEC 61000-3-2)	Class A	
Voltage fluctuation/flicker emissions (IEC 61000-3-3)	Complies	

The Nautilus VF ECMO Oxygenator is intended to be used in the electromagnetic environment specified below. The user should ensure that it is used in such an environment.

Immunity Test	IEC 60601-1-2 Test Level
Electrostatic Discharge (IEC 61000-4-2)	±8 kV contact
	±2, ±4, ±8, ±15 kV air
Electrical fast transients/bursts (IEC 61000-4-4)	±2 kV, 100 kHz repetition frequency for mains power
	±1 kV, 100 kHz repetition frequency for input/output
Surges (IEC 61000-4-5)	±0.5 kV, ±1 kV line to line
	±2 kV line to ground
Rated Power Frequency Magnetic Fields (IEC 61000-4-8)	30 A/m 50 or 60 Hz
Conducted disturbances induced by RF fields (IEC 61000-4-6)	3 V 0.15 MHz - 80 MHz 80% AM at 1 kHz
	6 V 0.15 MHz - 80 MHz within ISM and amateur radio bands 80% AM at 1 kHz
Radiated RF Electromagnetic Fields (IEC 61000-4-3)	3 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz
Proximity fields from RF Wireless communications Equipment (IEC 6100-4-3)	27 V/m 385 MHz 50% PM at 18 Hz
	28 V/m 450 MHz FM ± 5 kHz, 1 kHz sine
	9 V/m 710,745,780 MHz 50% PM at 217 Hz
	28 V/m 810,870,930 MHz 50% PM at 18 Hz
	28 V/m 1720, 1845, 1970 MHz 50% PM at 217 Hz
	28 V/m 2450 MHz PM 50%; 217 Hz
	9 V/m 5240, 5500, 5785 MHz PM 50%, 217 Hz

Voltage Dips (IEC 61000-4-11)	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
	0% UT; 1 cycle 70% UT; 25/30 cycles Single phase: at 0°
Voltage Interruptions (IEC 61000-4-11)	0% UT; 250/300 cycle
Note: UT is the AC supply voltage to be used for application of the test level.	

14 Benefit-risk assessment

Undesirable side effects must constitute an acceptable risk when weighed against the intended performance. The benefit of ECMO is life support for patients with acute respiratory failure or acute cardiopulmonary failure, where other available treatment options have failed, and continued clinical deterioration is expected or the risk of death is imminent. Therefore, it is concluded that any residual risks are acceptable because the benefit of the device outweighs each risk.

15 General information

ECLS General Information: This product is one component of an ECLS circuit. A typical ECLS circuit includes tubing, connectors, a gas exchange device and a blood pump.

Configuration of the circuit and selection of individual component types may vary according to each specific use and treatment center. Each component of the circuit is commercially available. The extracorporeal circuit is typically managed with a systemic anticoagulant (usually heparin). The anticoagulation is managed to titrate either the activated whole blood clotting time or activated partial thromboplastin time to a level above the normal ranges (typically 1-2 times normal). The principal pathway of the ECLS circuit consists of blood being drained from the drainage lumen of the catheter, passing through a blood pump which disperses the blood through a gas exchange device and back into the patient through the catheter's reinfusion lumen. The systemic temperature is controlled by a heat exchanger integrated into the gas exchange device. The components of the circuit are selected to meet the blood flow and gas exchange needs of each patient based on their metabolic requirements. The circuit has been widely described in the literature and is summarized in each iteration of the textbook published by the Extracorporeal Life Support Organization.² Maintenance of the ECLS circuit is managed by personnel who have been trained on the intricacies of extracorporeal support including physicians, perfusionists, nurses, respiratory therapists and other medical specialists. These individuals are qualified to manage the circuit based on training and continuing education per specific institutional guidelines.

16 Disclaimer of warranty

MC3, INC. MANUFACTURES ITS PRODUCTS (REFERRED TO AS THE "PRODUCTS") UNDER CONDITIONS THAT MEET APPLICABLE REGULATORY REQUIREMENTS, BUT HAS NO CONTROL OVER HOW THE PRODUCTS ARE ACTUALLY USED OR THE CONDITIONS UNDER WHICH THEY ARE USED. AS A RESULT, MC3, INC. DISCLAIMS ALL PRODUCT WARRANTIES INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MC3, INC. SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE.

Limitation on Liability

IN NO EVENT WILL MC3, INC.'S TOTAL LIABILITY TO PURCHASER OR ANY THIRD PARTY ARISING FROM THE SALE OF PRODUCTS EXCEED THE PURCHASE PRICE OF THE PRODUCTS SUBJECT TO THE CLAIM REGARDLESS OF THE THEORY OF LIABILITY OF THE CLAIM.

Limitation on Damages

MC3, INC. IS NOT RESPONSIBLE TO PURCHASER OR ANY THIRD PARTY FOR ANY DIRECT, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY CLAIMED BREACH OF WARRANTY, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT.

Applicable law in your jurisdiction may limit the ability of a manufacturer to disclaim some warranties or exclude or limit damages. While the above limitations are an integral part of the business relationship between MC3, Inc. and the purchasers of its Products, nothing in the above is intended to contravene these provisions of applicable law. If any part or term of this Disclaimer of Warranty or the Limitations on Liability or Damages is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of all remaining provisions shall not be affected and shall be construed and enforced as if the Disclaimer of Warranty or the Limitations on Liability or Damages did not contain the particular part or term held to be invalid.

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² Toomasian JM, Lawson DS, and Harris WH. ECLS Equipment. Chapter 8 in ECMO: Extracorporeal Cardiopulmonary Support in Critical Care, Fourth edition, Annich G, Bartlett RH, Lynch WR, MacLaren G, Wilson JM (eds). Extracorporeal Life Support Organization, Ann Arbor, 2012. pp 107-132.

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© A Nautilus VF ECMO oxigenátor nem MR-kompatibilis.
Medtronic and Medtronic logo are trademarks of Medtronic. TM* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. MC3 Cardiopulmonary is a Medtronic company.
MDSS GmbH
Medtronic B.V.
Medtronic, Inc.