
Bruksanvisning / Instructions for Use

Aliva® oral hydration system AV001

Bruksanvisning (NO)

4

Instructions for use (EN)

38

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Table of Contents

01. Preface	40	05. Maintenance	59
Support, Service, and Incident Reporting	40	Cleaning the Aliva pump	59
Retaining instructions	40	Cleaning the mouthpiece nozzle and tubing	59
Description of the user and intended use	40	Cleaning the Aliva bottle	59
Contraindications	41	Disinfection of the Aliva pump in healthcare settings	59
When to stop using the Aliva pump or the supplied accessories	41		
02. Explanation of safety warnings	42	06. Troubleshooting	60
Explanation of safety warnings used in the instructions for use:	42		
List of safety warnings	42	07. Disposal	62
03. Description of the product	44	Disposal of electronic components	62
Process overview – how it works	44	Disposal of batteries	62
Product illustrations and explanations	45	Disposal of packaging waste	62
		Disposal of the Aliva mouthpiece	63
04. Operation and how to use	52	How to dispose of the Aliva mouthpiece	63
How to unpack the Aliva oral hydration system for the first time	53	Disposal of the empty Aliva bottles	63
Setting up the Aliva pump	53		
Control panel LED lights	53	08. Appendixes	64
Installing the Aliva mouthpiece	54	Appendix 1 – Supplied accessories	64
How to install a new mouthpiece	54	Appendix 2 – Safety information symbols and explanations	64
Removing the Aliva mouthpiece	55	Appendix 3 – Technical data tables – Aliva pump and mouthpiece	66
Priming mode when installing a new Aliva mouthpiece	55	Appendix 4 – Electromagnetic Compatibility	67
Placing the Aliva mouthpiece nozzle in the mouth	56		
How to start and stop the flow of the water	56		
Battery level LED indications	57		
How to charge the internal battery	57		
How to operate the internal battery	58		
How to move and transport the Aliva system	58		
How to store the product safely	58		

01. Preface

This document is the English version of the Instructions for Use (IFU). The most recent version, including translations into official EU/EEA languages for countries where the product is marketed, is available at: www.aliva.no/EN/user-guide

These instructions apply to use within the European Economic Area (EEA), with national language versions provided for each country where the product is placed on the market.

Support, Service, and Incident Reporting

Any questions regarding the Aliva® oral hydration system may be directed to the legal manufacturer, Aliva AS. Please contact Aliva AS for support, technical assistance, or to report incidents.

Any incident related to the device or its use may be reported either to the manufacturer or to the competent authority of the relevant Member State.

Manufacturer:

Aliva AS
Rigedalen 41
4626 Kristiansand
Norway
Telephone: +47 488 48 838
E-mail: post@aliva.no
www.aliva.no

Retaining instructions

Read and understand the entire set of instructions before using this product. Failure to do so could result in serious injury or hazards to people or property.

Keep safety information and instructions for future reference and pass them on to subsequent users of the product.

The manufacturer is not liable for cases of material damage or personal injury caused by incorrect handling or non-compliance with the safety instructions. In such cases, the warranty will be voided.

Description of the user and intended use

This document is for adult (18 years and older) users of the Aliva oral hydration system, with oral dryness (xerostomia) who use the device at home, and for healthcare professionals who administer and support the use of the device.

The intended use of the Aliva oral hydration system is to continuously supply water to the oral cavity to alleviate oral dryness. The user adds clean, unthickened drinking water of choice (tap water or bottled).

The Aliva oral hydration system shall only be used in full accordance with the instructions described in this manual. Any use other than those described in this manual is considered as non-intended use. This will invalidate the warranty and liability.

Do not use the Aliva oral hydration system:

- In people under 18 years of age
- In people who have difficulty swallowing safely or who need thickened fluids
- With any liquids other than clean drinking water
- For longer than 12 hours in a single day
- If any part of the device has been changed or damaged

Contraindications

The Aliva oral hydration system must not be used by individuals who are unable to use the product safely in accordance with these instructions. Improper use could lead to serious health risks, such as aspiration, choking, or infection.

Do not use the Aliva oral hydration system in cases of:

- Impaired swallowing function with increased risk of aspiration, for instance due to neurological disorders like Parkinson's disease, or Amyotrophic lateral sclerosis (ALS).
- Swallowing difficulties (dysphagia).
- Prescribed fluid thickening agents due to impaired swallowing function.
- Cognitive impairment that prevents safe and independent use of the device, such as dementia.

If you are unsure whether the Aliva oral hydration system is suitable for your condition, consult your physician or qualified medical provider before use.

When to stop using the Aliva pump or the supplied accessories

If the power adapter is damaged, you should disconnect the adapter from the wall socket and the pump and stop using it.

If there is any visible damage to the Aliva oral hydration system, or any abnormal smell or sound coming from the Aliva pump, stop using the product and set the main power switch to the OFF position.

02. Explanation of safety warnings

Explanation of safety warnings used in the instructions for use:

⚠ WARNING

Warning indicates a hazard with a medium level of risk which, if not avoided, could result in death or serious injury.

⚠ CAUTION

Caution indicates a hazard with a low level of risk which, if not avoided, could result in minor or moderate injury.

NOTICE

Indicates information considered important, but not hazard-related.

List of safety warnings

⚠ WARNING

Read and understand this manual and its safety instructions before using this product. Failure to do so can result in serious or fatal injury.

Do not handle the device without good hand hygiene. Poor hand hygiene increases the risks of contamination and infection.

Do not use any other liquid than clean drinking water. Using other fluids may damage the device and increase the risk of contamination, infection or intoxication.

Do not add any substances to the water. This includes medications, flavourings, thickeners or similar. Additives can impair the function of the device and compromise user safety.

Do not use the device if you are on fluid restriction without consulting your medical doctor.

In conditions such as heart failure or kidney failure, excess fluid may lead to serious medical complications. Water used with the Aliva system must be included in your daily fluid balance.

Do not continue use if oral ulcers develop. Pause use until the area is healed. Watch for signs of infection (redness, swelling, discharge) and contact your doctor if symptoms worsen.

Do not use the device if you have severe swallowing difficulties without consulting your medical doctor. Risk of aspiration or chemical irritation may occur.

Do not use the device during MRI, HBO treatment or near a defibrillator in use. The device must be turned off before entering these environments (magnetic resonance imaging, hyperbaric chamber environment) to avoid safety risks and interference.

Do not stack the Aliva system with other equipment. Stacking or placing the device next to other systems can cause interference or malfunction. If this setup is necessary, monitor all devices closely for proper operation.

Do not use accessories or cables not provided or approved by the legal manufacturer. Using unapproved components like power adapter and cables can reduce safety, affect EMC performance, and damage the device.

Do not place RF equipment closer than 30 cm to the device. Portable Radio Frequency (RF) communications devices, including radios and external antennas, may interfere with the device. See details in Appendix.

Do not bite or chew on the nozzle of the mouthpiece. Chewing may cause parts to break off and pose a choking hazard.

Do not disassemble or modify any part of the device. Tampering may cause malfunction, loss of warranty and risk of harm.

Nozzle position. The nozzle must be placed underneath the upper lip as described in the *Operation and How To Use* section of this IFU. Misplacement can reduce the effectiveness of the system and lead to discomfort or accidental leakage.

Reposition after displacement. If the nozzle becomes dislodged during use (e.g. when eating), stop use and reposition it correctly before resuming operation.

Oral discomfort or soreness. If you experience soreness at the site of nozzle placement, consider limiting usage time to avoid blisters or ulcers.

Risk of contamination. Proper hygiene is essential for safe use. Always clean and handle the components as instructed to avoid contamination and reduce infection risk. Always use a sealed and undamaged package of mouthpiece when installing a new one.

Risk of damage. Do not pull on the Aliva mouthpiece. This could cause damage to the tubing or dislocating the nozzle.

Excessive amounts of liquid. The Aliva oral hydration system has built in safety measures to prevent overflow of water due to software malfunction. The safety does not rely on correct functioning of the software. The maximum flow of water is limited by hardware, so in case of a malfunction of the software the motor will be stopped by a hardwired circuit.

Nozzle position. Only place the nozzle of the mouthpiece in the oral cavity during use, according to instructions.

Misuse. Do not drop the Aliva pump to prevent damage to the device.

Cleanliness. This product is non-sterile, but manufactured and packaged under controlled conditions. Bioburden is tested on new and unused products and kept within safe limits.

03. Description of the product

The Aliva oral hydration system consists of the items in the table below.

Name	Article number
Aliva pump	07090057970022
Aliva mouthpiece	07090057970015
Power adapter	07090057970008
The Aliva bottle (accessory)	07090057970039

Process overview – how it works

The Aliva oral hydration system supplies the oral cavity with small amounts of water on a continuous basis during use. The nozzle of the Aliva mouthpiece is placed underneath the upper lip and is where the water is distributed. The nozzle is integrated with a tube, which is assembled to the Aliva pump using the Easy-bracket. The bracket is placed and connected in a designated area of the pump. The inlet of the tube is placed in the Aliva bottle.

The bottle represents a reservoir which is filled with clean, unthickened drinking water. The peristaltic Aliva pump enables transport of water from the bottle, through the tubing to the nozzle. The three options of preset flow rates mimic the levels of saliva that healthy salivary glands produce during rest (level 1), while unstimulated (level 2) and when stimulated (level 3). The device is used intermittently, depending on the user's need and level of oral dryness. The user of the device can set the flow rate to his or her need to a period of 60 minutes and up to a maximum of 12 hours per day/24 hours.

Product illustrations and explanations


Illustration 1: Aliva pump front and side

No.	Element	Function
A	Bottle holder	To hold the Aliva bottle.
B	Aliva peristaltic pump	Pumps the water from the bottle, through the tubing and out of the nozzle in the Aliva mouthpiece.
C	Lock	To securely connect the Easy-bracket to the Aliva pump.
D	Power port	To connect the power cable to the Aliva pump. The power port is a USB-C connector.


Illustration 2: Aliva pump front side and buttons panel

No.	Element	Function
E	RFID reader (Radio Frequency Identification)	Scans the RFID tag of the Aliva mouthpiece to activate the pump and to monitor the remaining usage time.
F	Buttons panel	Used to operate the device and view status from the LED indicators.


Illustration 3: Control panel

No.	Element	Function
G	Operation LED	Indicates the operating status of the Aliva pump: Green: Standby mode. Blue: operation mode. Yellow: alarm condition
H	Flow rate LEDs	Indicates the selected and active flow rate. The corresponding LED illuminates blue.
I	Flow level buttons 1, 2 and 3	Used to set the water flow rate. Level 1 is the lowest flow rate and level 3 is the highest.
J	Sleep mode symbol	Indicates sleep mode (flow level 1), recommended during sleep.
K	Braille knobs	Help visually impaired users identify and select the flow rate buttons. Level 1: A, level 2: B, and level 3: C.
L	Battery power indicator (LEDs)	Indicates the remaining charge of the internal rechargeable battery using five white LEDs. Each lit LED ≈ 20% battery capacity. All LEDs lit: battery fully charged. Bottom LED flashing: battery ≤10%. All LEDs off + yellow operation LED: the battery is depleted and pump stops.
M	Operation button	Used to switch the Aliva pump between Standby mode and operation mode.


Illustration 4: Aliva pump back side

No.	Element	Function
N	Main power switch	Used to turn the Aliva pump ON or OFF. Connects or disconnects the power adapter and the internal rechargeable battery.
O	Tube holder	Holds the Aliva mouthpiece when not in use.
P	Carrying grip	Enables a secure grip when carrying or moving the Aliva pump.


Illustration 5: Power adapter

No.	Element	Function
Q	Power adapter	Supplies power to the Aliva pump.


Illustration 6: Aliva mouthpiece

No.	Element	Function
R	Mouthpiece nozzle	Serves as the outlet of the Aliva mouthpiece. Is placed in the mouth during use and distributes water to the oral cavity.
S	Silicone tubing, nozzle side	Transports drinking water from the Aliva pump to the mouthpiece nozzle.
T	One-way check valve	Maintains water inside the tubing between uses.
U	Silicone bottle cap	Connects the Aliva mouthpiece to the Aliva bottle
V	Silicone tubing, bottle side	Transports drinking water from the Aliva bottle to the Aliva pump.
W	Easy-bracket®	Connects the Aliva mouthpiece to the Aliva pump.
X	RFID tag (embedded)	Communicates with the RFID reader to activate the Aliva pump and deactivate it once the allowed usage period has expired.


Illustration 7: Aliva bottle

No.	Element	Function
Y	Aliva bottle	Serves as a reservoir for filling and refilling with drinking water.
Z	Stainless steel cap	Keeps drinking water inside the Aliva bottle. Not used during active operation of the device.

04. Operation and how to use

Make sure to read the entire chapter 3 before first time assembling the product.

⚠ WARNING

- Note that the original packing of the Aliva oral hydration system is suitable for safe storage and/or travel, and we advise you to keep it for these purposes.
- Do not handle the device without good hand hygiene. Poor hand hygiene increases the risks of contamination and infection.
- Note that the Aliva pump shall only be powered using the supplied medical-grade power adapter or the integrated rechargeable battery.
- The mouthpiece can be used for 14 days starting from the moment it is installed in the Aliva pump. After 14 days, the Aliva pump stops working and a new Aliva mouthpiece needs to be installed.
- Wash or disinfect your hands before touching the mouthpiece. Unclean hands can introduce harmful microorganisms and cause infection.
- Make sure the packaging of the new, unused Aliva mouthpiece is sealed and intact before opening.
- The Aliva mouthpiece is for single-person use only.
- The Aliva mouthpiece needs to be replaced at a minimum every 14 days with a new, unused Aliva mouthpiece. The RFID tag of the Aliva mouthpiece is scanned by the RFID reader of the Aliva pump at the time of installation. Counting 14 days after first activation, the Aliva pump deactivates when the usage period ends. When a new mouthpiece is detected by the RFID reader, the pump is activated again.
- The device may be used for up to 12 hours per day, depending on individual needs. Check the tubing of the Aliva mouthpiece regularly for signs of damage or wear. If you see any damage, immediately stop the pump and replace the mouthpiece.
- Do not use any other liquid than clean, unthickened drinking water. Using other fluids may damage the device and increase the risk of contamination, infection or intoxication.
- **Do not add any substances to the water.** This includes medications, flavourings, thickeners or similar. Additives can impair the function of the device and compromise user safety.

How to unpack the Aliva oral hydration system for the first time

1. Unpack the Aliva oral hydration system.
2. Remove any packaging and/or restraints.
3. Make sure that the products are not damaged. Contact the manufacturer in case of damage.
4. Verify the contents of the Aliva oral hydration system:
 - Aliva pump
 - Power adapter
 - Aliva bottle (two Aliva bottles)
 - Aliva mouthpiece (two Aliva mouthpieces, in single-packs)


Setting up the Aliva pump

Before use, verify that the Aliva oral hydration system is placed on a clean, stable and horizontal surface.

1. Make sure the main power switch of the Aliva pump is in the OFF position.
2. Connect the USB-C plug of the power cable to the power port of the Aliva pump and connect the power adapter to the wall socket.
3. Set the main power switch of the Aliva pump to the ON position. The Operation LED lights up green to indicate the Standby mode.

The Aliva pump is now ready for use.

Control panel LED lights

- | | |
|---|--|
|  | Main Power OFF: No LEDs illuminated. |
|  | Main Power ON. Standby mode. Operation LED green. |
|  | Operation activated: Ready for use, Operation LED lights blue. |
|  | Priming mode: The LED above selected flow level (either 1, 2 or 3) is flashing blue. |
|  | Operation at level 3 Operation LED and level 3 LED light blue. |


-  Operation at level 2: Operation LED and level 2 LED light blue.
-  Operation at level 1: Operation LED and level 1 LED light blue.
-  Operation LED yellow. Easy-bracket may be misplaced or expired.

Illustration 8: Control panel and LED lights

Installing the Aliva mouthpiece

How to install a new mouthpiece


1. Make sure the Main power switch is in the OFF position.
2. Obtain a new, unused Aliva mouthpiece.
3. Fill a clean Aliva bottle with clean, unthickened drinking water.

4. Insert the end of the tubing into the Aliva bottle. Then press the integrated silicone cap firmly onto the bottle opening until it is securely in place.

5. Turn the lock into the open position.


6. Align the Easy-bracket with the receiving slot on the Aliva pump. Hook the top of the Easy- bracket® into place.

7. Gently push the lower part of the Easy-bracket into place.

8. Turn the lock to the closed position.
This will push the Easy-bracket upward, ensuring it is properly secured to the pump.


9. When not in use, you can wind the tubing around the Tube holder located on the back of the Aliva pump.

Removing the Aliva mouthpiece

1. Make sure the main power switch on the back of the Aliva pump is in the OFF position

2. Rotate the lock of the Aliva pump to the open position.

3. Gently lift and slide down the Easy-bracket out of the pocket in the Aliva pump.

4. Loosen the bottle cap from the Aliva bottle and remove the integrated bottle cap with tubing from the bottle.

Priming mode when installing a new Aliva mouthpiece

1. When the Aliva pump is in operation mode (blue LED), press and hold one of the three flow rate buttons for three seconds, to fill the tubing with water. The pump speeds up for approximately 30 seconds to fill the tubing. The selected flow level LED flashes blue. When the priming period is completed, the Aliva pump returns to the set flow rate.
2. To end the priming mode before the period has ended, press one of the flow rate buttons or the Operation button.

3. Priming Mode: The LED above selected flow level (either 1, 2 or 3) is flashing blue.

NOTICE


The priming mode can also be used to flush the tubing and the nozzle of the mouthpiece.

Placing the Aliva mouthpiece nozzle in the mouth


1. Ensure proper personal hygiene. Wash your hands thoroughly before handling the Aliva mouthpiece, to reduce risk of contamination.

NOTICE


It is recommended to moisture the nozzle of the Aliva mouthpiece before placing it, especially if you have a very dry mouth. Use clean water for this.

2. Hold the nozzle at the center, where it connects to the tubing. The nozzle is symmetrical with no front or back and can be placed either way.

3. Place the Aliva mouthpiece nozzle under your upper lip, resting against the upper teeth.

4. Use two fingers to gently push the nozzle further up under the upper lip, above the teeth, and against the gums.

5. Use your fingers to gently press the sides of the nozzle outward from center to improve contact between nozzle and gums.

How to start and stop the flow of the water

1. Set the main power switch in the ON position. The green LED above the Operation button indicates that the Standby mode is active.

2. Press the Operation button. The Operation LED lights up blue when a valid Easy-bracket is correctly installed.

3. Press the flow level button corresponding to the desired setting (1, 2 or 3). The LED for the selected flow level lights up blue. The Aliva pump starts, and when the water reaches the nozzle, it will begin to drip. To change the flow rate, press a different flow level button.

4. To stop the flow, press the Operation button. The LED above the Operation button lights up green to indicate Standby mode.

5. If the Operation LED lights up yellow, the RFID tag on the Easy-bracket is either not detected, or has expired.

NOTICE

A small amount of residual water in the nozzle may cause dripping after the pump has stopped.

If you plan to pause the use of Aliva oral hydration system while it is running on battery power, turn the Main power switch to the OFF position, to prevent the battery from draining.

Battery level LED indications

Battery level indicator: 5 white LEDs = 80-100% battery capacity.

Battery level indicator: 4 white LEDs = 60-80% battery capacity.

Battery level indicator: 3 white LEDs = 40-60% battery capacity.

Battery level indicator: 2 white LEDs = 20-40% battery capacity.

Battery level indicator: 1 white LEDs = 10-20% battery capacity.

Battery level indicator: 1 white LED light flashing, ≤ 10% battery capacity.

Operation LED lights up yellow and all white LED lights are off. Battery is out of charge and needs charging.

How to charge the internal battery

1. Set the main power switch to the OFF position.
2. Connect the power adapter to the power port and the wall socket.
3. Set the main power switch to the ON position.

The internal rechargeable battery is charging.

NOTICE

Battery charging time is approximately 6 hours. Battery operating time is up to 12 hours, depending on flow rate and ambient temperature.

How to operate the internal battery

1. Set the main power switch to the OFF position.
2. Disconnect the power adapter from the wall socket and the power port.
3. Set the main power switch of the Aliva pump to the ON position.

The Aliva pump is now powered by the internal battery and ready for use. The battery level indicator shows the remaining battery capacity.

How to move and transport the Aliva system

1. Set the main power switch to the OFF position.
2. Disconnect the power adapter from the wall socket.
3. Disconnect the power cable from the power port.
4. Wind the tube of the Aliva mouthpiece securely around the tube holder.
5. Use the carrying grip to lift the Aliva pump. Keep the Aliva pump level to prevent water from spilling from the Aliva bottle.

If you intend to transport the Aliva oral hydration system over longer distances, empty the bottle to prevent spills. Always use the original packaging or a suitable protective travel bag.

How to store the product safely

1. Set the main power switch to the OFF position.
2. Disconnect the power adapter from the wall socket and the power cable from the power port.
3. Remove the Aliva mouthpiece from the Aliva pump.
4. Remove the Aliva bottle from the Aliva pump.
5. Store the Aliva oral hydration system in a dry place, away from direct sunlight and heat sources.

NOTICE

The original packaging is recommended for safe and protective storage.

05. Maintenance

Cleaning the Aliva pump

Never submerge the Aliva pump. Doing so may cause serious damage to the device and its electronic components.

1. Ensure the main power switch is in the OFF position.
2. Disconnect the power adapter from the wall socket and the power port and clean the adapter separately.
3. Remove the Easy-bracket from the pump.
4. Remove the bottle from the bottle holder.
5. For home use, clean the outer surfaces of the Aliva pump with a clean, damp cloth using tap water and a food-safe detergent.

Cleaning the mouthpiece nozzle and tubing

1. Remove the Easy-bracket from the Aliva pump.
2. Rinse the nozzle with running tap water, then gently shake off excess water.
3. Wipe the tubing of the Aliva mouthpiece with clean, damp soft cloth using tap water and a food-safe detergent.

NOTICE

The priming function can be used to flush the tubing.

Cleaning the Aliva bottle

The Aliva bottle should be cleaned before first use and once daily thereafter.

1. The bottle is dishwasher safe. Place it upside down during the wash cycle.
2. Alternatively, clean the bottle using warm tap water and a small amount of food-safe detergent. Use a clean dish brush and rinse thoroughly thereafter.

Disinfection of the Aliva pump in healthcare settings

For disinfection in clinical environments, use only disinfectants that are approved on medical devices and suitable for healthcare settings (i.e. Virkon™, AntiBac 70%).

06. Troubleshooting

Error	Cause	Solution
The mouthpiece nozzle distributes too much liquid.	Flow rate is set too high.	Select a lower flow rate using the flow rate buttons.
The mouthpiece nozzle distributes too little liquid.	Flow rate is set too low.	Select a higher flow rate using the flow rate buttons.
	The Aliva mouthpiece is damaged and leaking to the surrounding environment.	Replace with a new Aliva mouthpiece.
The mouthpiece nozzle does not distribute any liquid.	A new Aliva mouthpiece has been installed and is not yet primed.	Use Priming mode to fill the tubing with liquid.
	The Easy-bracket is not installed, or installed incorrectly.	Check that the bracket is correctly installed.
	The inlet of tubing is incorrectly placed in the bottle.	Visually check and correct placement.
	The Aliva bottle is empty.	Refill the bottle with clean drinking water.
	The Aliva mouthpiece has expired.	Replace with a new, unused Aliva mouthpiece.
	The Aliva mouthpiece is damaged, or the flow is blocked.	Replace with a new Aliva mouthpiece.
	The Aliva mouthpiece is damaged and leaking.	Replace with a new Aliva mouthpiece.

Error	Cause	Solution
Operation LED lights yellow.	The Easy-bracket is not correctly installed.	Repeat the procedure for installing the Aliva mouthpiece.
	RFID usage time has expired.	Replace the old Aliva mouthpiece with a new.
Operations LED lights yellow. No other LEDs are lit.	Internal battery is empty.	Recharge the battery. Use the power adapter as alternative power source.
The Aliva pump does not enter Standby mode, even though the main power switch is ON.	The power adapter is not connected, and the internal battery is empty.	Connect the power adapter to supply power to the pump.
The internal rechargeable battery does not charge.	The main power switch is not in the ON position.	Check if the main power switch is set to ON .
	Power adapter is not properly connected to the wall socket.	Check that the power adapter is connected to the wall socket.
	Power cable is not connected to the pump.	Check that the power cable is securely connected to power port. If the internal rechargeable battery still cannot be charged, it may be defective and may need to be replaced. Contact Aliva AS or your supplier for support.

07. Disposal

Always check with your local waste management policy on how to dispose materials and device properly. Local regulations may differ from the instructions given in this section.

Disposal of electronic components

The symbol on the product, its accessory and packaging indicate that this device must not be treated as general waste but must be collected separately. This implies dispose of the device via a collection point for recycling of waste electrical and electronic equipment if you live within the EU and in other European countries. By disposing of the device in the proper manner, you help to avoid possible hazards for the environment and public health that could otherwise be caused by improper treatment of waste. Recycling of materials contributes to the conservation of natural resources.

Disposal of batteries

This product contains an internal, non-removable, rechargeable lithium-ion battery.

Batteries must not be disposed of as unsorted household waste. Improper disposal may pose risks to the environment and public health due to the presence of hazardous substances.

The battery is not intended to be removed by the user. Do not attempt to open the device or remove the battery.

At end of life, the product shall be disposed of via an authorised collection or recycling scheme for waste electrical and electronic equipment (WEEE) or waste batteries, in accordance with applicable local regulations.

By disposing of the product and battery correctly, you help prevent potential negative consequences for the environment and human health and support the conservation of natural resources.

Disposal of packaging waste

The packaging is made of recycled and environmentally friendly materials and should be disposed through your local recycling facilities. Follow the symbols on the packaging indicating what is recycled pulp, unbleached pulp, or pulp bleached without chlorine gas in accordance with the TCF or ECF method.

Disposal of the Aliva mouthpiece

The Aliva mouthpiece is made from silicone. Silicone cannot be treated as plastic waste but should preferably be deposited at your local recycling facility or disposed as general waste.

The Easy-bracket is made of plastic and can be disposed and recycled as plastic waste.

How to dispose of the Aliva mouthpiece

1. Separate the Easy-bracket from the silicone tubing of the Aliva mouthpiece.
2. Dispose of the silicone tube, bottle cap and nozzle as silicone waste at local collection point, or in general waste.
3. The Easy-bracket can be disposed of and recycled as plastic waste.

Disposal of the empty Aliva bottles

The Aliva bottles may be recycled as glass waste.

08. Appendixes

Appendix 1 – Supplied accessories

Supplied Accessories	Article number
Aliva bottle	07090057970039

Appendix 2 – Safety information symbols and explanations

Symbol	Meaning
	Aliva® logo
	Indicates the medical device legal manufacturer: Aliva AS, Rigedalen 41, 4626 Kristiansand, Norway, Telephone: +47 488 48 838, E-mail: post@aliva.no, www.aliva.no
	Indicates the date when the medical device was manufactured. Format: YYYY-MM-DD.
	Indicates the country of manufacture: China.
	Indicates the country of manufacture: Taiwan.
	Indicates that the device fulfils requirements from harmonized standards to apply the CE mark.
	Indicates the date after which the medical device is not to be used. Format: YYYY - MM - DD.
	Indicates the manufacturer's batch code so that the lot or batch can be identified.
	Indicates the item is a medical device.
	Indicates the catalogue number so that the device can be identified.
	Indicates the serial number so that a specific device can be identified.

Symbol	Meaning
	Indicates the quantity in the packaging.
	Indicates a carrier that contains unique device identifier (UDI) information.
	Indicates a medical device that may be used multiple times on a single patient.
	Indicates the need for the user to consult the instructions for use.
	Indicates the temperature limits to which the medical device can be safely exposed.
	Indicates Class II electrical equipment.
	Indicates a medical device that needs protection from sunlight.
	Indicates a medical device that needs to be protected from moisture.
	Indicates that the product is suitable for contact with food.
	Indicates not to use the product if the packaging is damaged.
	Indicates the medical device uses direct current.
	Indicates that the product should not be disposed with unsorted household waste. Must be collected separately for recycling.
IP22	Indicates protection class IP22: Protected against the ingress of solid objects ≥ 12.5 mm and against vertically falling water drops when tilted up to 15° .
	Indicates the humidity range to which the medical device can safely be exposed.
	Indicates that the user must consult the instructions before use.
	Indicates that the product is a Type BF applied part, with basic protection against electric shock
	Indicates the material is clear glass (GL 70) and suitable for recycling.

Appendix 3 – Technical data tables – Aliva pump and mouthpiece

Specification	Value
Product name	Aliva pump (component of Aliva oral hydration system)
Model /catalogue number	AV001
Energy consumption	15 W
Battery type	Type: Internal, non-replaceable lithium-ion battery Capacity: 7.3 V, 2450 mAh, 17.89 Wh Lifetime: >80% capacity after 500 charging cycles
Flow rates of water	Flow rate level 1: 0.1 ml/min Flow rate level 2: 0.5 ml/min Flow rate level 3: 1.5 ml/min
IP Code	IP22. Protected against the ingress of solid objects ≥ 12.5 mm and against vertically falling water drops when tilted up to 15°.
Operation conditions	Temperature: +5°C to +40°C Humidity: 15-93%, non-condensing Barometric pressure: 700-1060hPa
Storage/transportation conditions	Temperature: -20°C - +50°C Humidity: <93% non-condensing
Power adapter	Model: UES18LCP4-050300SPA Input: AC 100-240V, 50/60Hz, 500mA Output: DC 5.0V, 3.0A, 15.0W

Specification	Value
Product name	Aliva mouthpiece (component of Aliva oral hydration system)
Model/catalogue number	ALVM001
Material	Tubing, nozzle and one-way check valve: medical-grade silicone in accordance with ISO 10993 Bracket: Polypropylene Random Copolymer
Storage/transportation conditions	Temperature: +5 °C to +40 °C Avoid storage in direct sunlight
Sterility	Non-sterile

Specification	Value
Cleaning	Disinfection and sterilisation shall not be performed
Shelf life	36 months from date of manufacture
Use	May be used by the same patient for up to 14 days from first use

Appendix 4 – Electromagnetic Compatibility

⚠ WARNING

The Aliva Oral Hydration System is intended for use in hospitals or clinical settings. The device must not be used in magnetic resonance imaging (MRI) environments, hyperbaric chamber (HBO) environments, or in proximity to an operating defibrillator. Disconnect the device if any of these situations occur.

Use of Aliva oral hydration system adjacent to or stacked with other equipment should be avoided, as this may result in improper operation. If such use is necessary, the Aliva system and the other equipment must be observed to verify normal operation.

Use of accessories and cables other than those specified or provided by Aliva AS may degrade EMC performance.

Portable Radio Frequency (RF) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Aliva oral hydration system including manufacturer-specified cables. Otherwise, degradation of system performance may occur. See table below for additional details regarding recommended separation distances between portable and mobile RF communications equipment and the Aliva oral hydration system.

Guidance and manufacturer's declaration – electromagnetic emissions		
The Aliva oral hydration system is in compliance for each EMISSIONS test specified by the standard, e.g. EMISSIONS class and group.		
Emissions	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Aliva oral hydration system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF emissions CISPR 11	Class B	The Aliva oral hydration system is suitable for use in hospitals or clinical settings, including those directly connected to the public low voltage power supply network.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The **Aliva oral hydration system** is in compliance for each IMMUNITY test specified by the standard, e.g. IMMUNITY test level.

Immunity test	IEC 60601-1-2 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Radiated RF EM fields IEC 61000-4-3	10 V/m 80 MHz - 2.7GHz 80% AM at 1 kHz	10 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz
Electrical fast transient/burst IEC 61000-4-4	±2 kV 100 kHz repetition frequency	±2 kV 100 kHz repetition frequency
Surge IEC 61000-4-5	±0.5 kV, ±1 kV line-to-line; ±0.5 kV, ±1 kV and ±2 kV line-to-ground;	±0.5 kV, ±1 kV line-to-line; ±0.5 kV, ±1 kV and ±2 kV line-to-ground;
Conducted disturbances induced by RF fields IEC 61000-4-6	3V 0.15 MHz - 80 MHz 6V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	3V 0.15 MHz - 80 MHz 6V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT: 0.5 cycle ^{a)} At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°.	0% UT: 0.5 cycle ^{a)} At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°.
	0% UT: 1 cycle 70% UT: 25/30 cycles ^{b)} Single phase: at 0°	0% UT: 1 cycle 70% UT: 25/30 cycles ^{b)} Single phase: at 0°
	0% UT: 250/300 cycles ^{b)}	0% UT: 250/300 cycles ^{b)}

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz
NOTE	^{a)} UT is the a.c. mains voltage prior to application of the test level; ^{b)} E.g. 25/30 means 25 periods at 50 Hz or 30 periods at 60 Hz.	

Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test Frequency (MHz)	Band (MHz)	Service ^{a)}	Modulation ^{b)}	Max. Power (W)	Distance (m)	Test Level (V/m)	Compliance level
385	380-390	TETRA 400	Pulse Modulation 18 Hz	1.8	0.3	27	27
450	430-470	GMRS 460 FRS 460	FM ^{c)} ±5 kHz deviation 1 kHz sine	2	0.3	28	28
710	704-787	LTE-Brand 13, 17	Pulse Modulation 217 Hz	0.2	0.3	9	9
745							
780							
810	800-960	GSM 800/900, TETRA 800, Iden 820, CDMA 850, LTE Band 5	Pulse Modulation 18 Hz	2	0.3	28	28
870							
930							
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation 217 Hz	2	0.3	28	28
1845							
1970							
2450	2400-570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation 217 Hz	2	0.3	28	28
5240	5100-5800	WLAN 802.11 a/n	Pulse Modulation 217 Hz	0.2	0.3	9	9
5500							
5785							
NOTE	^{a)} For some services, only the uplink frequencies are included ^{b)} The carrier shall be modulated using a 50 % duty cycle square wave signal. ^{c)} As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case						

Recommended separation distances between portable and mobile RF communications equipment and the Aliva oral hydration system

The Aliva oral hydration system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Aliva oral hydration system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Aliva oral hydration system as recommended below, according to the maximum output power of the communications equipment.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Aliva oral hydration system.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = 0.58\sqrt{P}$	$d = 0.13\sqrt{P}$	$d = 0.25\sqrt{P}$
0.01	0.058	0.013	0.025
0.1	0.18	0.041	0.079
1	0.58	0.13	0.25
10	1.8	0.41	0.79
100	5.8	1.3	2.5

For transmitters rated at a maximum output power not listed above, the recommended separation distanced in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people