

The logo features a large, stylized red 'C' shape on the left. To its right, the word 'ThermaCor' is written in a bold, blue, sans-serif font with a registered trademark symbol. Above the 'Cor' part of 'ThermaCor', the number '1200' is written in a smaller, red, sans-serif font. Below 'ThermaCor', the tagline 'Infusing Warmth in a Heartbeat' is written in a blue, italicized, sans-serif font. At the bottom of the logo, the text 'SMISSON-CARTLEDGE BIOMEDICAL' is written in a black, bold, sans-serif font.

ThermaCor[®] 1200
Infusing Warmth in a Heartbeat
SMISSON-CARTLEDGE
BIOMEDICAL



RAPID THERMAL INFUSER



10 mL/HOUR to 1200 mL/MINUTE

Fast, Easy Set-Up

- Sterile, single patient disposable cassette
- Color-coded, mated, snap-in connectors for fast, correct set up
- Primed and ready in less than 60 seconds
- Targeted 40° C continuous heating at all flow rates

Automated Air Purge

- Patented ShurGard Vortex Air Trap
- Efficiently removes air at high flow rates
- Four separate air detection sensors
- Safety auto shutdown
- Fast, easy recovery from fluid out

Reduced Risk of Hypervolemia

- One touch bolus operation
- 4 Pre-set bolus selections - 100, 250, 500, 1000 mL
- Controlled rapid bolus delivery
- Auto-return to slow flow setting
- Precise, safe volume control



Advanced Safety Technology

- Continuous monitoring with 4 temp sensors, 4 air sensors and 3 pressure sensors
- Clear, concise, plain text directions speed recovery from alarms
- Ensures rapid return to continuous infusion

Pressure Controlled Flow

- Greater safety for pediatric and geriatric patients
- Automatically reduces flow rate based on pressure selected
- 300 or 100 mmHg pressure limits
- Facilitates connection to any catheter with dramatically reduced risk of hemolysis or catheter site damage

Portable, Battery Back-Up Power

- Up to 2-hours reserve power
- Minimum 30 minutes at maximum flow rate (1200mL/min)
- Full range of flow rates from 10 mL/hr to 1200 mL/min with Bolus capability

Product Listing

TIS-1200	ThermaCor® 1200 Rapid Infuser
DNC-1200	Cassette, Dual Patient Line
PNC-1200	Cassette, 6' patient Line
PTC-1200	Cassette, 6' Patient Line, 3 Spike Set Inflow Line
MR-1200	Rigid Fluid Reservoir
MH-1200	Rigid Fluid Reservoir Holder
TSS-1200	3 Spike Set Inflow Line
TH-1200	3 Spike Set Inflow Line Holder
PL-1200	Patient Line – 6'
EL-1200	Extension Line – 3'
DL-1200	Dual Patient Line – 6'
TD-1200	ThermaCor® 1200 Drape
SP-1200	Replacement Silpad
CG-1200	Replacement Component Guard
RC-1200	Replacement US Power Cord
RCB-1200	Replacement UK Power Cord
RCE-1200	Replacement European Power Cord



487 Cherry Street • Third Street Tower • Macon, GA 31201 • USA
 Phone: 478-744-9992 / Fax: 478-744-9996 • www.ThermaCor1200.com