

## Affinity™ AF100 Arterial Filter

Arterial Filter	
Min/Max Flow Rate	1-7 L/min blood flow
Prime Volume	100 mL
Inlet/Outlet Size	3/8" / 3/8"
Screen Size	30 µm
Biosurface Options	Cortiva™ or Balance™



### Arterial Filter

CFN	Product	Units per Case
CB851	Affinity AF100 Arterial Filter with Cortiva BioActive Surface	Tubing packs only
BB851	Affinity AF100 Arterial Filter with Balance Biosurface	4 per case