

# Centrifugal Pumps

AP40	
Style	Impeller
Prime Volume	40 mL
Inlet/Outlet Size	$\frac{3}{8}$ " inlet $\frac{3}{8}$ " outlet
Max Speed	4,000 RPM
Max Flow Rate	10 L/min

## Centrifugal Blood Pumps

CFN	Product	Units per Case
CBAP40	Affinity CP Centrifugal Blood Pump with Cortiva BioActive Surface	4 per case
BBAP40	Affinity CP Centrifugal Blood Pump Balance™ Biosurface	4 per case
AP40	Affinity CP Centrifugal Pump	4 per case

