

## Affinity™ Arterial Filter

### Arterial Filter

Min/Max Flow Rate 1-7 L/min blood flow

Prime Volume 212 mL

Inlet/Outlet Size 3/8" / 3/8"

	<b>351</b>	<b>353</b>
Screen Size	38 µm	20 µm

Biosurface Options Uncoated



### Adult Arterial Filter

CFN	Model #	Product	Units per Case
61399401124	351	Affinity Arterial filter with 38 µm screen	Tubing packs only
61399409434	353	Affinity Arterial Filter 20 µm Screen	Tubing packs only