

Affinity Fusion™ Oxygenation System*

Oxygenator	
Min/Max Flow Rate	1-7 L/min blood flow
Prime Volume	260 mL
Inlet/Outlet Size	3/8" / 3/8"
Membrane	Hollow Fiber Polypropylene
Membrane Surface Area	2.5 m ²
Biosurface Options	Cortiva™ or Balance™
Cardiotomy/Venous Reservoir for Open System	
Min/Max Flow Rate	1-7 L/min venous blood flow 1-6 L/min cardiomy blood flow
Volume	4,500 mL
Inlet/Outlet Size	1/2" or 3/8"* / 3/8"
Biosurface Options	Balance



Affinity Fusion Oxygenator with Integrated Arterial Filter

CFN	Product	Units per Case
CB841	Oxygenator with Cortiva BioActive Surface and Cardiotomy/Venous Reservoir with Balance Biosurface	1 per case
CB811	Oxygenator with Cortiva BioActive Surface	1 per case
BB841	Oxygenator and Cardiotomy/Venous Reservoir with Balance Biosurface	1 per case
BB811	Oxygenator with Balance Biosurface	1 per case
RCL841	Affinity Fusion Recirculation Line*	12 per case

Accessories and Holders for Affinity Fusion Oxygenation System

CFN	Product	Units per Case
ATP210	Affinity™ Temperature Probe	1 per case
AUH2093	Affinity Orbit™ Holder System	1 per case
AUH2070	Affinity Orbit Oxygenator Arm Replacement	1 per case
AMH2014	Affinity Manifold Holder	1 per case



*No recirculation line is included with Affinity Fusion product.

†Using a 3/8" to 1/2" adapter provided with product.

Products with Balance™ Biosurface Technology are licensed under agreement from BioInteractions, Limited, United Kingdom.