

## Affinity Pixie™ Oxygenation System

### Oxygenator

Min/Max Flow Rate	0.1-2.0 L/min blood flow
Prime Volume	48 mL
Inlet/Outlet Size	1/4" / 1/4"
Membrane	Hollow Fiber Polypropylene
Membrane Surface Area	0.67 m <sup>2</sup>
Biosurface Options	Cortiva™ or Balance™
Included Accessories	2 – Tubing adapter – 3/16" 1 – Recirculation line – 3/16"

### Cardiotomy/Venous Reservoir

Min/Max Flow Rate	0.1-2.0 L/min venous blood flow
Max Cardiotomy Flow Rate	2.0 L/min cardiotomy blood flow
Volume	1,200 mL
Inlet/Outlet Size	1/4" / 1/4"
Biosurface Options	Cortiva or Balance
Included Accessories	2 – Tubing adapters – 3/16" 2 – Tubing adapters – 3/8" 2 – Flexible Luer Lock adapters

1 – Sampling manifold assembly



## Affinity Pixie Oxygenation System

CFN	Product	Units per Case
CBP241	Affinity Pixie Oxygenator and Cardiotomy/Venous Reservoir with Cortiva BioActive Surface	Tubing packs only
CBP211	Affinity Pixie Oxygenator with Cortiva BioActive Surface	Tubing packs only
BBP241	Affinity Pixie Oxygenator and Cardiotomy/Venous Reservoir with Balance Biosurface	4 per case
BBP211	Affinity Pixie Oxygenator with Balance Biosurface	Tubing packs only

Products with Balance™ Biosurface Technology are licensed under agreement from BioInteractions, Limited, United Kingdom.