


## Cardiotomy Reservoir

	1351	EL2120	EL240	EL400	EL402	EL404
Style	Hard Shell	Hard Shell	Hard Shell	Hard Shell	Hard Shell	Hard Shell
Volume	2,600 mL	4,000 mL	4,000 mL	4,000 mL	4,000 mL	4,000 mL
# Ports	4	2	2	4	4	4
Screen Size	20 µm	120 µm	40 µm	120 µm	20 µm	40 µm

### Intersept™ Cardiotomy Reservoir

Used for Cardiopulmonary Bypass and Post-op Chest Drainage

CFN	Product	Units per Case
1351	Intersept Cardiotomy Reservoir with Filter	6 per case



### EL Series Blood Collection Reservoir — 2 Port

Used for Autotransfusion

CFN	Product	Units per Case
EL2120	Blood Collection Reservoir	6 per case
EL240	Blood Collection Reservoir	6 per case



### EL Series Blood Collection Reservoir — 4 Port

Used for Cardiopulmonary Bypass and Autotransfusion

CFN	Product	Units per Case
EL400	Cardiotomy Reservoir	6 per case
EL402	Cardiotomy Reservoir	6 per case
EL404	Cardiotomy Reservoir	6 per case

