

Soft Shell Venous Reservoir

	MVR800	MVR1600
Style	Soft Bag	Soft Bag
Volume	800 mL	1,600 mL
Inlet/Outlet Size	1/2" / 3/8"	1/2" / 3/8"
Screen Size	105 µm	105 µm
Biosurface Options	as listed below	

Venous Reservoirs — Uncoated

CFN	Product	Units per Case
MVR800	MVR™ Venous Reservoir Bag 800 mL	6 per case
MVR1600	MVR Venous Reservoir Bag 1,600 mL	6 per case

