

MYOthem XP™ with 4:1 Blood to Crystalloid

Cardioplegia Delivery System

Recommended Flow Rate	0.5 L/min blood flow
Prime Volume	44 mL
Inlet/Outlet Size	¼" / ⅜"
Heat Exchanger Material	Stainless steel
Blood Contact Surface Area	0.06 m ²
Outer Filter Screen Size	150 µm
Pressure Relief Valve Opening Pressure	773 mmHg
Biosurface Options	Trillium™, Cortiva™, or Uncoated



MYOthem XP Blood Cardioplegia Delivery Systems

CFN	Model #	Product	Units per Case
61399404996	XP41	MYOthem XP Cardioplegia Delivery System	4 per case
61399405331	XP41B	MYOthem XP Cardioplegia Delivery System with Bridge Between Blood and Crystalloid Lines	4 per case
61399401032	61399401032	MYOthem XP Delivery Line Single Table Delivery System	4 per case
XP41T	XP41T	MYOthem XP Cardioplegia Delivery System with Trillium Biosurface	4 per case
XP41BT	XP41BT	MYOthem XP Cardioplegia Delivery System with Trillium Biosurface with Bridge Between Blood and Crystalloid Lines	4 per case
CBXP41	CBXP41	MYOthem XP Cardioplegia Delivery System with Cortiva BioActive Surface	4 per case
CBXP41B	CBXP41B	MYOthem XP Cardioplegia Delivery System with Cortiva BioActive Surface with Bridge Between Blood and Crystalloid Lines	4 per case
CB4Z25R	CB4Z25R	MYOthem XP Delivery Line with Cortiva BioActive Surface Single Table Delivery System	4 per case