

BLOOD & INFUSION WARMER

QW618

Dry heating blood and infusion warmer for routine and massive clinical applications:

- Blood transfusion / infusion for the patients during or after operation
- Massive blood transfusion / infusion for emergency care
- Routine blood transfusion / infusion in the clinical settings
- Enteral / parenteral nutrition and rinsing solutions



Extension tube for medical fluids or blood

Grooved cylinder heating exchanger

Display set/heating temperature

Temperature set freely

Microprocessor control

BLOOD & INFUSION WARMER

QW618

FEATURES

Technology of control:

- Microprocessor system integrated with sensors provides high precise control, permanent running self-tests and malfunction alarms

Low cost disposables:

- Grooved dry heating cylinder enables to use standard extension tube or routine IV tube at slow fluid rate

Accurate temperature monitoring:

- Three sets of independent temperature sensor and numerical display for real-time temperature or set temperature in 0.1°C

Multiple protections for high safety:

- High temperature cut-off, low temperature alarm and sensor fault alarm

Unique high temperature alarm testing function:

- Extra safe in use and easy for maintenance

Warming capacity:

- Up to 6000ml/hour (flow rate kvo – 100ml/min)

Good for simultaneous warming multiple infusions

Quick warming up & 24 hours continuous operating

TECHNICAL SPECIFICATIONS

Model:	QW618
Temperature setting:	37°C - 42°C
Power supply:	AC230V±10%/50~60Hz
Power consumption:	Max. 250VA
Type of protection against electric shock:	Class I
Degree of protection against electric shock:	Type B
Degree of protection against ingress of liquids:	IPX4 splash proof
Temperature accuracy:	±0.5°C
Overheat protection:	43°C/45°C/48±3°C
Low temperature alarm:	36°C
Warming up time:	< 45sec (20°C - 36°C)
Operating mode:	Continuous
Dimension (W*D*H):	100×170×210mm
Net Weight:	1.3kg

* Subject to technical modifications without notice



Infusion stand not included in delivery

keewell
Medical Technology

Keewell Medical Technology Co.,Ltd.

Xiananyi Industrial Zone, Pingzhou, Nanhai, Foshan,
Guangdong 528251, P.R.China

T.: +86 757 86260855
marketing@keewell.cn

F.: +86 757 86260100
www.keewell.com.cn

CE 0482

Distributor:



BITMEDICAL

MEDIZINUMFASSEND

www.bitmedical.com